

3050 Spruce Street Saint Louis, Missouri 63103 USA Telephone 800-325-5832 • (314) 771-5765 Fax 800-325-5052 • (314) 286-3828 email: sigma-techserv@sial.com Internet: http://www.sigmaaldrich.com

ProductInformation

PHOSPHODIESTERASE 3':5'-CYCLIC NUCLEOTIDE

ACTIVATOR

from Bovine Brain

Product Number P2277

Cas #: 73298-54-1

Synonyms: calmodulin, calcium dependent regulator protein (CDR), and troponin C - Like peptide.

Product Description

Appearance: White powder. Molecular weight:

- 16,790 by amino acid sequence¹
- 18,700 by sedimentation equilibrium²
- 19,000 by SDS gel (in presence of EGTA)

NOTE: Migration rate in SDS is faster when Ca+2 is present and slower when EGTA removes the Ca+2.)³ Stokes Radius: 20.9 angstroms as calculated from above values²

 $E^{1\%}(278nm) = 1.8$ Solvent:0.1M Imidazole-HCl - 0.001M EDTA, pH 7.0⁴

If calcium is removed, E1% increases to 2.0⁵ Isoelectric point: The literature reports various pl values ranging from 3.9-4.3^{6,7}

Purity: Greater than 98% by SDS PAGE using brilliant blue stain. Supplied as an essentially salt-free powder. Activity:>40,000 units per mg protein (Lowry) Structure: Calmodulin has a highly conserved structure and is ubiquitous to eukaryotes.^{8,9,10} A Posttranslational modification of most species is a trimethylated lysine at the 115 position,^{10,11} and acetylation of the terminal alanine residue.¹¹

Calmodulin contains many acidic amino acids. The protein lacks cysteine, hydroxyproline and tryptophan. The abundance of acidic amino acid carboxylate groups allows reversibly binding to calcium 2+, while the absence of cysteine and hydroxyproline allows for a very flexible tertiary structure for interaction with various calmodulin regulated proteins.

Calmodulin also has a high ratio of phenylalanine (8 residues) to tyrosine (2 residues) and has a distinctive ultraviolet spectrum with five peaks at 253 nm, 259 nm, 365 nm, 269 nm and 277 nm, and a shoulder at 282 nm.¹¹ Phosphorylation of Calmodulin in vivo has been reported when cells are stimulated with insulin,^{12,13,14} and in vitro by various protein kinases.^{15,16,17,18,19,20}

Xray crystallography of calmodulin, in the presence of calcium, indicates that it exists as a long dumbbell shaped molecule 65 anstroms long. Each globular end contains two calcium binding domains. These domains are common among many calcium binding proteins, and are described as helix-loop-helix "EF-Hand" regions. The calcium binding regions are connected by an extended 40 angstrom, 28 amino acid alpha helical region.^{21,22,23} Upon binding to calcium, calmodulin undergoes a conformational change in which the hydrophobic regions are believed to be involved with enzyme binding.^{24,25}

The four calcium binding sites of calmodulin are designated as I, II, III and IV, starting from the site closest to the amino terminus. The order of calcium to calmodulin is believed to occur in the order of III, IV, I and II. Sites III and IV have an affinity for calcium 10-20 times higher than sites I and II.

Composition: Calmodulin P2277 is supplied as a saltfree lyophilized powder. The calcium content of a representative lot was found to be less than 0.05%

Applications:

Calmodulin has been found to activate several enzymes.

Cyclic Nucleotide Phosphodiesterase^{29, 30, 31} Nitric Oxide Synthetase^{32, 33} ATPase^{34, 35, 36, 37} Phosphorylase Kinase^{38, 39, 40} Myosin Light-Chain Kinase⁴¹ Protein Phosphatase (Calcineurin)^{42, 43} Adenylate Cyclase⁴⁴ NADPH Oxidase⁴⁵ CaM Kinases I, II and III^{46, 47, 48, 49}

Calmodulin is believed to be involved in intracellular Ca 2+ homeostasis, cell proliferation, smooth muscle contraction, microtubular function, exocytotic secretion of cellular products, and cell motility. Several review articles on calmodulin have been published.^{50,51}

Preparation Instructions

Sigma's method of preparation is based on a modification of the procedure of Gopalakrishna.⁵² Solutions at 1 mg/mL or higher in water or saline (no phosphate or other calcium binding buffers)at refrigerator temperature are stable for 1-2 days or frozen at 1 mg/mL or higher for approximately one week.

Storage/Stability

Store at less than 0°C. Stored properly, as supplied in the powder form, calmodulin should be stable for a minimum of two to three years.

Product Profile

One unit will stimulate 0.016 activated unit of phosphodiesterase, 3':5' - cyclic nucleotide, P 0520, in a 3 ml reaction volume at pH 7.5 and 30°C, to 50% of the maximum activity of the enzyme when saturated with activator, in the presence of 0.01 mM Ca⁺⁺.

References

- Watterson, D. M., *J. Biol. Chem.*, vol. 255, 962-975 (1980).
- 2. Ann. Rev. Biochem., vol 49, 493 (1980).
- 3. Klee, C.B. et al., PNAS (USA), vol 76, 6270 (1979).
- 4. Watterson, D.M. et al., *J. Biol. Chem.*, vol. 251, 4501-4513 (1976).
- 5. Klee, C.B., *Biochem.*, vol. 16,1017 (1977).
- Lin, Y.M. et al., *J. Biol. Chem.*, vol. 249, 4943 (1974).
- Crouch, T. H. et al., *Biochem.*, vol 19, 3692-3698 (1980).
- Klee, C. and Vanaman, T., *Adv. Protein Chem.*, 35, 213-321 (1982).
- 9. Means, A. et al., *Physiol. Reviews*, 62, 1-39 (1982)
- 10. Wylie, D. and Vanaman, T., in *Calmodulin*, (Cohen, P. and Klee, C. eds) Elsevier, New York, pp. 1-15
- 11. Cheung, W.Y., *Science*, 207, January 4, 1980, p. 19.
- 12. Sacks, D. et al., Biochem. J., 286, 211-216 (1992)
- 13. Fukami, Y., et al., *Proc. Natl. Acad. Sci.* USA, 83 4190-4193 (1986).
- 14. Colca, J. et al., *J. Biol. Chem.*, 262, 11399-11402 (1987).
- 15. Graves, C. et al., *J. Biol. Chem.*, 261, 10429-10438 (1986).
- 16. Sacks, D., and McDonald, J., *J. Biol. Chem.*, 263, 2377-2383 (1988).
- 17. Sacks, D. et al., Biochem. J., 262, 803-812 (1989).
- 18. Meggio, F. et al., FEBS Lett. 215, 241-246 (1987).
- 19. Sacks, D. et al., Biochem. J., 283-21-24 (1992).
- Sacks, D. et al., *Biochem. Biophys.* Res. Comm., 188, 754-759 (1992).
- 21. Babu, S. et al., Nature, 315, 37-40 (1985).

- 22. Babu, S. et al., J. Mol. Biol., 204, 191-204 (1988).
- Krestinger, R. in Calcium Transport in Contraction and Secretion (Carafoli, E., et al ed.), pp. 46-478 (1975).
- 24. Newton, D. et al., *J. Biol. Chem.*, 259, 4419-4426 (1984).
- 25. Newton, D. et al., *Biochim. Biophys. Acta.*, 845, 533-539 (1985)
- 26. Olwin, B. et al., *J. Biol. Chem.*, 259, 10949-10955 (1984).
- Haiech, J. et al., *J. Biol. Chem.*, 266, 3427-3431 (1991).
- 28. Starovasnik, M. et al., *Protein Sience*, 1, 245-253 (1992).
- Beavo, J. et al., *Mol. Cell. Endocrinol.*, 28, 387-410 (1982).
- 30. Dedman, J. et al., *J. Biol. Chem.*, 252, 8415-8422 (1977).
- 31. Cheung, W., *Biochem. Biophys. Res. Commun.*, 38, 533-538 (1970).
- 32. Bredt, D. and Snyder, S., *Proc. Nat, Acad. Sci.* USA, 87, 682-685 (1990).
- Busse, R. and Mulsch, A. *Febs Lett.*, 265, 133-136 (1990).
- 34. Watterson, D. et al., *J. Biol. Chem.*, 251, 4501-4513 (1976).
- 35. Blum, J. et al., J. Cell Biol., 87, 386-397 (1980).
- Carafoli, E. et al., Ann. N.Y. Acad. Sci., 402, 304-328 (1982).
- Caroni, P. and Carafoli, E., J. Biol. Chem., 256, 3263-3270(1981).
- Newsholme, P. et al., J. Biol. Chem., 267, 810-818 (1992).
- Chan, K and Graves, D., J. Biol. Chem., 257 5948-5955 (1982).
- Cohen, P., *Mol. Aspects Cell. Regul.*, 5, 123-144 (1988).
- 41. Stull, J., *Mol. Aspects Cell. Regul.*, 5, 91-122 (1988).
- 42. Klee, C., Methods Enzymol., 102, 227-244 (1983).
- 43. Stewart, A. et al. *Eur. J. Biochem.*, 132, 289-295 (1983).
- 44. Bronstrom, C. et al., *Proc. Nat. Acad. Sci.* USA, 72, 64-68 (1976).
- 45. Jones, H. et al., *Biochim. Biophys. Acta.*, 714, 152-156 (1982).
- 46. Nairn, A. and Greengard, P., *J. Biol. Chem.*, 262, 7273-7281 (1987).
- 47. Colbran, R and Soderling, T., Curr. *Top. Cell.* Regul., 31, 181-221 (1990).

- 48. Nairn, A. and Palfrey, H., *J. Biol. Chem.*, 262, 17299-17303 (1987)
- 49. Nairn, A. et al., *Proc. Natl. Acad. Sci.* USA, 82, 7939-7943 (1985).
- 50. Cheung, W., Science, 207, 19-27 (1980).
- 51. Hinrichsen, R., *Biochim. Biophys. Acta.*, 1155, 277-293 (1993).
- 52. Gopalakrishna, R. and Anderson, W., *Biochem. Biophys. Res. Commun.*, 104, 830-836 (1982).

rbg 10/26/98

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.