

Product Information

S-Nitroso-N-acetyl-DL-penicillamine

Catalog Number **N3398**Storage Temperature -20°C


CAS RN: 67776-06-1

Synonym: SNAP

Product Description

Molecular Formula: $\text{C}_7\text{H}_{12}\text{N}_2\text{O}_4\text{S}$

Molecular Weight: 220.25

Purity: $\geq 97\%$ Molar extinction coefficients:¹(in 0.1 M citrate/HCl Buffer, pH 2.0, or 0.5 M HCl) $E_{335} = 519 \text{ M}^{-1} \text{ cm}^{-1}$;
 $E_{591} = 7.0 \text{ M}^{-1} \text{ cm}^{-1}$

S-Nitroso-N-acetyl-DL-penicillamine is synthetically prepared and is a racemic mixture of the D and L isomers. Methods of preparation¹⁻³ and quantitative methods of analysis of SNAP¹, including biological samples⁴, have been reported.

SNAP, a nitrosothiol derivative, releases nitric oxide (NO) under physiological conditions, thus making it a useful tool for studying pharmacological and physiological actions of NO. Because of its NO releasing properties, SNAP is a potent vasodilator *in vitro*^{5,6} and *in vivo*^{7,8} and is less prone to produce pharmacological tolerance.⁵⁻⁸ SNAP inhibited cultured rat vascular smooth muscle cell mitogenesis and proliferation at about 0.1 mM.⁹ It inhibited the dioxygenase activity of lipoxygenase (lipoxidase) in human platelets and in human CHP100 neuroblastoma cells. The effect of the inhibition on the dioxygenase activity was also studied using the model soybean lipoxygenase Type II (LOX-2). SNAP is a competitive inhibitor of the enzyme ($K_i = 710 \mu\text{M}$) and inactivated the enzyme by the reduction of the catalytic iron to the inactive Fe(II) form and inhibited the hydrogen peroxide mediated activation of the lipoxygenase catalyzed dioxygenase reaction.¹⁰

SNAP use in tissue samples resulted in nitric oxide mediated inhibition of lipoxygenase, and inhibition of superoxide anion- and peroxynitrite-dependent liposomal, lipoprotein, and fatty acid oxidation.¹¹ It was shown to have a reproducible vasorelaxant effect on isolated canine coronary arteries ($\text{IC}_{50} = 113 \text{ nM}$).¹² Denitrosation of SNAP and the subsequent release of NO may not be spontaneous but may be due to catalysis by specific enzymes located at external vascular cell membranes.⁶

SNAP (10^{-5} M) caused a decrease in the calcium ion concentration in rat vascular smooth muscle cells.¹³ It inhibited *in vitro* clot formation and formation of γ dimers of fibrin at 10 mM, both of which are catalyzed by Factor 13 (F13), a transglutaminase. SNAP inhibited F13 activity ($\text{IC}_{50} = 230 \mu\text{M}$) in the presence of dithiothreitol. These effects are presumably due to the S-nitrosylation of the highly reactive cysteine residue of the F13 enzyme.¹⁴ It decreased α -Tumor Necrosis factor induced apoptosis (which is partly mediated through the cGMP pathway) maximally (60%) at 250 μM in endothelial cells. In contrast, high concentrations (1 mM) induced endothelial apoptosis via cGMP-independent pathways.¹⁵ SNAP (10^{-4} M) stimulated human colonic ion transport *in vitro*¹⁶; inhibited (5-500 μM) aggregation of thrombin induced human platelets¹⁷ and increased ($\text{EC}_{50} \sim 50 \mu\text{M}$) cGMP levels in rat cerebellar slices up to 300 fold.¹⁸

In resting human peripheral blood mononuclear cells, SNAP enhanced the rate of glucose transport, induced TNF- α secretion and NF- κB binding activity and enhanced activity (10 μM) of membrane-associated protein tyrosine phosphatase.¹⁹ It induced a loss of about 90% cell viability of cultivated endothelial cells under hypoxic conditions compared to 45% under normoxic conditions at 5 mM.²⁰ It induced manganese superoxide dismutase (MnSOD) mRNA in rat vascular smooth muscle cells²¹ and MnSOD and copper zinc SOD in a murine macrophage cell line at 100 μM .²² Through an S-nitrosylation mediated inhibition of phosphodiesterase, SNAP increased basal lipolysis in white adipose tissue.²³

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

Storage/Stability

Stable for at least one year if stored desiccated at -20°C and protected from exposure to light.²⁴ Soluble in DMSO (57.5 mg/ml) and in water (2.1 mg/ml). Solutions of SNAP (10 mg/ml) have been prepared in methanol. The half life of SNAP in aqueous media is about 5 hours.²⁵ Stock solutions should be made fresh using deoxygenated citrate/HCl buffer, pH 2.0 or 0.5-1 M HCl. Solutions can be kept on ice protected from light for several hours.^{6,24} Stock and working solutions were freshly prepared before use.^{6,26}

Controlled studies of solutions in oxygenated Krebs buffer showed degradation of SNAP at a concentration of 4.4 μM to about 3 μM on incubation at 37°C for 10 min.⁶ For application of SNAP to thymic cells, a 25 mM stock solution (1.10 mg/194 μL of PBS + 6 μL of 1 M sodium hydroxide) on ice was prepared and diluted right before use.² The decomposition of 1 mM SNAP (absorbance at 340 nm) in PBS, pH 7.2, at 37°C in the absence and presence of 0.1-1 mM cysteine (which enhances decomposition) and the effect of cysteamine on the rate of NO-release from SNAP have been reported.^{2,27} NO formation from SNAP is high (100 μM gives about 1.4 μM NO/minute at 37°C) and is linear over a wide concentration range²⁸. Metal chelators, such as EDTA, may stabilize SNAP solutions.^{1,24} The effects of trace metals and EDTA on the decomposition of SNAP was reported.²⁸ It is an effector releaser of NO in the pH range studies (pH 5-8).²⁶

References

1. Stamler, J.S. and Feelisch, M. "Preparation and Detection of S-Nitrosothiols" in *Methods in Nitric Oxide Research*, Chapter, 36, 521, Feelisch, M. and Stamler, J.S. Eds., John Wiley & Sons, Ltd, New York, 1996.
2. Fehsel, K. et al. *Methods in Enzymology*, **269**, 426, (1996).
3. Field, L. et al., *J. Chem. Soc. Perkin Trans. I*, 249, (1978).
4. Samouilov, A. and Zweier, J.L., *Anal. Biochem.* **258**, 322, (1998).
5. Kowaluk, E.A. et al., *Eur. J. Pharmacol.* **144**, 379, (1987).
6. Kowaluk, E.A. and Fung, H.-L., *J. Pharmacol. Exp. Ther.* **255**, 1256, (1990).
7. Shaffer, J.E. et al., *J. Pharmacol. Exp. Ther.* **260**, 286, (1992).
8. Bauer, J.A. and Fung, H.-L., *J. Pharmacol. Exper. Ther.* **256**, 249, (1991).
9. Garg, U.C. and Hassid, A., *J. Clin. Invest.* **83**, 1774, (1989).
10. Maccarrone, M. et al., *Biochem. Biophys. Res. Commun.* **219**, 128, (1996).
11. Rubbo, H. and Freeman, B.A., *Methods in Enzymology*, **269**, 385, (1996).
12. Ku, D.D., *Methods in Enzymology*, **269**, 107, (1996).
13. Okishio, M. et al., *Biochem. Biophys. Res. Commun.* **183**, 849, (1992).
14. Catani, M.V. et al., *Biochem. Biophys. Res. Commun.* **249**, 275, (1998).
15. Shen, Y. H. et al., *FEBS Letters*. **433**, 125, (1998).
16. Stack, W.A. et al., *Gut*, 39, 93, (1996).
17. Azula, F.J. et al., *Mol. Pharmacology* **50**, 367, (1996).
18. Southam, E. and Garthwaite, J., *Neuroscience Letters*, **130**, 107, (1991).
19. Lander, H.M. et al., *J. Immunology*, **150**, 1509, (1993).
20. Ioannidis, I. et al., *Biochem. J.* **318**, 789, (1996).
21. Sano, H. et al., *J. Cellular Biochem.* **62**, 50, (1996).
22. Sano, H. et al., *Immunology*, **92**, 118, (1997).
23. Gaudiot, N. et al., *J. Biol. Chem.* **273**, 13475, (1998).
24. Feelisch, M. and Stamler, J.S. "Donors of Nitrogen Oxides" *ibid*, Chapter 7, 71.
25. Ignarro, L.J. et al., *J. Pharmacol. Exp. Ther.* **218**, 739, (1981).
26. Gryglewski, R.J. et al., *Int. J. Tiss. Reac.* **11**, **269**, (1989).
27. Noack, E. and Feelisch, M. *Basic Research in Cardiology*, **86**, 37, (1991).
28. Feelisch, M., *J. Cardiovascular Pharmacology*, **17**, S25, (1991).

ES-S,PHC 01/07-1

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.