cell signaling

The Effective Use of Protein Kinase Inhibitors

By Philip Cohen, Ph.D.

Royal Society Research Professor, University of Dundee, Scotland

Introduction

The advent of relatively specific cell-permeable inhibitors of protein kinases in the mid-1990s has had a major impact on the study of signal transduction. The ability to rapidly suppress the cellular activity of a particular protein kinase has proved to be a powerful method for identifying the physiological substrates of these enzymes and the roles of the signaling pathways in which they participate. These compounds enter cells within minutes, so that indirect effects caused, for example, by changes in gene expression (a potential hazard when using cells deficient in a particular protein kinase), are excluded. Moreover, the use of protein kinase inhibitors avoids the need for transfectionbased approaches that have the potential to give misleading results, since the fidelity of signaling can break down when components are overexpressed. Nevertheless, in order to use protein kinase inhibitors effectively it is important to realize their limitations, as well as their strengths.

An appreciation of the degree of specificity of any particular inhibitor is clearly a critical issue. There are just over 500 protein kinases encoded by the human genome, most of which belong to the same superfamily. It is therefore a challenging and difficult task to develop compounds that inhibit one particular protein kinase, without inhibiting several related enzymes. Table 2 provides current information about the specificities of 41 protein kinase inhibitors, many of which are available from Sigma-RBI, on a wide range of protein kinases.

Protein kinases in cell-based assays

Most inhibitors of protein kinases target more than one enzyme. There is, therefore, a danger that, in cell-based assays, the observed effects do not result from inhibition of the kinase of interest, but rather from inhibition of another protein kinase. In order to exclude this possibility, it is necessary to show that the effects of an inhibitor disappear in cells that express an inhibitor-resistant mutant of the kinase of interest. However, at present, the availability of such cells is very limited.^{1,2} In order to reduce this risk, it is important to examine, wherever possible, the effects of at least two structurally unrelated inhibitors of the same protein kinase. For example, kenpaullone (Product Code K 3888) and roscovitine (Product Code R 7772), which are relatively specific inhibitors of cyclin-dependent protein kinases (CDKs), also inhibit a few other protein kinases. However, the other enzymes inhibited by roscovitine are not the same as those inhibited by kenpaullone (Table 2). Thus, if identical effects are observed with roscovitine or kenpaullone, one can have greater confidence that the effects are mediated by a CDK.

For similar reasons, it is advisable to use both LiCl (Product Code <u>L 0505</u>) and kenpaullone to study glycogen synthase kinase-3 (GSK-3). Wortmannin (Product Code <u>W 1628</u>) and LY 294002 (Product Code <u>L 9908</u>) are used to identify potential roles of phosphoinositide 3-kinase (PI3K), PP1 or PP2 and SU6656 (Product Code <u>S 9692</u>) for Src family kinases. Y 27632 and HA1077 (Product Code <u>H-139</u>) are used for Rho kinase (ROCK) and protein kinase C-related kinase 2 (PRK2) (Tables 1 and 2).^{3,4}

Even compounds that inhibit a number of protein kinases can sometimes be useful in excluding the involvement of one or more protein kinases in the control of a particular process. For example, H89 (Product Code B 1427), which inhibits isoforms of mitogen- and stress-activated kinase (MSK), but not the structurally related isoforms of ribosomal S6 kinase (RSK), has been used to provide evidence that RSKs do not mediate the growth factor-induced phosphorylation of the transcription factor cAMP-responseelement-binding protein (CREB).⁵ MSK isoforms were later shown to be the physiologically relevant protein kinases using cells deficient in these kinases.⁶ Similarly, UCN01, which inhibits checkpoint kinase 1 (CHK1), but not CHK2, can be used to exclude the involvement of CHK2 in the control of responses to DNA damage or cell cycle checkpoints (Tables 1 and 2).³

It is also possible to vary the concentrations of inhibitors in the culture medium to differentially inhibit particular protein kinases. For example, at low concentrations PD184352 inhibits the classical mitogen-activated protein kinase (MAPK) cascade specifically, but at higher concentrations it also blocks the mitogen-activated protein kinase 5 (MKK5/ERK5) pathway.⁷ However, the precise concentrations needed may vary from cell to cell. For this reason, it is essential to define the minimum concentration of an inhibitor required to suppress activity by 80-90% by examining the phosphorylation of a validated substrate of the protein kinase that is under investigation.

Application Notes

Ε

ma-aldrich.co

s i g

The vast majority of protein kinase inhibitors target the adenosine 5'-triphosphate (ATP)-binding site of a protein kinase. For this reason, much higher concentrations of inhibitors are generally needed to suppress the activity of a protein kinase in cells (where the ATP concentration is in the millimolar range), compared to the amounts required for inhibition in vitro (where assays are performed at much lower ATP concentrations, typically 0.01-0.1 mM). There are, however, a few inhibitors that are actually more potent in cell-based assays than they are in vitro. For example, PD 98059 (Product Code P-215) and U0126 (Product Code U-120), which are non-competitive inhibitors of mitogen-activated protein kinase kinase 1 (MKK1), bind much more strongly to the dephosphorylated, inactive form of this protein kinase than the phosphorylated, active enzyme.

These compounds prevent the conformational change required for activation of MKK1 and therefore suppress the classical MAPK cascade at much lower concentrations than those needed to inhibit activated MKK1 *in vitro*.^{3,8} Similarly, lithium ions, a relatively specific inhibitor of GSK-3, compete for binding with magnesium ions. The free concentration of magnesium ions in cells is less than 0.5 mM, much lower than the concentration used to assay

GSK-3 routinely (10 mM, see Table 1). For this reason, lithium ions inhibit GSK-3 more potently in cells than *in vitro*.³

In recent years, potent and highly specific inhibitors of a variety of protein kinases have been developed by several pharmaceutical companies. Many have entered human clinical trials and, in two cases (Glivec and Iressa), have been approved for the treatment of different types of cancer.⁹ Over the next ten years one can therefore expect that many more protein kinase inhibitors will become available to the scientific community, which should advance at an even faster pace our understanding of the function of these enzymes.

References

- 1. Eyers, P. et al. FEBS Lett., 451, 191-196 (1999).
- 2. Brown, E.J. et al. Nature, **377**, 441-446 (1995).
- 3. Davies, S.P. et al. Biochem. J., **351**, 95-105 (2000). 4. Bain, J. et al. Biochem. J., **371**, 199-204 (2003).
- 5. Thomson, S. et al. EMBO J., **18**, 4,779-4,793 (1999).
- Wiggin, G.R. et al. Mol. Cell Biol., 22, 2,871-2,881 (2002).
- 7. Mody, N. et al. FEBS Lett., **502**, 21-24 (2001).
- 8. Alessi, D.R. et al. J. Biol. Chem., 270, 27,489-27,494 (1995).
- 9. Cohen, P. Nature Rev. Drug Discovery, 1, 309-316 (2002).

| Table 1. How to use the more specific inhibitors of | protein kinases in cell-based assays |
|---|--------------------------------------|
| | |

| Inhibitor | Specificity | Target Kinase(s)*** | Concentration to use in Culture Medium (µM)* |
|-----------------------------|-------------|-----------------------------|---|
| Rapamycin (R 0395) | Very high | mTOR | 0.1 |
| PD 98059 (P-215) | High | MKK1 | 50 |
| PD 184352** | High | MKK1 | 1-2 |
| PD 184352** | High | MKK1, MKK5 | 10-20 |
| U0126 (U-120) | High | MKK1, MKK5 | 5-10 |
| SB-203580 (5 8307) | High | SAPK2a/p38α, SAPK2b/p38β2 | 1-10 |
| SB-202190 (S 7067) | High | SAPK2a/p38α, SAPK2b/p38β2 | 1-5 |
| KN62 (I 2142) | High | CaM-KII, other CaM-Ks | 10 |
| Wortmannin (W 1628) | High | PI3K | 0.1 |
| LY 294002 (L 9908) | Quite high | PI3K | 50-100 |
| Y27632 | Quite high | ROCK, PRK2 | 10-20 |
| HA1077 (H-139) | Medium | ROCK, PRK2 | 10-100 |
| LiCI (L 0505) | Quite high | GSK-3 | 10 |
| Kenpaullone (K 3888) | Quite high | GSK-3, CDKs | 10 |
| Roscovitine (R 7772) | High | CDKs | 10-100 |
| PP1 | Quite high | Src, Fyn, Lck | 0.1-1.0 |
| PP2 | Quite high | Src, Fyn, Lck | 0.1-1.0 |
| SU6656 (S 9692) | Medium | Src, Fyn, Lck | 10-50 |
| ML7 (I 2764) | Quite high | Sm-MLCK | 50-100 |
| H89** (B 1427) | Medium | РКА | 5-10 |
| H89** (B 1427) | Medium | PKA, MSKs, but not RSKs | 10-25 |
| Ro 31-8220** (R-136) | Medium | Conventional PKCs | 1 |
| Ro 31-8220** (R-136) | Medium | PKCs, MSKs, RSKs, etc | 5 |
| UCN01 | Quite low | PDK1 and CHK1, but not CHK2 | 0.3-1 |

* The suggested concentrations are only guidelines. The optimal concentrations can vary and need to be defined for each cell used, as discussed in the text. Sigma-RBI product numbers are shown in red.

Depending on the concentration range at which they are used, these kinase inhibitors can be used to target different groups of protein kinases. *Kinase names are provided in the Table 1 legend. m ю . ص \sim ω 0 0 ø Ø υ r v i Ð ŝ i c a ⊆ ч о Ч е 0 0 m 25. m rder: 1.800.

 \sim

Table 2. Inhibition of protein kinases by various inhibitors. Results indicate percent activity observed in the presence of inhibitor expressed as a percentage of control incubations in the absence of inhibitor. Data are the means of duplicate determinations. Data highlighted in boxes indicate instances when a given inhibitor reduced kinase activity to \leq 25% of control values. Assays were carried out at a magnesium ion concentration of 10 mM and an ATP concentration of 0.1 mM. Column headers indicate the protein kinase inhibitors tested, together with the concentrations at which they were used and their Sigma-RBI product numbers in red.

| | | (U17) | - | (NI) | OLINA | 1. m. | (11) | 2 | , | 0 0 LIMAI | 0 m | nin Time | (INI) CINN) | 20 LMAI | in 1 | 10VILIO | in low | 20 1/1 | (hvi) | (Mr) |
|--|--|--|--|----------------------------------|----------------|-----------------------------|--------------|-----------|---|--|--|--------------------------------|---|------------|-------------|-------------|----------------------------|-------------|-------------------------------------|------------|
| Protein kinase | 489 1489 | 1 2 2 0 LIN | (111) HA 10>> | Rotter, 139 | 17 56481 COUNT | (1) 25 (10) (10) 25 (10) | (1) 180, LIA | PO 000 | (P. 2, 3059 58, 50 58, 2035, 111) | (S & 20) (10 (S & 20) (10 (S & 20) (10 | 1000, 100 1000, 100 1000, 100 1000, 100 1000, 100 | (W (man) 1628) 11 290 LW | (1) 3000 (1) (1) 3000 (50 (1) 50 (50) (1) 50) (1) 50 (50) (1) 50 (50) (1) 50) (1) 50 (50) (1) 50 (50) (1) 50) (1) 50) | Rapamin MA | Lice (1 LIL | KCI (10 mm) | Ro 37.00 mm | Bis. 1 LINA | (8 5.9) Bis 3 (1 11) (8 3 78) | (evr) (Co. |
| (A) Core panel | · - | | | • | * | • | ÷ | 2 | | | | ÷ | 2 | | | | | | ÷ | |
| VIKK1 | 90 | 103 | 89 | 106 | 101 | 56 | 5 | 89 | 99 | 93 | 96 | 101 | 94 | 99 | 116 | 109 | 94 | 94 | 90 | |
| MAPK2/ERK2 | 87 | 94 | 94 | 139 | 92 | 92 | 107 | 85 | 85 | 89 | 90 | 114 | 113 | 90 | 107 | 102 | 97 | 98 | 107 | |
| JNK1α1/SAPK1c | . 97 | 98 | 96 | 49 | 104 | 96 | 102 | 111 | 101 | 93 | 97 | 108 | 101 | 98 | 92 | 89 | 95 | 91 | 100 | |
| JNK/SAPK1c | | | | | | | | | | | | | | | | | | | | |
| SAPK2a/p38 | 99 | 94 | 93 | 111 | 95 | 75 | 100 | 85 | 2 | 0 | 86 | 98 | 138 | 93 | 108 | 102 | 84 | 88 | 104 | |
| | 97 | 107 | 97 | 127 | 98 | 90 | 119 | 95 | 10 | 3 | 74 | 98 | 150 | 88 | 96 | 104 | 97 | 92 | 115 | |
| | 106 | 100 | 87 | 146 | 95 | 100 | 100 | 96 | 96 | 80 | 75 | 97 | 132 | 100 | 99 | 108 | 92 | 89 | 116 | |
| • | 105 | 95 | 103 | 130 | 110 | 111 | 98 | 94 | 93 | 87 | 79 | 94 | 103 | 82 | 84 | 99 | 104 | 113 | 133 | |
| МАРКАР-К1а | | | | | | | | | | | | | | | | | _ | _ | _ | |
| | 16 | 72 | 37 | 79 | 89 | 88 | 86 | 93 | 83 | 95 | 92 | 70 | 20 | 95 | 95 | 78 | 2 | 9 | 2 | |
| МАРКАР-К2 | 99 | 99 | 90 | 5 | 59 | 102 | 98 | 95 | 93 | 97 | 102 | 74 | 90 | 125 | 72 | 98 | 103 | 90 | 97 | |
| | 3 | 57 | 19 | 38 | 81 | 104 | 118 | 86 | 86 | 88 | 99 95 | 83 | 37 | 104 | 104 | 105 | 2 | 21 | 9 | |
| | 81 | 104 | 91 25 | 6 | 36 | 93 05 | 71 | 108 | 112 | 88 | 85 | 68 | 51 | 74 | 76 | 104 | 96 | 98 | 89 | |
| | 2 | 91 98 | 35 | 17 | 94 95 | 95 | 105 | 106 | 96 89 | 66 | 97 100 | 91 01 | 104 | 104 00 | 96 98 | 96 97 | 70 | 99 | 87 | |
| | 79 104 | 98 115 | 86 92 | 95 36 | 95 70 | 92 99 | 99 85 | 93 86 | 89 89 | 92 87 | 100 88 | 91 76 | 70 81 | 99 110 | 98 105 | 97 98 | 84 | 4 | 7 | |
| | 104 | 115 90 | 92 88 | 36 27 | 70 67 | 99 79 | 85 89 | 86 82 | 89 62 | 87 53 | 88 96 | 76 60 | 81 99 | 110 91 | 105 95 | 98 96 | 84 73 | 85 77 | 69 99 | |
| | 17 25 | 90 109 | 88 92 | 27 81 | 67 78 | 79 91 | 89 111 | 82 90 | 62 83 | 53 98 | 96 101 | 60 72 | 99 35 | 91 108 | 95 99 | 96 100 | 73 21 | 77 63 | 99 26 | |
| | 25 | 109 94 | 92 32 | 81 98 | 78 93 | 91 92 | 111 86 | 90 100 | 83 87 | 98 75 | 101 | 72 81 | 35 25 | 108 | 99 95 | 100 | 6 | 63 32 | 26 18 | |
| | 107 | 94 92 | 32 90 | 98 | 93 38 | 92 105 | 86 | 100 | 87 66 | 75 61 | 85 | 53 | 30 | 109 89 | 95 58 | 99 | 5 | 32 50 | 18 46 | |
| GSK-3p ROCK-II | 107 | 92 | 90 | 13 88 | 38 88 | 105 94 | 83 107 | 80 | 66 77 | 61 | 85 91 | 53 104 | 30 55 | 89 92 | 58 101 | 99 102 | 92 | 50 90 | 46 89 | |
| | 19 | 95 | 77 | 88 98 | 88 97 | 94 85 | 89 | 80 97 | 96 | 94 | 106 | 104 | 16 | 92 106 | 101 | 102 | 92 42 | 90 54 | 23 | |
| CK1 | | 55 | 11 | 50 | 51 | 0.5 | 05 | 5. | 50 | 51 | 100 | 102 | <u> </u> | 100 | 100 | 105 | "TL | 5.1 | | |
| | 104 | 98 | 102 | 103 | 103 | 107 | 96 | 87 | 97 | 93 | 98 | 18 | 19 | 104 | 73 | 112 | 104 | 101 | 106 | |
| РНК | 51 | 81 | 58 | 63 | 105 | 107 | 117 | 87 | 104 | 91 | 100 | 44 | 32 | 104 | 96 | 93 | 57 | 54 | 20 | |
| LCK | 76 | 109 | 94 | 70 | 92 | 87 | 99 | 85 | 32 | 37 | 95 | 85 | 83 | 102 | 99 | 105 | 79 | 86 | 86 | |
| | 21 | 99 | 82 | 107 | 104 | 95 | 104 | 99 | 95 | 95 | 99 | 90 | 56 | 102 | 96 | 97 | 42 | 60 | 40 | |
| снк2 | | | Ż | · | | | Ż | | Ż | ġ | | | | | ġ, | ġ | | | | |
| CSK | | | | | | | | | | | | | | | | | | | | |
| CDK2/Cyclin A | | | | | | | | | | | | | | | | | | | | |
| DYRK1A | | | | | | | | | | | | | | | | | | | | |
| PKG | | | | | | | | | | | | | | | | | | | | |
| (B) Other kinases | 35 | | | | | | | | | | | | | | | | | | | |
| CAM-KII | | | | | 0 | | | | | | | | | | | | | | | |
| SkMLCK | | | 93 | | | | | | | | 104 | | | | | | | | | |
| SmMLCK | | | 93 | | 96 | | | | | | 4 | | | | | | | | | |
| ΡΚϹδ | | | | 101 | | | | | | | | | | | | | | | | |
| | 95 | | | | | 114 | | | | | | | | | 94 | 109 | | | | |
| МКК4 | 80 | | | | | 81 | | | | | | | | | 87 | 94 | | | | |
| МКК6 | 86 | | | | | 79 | | | | | | | | | 108 | 113 | | | | |
| MKK7 | 91 | | | | | 89 | | | | | | | | | 100 | 102 | | | | |
| PI3K PRK2 | | 6 | 15 | | | | | | | | 0 | 13 | 18 | | | | | | | |
| Abbreviations AMPK: CaM-KII: CDK2/Cyclin A: CHK1: CHK2: CK1: CK2: CSK: | AMP-aci Calcium Cyclin-d Checkpo Checkpo Casein k Casein k COOH-t | citvated pro m/Calmodu dependent point kinase kinase 1 kinase 2 -terminal Sr pecificity typ | protein kina lulin protei 1t kinase 2, se 1 se 2 Src kinase | ein kinase I 2/Cyclin A se | A complex | | inase 1A | _ | _ | JNK/SA JNK1α LCK: MAPKA MAPKA MAPKA MAPKA MAPKA | GSK-3β: Glycogen synthase kinase-3β JNK/SAPK1c:: c-jun N-terminal kinase JNK1/x1/SAPK1c:: c-jun N-terminal kinase LCK: T-cell specific kinase; lymphocyte-specific kinase MAPK2/ERK2: Mitogen-activated protein kinase-activated protein kinase kinase 1 MKK3: Mitogen-activated protein kinase kinase 3 | | | | | | tein kinase tein kinase | se-1b | | |

22

6

sigma-aldrich.com

Biochem. J., **371**, 199-204. © The Biochemical Society.

For more information on Sigma-RBI's extensive range of protein kinase inhibitors, natius, maine in a ma visit our website at sigma-aldrich.com/cellsignaling.

| visit our website at sigma-aldrich.com/cellsignaling. | | | | | | | | | | | | | | | | | | | | | | |
|---|------------|------------|------------|--------------|---------------|-----------------|------------|-----------|--------------|----------|--------------|------------------------------------|--------------------|-----------------|-----------|----------------|-------------|---|------------|-----------------|-----------------|-------------------|
| | | | | | | | - | | | | oth | | | | | | | | | | | |
| Bisca | 5 (1 MJ) | Bis. 8 LIN | Bis 70 LIN | 60 69 (1 km) | 125 10 LIN | 47 4 3,52 | UCULO LIVI | Roscow. | Purvelon und | Indirubi | Kenoguisting | 415terno 115terno | PP7 (10 MG) | £ •. | 2 . | SU SU Se | 50 6001 (10 | 1000 100 100 100 100 100 100 100 100 10 | EGC 10 110 | (111 01) (-111) | 100 (C, 30 LIN) | (MN) 001/12/10 |
| | | | _ | | | | | | | | | | | | | | | | | | | 2 |
| 84 | 99 | 99 | 9 | 64 | 59 | 20 | 90 | 93 | 80 | 94 | 103 | 96 | 52 | 55 | 89 | 82 | 89 | 99 | 94 | 102 | 95 | |
| 86 | 101 | 105 | 106 | 34 | 37 | 57 | 87 | 81 | 26 | 104 | 70 | 64 | 75 | 61 | 92 | 102 | 55 | 106 | 53 | 85 | 101 | |
| 95 | 100 | 95 | 103 | 88 | 82 | 81 | 98 | 0.1 | 0.1 | | 120 | 00 | 00 | 00 | 00 | 100 | 20 | 00 | 76 | 07 | 4.00 | |
| 00 | 4.0.4 | 445 | 447 | 0.2 | 400 | 102 | 0.2 | 91 | 84 | 21 | 128 | 89 | 98 | 98 | 99 | 103 | 38 | 99 | 76 | 97 | 102 | |
| 90 | 101 104 | 115 | 117 | 92 81 | 100 | 103 | 92 91 | 96 | 101 | 83 | 84 | 79 | 13 | 21 26 | 100 84 | 70 90 | 86 | 102 | 105 | 100 89 | 94 | |
| 88 83 | 99 | 104 107 | 102 116 | 72 | 88 64 | 98 95 | 76 | 104 94 | 117 105 | 45 74 | 108 87 | 64 55 | 22 92 | 26 90 | 84 83 | 90 75 | 75 82 | 100 103 | 87 85 | 89 95 | 98 98 | |
| 90 | 88 | 129 | 112 | 96 | 93 | 95 | 91 | 106 | 98 | 65 | 104 | 73 | 102 | 96 | 95 | 78 | 98 | 105 | 60 | 75 | 86 | |
| 50 | 00 | 125 | 112 | 50 | 55 | 55 | 51 | 69 | 15 | 20 | 86 | 82 | 69 | 114 | 115 | 27 | 46 | 105 | 109 | 99 | 44 | |
| 54 | 92 | 4 | 5 | 5 | 29 | 23 | 5 | 05 | | | 00 | 02 | 05 | 114 | 115 | 27 | -0 | 105 | 105 | 55 | | |
| 113 | 124 | 104 | 105 | 80 | 83 | 62 | 90 | 85 | 97 | 69 | 110 | 40 | 96 | 106 | 98 | 96 | 72 | 71 | 79 | 93 | 106 | |
| 56 | 103 | 13 | 5 | 2 | 62 | 39 | 19 | 92 | 88 | 40 | 92 | 98 | 55 | 57 | 93 | 94 | 51 | 74 | 88 | 59 | 14 | |
| 96 | 85 | 102 | 101 | 76 | 50 | 58 | 77 | 87 | 76 | 38 | 83 | 78 | 85 | 85 | 97 | 95 | 39 | 83 | 10 | 78 | 91 | |
| 93 | 95 | 94 | 90 | 85 | 99 | 39 | 39 | 93 | 82 | 80 | 109 | 79 | 76 | 74 | 89 | 73 | 82 | 102 | 109 | 85 | 71 | |
| 57 | 104 | 14 | 1 | 7 | 60 | 87 | 1 | 101 | 97 | 68 | 120 | 93 | 80 | 66 | 100 | 90 | 79 | 102 | 76 | 85 | 80 | |
| 95 | 109 | 107 | 122 | 18 | 6 | 7 | 0 | 104 | 103 | 27 | 52 | 41 | 97 | 99 | 82 | 94 | 62 | 114 | 50 | 107 | 113 | |
| 101 | 91 | 98 | 104 | 87 | 65 | 36 | 31 | 89 | 88 | 62 | 86 | 61 | 97 | 77 | 108 | 100 | 95 | 83 | 87 | 88 | 76 | |
| 101 | 96 | 28 | 34 | 77 | 28 | 29 | 63 | 104 | 91 | 9 | 92 | 42 | 94 | 111 | 94 | 70 | 22 | 72 | 104 | 101 | 84 | |
| 87 | 110 | 29 | 12 | 12 | 49 | 66 | 28 | 94 | 79 | 20 | 117 | 64 | 43 | 70 | 93 | 71 | 22 | 90 | 87 | 74 | 27 | |
| 49 | 109 | 21 | 42 | 33 | 84 | 39 | 83 | 106 | 84 | 6 | 6 | 4 | 86 | 113 | 99 | 83 | 60 | 52 | 78 | 107 | 107 | |
| 100 | 107 | 106 | 103 | 41 | 19 | 84 | 22 | 94 | 57 | 15 | 53 | 45 | 65 | 75 | 93 | 43 | 59 | 97 | 53 | 74 | 23 | с с |
| 96 | 102 | 60 | 71 | 50 | 30 | 22 | 1 | 95 | 59 | 6 | 72 | 36 | 93 | 84 | 94 | 14 | 26 | 90 | 101 | 82 | 61 | c |
| | | | | | | | | 48 | 75 | 38 | 99 | 86 | 8 | 6 | 54 | 66 | 10 | 83 | 105 | 106 | 103 | |
| 99 | 98 | 103 | 104 | 101 | 67 | 95 | 77 | 101 | 95 | 82 | 99 | 67 | 96 | 90 | 76 | 93 | 63 | 98 | 105 | 83 | 69 | |
| 90 | 100 | 58 | 60 | 1 | 0 | 0 | 0 | 80 | 36 | 102 | 110 | 77 | 70 | 93 | 99 | 21 | 34 | 99 | 78 | 85 | 79 | ſ |
| 91 | 96 | 86 | 98 | 58 | 68 | 46 | 7 | 97 | 22 | 11 | 15 | 19 | 0 | 1 | 77 | 8 | 53 | 94 | 89 | 87 | 116 | c c |
| 91 | 104 | 71 | 50 | 3 | 4 | 13 | 3 | 95 | 78 | 38 | 65 | 34 | 91 | 93 | 96 | 66 | 39 | 93 | 79 | 94 | 103 | c |
| | | | | | | | 69 | | | | | | | | | | | | | | | - |
| | | | | | | | | 109 | 20 | 105 | 86 | 90 | 3 | 3 | 96 | 89 | 71 | 93 | 76 | 94 | 96 | - |
| | | | | | | | | 2 | 2 | 18 | 22 | 5 | 75 | 68 | 84 | 54 | 20 | 94 | 97 | 73 | 75 | |
| | | | | | | | | 15 | 6 | 62 | 116 | 102 | 101 | 94 | 60 | 23 | 16 | 112 | 9 | 66 | 38 | 2 |
| | | | | | | | | | | | | | | | | | | 109 | | | | |

63 54 10 25

MKK4: MKK6: MKK7: MSK1: mTOR: PDK1: PHK: PI3K: PKA: PKA: PKCα: PKCα: PKCδ: PKG:

Mitogen-activated protein kinase kinase 4 Mitogen-activated protein kinase kinase 6 Mitogen-activated protein kinase kinase 7 Mitogen- and Stress-activated protein kinase-1 Mammalian target of rapamycin 3-Phosphoripase kinase Phosphoripase kinase Phosphoripase kinase Protein kinase A Protein kinase Bα Protein kinase Cα Protein kinase Cα Protein kinase G

PRAK: PRK2: ROCK: ROCK-II: RSK: S6K1: SAPK2a/p38; SAPK2b/p38j2: SAPK3/p38; SAPK4/p383; SGK: SGK: SK-MLCK: Sm-MLCK:

p38-regulated/activated kinase PKC-related kinase 2 Rho-associated protein kinase Rho-associated coiled-coil forming protein kinase-II Ribosomal S6 kinase p70 S6 kinase p38 kinase p38 kinase p38 kinase p38 kinase p38 kinase serum and glucocorticoid-induced kinase Skeletal myosin light chain kinase Smooth myosin light chain kinase

 \sim Μ ø ß . م