# CCMCL1 Mantle Cell Lymphoma Cell Line

Cancer Cell Line Cat. # SCC132

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES. NOT FOR HUMAN OR ANIMAL CONSUMPTION. Pack size: ≥1X10^6 viable cells/vial

Store in liquid nitrogen



page 1 of 3

## Background

Mantle Cell Lymphoma (MCL) is a rare B-cell non-Hodgkin lymphoma (NHL) that primarily affects men over the age of 50. In MCL, B-lymphocytes located in the "mantle zone" of the lymph nodes undergo malignant transformation resulting in uncontrolled proliferation that causes enlargement of the lymph nodes. MCL cells can spread to other tissues such as bone marrow, liver and the gastrointestinal tract. Over 90% of MCL patients overexpress cyclin D1, a cell cycle protein that contributes to abnormal proliferation of the malignant cells. Although cyclin D1 overexpression is believed to directly contribute to the tumorigenesis of MCL, the pathogenesis of the disease is complex and is not fully understood<sup>1, 2</sup>.

CCMCL1 cell line is a new cell culture model of aggressive mantle cell lymphoma<sup>3</sup>. Primary MCL cells of a 58 year old man with progressive MCL and hyperleukocytosis were engrafted in NOD-*SCID-* $\gamma$  mice. Splenocytes obtained from the engrafted mouse were cultured in vitro and yielded the CCMCL1 cell line. CCMCL1 cells possess a blastoid nuclear morphology and proliferate readily in suspension culture<sup>3</sup>. CCMCL1 cells express CD19, CD20, and CD45 and are negative for CD3 and CD25.

## Short Tandem Repeat (STR) Profile

D16S539: 12, 13 CSF1PO: 10, 13 Penta D: 12, 13
vWA: 14, 18
D8S1179: 13, 14
TPOX: 8, 11
FGA: 19, 22.2
Amelogenin: X, Y

Cancer cell lines are inherently genetically unstable. Genetic instability may arise in the form of loss of heterozygosity of alleles at one or more genetic sites with increased passages.

## Storage and Handling

CCMCL1 Mantle Cell Lymphoma Cell Line should be stored in liquid nitrogen. The cells can be cultured for at least 10 passages after initial thawing without significantly affecting the cell marker expression and functionality.

# **Quality Control Testing**

- Each vial contains ≥ 1X10<sup>6</sup> viable cells.
- Cells are tested negative for Epstein-Barr virus, HPV-16, HPV-18, Hepatitis A, C, Herpesvirus type 6, 7, 8 and HIV-1 & 2 viruses by PCR
- Cells are negative for mycoplasma contamination.
- Each lot of cells is genotyped by STR analysis to verify the unique identity of the cell line.

## **Representative Data**



Figure 1. Day 1 after thaw (A). Flow analysis of cell surface molecules. CCMCL1 cells express high levels of CD19 (B), CD20 (C) and CD45 (D) and are negative for CD3 (E) and CD25 (F).

SPECIES LEGEND: H Human Ca Canine M Mouse R Rat Rb Rabbit B Bovine P Porcine WR Most Common Vertebrates

Please visit www.millipore.com for additional product information and references.

Submit your published journal article, and earn credit toward future purchases. Visit www.millipore.com/publicationrewards to learn more!

# Protocols

CCMCL1 are suspension cells that are very small in size. Passage when the cell density reaches 1–1.5 million cells/mL. Optimal plating density should be ~250,000 cells/mL.

#### Thawing Cells

 Do not thaw the cells until the recommended medium is on hand. Cells can grow on normal tissue culture ware surfaces without any additional coating.

Cells are thawed and expanded in RPMI-1640 (Sigma Cat. No. R0883), 10% FBS (Cat. No. ES-009-B), 1X L-Glutamine (Cat. No. TMS-002-C) and 1X Penicillin-Streptomycin Solution (Cat. No. TMS-AB2-C) (optional).

 Remove the vial of frozen CCMCL1 cells from liquid nitrogen and incubate in a 37°C water bath. Closely monitor until the cells are completely thawed. Maximum cell viability is dependent on the rapid and complete thawing of frozen cells.

#### IMPORTANT: Do not vortex the cells.

- 3. As soon as the cells are completely thawed, disinfect the outside of the vial with 70% ethanol. Proceed immediately to the next step.
- In a laminar flow hood, use a 1 or 2 mL pipette to transfer the cells to a sterile 15 mL conical tube. Be careful not to introduce any bubbles during the transfer process.
- 5. Using a 10 mL pipette, slowly add dropwise 9 mL of CCMCL1 Expansion Medium (Step 1 above) to the 15 mL conical tube.

IMPORTANT: Do not add the entire volume of media all at once to the cells. This may result in decreased cell viability due to osmotic shock.

6. Gently mix the cell suspension by slowly pipetting up and down twice. Be careful not to introduce any bubbles.

#### IMPORTANT: Do not vortex the cells.

- 7. Centrifuge the tube at 300 x g for 2-3 minutes to pellet the cells.
- Decant as much of the supernatant as possible. Steps 5-8 are necessary to remove residual cryopreservative (DMSO).
- 9. Resuspend the cells in 10-15 mL of CCMCL1 Expansion Medium.
- 10. Transfer the cell suspension to a T25 flask.
- 11. Incubate the cells at 37°C in a humidified incubator with 5%  $\mbox{CO}_2.$

## Subculturing Cells

CCMCL1 suspension cells require media replenishment every 2-3 days. Passage cells when the cell density is at 1 -1.5 million cells/mL.

- 1. Remove flask from incubator, tighten cap and place in tissue culture hood.
- Dislodge any cells that may adhere to the flask by firmly rapping the side of the flask with the palm of the hand and gently swirl the medium over the cells to mix. Visually inspect flask to ensure the cells have been dislodged and the suspension is free of contaminants.
- Determine cell count and viability using a hemocytometer or automated cell counter.
- 4. Cells are typically plated at a density of 250,000 cells/mL

### Cryopreservation of Cells

CCMCL1 Human Mantle Cell Lymphoma Cell Line may be frozen in RPMI-1640 medium containing 30-40% FBS and 10% DMSO using a Nalgene slow freeze Mr. Frosty container.

## References

- Pérez-Galán P, Dreyling M, Wiestner A (2011) Mantle cell lymphoma: biology pathogenesis and the molecular basis of treatment in the genomic era. *Blood* 117(1): 26-38.
- Inamdar AA, Goy A, Ayoub NM, Attia C, Oton L, Taruvai V, Costales M, Lin YT, Pecora A, Suh KS (2016) Mantle cell lymphoma in the era of precision medicine-diagnosis, biomarkers and therapeutic agents. *Oncotarget* 7(30: 48692-48731.
- Zhao X, Chen-Kiang S, Shetty S, Di Liberto M, Bodo J, Durkin L, Eng K, Elemento O, Smith MR, Hsi ED (2015) CCMCL1: a new model of aggressive mantle cell lymphoma. *Blood* 125(17): 2730-2732.

📕 antibodies 📕 Multiplex products 📕 biotools 📕 cell culture 🔳 enzymes 📕 kits 📕 proteins/peptides 📒 siRNA/cDNA products

#### Please visit www.millipore.com for additional product information, test data and references

FOR RESEARCH USE ONLY. Not for use in diagnostic or therapeutic applications. Purchase of this Product does not include any right to resell or transfer, either as a stand-alone product or as a component of another product. Any use of this Product for purposes other than research is strictly prohibited without prior written authorization from an authorized officer of Millipore Corporation. Upstate®, Chemicon® and all other trademarks are owned by Millipore Corporation. Copyright ©2007 Millipore Corporation. All rights reserved.

# ACADEMIC USE AGREEMENT (subject to local law)

THIS PRODUCT MAY ONLY BE USED BY INDIVIDUALS EMPLOYED BY AN ACADEMIC INSTITUTION AND IS INTENDED SOLELY TO BE USED FOR ACADEMIC RESEARCH, WHICH IS FURTHER DEFINED BELOW. BY OPENING THIS PRODUCT, YOU ("PURCHASER") HEREBY REPRESENT THAT YOU HAVE THE RIGHT AND AUTHORITY TO LEGALLY BIND YOURSELF AND/OR YOUR EMPLOYER INSTITUTION, AS APPLICABLE, AND CONSENT TO BE LEGALLY BOUND BY THE TERMS OF THIS ACADEMIC USE AGREEMENT. IF YOU DO NOT AGREE TO COMPLY WITH THESE TERMS, YOU MAY NOT OPEN OR USE THE PRODUCT AND YOU MUST CALL MILLIPORESIGMA ("SELLER") CUSTOMER SERVICE (1-800-645-5476) TO ARRANGE TO RETURN THE PRODUCT FOR A REFUND.

"Product" means CCMCL1 Mantle Cell Lymphoma Cell Line (SCC132)

"Academic Research" means any internal *in vitro* research use by individuals employed by an academic institution. Academic Research specifically excludes the following uses of whatever kind or nature:

- Re-engineering or copying the Product
- Making derivatives, modifications, or functional equivalents of the Product
- Obtaining patents or other intellectual property rights claiming use of the Product
- Using the Product in the development, testing, or manufacture of a Commercial Product
- Using the Product as a component of a Commercial Product
- Reselling or licensing the Product
- Using the Product in clinical or therapeutic applications including producing materials for clinical trials
- Administering the Product to humans
- Using the Product in collaboration with a commercial or non-academic entity

"Commercial Product" means any product intended for: (i) current or future sale; (ii) use in a fee-for-service; or (iii) any diagnostic, clinical, or therapeutic use.

Access to the Product is limited solely to those officers, employees, and students of PURCHASER's academic institution who need access to the Product to perform Academic Research. PURCHASER shall comply with all applicable laws in its use and handling of the Product and shall keep it under reasonably safe and secure conditions to prevent unauthorized use or access.

These use restrictions will remain in effect for as long as PURCHASER possesses the Product.

COMMERCIAL OR NON-ACADEMIC ENTITIES INTERESTED IN PURCHASING OR USING THE PRODUCT MUST CONTACT licensing@emdmillipore.com AND AGREE TO SEPARATE TERMS OF USE PRIOR TO USE OR PURCHASE.

📕 antibodies 📕 Multiplex products 📕 biotools 📕 cell culture 📕 enzymes 📕 kits 📕 proteins/peptides 📙 siRNA/cDNA products

We Buy 100% Certifier Renewable Energy

Please visit www.millipore.com for additional product information, test data and references EMD Millipore Corporation, 28820 Single Oak Drive, Temecula, CA 92590, USA 1-800-437-7500 Technical Support: T: 1-800-MILLIPORE (1-800-645-5476) • F: 1-800-437-7502 FOR RESEARCH USE ONLY. Not for use in diagnostic procedures. Not for human or animal consumption. Purchase of this Product does not include any right to resell or transfer, either as a stand-alone product or as a component of another product. Any use of this Product for purposes other than research is strictly prohibited. Millipore®, the M mark, Upstate®, Chemicon®, Linco® and all other trademarks, unless specifically identified above in the text as belonging to a third party, are owned by Merck KGaA, Darmstadt. Copyright ©2008-2017 Merck KGaA, Darmstadt. All rights reserved.