



Product Information

2,3-Dimercapto-1-propanol

Product Number **D 1129**
Storage Temperature 2-8 °C

Product Description

Molecular Formula: $C_3H_8OS_2$
Molecular Weight: 124.2
CAS Number: 59-52-9
Boiling Point: 140 °C (40 mm Hg)¹
Density: 1.2385 g/ml¹
Synonym: dimercaprol; 2,3-dimercaptopropanol;
1,2-dithioglycerol; British Anti-Lewisite; BAL

2,3-Dimercapto-1-propanol is a reducing agent that serves a similar function to dithiothreitol. It is also a chelating agent for heavy metals such as lead, mercury, and arsenic, and has been studied in toxicology research.² Dimercaprol can maintain sulfhydryl groups in the reduced state and can reverse arsenic-related inhibition of various enzymes, such as pyruvic and α -ketoglutaric oxidases.³

Dimercaprol has been used to chelate labile zinc and subsequently alter eotaxin, RANTES, and MCP-1 production in stimulated human airway epithelium and fibroblasts.⁴ Dimercaprol (3 mM) has been utilized to probe the inhibition of activity of voltage-gated Ca^{2+} channels in *Xenopus* oocytes by syntaxin 1A.⁵ A comparative study has indicated that dimercaprol can inhibit Ca^{2+} transport by Ca^{2+} -ATPases in the microsome sarco/endoplasmic reticulum from brain, but not from fast-skeletal muscle.⁶

Precautions and Disclaimer

For Laboratory Use Only. Not for drug, household or other uses.

Preparation Instructions

This product is miscible in ethanol (0.1 ml/ml, v/v), yielding a clear, colorless solution. It is also soluble in chloroform (100 mg/mL, w/v). This product is soluble in ether and in vegetable oils, and in water (68 mg/ml, w/v) with slow decomposition.²

Storage/Stability

This product oxidizes in aqueous solution to a disulfide, in the presence of air. At pH 5, the rate is negligible, but in alkaline solution, the rate is rapid.²

References

1. The Merck Index, 12th ed., Entry# 3255.
2. Aposhian, H. V., et al., Mobilization of heavy metals by newer, therapeutically useful chelating agents. *Toxicology*, **97(1-3)**, 23-38 (1995).
3. Data for Biochemical Research, 3rd ed., Dawson, R. M. C., et al., Oxford University Press (New York, NY: 1986), p. 380-381.
4. Richter, M., et al., Zinc chelators inhibit eotaxin, RANTES, and MCP-1 production in stimulated human airway epithelium and fibroblasts. *Am. J. Physiol. Lung Cell. Mol. Physiol.*, **285(3)**, L719-729 (2003).
5. Arien, H., et al., Syntaxin 1A modulates the voltage-gated L-type calcium channel ($Ca_v1.2$) in a cooperative manner. *J. Biol. Chem.*, **278(31)**, 29231-29239 (2003).
6. Quinhones, E. B., et al., 2,3-Dimercaptopropanol inhibits Ca^{2+} transport in microsomes from brain but not from fast-skeletal muscle. *Neurochem Res.*, **26(3)**, 251-256 (2001).

GCY/CRF 1/04

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.