

For life science research only.
Not for use in diagnostic procedures.



Biotin-High Prime

 **Version: 09**

Content Version: November 2020

For the nonradioactive labeling of DNA with Biotin-16-dUTP, using random oligonucleotides as primers. Premixed solution for 25 labeling assays.

Cat. No. 11 585 649 910 100 µl
25 labeling assays

Store product at –15 to –25°C.

1.	General Information	3
1.1.	Contents	3
1.2.	Storage and Stability	3
	Storage Conditions (Product)	3
1.3.	Additional Equipment and Reagent required	3
1.4.	Application	3
2.	How to Use this Product	4
2.1.	Before you Begin	4
	Sample Materials	4
	Templates for labeling reaction	4
	Safety Information	4
	For customers in the European Economic Area	4
	Working Solution	4
2.2.	Protocols	4
	Random primed DNA labeling	4
	Labeling DNA isolated from low-melting point agarose	5
	Labeling reaction yield	5
	Detection of biotin-labeled DNA	5
2.3.	Parameters	6
	Chemical Name	6
	Structural formula	6
3.	Additional Information on this Product	6
3.1.	Test Principle	6
4.	Supplementary Information	7
4.1.	Conventions	7
4.2.	Changes to previous version	7
4.3.	Ordering Information	7
4.4.	Trademarks	8
4.5.	License Disclaimer	8
4.6.	Regulatory Disclaimer	8
4.7.	Safety Data Sheet	8
4.8.	Contact and Support	8

1. General Information

1.1. Contents

Vial / Bottle	Label	Function / Description	Content
1	Biotin-High Prime, 5x conc.	Random prime labeling mixture: Premixed solution of 1 U/ μ l Klenow polymerase, labeling grade, 1 mM dATP, 1 mM dCTP, 1 mM dGTP each, 0.65 mM dTTP, 0.35 mM Biotin-16-dUTP, and 5x stabilized reaction buffer in 50% (v/v) glycerol.	1 vial, 100 μ l

1.2. Storage and Stability

Storage Conditions (Product)

When stored at -15 to -25°C , the product is stable through the expiry date printed on the label.

Vial / Bottle	Label	Storage
1	Biotin-High Prime, 5x conc.	Store at -15 to -25°C . ⚠ Avoid repeated freezing and thawing. ⚠ To avoid contamination, aliquot and store the solution in 2 to 3 vials.

1.3. Additional Equipment and Reagent required

For random primed DNA labeling

- 0.2 M EDTA, pH 8.0
- Autoclaved, double-distilled water
- Water bath
- Ice bath

For labeling DNA isolated from agarose

- 0.2 M EDTA, pH 8.0
- Autoclaved, double-distilled water
- Water bath

For detection of biotin-labeled DNA

- CSPD*
- CDP-*Star**
- NBT/BCIP*

1.4. Application

Biotin-High Prime-labeled probes are used in a variety of hybridization reactions:

- Southern blots
- Northern blots
- Dot/slot blots
- Screening of gene libraries
- *In situ* hybridizations

2. How to Use this Product

2.1. Before you Begin

Sample Materials

Templates for labeling reaction

- DNA fragments of at least 100 bp.
- Linearized plasmid, cosmid, or λ DNA.
- Supercoiled DNA
- Minimal amounts of DNA (10 ng), such as DNA restriction fragments isolated from gels or in molten agarose.

i For the Random primed DNA labeling protocol, 0.01 to 3 μ g of template can be labeled.

Safety Information

For customers in the European Economic Area

Contains SVHC: octyl/nonylphenol ethoxylates. For use in research and under controlled conditions only – acc. to Art. 56.3 and 3.23 REACH Regulation.

Working Solution

Solution	Composition	Use	Storage and Stability
Water	Autoclaved, double-distilled water	Dilution of DNA.	Store at +15 to +25°C.
EDTA	0.2 M ethylenediaminetetraacetic acid, pH 8.0	Stops the reaction.	Store at +15 to +25°C.

2.2. Protocols

Random primed DNA labeling

Perform the standard random primed DNA labeling according to the following steps.

- 1** To a reaction vial, add 1 μ g template DNA (linear or supercoiled) and autoclaved, double-distilled water to a final volume of 16 μ l.
- 2** Denature the DNA by heating in a boiling water bath for 10 minutes; quickly chill in an ice/water bath.
i Full denaturation is essential for efficient labeling.
- 3** Mix Biotin-High Prime thoroughly and add 4 μ l to the denatured DNA; mix, and centrifuge briefly.
– Incubate for 1 hour or overnight at +37°C.
i Longer incubations up to 20 hours increase the yield of Biotin-labeled DNA, see **Table, Labeling reaction yield**.
- 4** Stop the reaction by adding 2 μ l 0.2 M EDTA, pH 8.0, and/or by heating to +65°C for 10 minutes.
i The length of the Biotin-labeled fragments obtained with Biotin-High Prime range from 200 bp to $\geq 1,000$ bp, depending on the length of the original template.

Labeling DNA isolated from low-melting point agarose

- 1 Excise the DNA fragment to be labeled cleanly from a low-melting point agarose gel and transfer it to a 1.5 ml microfuge tube.

- 2 Add sterile, double-distilled water to a ratio of 3 ml/g gel and heat the tube for 7 minutes at +100°C to melt the gel and denature the DNA.
– After cooling to +37°C, the DNA/agarose mixture can be used directly for labeling.

- 3 Mix Biotin-High Prime thoroughly; add 4 µl, mix, and centrifuge briefly.

- 4 Incubate overnight at +37°C.

- 5 Stop the reaction by adding 2 µl 0.2 M EDTA, pH 8.0 and/or by heating to +65°C for 10 minutes.

Labeling reaction yield

The labeling efficiency is shown in the following table and in Figure 1.

- i* Reactions were performed with increasing amounts of different template DNAs for 1 hour and 20 hours. The yield of biotin-labeled DNA was determined by the incorporation of a radioactive tracer and confirmed by a dot blot. Numbers shown are the average of independent labeling assays.




Fig. 1: Yield of biotin-labeled DNA from different amounts of template DNAs for 1 hour and 20 hours incubation of the Biotin-High Prime reaction at +37°C.

Detection of biotin-labeled DNA

Biotin-labeled DNA is detected by streptavidin conjugated to the enzyme alkaline phosphatase which catalyzes a color reaction with 5-bromo-4-chloro-3-indolyl-phosphate and 4-nitro blue tetrazolium chloride (NBT/BCIP*) or by a chemiluminescent reaction with CSPD* or CDP-Star*. Alternatively, especially for *in situ* applications, biotin-labeled hybrids can also be detected by streptavidin conjugated to different fluorochromes.

2.3. Parameters

Chemical Name

Structural formula




Fig. 2: Structure of Biotin-16-dUTP.

3. Additional Information on this Product

3.1. Test Principle

Biotin-labeled DNA probes are generated with Biotin-High Prime according to the random primed labeling technique.

- ① The complementary DNA strand is synthesized by Klenow polymerase using the 3'OH termini of the random oligonucleotides as primers.

 - ② Biotin-16-dUTP is incorporated into the newly synthesized complementary DNA strand.
 - The premixed Biotin-High Prime is ideal for for convenient and efficient nonradioactive labeling of DNA with biotin.
-

4. Supplementary Information

4.1. Conventions

To make information consistent and easier to read, the following text conventions and symbols are used in this document to highlight important information:

Text convention and symbols

 *Information Note: Additional information about the current topic or procedure.*

 **Important Note: Information critical to the success of the current procedure or use of the product.**

① ② ③ etc. Stages in a process that usually occur in the order listed.

① ② ③ etc. Steps in a procedure that must be performed in the order listed.

* (Asterisk) The Asterisk denotes a product available from Roche Diagnostics.

4.2. Changes to previous version

Layout changes.

Editorial changes.

New information added related to the REACH Annex XIV.

4.3. Ordering Information

Product	Pack Size	Cat. No.
Reagents, kits		
CSPD, ready-to-use	2 x 50 ml	11 755 633 001
CDP-Star, ready-to-use	2 x 50 ml	12 041 677 001
NBT/BCIP Stock Solution	8 ml	11 681 451 001

4. Supplementary Information

4.4. Trademarks

All product names and trademarks are the property of their respective owners.

4.5. License Disclaimer

For patent license limitations for individual products please refer to:

List of biochemical reagent products.

4.6. Regulatory Disclaimer

For life science research only. Not for use in diagnostic procedures.

4.7. Safety Data Sheet

Please follow the instructions in the Safety Data Sheet (SDS).

4.8. Contact and Support

To ask questions, solve problems, suggest enhancements or report new applications, please visit our **Online Technical Support Site.**

To call, write, fax, or email us, visit **sigma-aldrich.com**, and select your home country. Country-specific contact information will be displayed.

