


Product Information



15180 hemi-Calcium bis[4-(1,1,3,3-tetramethylbutyl)phenyl]phosphate

(Bis[4-(1,1,3,3-tetramethylbutyl)phenyl]phosphate Calcium salt)
Selectophore®, function tested

Electrochemical Transduction

Ion-Selective Electrodes

Application 1 and Sensor Type¹

Assay of Ca^{2+} activity with solvent polymeric membrane electrodes based on the liquid ion-exchanger Bis[4-(1,1,3,3-tetramethylbutyl)phenyl]phosphate Calcium salt.

Recommended Membrane Composition

- 7.00 wt% Bis[4-(1,1,3,3-tetramethylbutyl)phenyl]phosphate Calcium salt ([15180](#))
- 29.86 wt% Poly(vinyl chloride) high molecular weight ([81392](#))
- 63.14 wt% Di-n-octylphenylphosphonate ([12584](#))

Recommended Cell Assembly

Reference || sample solution || ion-selective electrode | 0.1 M CaCl_2 | AgCl, Ag

Electrode Characteristics and Function

Selectivity coefficients $\log K_{Ca,M}^{Pot}$ as obtained by the separate solution method in (0.1 M solutions of the chlorides).

$\log K_{Ca,Na}^{Pot}$	-4.4
$\log K_{Ca,K}^{Pot}$	-4.5
$\log K_{Ca,Mg}^{Pot}$	-4.9

Slope of linear regression:	30.5 at 25°C (10^{-5} to 10^{-0} Ca^{2+})
Detection limit:	$3.2 \cdot 10^{-6}$ M Ca^{2+}
pH range for 10^{-3} M CaCl_2 :	4.8 to 8.8
Response time:	<10 s
Operational lifetime:	3 months
Membrane resistance:	3 MΩ



Application 2 and Sensor Type²

Assay of Ca^{2+} activity with solvent polymeric membrane electrodes based on the liquid ion-exchanger Bis[4-(1,1,3,3-tetramethylbutyl)phenyl]phosphate Calcium salt.

Recommended Membrane Composition

0.10 wt% Bis[4-(1,1,3,3-tetramethylbutyl)phenyl]phosphate Calcium salt ([15180](#))
33.43 wt% Poly(vinyl chloride) high molecular weight ([81392](#))
66.43 wt% Bis(2-ethylhexyl)sebacate ([84818](#))

Recommended Cell Assembly

Reference || sample solution || ion-selective electrode | 0.1 M CaCl_2 | AgCl, Ag

Electrode Characteristics and Function

Selectivity coefficients $\log K_{Ca,M}^{Pot}$ as obtained by the separate solution method in (0.1 M solutions of the chlorides).

$\log K_{Ca,NH_4}^{Pot}$	-3.3	$\log K_{Ca,H}^{Pot}$	-2.4
$\log K_{Ca,Li}^{Pot}$	-4.1	$\log K_{Ca,Mg}^{Pot}$	-6.2
$\log K_{Ca,Na}^{Pot}$	-4.0	$\log K_{Ca,Ba}^{Pot}$	-3.0
$\log K_{Ca,K}^{Pot}$	-3.0		

Slope of linear regression: 31.9 ± 2.3 mV (10^{-4} to 10^{-1} CaCl_2)
Membrane resistance: 9.9 ± 1.0 M Ω

¹ The role of polymeric materials in the fabrication of ion-selective electrodes and biosensors. G.J. Moody, Polym. Mat. Sci. Eng. 64, 362, (1991) and ref. cited therein.

² Carrier mechanism of acidic ionophores in solvent polymeric membrane ion-selective electrodes. U. Schaller, E. Bakker, E. Pretsch, Anal. Chem. 67, 3123 (1995).

