

Product Information

4-CHLORO-1-NAPHTHOL SOLUTION

Product No. **C 8302**

Store at 2-8 °C

Product Description

Solution contains 0.48 mM 4-Chloro-1-Naphthol, 50 mM Tris-HCl and 0.2 M NaCl in 17% methanol. This is a substrate solution designed for visualizing horseradish peroxidase conjugates in Western Blotting.

Items not provided:

10X TBS (Product No. T 5912)

20% Tween 20 (Product No. P 9203)

BSA Fraction V powder (Product No. A 9647)

Nonfat dry milk (Product No. M 7409)

Instructions for Use:

1. After the gel is transferred onto a blotting membrane, wash the membrane for 5 minutes with the washing solution (See Table A for selection of appropriate solutions).
2. Incubate the rinsed membrane with primary antibody diluted in blocking solution (See Table A for selection of appropriate solutions) for 2 hours at room temperature with gentle agitation. Use a dilution based on the manufacturer's recommendation. A blocking step prior to step 2 is usually not necessary.
3. Wash the membrane for 5 minutes with the washing solution.
4. Incubate the washed membrane with secondary antibody peroxidase conjugate in blocking solution for 2 hours at room temperature with gentle agitation. (A 1:1000 dilution of antibody is recommended).
5. Wash the membrane 3 times for 5 minutes each in washing solution.
6. Add hydrogen peroxide (Product No. H 1009) to 4-Chloro-1-Naphthol Solution (Product No. C 8302) to obtain a final concentration of 0.01% hydrogen peroxide (v/v). Prepare immediately before use.
7. Cover the membrane with the substrate solution for 1-5 minutes at room temperature until the desired color is obtained. Use 10-20 ml for a 8x10 cm membrane. **MAKE SURE THE SURFACE OF THE MEMBRANE IS COMPLETELY COVERED WITH SUBSTRATE SOLUTION.**
8. The color development can be stopped by extensive washing with water. **WATCH CAREFULLY DURING COLOR DEVELOPMENT TO AVOID OVER DEVELOPING.**

TABLE A: Several commonly used washing and blocking solutions are listed. Selection depends on the stringency required. Condition 1 is the least stringent and condition 3 is most stringent. Washing solution 3 and blocking solution 3 with 2% nonfat milk is recommended for most applications.

Condition	Washing Solution	Blocking Solution
1	1X TBS*	Washing solution + 0.5% Tween 20
2	1X TBS + 0.05% Tween 20	Washing solution + 1-5% BSA
3	1X TBS + 0.05% Tween 20	Washing solution + 2-4% Nonfat dry milk

*1X TBS = 20 mM Tris-HCl, pH 7.5 and 0.9% NaCl

References

1. B. Batteiger, Journal of Immunological Methods, **55** (1982) 297-307.
2. D.I. Scott, Journal of Immunological Methods, **119** (1989) 153-187.

3. F. Miescher, Analytical Biochemistry, **119** (1982) 142-147.

Pcs2/01

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications.

Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply.

Please see reverse side of the invoice or packing slip.