sigma-aldrich.com



3050 Spruce Street, St. Louis, MO 63103 USA Tel: (800) 521-8956 (314) 771-5765 Fax: (800) 325-5052 (314) 771-5757 email: techservice@sial.com sigma-aldrich.com

Product Information

18723 Mant-ATP triethylammonium salt solution

Molecular Formula: Molecular Weight: Form: Concentration: pH: Spectroscopic Properties:

Package quantity:

Unit Definition:

 $\begin{array}{l} C_{18}H_{23}N_6O_{14}P_3 \mbox{ (free acid)} \\ 640.33 \mbox{ g/mol (free acid)} \\ clear aqueous solution \\ 10 \mbox{ mM} \\ 7.5 \mbox{ +/-}0.5 \\ \lambda_{max} \mbox{ 255/355 nm}; \mbox{ } \lambda_{exc} \mbox{ 355 nm}; \mbox{ } \lambda_{em} \mbox{ 448 nm}; \\ \epsilon \mbox{ 23.3/5.8 L mmol}^{-1} \mbox{ cm}^{-1} \mbox{ (Tris-HCl pH 7.5)} \\ 150 \mbox{ Units} \\ 1 \mbox{ unit = 1 ul of a 10 mM solution} \end{array}$



Structural formula of Mant-ATP

Applications:

Conformational dynamic: DnaB/C-protein^[1], Csk^[2] Inhibition of AC-isoforms^[3] Fluorescence stop-flow kinetics: PKA^[4] FRET: AC^[5], myosin V^[6]

Specific Ligands:

DnaC-protein^[7] Factor Rho^[8]

Selected References:

[1] Roychowdhury *et al.* (2009) Interactions of Escherichia coli DnaB-DnaC protein complex with nucleotide cofactors. 1. Allosteric conformational transitions of the complex. *Biochemistry* **48**:6712.

[2] Shaffer *et al.* (2001) Nucleotide release and associated conformational changes regulate function in the COOH-terminal Src kinase, Csk. *Biochemistry* **40**:11149.

[3] Gille *et al.* (2004) Differential inhibition of adenylyl cyclase isoforms and soluble guanylyl cyclase by purine and pyrimidine nucleotides. *J. Biol. Chem.* **279**:19955.

[4] Ni *et al.* (2000) Insights into nucleotide binding in protein kinase A using fluorescent adenosine derivatives. *Protein Science* **9**:1818.

[5] Goettle *et al.* (2007) Molecular analysis of the interaction of Bordetella pertussis adenylyl cyclase with fluorescent nucleotides. *Molecular Pharmacology* **72** (3):526.

[6] Sun *et al.* (2006) Dynamics of the upper 50-kDa domain of myosin V examined with fluorescence resonance energy transfer. *J. Biol. Chem.* **281**:5711.

[7] Galletto *et al.* (2005) The nucleotide-binding site of the Escherichia coli DnaC protein: Molecular topography of DnaC protein-nucleotide cofactor complex. *Cell Biochem. And Biophys.* **43**:331.

[8] Jeong *et al.* (2004) Nucleotide binding induces conformational changes in Escherichia coli transcription termination factor Rho. *J. Biol. Chem.* **279**:1837.

Seifert et al. (2012) Inhibitors of membranous adenylyl cyclases. Trends Pharmacol. Sci. 33(2):64.

Wang *et al.* (2011) Charge isomers of myelin basic protein: structure and interactions with membranes, nucleotide analogues, and calmodulin. *PLoS One.* **6(5)**:e19915.

Pinto *et al.* (2011) Structure-activity relationships for the interactions of 2'- and 3'-(O)-(N-methyl)anthraniloyl-substituted purine and pyrimidine nucleotides with mammalian adenylyl cyclases. *Biochem. Pharmacol.* **82(4)**:358.

Hübner *et al.* (2011) Effect of MANT-nucleotides on L-type calcium currents in murine cardiomyocytes. *J. Naunyn Schmiedebergs Arch.Pharmacol.* **383(4)**:573.

Spangler *et al.* (2011) Interaction of the diguanylate cyclase YdeH of Escherichia coli with 2',(3')-substituted purine and pyrimidine nucleotides. *J. Pharmacol. Exp. Ther.* **336(1)**:234.

Erdorf *et al.* (2011) Pharmacological characterization of adenylyl cyclase isoforms in rabbit kidney membranes. *J. Naunyn Schmiedebergs Arch. Pharmacol.* **383(4)**:357.

Goettle *et al.* (2010) Cytidylyl and uridylyl cyclase activity of bacillus anthracis edema factor and Bordetella pertussis CyaA. *Biochemistry***49(26)**:5494.

Chen *et al.* (2009) Structural Basis of Mechanochemical Coupling in a Hexameric Molecular Motor. *J. Biol. Chem.* **283** (6):3607.

Kainov *et al.* (2008) ADP but Not Pi Dissociation Contributes to Rate Limitation for *Escherichia coli* Rho. *J. Biol. Chem.* **284** (**49**):33773.

Kambach *et al.* (2007) Human OLA1 Defines an ATPase subfamily in the Obg Family of GTP-binding proteins. *J. Biol. Chem.* **282** (27):19928.

Scott *et al.* (2006) Chain-Terminating Dinucleoside Tetraphosphates are substrates for DNA Polymerization by Human Immunodefiency Virus Type 1 Reverse Transcriptase with Increased Activity against Thymidine Analogue-Resistant Mutants. *AAC* **50** (11):3607.

Lísal et al. (2005) Cooperative Mechanism of RNA Packaging Motor. J. Biol. Chem. 280 (24):23157.

Tung-Chung Mou *et al.* (2005) Structural Basis for the Inhibition of Mammalian Membrane Adenylyl Cyclase by 2'(3')-O-(N-Methylanthraniloyl)-guanosine 5'-Triphosphate. *J. Biol.Chem.* **280** (8):7253.

Fujita *et al.* (1999) Fluorescence changes of a label attached near the myosin active site on nucleotide binding in rat skeletal muscle fibres. *J. Physiol.* **515**:869.

Thoenges *et al.* (1999) Tight binding of bulky fluorescent derivatives of adenosine to the low affinity E2ATP site leads to inhibition of Na+/K+-ATPase. Analysis of structural requirements of fluorescent ATP derivatives with a Koshland-Nemethy-Filmer model of two interacting ATP sites. *J. Biol. Chem.* **274**:1971.

Cheng et al. (1998) Interaction of mant-adenosine nucleotides and magnesium with kinesin. Biochemistry 37:5288.

Mocz *et al.* (1998) Probing the nucleotide binding sites of axonemal dynein with the fluorescent nucleotide analogue 2'(3')-O-(-N-Methylanthraniloyl)-adenosine 5'-triphosphate. *Biochemistry* **37**:9862.

Jameson et al. (1997) Fluorescent nucleotide analogs: synthesis and applications. Methods Enzymol. 278:363.

Hiratsuka (1983) New ribose-modified fluorescent analogs of adenine and guanine nucleotides available as substrates for various enzymes. *Biochim. Biophys. Acta* **742**:496.

Storage Conditions: store at -20 °C, protect from light

Precautions and Disclaimer:

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

©2016 Sigma-Aldrich Co. LLC. All rights reserved. SIGMA-ALDRICH is a trademark of Sigma-Aldrich Co. LLC, registered in the US and other countries. Sigma brand products are sold through Sigma-Aldrich, Inc. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see product information on the Sigma-Aldrich website at www.sigmaaldrich.com and/or on the reverse side of the invoice or packing slip.