



**SIGMA-ALDRICH**

3050 Spruce Street  
Saint Louis, Missouri 63103 USA  
Telephone 800-325-5832 • (314) 771-5765  
Fax (314) 286-7828  
email: techserv@sial.com  
sigma-aldrich.com

## Product Information

### Tunicamycin from *Streptomyces sp.*

Product Numbers **T7765, T1152, T0527**

#### Product Description

Tunicamycin is a mixture of homologous antibiotics which contain the following common building blocks: uracil, N-acetyl glycosamine, an 11-carbon aminodialdose called tunicamine and a fatty acid linked to the amino group. Tunicamycin is an inhibitor of bacterial and eukaryote N-acetylglucosamine transferases; preventing formation of N-acetylglucosamine lipid intermediates and glycosylation of newly synthesized glycoproteins.<sup>1</sup> It blocks the formation of protein N-glycosidic linkages by inhibiting the transfer of N-acetylglucosamine 1-phosphate to dolichylmonophosphate.<sup>2</sup> It enhanced the activity of interferon<sup>3</sup> and has been shown to inhibit 99% of replication of Vesticular Stomatitis virus in mammalian cells at 0.5 µg/ml. Tunicamycin is also an inhibitor of ceruloplasmin, of α2-macroglobulin, and of the secretion of α1-protease inhibitor.<sup>4</sup> It has been shown that tunicamycin induced a programmed cell death in plant cells<sup>5</sup> and in mammalian cells via stimuli of ER stress.<sup>6</sup>

There are at least 10 homologs, the main components being A, B, C, and D.<sup>7</sup> The homologs differ in their fatty acid components, which vary in chain length.<sup>7,8,9,10</sup> Homologs A1, A2, B1, and C2 inhibit N-linked glycosylation of lipids at 1-2 µg/ml, while homologs B2, C1, D1, and D2 require 5 µg/ml to achieve the same amount of inhibition. Homologs A1, B2, and C2 do not inhibit protein synthesis, at up to 100 ng/ml, whereas homologs A2, B1, C1, D1, and D2 markedly inhibit tyrosine incorporation at less than 50 ng/ml.<sup>11</sup>

#### Precautions and Disclaimer

These products are for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

#### Preparation Instructions

Tunicamycin is soluble in excess of 10 mg/ml in DMF, DMSO, and pyridine. It is soluble at less than 5 mg/ml in water at pH = 9.0, warm methanol, and warm 95%

ethanol. It is soluble at less than 1 mg/ml in dioxane and THF. It is insoluble in other organic solvents such as acetone, chloroform, and ethyl acetate and in aqueous solutions with pH less than 6. Aqueous solutions can be prepared from stock solutions by diluting with water at pH 8-10 or with buffers with pH above 7, preferably greater than 8. It will not dissolve directly in phosphate buffer, pH 8, at 1 mg/ml, even with heating, but solubility can be achieved by raising the pH to 9 and back titrating to pH 7 to 8.

#### Storage/Stability

The products are stable for 3 years if stored at the recommended temperature. Tunicamycin is unstable in acid solutions, but stable at alkaline pH.<sup>10</sup>

#### T7765 Tunicamycin

This product is the mixture of homologous antibiotics, A, B, C and D  
Molecular weights: A = 817, B = 831, C = 845, D = 859<sup>7</sup>  
Storage Temperature 2-8 °C

CAS RN: (for the complex of homologs): 11089-65-9

Melting Point: 234 - 235 °C with decomposition<sup>12</sup>  
 $\lambda_{\text{max}}$ : 205 nm, 260 nm; (methanol)  
( $E_{1\%}^{1\text{cm}} = 230, 110$ )<sup>12</sup>  
 $[\alpha]_{\text{D}}^{20} : +52^{\circ}$  ( $c = 0.5$  in pyridine)<sup>12</sup>

Sigma determines the amounts of A, B, C, and D homologs in the product. Composition may vary from lot to lot and actual content is given on the label.

#### T1152 C2 homolog

Purity 95% (HPLC)  
Storage Temperature -20 °C

CAS RN: 73942-07-1

Molecular formula: C<sub>39</sub>H<sub>64</sub>N<sub>4</sub>O<sub>16</sub>  
Molecular weight: 844.94

**T0527 A1 homolog**

Purity 95% (HPLC)

Storage Temperature -20 °C

CAS RN: 66081-37-6

Molecular formula: C<sub>37</sub>H<sub>60</sub>N<sub>4</sub>O<sub>16</sub>

Molecular weight: 816.89

**References**

1. Dawson, R.M.C., et al., Data for Biochemical Research, 3rd ed., p. 335., Oxford University Press, New York, (1986).
2. Heifetz A., et al., Mechanism of action of tunicamycin on the UDP-GlcNAc: dolichyl-phosphate transferase., *Biochem.*, **18**, 2186-92 (1979).
3. Maheshwari, R.K., et al., Tunicamycin enhances the antiviral and anticellular activity of interferon., *Science*, **219**, 1339-1341 (1983).
4. Bauer, H.C., et al., Role of carbohydrate in glycoprotein secretion by human hepatoma cells. *Biochem. Biophys. Res. Comm.*, **128**, 368-375 (1985).
5. Crosti P., et al., Tunicamycin and Brefeldin A induce in plant cells a programmed cell death showing apoptotic features. *Protoplasma*, **216**, 31-8 (2001).
6. Fujita E, et al., Caspase-12 processing and fragment translocation into nuclei of tunicamycin-treated cells., *Cell Death Differ.*, **9**, 1108-14 (2002).
7. Dictionary of Antibiotics, 1, 714 (1988).
8. Takatsuki A., et al., Tunicamycin, a new antibiotic. I. Isolation and characterization of tunicamycin.: *J. Antibiotics*, **24**, 215-23 (1971).
9. Mahoney, W.C. and Duksin, D., Separation of tunicamycin homologues by reversed-phase high-performance liquid chromatography., *J. Chromatogr.*, **198**, 506-10 (1980).
10. Ito, T. et al.; Isolation and structure of tunicamycin components., *Agric. Biol. Chem.* **44**, 695 (1980).
11. Duksin, D. and Mahoney, W.C., Relationship of the structure and biological activity of the natural homologues of tunicamycin., *J. Biol. Chem.*, **257**, 3105-3109 (1982).
12. The Merck Index, 12th ed., No. 9949 (1996).

NDH/PHC 12/04

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.