

## Product Information

### Trypsin from porcine pancreas

Catalog Number **T0303**  
Storage Temperature  $-20\text{ }^{\circ}\text{C}$

CAS RN 9002-07-05  
EC 3.4.21.4  
Molecular mass:<sup>1</sup> 23.4 kDa  
Extinction Coefficient:<sup>2</sup>  $E^{1\%} = 15.0$  (280 nm)  
pI:<sup>1,2</sup> 10.2–10.8  
pH optimum:<sup>2</sup> 7–9  
Synonyms: Tryptase, Tryptar, Cocoonase, Parenzyme, Parenzymol

#### Product Description

Trypsin is a member of the serine protease family. The active site amino acid residues of trypsin include His<sup>46</sup> and Ser<sup>183</sup>.<sup>1,3</sup> Trypsin consists of a single chain polypeptide of 223 amino acid residues. Trypsin is produced by the cleavage of the N-terminal hexapeptide from its precursor, trypsinogen, at the Lys<sup>6</sup>–Ile<sup>7</sup> bond. The amino acid sequence of trypsin is crosslinked by 6 disulfide bridges. This native form of trypsin is referred to as  $\beta$ -trypsin. Autolysis of  $\beta$ -trypsin by cleavage at its Lys<sup>125</sup>–Ser<sup>126</sup> bond<sup>4</sup> results in  $\alpha$ -trypsin, which is held together by disulfide bridges.

Trypsin will cleave peptides on the C-terminal side of lysine and arginine amino acid residues. The rate of hydrolysis is slower if an acidic residue is on either side of the cleavage site and no cleavage occurs if a proline residue is on the carboxyl side of the cleavage site. Trypsin will also hydrolyze ester and amide linkages of synthetic derivatives of amino acids such as: benzoyl-L-arginine ethyl ester (BAEE), *p*-toluenesulfonyl-L-arginine methyl ester (TAME), tosyl-L-arginine methyl ester, *N* $\alpha$ -benzoyl-L-arginine *p*-nitroanilide (BAPNA), L-lysyl-*p*-nitroanilide, and benzoyl-L-arginamide.<sup>2,3,5,6</sup>

Assuming that the pH and temperature are the same and using a molar extinction coefficient of 808 at 254 nm for BAEE, the following conversions are valid:

1 BAEE  $\mu\text{M}$  Unit = 200  $A_{253}$  BAEE Units  
1 TAME  $\mu\text{M}$  Unit = 0.27 BAEE  $\mu\text{M}$  Units  
1 BAEE  $\mu\text{M}$  Unit = 3.64 TAME Units  
1 TAME  $\mu\text{M}$  Unit = 55 BAEE  $A_{253}$  Units  
1 BAEE  $A_{253}$  Unit = 0.018 TAME  $\mu\text{M}$  Unit  
1 TAME  $\mu\text{M}$  Unit = 180 TAME  $A_{247}$  Units  
1 TAME  $A_{247}$  Unit = 0.33 BAEE Units  
1 USP Unit =  $\Delta A_{253}$  of 0.003 per minute  
1 NF Unit = 3.3  $A_{253}$  BAEE Units.<sup>7</sup>

**Note:** These activity conversions were determined using bovine trypsin. However, they are thought to be similar for porcine trypsin.

The oxidized B chain of insulin is often used as a substrate to determine the suitability of trypsin for use in protein sequencing. The presence of two peptide bonds (Arg<sup>22</sup>–Gly<sup>23</sup> and Lys<sup>29</sup>–Ala<sup>30</sup>) makes it an ideal peptide for use in this kind of application.<sup>8</sup>

Serine protease inhibitors that will inhibit trypsin include DFP (diisopropyl fluorophosphate), TLCK (*N* $\alpha$ -*p*-tosyl-L-lysine chloromethyl ketone), PMSF (phenylmethanesulfonyl fluoride), APMSF (4-amidinophenylmethanesulfonyl fluoride), AEBSEF (4-(2-aminoethyl)benzenesulfonyl fluoride), aprotinin, leupeptin,  $\alpha_2$ -macroglobulin,  $\alpha_1$ -antitrypsin, *p*-aminobenzamidine, benzamidine (reversible), soybean trypsin inhibitor, lima bean inhibitor, bovine pancreas trypsin inhibitor, chicken egg white inhibitor, and turkey egg white inhibitor.<sup>1,9</sup>

#### Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

#### Preparation Instructions

This enzyme is soluble in 1 mM HCl (1 mg/ml).

### Storage/Stability

Solutions in 1 mM HCl (pH 3) remain active for ~1 year when aliquoted and stored at  $-20^{\circ}\text{C}$ . The presence of  $\text{Ca}^{2+}$  (20 mM) will also retard the autolysis of trypsin and will maintain the stability of the trypsin in solution.<sup>1,11</sup>

Trypsin retains most of its activity in 2.0 M urea, 2.0 M guanidine HCl, or 0.1% (w/v) SDS.<sup>12</sup> Trypsin is reversibly denatured at high pH (above 11), by precipitation with TCA, or by high concentrations of urea (greater than 6.5 M).<sup>3</sup> In order to abolish all trypsin activity, heating at  $100^{\circ}\text{C}$  in 1% (w/v) SDS for 5 minutes is required.<sup>13</sup>

### Procedure

For trypsin digestion of proteins, use a ratio (w:w) of 1:100 to 1:20 for trypsin:protein. Trypsin preparations usually contain some contaminating chymotrypsin, and this should be inhibited with N-tosyl-L-phenylalanyl chloromethyl ketone (TPCK).<sup>10</sup>

### References

1. Walsh, K.A., Trypsinogens and trypsins of various species. *Meth. Enzymol.*, **19**, 41-63 (1970)
2. Buck, F.F., *et al.*, On the mechanism of enzyme action. LXXII. Studies on trypsins from beef, sheep, and pig pancreas. *Arch. Biochem. Biophys.*, **97**, 417-424 (1962).
3. Keil, B., in *The Enzymes*, 3rd ed., Vol. III, Boyer, P. D., Academic Press (New York, NY: 1971), pp. 250-275.
4. Hermodson, M.A., *et al.*, Determination of the Amino Acid Sequence of Porcine Trypsin by Sequenator Analysis. *Biochemistry*, **12(17)**, 3146-3153 (1973).
5. Burdon, R.H., *et al.*, in *Laboratory Techniques in Biochemistry and Molecular Biology*, Vol. 9, 2nd ed. (Allen, G., ed.), Elsevier/North (New York, NY: 1989), pp. 73-104.
6. *Enzyme Handbook, Vol. II*, Barman, T.E., Springer-Verlag (New York, NY: 1969), pp. 618-619.
7. USP, Vol. 23, p. 1611.
8. Wang, S.-S., and Carpenter, F.H., Kinetics of the tryptic hydrolysis of the oxidized B chain of bovine insulin. *Biochemistry*, **6(1)**, 215-224 (1967).
9. *Proteolytic Enzymes, A Practical Approach*, Beynon, R.J., ed., IRL Press (New York, NY: 1989), p. 240.
10. Carpenter, F.H., Treatment of trypsin with TPCK. *Methods Enzymol.*, **11**, 237 (1967).
11. Sipos, T., and Merkel, J.R., An effect of calcium ions on the activity, heat stability, and structure of trypsin. *Biochemistry*, **9(14)**, 2766-2775 (1970).
12. *Methods of Molecular Biology*, Vol. 3, Smith, B.J., Humana Press, (New Jersey, 1988), pp. 57-69.
13. Porter, W.H., and Preston, J.L., Retention of trypsin and chymotrypsin proteolytic activity in sodium dodecyl sulfate solutions. *Anal. Biochem.*, **66(1)**, 69-77 (1975).

TMG,GCY,RXR,MAM 05/16-1