

Data Sheet

# BioTracker™ Polyamine Spermine Live Cell Probe

Live Cell Dye

**SCT250****Pack Size: 500 µg****Store at -20 °C**FOR RESEARCH USE ONLY**Not for use in diagnostic procedures. Not for human or animal consumption.**

## Background

Polyamines play a variety of critical roles in cellular processes. Some of these include gene expression, cell proliferation and differentiation, and cellular stress responses. Polyamine dysregulation has been shown to be involved in various cancers. In colorectal cancer, for example, polyamine content and biosynthesis is significantly higher when compared to normal colorectal tissue. A relationship between polyamines and their metabolites has been used to observe progression in breast, lung, colorectal, ovarian, prostate, and pancreatic cancers. In normal cells, there is close regulation of polyamines via usual cellular mechanisms. One of these is known as the polyamine transport system (PTS), the characteristics of which are not fully understood. Studies of the PTS have shown that cancer cells have high PTS activity, suggesting this system as a candidate for cancer research.

The relationship between polyamine systems and cancer led to the creation of clickable polyamine derivatives which can functionally be used to study the PTS. The BioTracker™ spermine live cell dye is a spermine moiety attached to a BODIPY fluorophore through click chemistry. The BioTracker™ spermine live cell dye permits evaluation of polyamine.

## Source

SCT250 does not contain genetically modified organisms.

## Spectral Properties

Absorption: 494 nm

Emission: 500-510 nm

## Quality Control Testing

Purity: ≥ 80% confirmed by HNMR, LC-MS and HPLC and elemental analysis

Molar Mass: 743.35 g/mol

## Storage and Handling


Store BioTracker™ Spermine Live Cell Dye at -20 °C, desiccate and protect from light.

**Note:** Centrifuge vial briefly to collect contents at bottom of vial before opening.


## Presentation

Lyophilized. Red-brown solid.


## Representative Data



**Figure 1:** (A) Confocal imaging of SCT250 spermine probe (green) in MCF7 cells and (B) MCF10A cells. Cells were treated with 10  $\mu$ M spermine probe for 3 hours at 37  $^{\circ}$ C and counterstained with Hoechst 33342 nuclear dye (cyan) for 10 minutes. The highly metastatic MCF7 breast cancer cells are characterized by high polyamine uptake, while MCF10A cells demonstrate low polyamine uptake.



**Figure 2:** Probe excitation (blue) and emission (orange) spectral data. 7  $\mu$ L of probe at stock concentration (10 mM) was diluted in 1 mL of Tris HCl buffer (10  $\mu$ M, pH 7.0). Spectral scans were conducted using a PerkinElmer FL8500 Fluorescence Spectrophotometer.



**Figure 3:** Probe absorbance data. 7  $\mu$ L of probe at stock concentration (10 mM) was diluted in 1 mL of Tris HCl buffer (10  $\mu$ M, pH 7.0). Spectral scans were conducted using a PerkinElmer FL8500 Fluorescence Spectrophotometer.

## Protocols

### Reagent Preparation

1. Before opening the vial, spin down the solid to the bottom by a microcentrifuge or by a desktop centrifuge.
2. Warm the vial to the room temperature and add DMSO to make a 1000X stock solution of 10 mM (freeze aliquots at  $-20^{\circ}\text{C}$ ).
3. Dilute in cell culture media at a final concentration of 10  $\mu$ M and add to cells in culture. Incubate at  $37^{\circ}\text{C}$  for 3-4 hours.
4. If desired, counterstain with nuclear stain (for example, Hoechst 33342 at 0.1 mg/mL) for 10 minutes.
5. Wash cells with PBS buffer before imaging.

**Note:** Optimal concentration must be determined by end user.

## References

1. Vanhoutte, R., Kahler, J. P., Martin, S., van Veen, S., & Verhelst, S. H. (2018). Clickable polyamine derivatives as chemical probes for the polyamine transport system. *ChemBioChem*, 19(9), 907-911.
2. Li, J., Meng, Y., Wu, X., & Sun, Y. (2020). Polyamines and related signaling pathways in cancer. *Cancer cell international*, 20, 1-16.

---

## Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

## Technical Assistance

Visit the tech service page at [SigmaAldrich.com/techservice](https://SigmaAldrich.com/techservice).

## Terms and Conditions of Sale

Warranty, use restrictions, and other conditions of sale may be found at [SigmaAldrich.com/terms](https://SigmaAldrich.com/terms).

## Contact Information

For the location of the office nearest you, go to [SigmaAldrich.com/offices](https://SigmaAldrich.com/offices).

The life science business of Merck KGaA, Darmstadt, Germany  
operates as MilliporeSigma in the U.S. and Canada.

MilliporeSigma, BioTracker and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2024 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

Document Template 20306518 Ver 6.0

23081502 Ver 1.0, Rev 08JUL2024, LB

**Millipore  
Sigma**