

For life science research only.  
Not for use in diagnostic procedures.



# Endoproteinase Asp-N Sequencing Grade from *Pseudomonas fragi* mutant

 **Version: 18**  
Content Version: June 2021

Lyophilized

**Cat. No. 11 420 488 001**    2 µg  
**Cat. No. 11 054 589 001**    3 x 2 µg

**Store the product at +2 to +8°C.**

<b>1.</b>	<b>General Information .....</b>	<b>3</b>
1.1.	Contents .....	3
1.2.	Storage and Stability .....	3
	Storage Conditions (Product) .....	3
1.3.	Additional Equipment and Reagent required .....	3
	For preparation of digestion buffer .....	3
	For solubilization of proteins .....	3
1.4.	Application .....	3
<b>2.</b>	<b>How to Use this Product .....</b>	<b>4</b>
2.1.	Before you Begin .....	4
	General Considerations .....	4
	General handling recommendations .....	4
	Activity determination .....	4
	Working Solution .....	4
2.2.	Protocols .....	5
	Digestion of proteins in solution .....	5
2.3.	Parameters .....	5
	Molecular Weight .....	5
<b>3.</b>	<b>Results .....</b>	<b>6</b>
	Verification of specificity and nonspecificity of Endoproteinase Asp-N .....	6
<b>4.</b>	<b>Additional Information on this Product .....</b>	<b>8</b>
4.1.	Quality Control .....	8
<b>5.</b>	<b>Supplementary Information .....</b>	<b>9</b>
5.1.	Conventions .....	9
5.2.	Changes to previous version .....	9
5.3.	Ordering Information .....	9
5.4.	Trademarks .....	10
5.5.	License Disclaimer .....	10
5.6.	Regulatory Disclaimer .....	10
5.7.	Safety Data Sheet .....	10
5.8.	Contact and Support .....	10

# 1. General Information


## 1.1. Contents

Vial / Bottle	Label	Function / Description	Catalog Number	Content
1	Endoproteinase Asp-N Sequencing Grade	Highly purified and specific protease.	11 420 488 001	1 vial, 2 µg
			11 054 589 001	3 vials, 2 µg each

## 1.2. Storage and Stability

### Storage Conditions (Product)

When stored at +2 to +8°C, the product is stable through the expiry date printed on the label.

Vial / Bottle	Label	Storage
1	Endoproteinase Asp-N Sequencing Grade	Store at +2 to +8°C.  <b>Store dry.</b>

## 1.3. Additional Equipment and Reagent required

### For preparation of digestion buffer

 See section, **Working Solution** for additional information on preparing solutions.

- Sodium phosphate buffer

### For solubilization of proteins

- Sodium dodecyl sulfate (SDS\*)
- Urea
- Methylamine
- Guanidine hydrochloride
- Acetonitrile

## 1.4. Application

Use Endoproteinase Asp-N for the specific cleavage of proteins and peptides for:

- Protein structure
- Sequence analysis

## 2. How to Use this Product

### 2.1. Before you Begin

#### General Considerations

##### General handling recommendations

The content of one vial may be used for several simultaneous digests.

**⚠ Take a new vial when repeating a digest in order to minimize the risk of contamination or autolysis.**

##### Activity determination

Activity determination of Endoproteinase Asp-N, with azocoll as substrate in the presence of stated concentrations of denaturing agents. Incubation of Endoproteinase Asp-N 200 µg/ml with denaturing agent for 6 hours at +25°C in 25 mM sodium phosphate buffer, pH 7.8.

**i** Add 20 mM methylamine when applying urea.

Denaturing agent	Concentration	Enzyme activity [%]
without addition (control)	–	100
SDS	0.001% (w/v)	113
	0.01% (w/v)	122
	0.1% (w/v)	10
Urea (+ methylamine)	0.1 M	100
	0.5 M	108
	1 M	105
Guanidine hydrochloride	0.1 M	100
	0.5 M	85
	1 M	80
Acetonitrile	1% (v/v)	90
	5% (v/v)	115
	10% (v/v)	125

#### Working Solution

Solution	Preparation/Composition	Storage and Stability	For use in...
Endoproteinase Asp-N Sequencing Grade	<ul style="list-style-type: none"> <li>Add 50 µl double-distilled water to the lyophilizate to a final concentration of 10 mM Tris-HCl, pH 7.5.</li> </ul> <p><b>⚠ To avoid autolysis, the incubation temperature must not exceed +37°C.</b></p>	Store 1 week at +2 to +8°C.	Digestion mixture
Digestion buffer	50 mM sodium phosphate buffer, pH 8.0.	–	Dissolution of the proteins to be sequenced.

## 2.2. Protocols

### Digestion of proteins in solution

- i* See section, **Working Solution** for information on preparing solutions.
- 1 Dissolve the proteins to be sequenced in Digestion buffer.
    - i* For proteins that are hard to solubilize, add urea, SDS, or guanidine hydrochloride to the Digestion buffer prior to solubilizing the protein. When applying urea, also add 20 mM methylamine.
- 
- 2 Dilute protein solution with buffer, see section, **General Considerations** to achieve a suitable concentration of the denaturing agent in the digest.
    - ⚠ Avoid complexing agents, such as EDTA or 2-phenanthroline in the samples, since Endoproteinase Asp-N is a metalloprotease.**
    - i* The recommended amount of enzyme is 1/200 to 1/20 of the protein by weight.
- 
- 3 Choose an incubation time between 2 and 18 hours at +37°C, depending on the amount of enzyme.
    - i* Under these conditions, additional cleavage at glutamyl residues can be observed. The cleavage velocity of Endoproteinase Asp-N revealed that the aspartyl specific cleavage is at least 2,000-fold faster than the glutamyl side activity of the Endoproteinase Asp-N.
    - ⚠ Reduce enzyme concentration (enzyme-substrate ratio of 1:1,000 [w/w]) at an incubation time of 2 to 6 hours to prevent the additional cleavage at glutamyl residues.**
- 

## 2.3. Parameters

### Molecular Weight


27 kDa

## 3. Results

### Verification of specificity and nonspecificity of Endoproteinase Asp-N


Endoproteinase Asp-N is a metalloprotease that specifically cleaves peptide bonds N-terminally at aspartic and cysteic acids in phosphate, acetate, or Tris buffers at pH 6.0 to 8.5.

The specificity and nonspecificity of Endoproteinase Asp-N is verified with glucagon or melittin as substrate.



**Fig 1:** Specificity of Endoproteinase Asp-N in reversed phase HPLC. High concentrations of Endoproteinase Asp-N (1 part by weight enzyme with 100 parts by weight glucagon) are incubated for 1 hour to detect the fragments of the specific digested substrate.

Digest	40 µg glucagon in 200 µl 45 mM NaH <sub>2</sub> PO <sub>4</sub> buffer, 0.09% TFA at pH 8.0 + 0.4 µg Endoproteinase Asp-N dissolved in 10 µl double-distilled water; 1 hour at +37°C; reversed phase HPLC: undiluted.
Column	Nucleosil 100-5-C18 4 × 100 mm, 5 µm
Solvent A	0.1% TFA (v/v) in double-distilled water
Solvent B	0.1% TFA (v/v) in double-distilled water; 70% acetonitrile (v/v)
Gradient	40 minutes linearly 0 to 100% B
Flow rate	1 ml/minute
Wavelength	215 nm
Fragments	4.49 min Asp (15) – Gln (20) 8.89 min His (1) – Ser (8) 13.78 min Asp (9) – Leu (14) 22.73 min Asp (21) – Thr (29)



**Fig 2:** Nonspecificity of Endoproteinase Asp-N in reversed phase HPLC. High concentrations of Endoproteinase Asp-N (1 part by weight enzyme with 20 parts by weight melittin) are incubated for 4 hours to detect traces of impurities.

Digest	40 µg melittin in 200 µl 50 mM NaH <sub>2</sub> PO <sub>4</sub> buffer at pH 8.0 + 2 µg Endoproteinase Asp-N dissolved in 50 µl double-distilled water; 4 hours at +37°C; reversed phase HPLC: undiluted.
Column	Nucleosil 100-5-C18 4 × 100 mm, 5 µm
Solvent A	0.1% TFA (v/v) in double-distilled water
Solvent B	0.1% TFA (v/v) in double-distilled water; 70% acetonitrile (v/v)
Gradient	40 minutes linearly 0 to 100% B
Flow rate	1 ml/minute
Wavelength	215 nm

## 4. Additional Information on this Product

### 4.1. Quality Control

For lot-specific certificates of analysis, see section, **Contact and Support**.



## 5. Supplementary Information

### 5.1. Conventions

To make information consistent and easier to read, the following text conventions and symbols are used in this document to highlight important information:

#### Text convention and symbols

 *Information Note: Additional information about the current topic or procedure.*

 **Important Note: Information critical to the success of the current procedure or use of the product.**

① ② ③ etc. Stages in a process that usually occur in the order listed.

① ② ③ etc. Steps in a procedure that must be performed in the order listed.

\* (Asterisk) The Asterisk denotes a product available from Roche Diagnostics.

### 5.2. Changes to previous version

Layout changes.

Editorial changes.

### 5.3. Ordering Information

Product	Pack Size	Cat. No.
Reagents, kits		
Sodium Dodecyl Sulfate (SDS)	1 kg	11 667 289 001

## 5. Supplementary Information

### 5.4. Trademarks

All product names and trademarks are the property of their respective owners.

### 5.5. License Disclaimer

For patent license limitations for individual products please refer to:

**List of biochemical reagent products.**

### 5.6. Regulatory Disclaimer

For life science research only. Not for use in diagnostic procedures.

### 5.7. Safety Data Sheet

Please follow the instructions in the Safety Data Sheet (SDS).

### 5.8. Contact and Support

To ask questions, solve problems, suggest enhancements or report new applications, please visit our **Online Technical Support Site.**

To call, write, fax, or email us, visit **sigma-aldrich.com**, and select your home country. Country-specific contact information will be displayed.

