

24-Hour Emergency Number (ChemTel) 888.255.3924 in U.S. cleaning@alconox.com

alconox.com

- distributors
- technical information
- free samples
- new developments
- white papers
- prices



DETONOX®

Ultimate Precision Cleaner

- Concentrated to save you money
- Exceptional cosmetic and pharmaceutical polymer, cream and lotion removal
- Non-caustic for greater worker safety and greater material compatibility
- Free rinsing to provide reliable results and no interfering residues

Used to remove: Cosmetic bases and formulations, organics, grease, tars, resins, waxes, fats, oils, titanium and zinc oxides, ointments and lotions.

Surfaces cleaned: Recommended for glass, stainless steel, PTFE, ceramic, cement and solvent resistant plastics and metals. Corrosion testing is advisable.

Cleaning method: Soak, circulate under immersions, ultrasonic cleaners, and manual cleaning.

Directions: Make a fresh 1–3% solution (1 1/4–3 3/4 oz per gal, or 10–30 ml per L) in cold, warm or hot water. If available, use very hot water to clean cosmetic residues. For difficult soils, raise water temperature and use more detergent. Elevation of temperature is true under most circumstances. Check with Alconox, Inc. for use of hot water if products contain titanium oxide (TiO) or zinc oxide (ZnO) as they are notable exceptions to heat-expedited detergency. Clean by soak, circulate, wipe or ultrasonic method. Not for spray machines: will foam. RINSE THOROUGHLY—preferably with running water. For critical cleaning, do final or all rinsing in purified water.

Caution: Contains anionic and nonionic surfactants with phosphate chelating agents. May cause skin and eye irritation. Wear eye protection and chemical resistant gloves. Material compatibility guidelines available on request.

					0		14
Convenient Sizes	Alco	nox Cat. #		-		U	1
1 gal bottle (3.8 L)		2301-1					
Case of 4 x 1 gal (4 x 3.8L)		2301	ET	-		ELE	
5 gal (19 L) drum		2305					
55 gal (208L) drum		2355	-				

1 gal makes 100 gal (1 L makes 100L) solution

DETONOX is available from leading laboratory, hospital, clinical and industrial suppliers. To find a distributor for Alconox, Inc. detergents, use the Locate Dealers search tool in the Support section of alconox. com. To request FREE samples, use the Get Samples form at alconox.com. Or call Alconox, Inc. at +914.948.4040, today.

DXTB.1



24-Hour Emergency Number (ChemTel) 888.255.3924 in U.S. cleaning@alconox.com

alconox.com

- distributors
- technical information
- free samples
- new developments
- white papers
- prices

Physical Data		Typical Value		
pH (1%)		10		
Surface Tension 1	l% Sol'n (dyne/cm)	32		
Conductivity 1%	Sol'n (mS/cm at 21 degrees (0.7		
Form		Liquid		
Color		Pale yellow		
Solubility in Wate	c Completely soluble	e in all proportions		
Hard Water Effect	riveness	Highly effective		
Biodegradability	Yes, m	eets EC648/2004		
Foam Tendency		High foaming		
Shelf Life	24 months from end of mon	th of manufacture		
Fragrance and Dy	re Content	0%		
Phosphate Conte	nt (as phosphorus)	0.9% as P		
Flash Point (degree	es F)	None		
Organic Carbon (% calculated w/w)		19%		
Density (g/mL − lb/	'gal)	9.0		

Chemical Description

DETONOX contains anionic and nonionic detergent with wetting agents and emulsifiers with phosphate chelating and sequestering agents.

Cleaning Validation Methods

To quantify trace levels of Citranox for cleaning validation purposes, use conductivity, UV, TOC, HPLC or phosphate analysis. For details see the Cleaning Validation page in Support at alconox.com. For lab accreditation tests to show the detergent is rinsed away, see the Accreditation Residue Detection in the Support section of alconox. com.

Health Safety Information

OSHA Hazardous Ingredients	None	
RCRA Hazard Class	Nonhazardous	
EPA Priority Pollutants	None	
DOT Hazard Class	Not regulated	
Flammability	Nonflammable	
Latex Content	None in product or packaging.	
Irritation	Cat. 2 skin, Cat. 2A eye irritant	
Oral Toxicity	$LD_{50} > 5000$ mg/kg oral rat.	
Inhalation Toxicity	Nonirritating	
VOC Content	0%	
GMO Content	None	
BSE/TSE Source	None	
2007/68/EC Annex IIIa Allergens		
REACH Art. 57 (SVHC)	None	
RoHS	No RoHS hazardous ingredients.	
Carcinogenicity	NTP = No IARC = No OSHA = No	
Indgredients on TSCA Inventory		
	ains no enzymes, optical brighteners, c fragrances, or preservative agents.	

Precautions

Avoid contact with skin, eyes and clothing. Do not take internally. In case of contact immediately flush skin with plenty of water. For eyes, flush with plenty of water for at least 15 minutes and get medical attention immediately. A Safety Data Sheet is available at alconox.com.

Contact Alconox, Inc. for purchase specifications. Information presented is typical. Not to be taken as specifications. Typical data is not a specification.

While the information in this report should not be considered to be a product warranty, we urge you to investigate, test and verify the suitability of Alconox, Inc. detergents for your specific application. We, of course, can not give permission to use, or recommend the use of, our detergents where they infringe patents. No representation or warranty is made as to the safety of products or materials mentioned under the Federal Food Additives Amendment of 1958.

DXTB.1