

Technical Data Sheet

IsoBag™ TSA + LTHTh Settle_{plus}

Ordering number: 1.46755.0080

The IsoBag™ TSA + LTHTh Settle_{plus} is designed for quick and convenient transfer of culture media for air monitoring into isolators or RABS systems containing Getinge La Calhène DPTE® Alpha Ports of 190 mm diameters. The plate can be transferred into the isolator without the need for additional decontamination of the outside packaging of plates. The product allows increased production time by decreased requirements for decontamination cycles, e.g. in case of campaign productions. Furthermore this procedure will save space for storage of culture media in the isolator.

The IsoBag™ TSA + LTHTh Settle_{plus} is a sterile polyurethane DPTE® Beta Bag (DPTE is a trademark by Getinge La Calhène) compatible with 190 mm DPTE® Alpha Ports and filled with:

- 8 packs of 10 single bagged, lockable 90 mm settle plates (article 146683 – TSA + LTHTh ICR_{plus}) and
- 4 x 5 single bagged, pre-sterilized zip bags for safe transport of used plates into non-controlled areas for incubation and colony counting.

The filled IsoBag™ is double bagged and irradiated at a dose of 10-20 kGy. Further IsoBag designs are available with the identical media formulation:

- IsoBag™ TSA + LTHTh Settle (article number 1467560080): 190 PU DPTE® Beta Bag filled with 8 x 10 single bagged settle plates (article number 146069) and 4 x 5 single bagged transparent, zip-lock transport bags. The product is intended for microbial monitoring of air (passive and active) in isolators and RABS. The plate design allows aerobic incubation only.
- IsoBag™ TSA + LTHTh Contact (article number 1467540100): 190 PU DPTE® Beta Bag filled with 10 x 10 single bagged contact plates (article number 146231) and 4 x 5 single bagged transparent, zip-lock transport bags. The product is intended for microbial monitoring of surfaces in isolators and RABS. The plate design allows aerobic incubation only.



Typical Composition of included lockable settle plates

Casein Peptone	15 g/l
Soy Peptone	5 g/l
NaCl	5 g/l
Polysorbate (Tween®) 80	5 ml/l
Lecithin	0.7 g/l
Histidine	0.5 g/l
Sodium Thiosulfate	0.5 g/l
Agar	15 g/l

The appearance of the medium is clear and yellowish. The pH value is in the range of 7.1-7.5. The medium can be adjusted and/or supplemented according to the performance criteria required. (Tween is a trademark of ICI Americas Inc.)

Application and Interpretation

The plates are introduced into isolators or RABS by connecting the IsoBag™ to a 190 mm DPTE® Alpha Port. Therefore please follow the instructions for use of Getinge La Calhène for the DPTE bidirectional transfer system. The integrity of the bag remains stable up to 14 connections to the Alpha Port.

Single packs of plates may be transferred into the isolator by the possibility of multiple connections. Up to 14 connections of the IsoBag™ with the Alpha Port are possible. Please do not open the IsoBag™ outside of the Isolator or RABS, if further connections to the Alpha Port will be performed.

There is no need to decontaminate the plates, which were transferred into the isolator. As long as the packaging of plate packs is not damaged a VHP decontamination is allowed, however.

For the use of the included settle plates for active or passive air monitoring within the isolator please refer to the Technical Data Sheet of Tryptic Soy Agar + LTHTh – ICR_{plus} in 90 mm, lockable settle plates (article 146683).

For transfer of used plates to the incubator 4 settle plates may be packed into the zip bags included in the IsoBag™. The outtake can be performed by a new DPTE® Beta Bag from Getinge La Calhène or by using an empty IsoBag, if the maximum connection number of 14 has not been applied before.

Storage and Shelf Life

The product can be used for sampling until the expiry date if stored upright, protected from light and properly sealed at +15°C to +25°C. The testing procedures as described on the CoA can be started up to the expiry date printed on the label.

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose. The M mark is a registered trademark of Merck KGaA, Darmstadt, Germany.

Lit. No. TN6481FN00

Condensation can be prevented by avoiding quick temperature shifts and mechanical stress.

Disposal

Please mind the respective regulations for the disposal of used culture medium (e.g. autoclave for 20 min at 121 °C, disinfect, incinerate etc.).

Quality Control

Please refer to the actual batch related Certificate of Analysis as well as the batch related Certificate of Analysis of the included batch of Tryptic Soy Agar + LTHTh – ICR_{plus} in lockable, 90 mm settle plates (article 146683 download from webpage using either 1466830020 or 1466830120 combined with the concerned batch number of included plates).

Merck KGaA
64271 Darmstadt, Germany
Fax: +49 (0) 61 51 / 72-60 80
mibio@merckgroup.com

Find contact information for your
country at:
www.merckmillipore.com/offices

For Technical Service, please visit:
www.merckmillipore.com/techservice