

71725 / 71729 Sodium dodecyl sulfate

Synonyms: Dodecyl sulfate sodium salt

SDS

Sodium lauryl sulfate Lauryl sulfate sodium salt Dodecyl sodium sulfate CH₃(CH₂)₁₀CH₂O-S-ONa

CAS number: 151-21-3

Product Description:

Molecular Formula: $C_{12}H_{25}NaO_4S$ Molecular Weight: 288.4 g/mol Melting Point: 204-207 °C

The BioUltra grade SDS has trace element analysis performed. Trace element analysis is often important for experiments where tight control of elemental contamination is required (i.e. NMR, ICP atomic emission and purification steps). The molecular biology quality is additionally tested on the absence of DNases, RNases, Phosphatases and Proteases. For common biochemistry work the BioUltra quality is recommended.

SDS is an anionic detergent and wetting agent that is effective in both acid and alkaline solutions. SDS has a wide variety of applications, but is most often used as a protein and lipid solubilization reagent. As a general rule for the solubilization of proteins, SDS should be used at its critical micelle concentration. ²

SDS is also a powerful protein denaturant. The effects of SDS on protein conformation has been published.^{3,4}

Comparisons between SDS and other detergents for solubilization of lipids, proteins, and its effect on enzymes activity has been published.^{2,5}

To remove SDS from protein samples, it is recommended to use Product No. 428698, Dowex[®] Retardation 11A8 Ion-Exchange Resin. This resin contains paired anion and cation exchange sites. Methods for SDS removal by ion exchange chromatography have been published.⁶ Methylene blue can be used to determine the amount of SDS remaining following removal of SDS by ion exchange chromatography.⁷

Preparation Instructions:

This product can be dissolved in water (200 mg/ml), yielding a clear, colorless solution.

Storage/Stability:

SDS undergoes hydrolysis at elevated temperatures especially in acidic medium. Prolonged heating at 40°C or greater causes decomposition of alkyl sulfates into fatty alcohols and sodium sulfate.⁸



References:

- 1. Martindale The Extra Pharmacopoeia, 29th ed., Reynolds, J. E. F., ed., The Pharmaceutical Press (London, England: 1989), p. 1417.
- 2. Womack, M. D., et al., Detergent Effects on Enzyme Activity and Solubilization of Lipid Bilayer Membranes. Biochim. Biophys. Acta, 733(2), 210-215 (1983).
- 3. Protein Structure: A Practical Approach, 2nd ed., Creighton, T. E., ed., IRL Press at Oxford University Press (New York, NY: 1997), pp. 1-19.
- 4. Protein Structure: A Practical Approach, 2nd ed., Creighton, T. E., ed., IRL Press at Oxford University Press (New York, NY: 1997), pp. 198-199.
- 5. Weber, K., and Kuter, D. J., Reversible Denaturation of Enzymes by Sodium Dodecyl Sulfate. J. Biol. Chem., 246(14), 4504-4509 (1971).
- 6. Kapp, O. H., and Vinogradov, S. N., Removal of Sodium Dodecyl Sulfate from Proteins. Anal. Biochem., 91(1), 230-235 (1978).
- 7. Hayashi, K., A Rapid Determination of Sodium Dodecyl Sulfate with Methylene Blue. Anal. Biochem., 67(2), 503-506 (1975).
- 8. Chem. Abstract, vol. 72:22811d (1970). Dowex is a registered trademark of the Dow Chemical Co.
- 9. CRC Handbook of Electrophoresis Boca Raton, Florida 1, (1980)
- 10. R.A. Grohs, Analysis of drugs in the presence of serum albumin by LC with eluants containing surfactants Anal. Chem. 63, 384 (1991)
- 11. U.K. Laemmli, PAGE a SDS-discontinuous system Nature 227, 680 (1970)
- 12. T.B. Nelsen, J.A. Reynolds, Measurement of molecular weight by gel electrophoresis Meth. Enzymol. 48, 3 (1978)
- 13. K. Harewood, J.S. Wolff, Anionic detergent for protein solubilization Anal. Biochem. 55, 573 (1973)
- 14. eds. L.A. Lewis, J.J. Oppelt, CRC Handbook of Electrophoresis Boca Raton, Florida: 1, (1980)
- 15. Sigma-Aldrich quality control

Dowex is a registered trademark of the Dow Chemical Co.

Precautions and Disclaimer:

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

The vibrant M and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. Detailed information on trademarks is available via publicly accessible resources.

© 2018 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

