

3050 Spruce Street
Saint Louis, Missouri 63103 USA
Telephone 800-325-5832 • (314) 771-5765
Fax (314) 286-7828
email: techserv@sial.com
sigma-aldrich.com

ProductInformation

N,O-Bis(trimethylsilyl)trifluoroacetamide with 1% trimethylchlorosilane

Product Number **T 6381** Storage Temperature 2-8 °C

Product Description

Molecular Formula (BSFTA): C₈H₁₈F₃NOSi₂ Molecular Weight (BSTFA): 257.4

CAS Number (BSTFA): 25661-30-2

Density: 0.97 g/ml Synonym: BSFTA

The position of the silyl groups has not been unequivocally established. BSTFA has been referred to as both N,N and N,O.

Bis(trimethylsilyl)trifluoroacetamide BSTFA is a powerful silylating reagent for the derivatization of a wide range of compounds. BSTFA replaces labile hydrogens on compounds with trimethylsilyl groups, to form volatile and thermally stable derivatives for GC and GC-MS analysis. BSTFA is similar in strength to N,O-bis(trimethylsilyl)acetamide (BSA) and possesses advantages over BSA:

 BSFTA and its by products are more volatile, resulting in less chromatographic interference.
 The presence of fluorine reduces detector (FID) contamination.

The silyl donor properties of silylating reagents are distinctly influenced by the solvent systems and any added catalysts. Thus, trimethylchlorosilane (TMCS) is often added to the BSTFA to increase the silyl donor strength of the BSTFA.² A comparison of the relative efficiency of derivatization of carbonyl compounds by BSTFA and BSTFA/TMCS has been reported.³

BSTFA is a liquid at room temperature. It is very miscible with most commonly used silylation solvents, and BSTFA itself has good solvent properties. Thus, BSTFA can be used as a silylation reagent, without additional solvents. In general, BSTFA can effectively derivatize compounds at room temperature.

BSTFA is readily susceptible to hydrolysis, so it is strongly advised to exclude water from the samples before derivatization. Because of the reactivity of BSTFA with labile hydrogens, solvents which contain labile hydrogen atoms should not be used in the derivatization process.

BSTFA has been applied to the analysis of compounds from plants,⁴ human biological fluids,⁵⁻⁷ and the ambient atmosphere.⁸

Precautions and Disclaimer

For Laboratory Use Only. Not for drug, household or other uses.

Preparation Instructions

BSFTA is miscible in hexane (100 mg/ml, w/v), yielding a clear to slightly hazy, colorless solution.

References

- Silylating Agents, 2nd ed., Heberle, J., and Simchen, G., Fluka Chemie AG (Buchs. 1995).
- 2. Knapp, D., Handbook of Analytical Derivatization Reactions, 1st ed., Knapp, D., John Wiley & Sons (New York, NY: 1979), p. 9.
- Spaulding, R.S., and Charles, M.J., Comparison of methods for extraction, storage, and silylation of pentafluorobenzyl derivatives of carbonyl compounds and multi-functional carbonyl compounds. Anal. Bioanal. Chem., 372(7-8), 808-16 (2002).
- Greenway, W., and Whatley, F.R., Resolution of complex mixtures of phenolics in poplar bud exudate by analysis of gas chromatography-mass spectrometry data. J. Chromatogr. A, 519(1), 145-158 (1990).

- Hammargren, W.M., et al., Identification of a novel nucleoside, 1,N⁶-dimethyladenosine, in human cancer urine. Anal. Chim. Acta, 247(2), 201-209 (1991).
- 6. Schaefer, W.R., et al., Cytochrome P450 metabolites of arachidonic acid in human placenta. Prostaglandins, **54(4)**, 677-687 (1997).
- Soleas, G.J., et al., Ultrasensitive assay for three polyphenols (catechin, quercetin and resveratrol) and their conjugates in biological fluids utilizing gas chromatography with mass selective detection. J. Chromatogr. B Biomed. Sci. Appl., 757(1), 161-172 (2001).
- 8. Destaillats, H., et al., Ambient air measurement of acrolein and other carbonyls at the Oakland-San Francisco Bay Bridge toll plaza. Environ. Sci. Technol., **36(10)**, 2227-2235 (2002).

GCY/RXR 2/03