

Technical Bulletin

Cytotoxicity Assay Kit

Catalogue number MAK536


Product Description

Adenosine 5'-triphosphate (ATP) is the chemical energy for cellular metabolism and is often referred to as "energy currency" of the cell. ATP is produced only in living cells during photosynthesis and cellular respiration and consumed in cellular processes including biosynthetic reactions, motility and cell division. It is a key indicator of cellular activity and has been utilized as a measure of cell viability and cytotoxicity in research and drug discovery.

The Cytotoxicity Assay Kit provides a rapid method to measure intracellular ATP, cell viability and cytotoxicity.

The single working reagent lyses cells to release ATP, which, in the presence of luciferase, immediately reacts with the Substrate D-luciferin to produce light. The light intensity is a direct measure of intracellular ATP concentration and hence number of living cells.

This non-radioactive, homogeneous cell-based assay can be conveniently performed in microplates. The reagent is compatible with all culture media and liquid handling systems for high-throughput screening applications in 96-well and 384-well plates. This kit can detect as low as 50 cells in a given sample. The kit is used to study effects of cytokines, growth factor, nutrients on cells.



Components

The kit is sufficient for 100 bioluminescent assays in 96-well plates.

- | | |
|--------------------------|--------|
| • Assay Buffer | 10 mL |
| Catalogue Number MAK536A | |
| • Substrate | 120 µL |
| Catalogue Number MAK536B | |
| • ATP Enzyme | 120 µL |
| Catalogue Number MAK536C | |

Reagents and Equipment Required but Not Provided

- Pipetting devices and accessories (example., multichannel pipettor)
- Luminometer.
- White opaque tissue culture plates.
- 1.5 mL microcentrifuge tubes
- PBS (Catalogue number P4474 or equivalent)
- Cell Culture Medium.

Precautions and Disclaimer

For Research Use Only. Not for use in diagnostic procedures. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

Storage/Stability

The kit is shipped on wet ice. Store components at -20 °C.

Preparation Instructions

Briefly centrifuge small vials prior to opening. Equilibrate all components to room temperature prior to use.

Keep thawed ATP Enzyme on ice or 4 °C.

Note: After adding the Reconstituted Reagent, the luminescence signal is stable for about 15 min and decreases slow thereafter. Reading is best performed within 30 min.

Procedure

All Samples and Standards should be run in duplicate.

Assay Procedure In 96-Well Plates

Cell Culture

1. Plate cells at 100 μL /well in white opaque tissue culture plates.
2. If desired, add 5 μL test compounds and controls dissolved in PBS or culture medium per well.
3. Rock plate lightly to mix and incubate for desired period of time (example, overnight).

Assay Reaction

Prepare Working Reagent as per Table 1. Each well requires 90 μL of Working Reagent. Prepare enough for the assay.

Table 1.
96-well Working Reagent Preparation

Reagent	Volume (μL)
Assay Buffer	95
Substrate	1
ATP Enzyme	1

1. Add 90 μL Reconstituted Reagent to each test well.
2. Tap the plate to mix well.
3. Incubate for 2 minutes at room temperature.
4. Read luminescence on a luminometer.

Note: For most luminometers (Berthold Luminometer, LJJ Analyst, Top Count, MicroBeta Counters, CLIPR and LeadSeeker), integration time of 0.1 to 5 sec is appropriate.

Assay Procedure In 384-Well Plates

Cell Culture

1. Plate cells at 25 μL /well in white opaque tissue culture plates.
2. If desired, add 5 μL test compounds and controls dissolved in PBS or culture medium per well.
3. Rock plate lightly to mix and incubate for desired period of time (example, overnight).

Assay Reaction

Prepare Working Reagent as per Table 2. Each well requires 25 μL of Working Reagent. Prepare enough for the assay.


Table 2.
384-well Working Reagent Preparation

Reagent	Volume (μL)
Assay Buffer	30
Substrate	0.3
ATP Enzyme	0.3

1. Add 25 μL Reconstituted Reagent to each test well.
2. Tap the plate to mix well.
3. Incubate for 2 minutes at room temperature.
4. Read luminescence on a luminometer.

Note: For most luminometers (Berthold Luminometer, LJJ Analyst, Top Count, MicroBeta Counters, CLIPR and LeadSeeker), integration time of 0.1 to 5 sec is appropriate.

Figure 1.
Typical Linearity of Luminescence in 96 well plate



References

1. Li W., *et al.*, Human primary renal cells as a model for toxicity assessment of chemo-therapeutic drugs. *Toxicol In Vitro.*, **20(5)**, 669-76 (2006).
2. Zhelev Z., *et al.*, Phenothiazines suppress proliferation and induce apoptosis in cultured leukemic cells without any influence on the viability of normal lymphocytes. Phenothiazines and leukemia. *Cancer Chemother Pharmacol.*, **53(3)**, 267-75 (2004).
3. Ingram PR., *et al.*, A comparison of the effects of ocular preservatives on mammalian and microbial ATP and glutathione levels. *Free Radic Res.*, **38(7)**, 739-50 (2004).

Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

Technical Assistance

Visit the tech service page at SigmaAldrich.com/techservice.

Terms and Conditions of Sale

Warranty, use restrictions, and other conditions of sale may be found at SigmaAldrich.com/terms.

Contact Information

For the location of the office nearest you, go to SigmaAldrich.com/offices.

The life science business of Merck operates
as MilliporeSigma in the U.S. and Canada.

Merck and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates.
All other trademarks are the property of their respective owners. Detailed information on
trademarks is available via publicly accessible resources.

© 2021 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

mak536 Rev 09/23

