



3050 Spruce Street
Saint Louis, Missouri 63103 USA
Telephone 800-325-5832 • (314) 771-5765
Fax (314) 286-7828
email: techserv@sial.com
sigma-aldrich.com

Product Information

Pyrophosphatase, inorganic from baker's yeast (*S. cerevisiae*)


Product Number **I 1643**
Storage Temperature -0 °C

Product Description

Enzyme Commission (EC) Number: 3.6.1.1
CAS Number: 9024-82-2
Molecular Weight: 71 kDa¹
Extinction Coefficient: $E^{1\%} = 14.5$ (280 nm)²
pI: 4.75²
Synonyms: Inorganic Pyrophosphatase, PPI,
Pyrophosphate phosphohydrolase

Inorganic pyrophosphatase from baker's yeast is a homodimer consisting of two equal subunits of molecular weight 32-35 kDa.^{2,3}

Inorganic pyrophosphatase catalyzes the following reaction:



This ubiquitous enzyme serves to drive metabolic reactions that produce pyrophosphate, since these reactions typically have equilibrium constants near unity. The catalytic mechanism has been described in the literature.³ Inorganic pyrophosphatase is a metalloprotease that requires Mg^{2+} for maximal activity. Although the hydrolysis of inorganic pyrophosphate is specific in the presence of Mg^{2+} , both ADP and ATP can be hydrolyzed if zinc is present. The following metals can act as activators: $\text{Mg}^{2+} > \text{Zn}^{2+} > \text{Co}^{2+} > \text{Mn}^{2+} > \text{Ca}^{2+}$.^{2,3,4}

Inorganic pyrophosphatase from yeast is strongly inhibited by EDTA.²

Precautions and Disclaimer

For Laboratory Use Only. Not for drug, household or other uses.

Preparation Instructions

The lyophilized powder contains approximately 85% buffer salts (Tris, citric acid and magnesium chloride). This product is soluble in deionized water (2 mg/ml), yielding a clear and colorless solution.

Storage/Stability

A frozen solution containing Mg^{2+} is stable for at least one year at -20 °C.

References

1. Specifications and Criteria for Biochemical Compounds, 3rd ed., National Academy of Sciences (Washington, DC: 1972), p. 108.
2. Butler, L. G., in The Enzymes, 3rd ed., Vol. IV, Boyer, P. D., ed., Academic Press (New York, NY: 1971), pp. 529-540.
3. Cooperman, B. S., The mechanism of action of yeast inorganic pyrophosphatase. *Methods Enzymol.*, **87(pt. C)**, 526-548 (1982).
4. Knight, W. B., et al., Investigations of the metal ion binding sites of yeast inorganic pyrophosphatase. *J. Biol. Chem.*, **259(5)**, 2886-2895 (1984).

TMG/RXR 12/03

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.