



pEx/Bac™-3 Expression Vector

TB470 0807

	Cat No.
pEx/Bac-3 DNA	71726-3
pEx/Bac-3 sequence landmarks	
hr5 enhancer	1196-1677
IE1 promoter	1682-2228
IE1 transcription start	2212
p10 promoter	2279-2391
p10 transcription starts	2323, 2324
His•Tag® coding sequence	2400-2429
Strep•Tag® II coding sequence	2529-2552
Multiple cloning sites (Nco I–Bsu36 I)	2393-2566
IE1 terminator	2571-2877
pUC ori	4527
<i>bla</i> (Amp ^R)	5288-6145
f1 origin	6277-6732

The dual-purpose pEx/Bac™ vectors are designed for cloning and high-level expression of proteins by transiently transfecting *Spodoptera*-derived insect cells or by generating baculovirus recombinants. Transient transfection and early baculovirus expression is driven by a promoter/enhancer combination that recruits endogenous insect cell transcription machinery, the AcNPV derived hr5 enhancer and ie1 promoter. Late/very late expression in the baculovirus mode is driven by the strong p10 promoter. The pEx/Bac-3 vector carries an N-terminal His•Tag® coding sequence followed by a recognition site for HRV 3C protease. The multiple cloning region is followed by an optional C-terminal Strep•Tag® II coding sequence (1). The presence of two “gentle elution” tags at both the N- and C-terminus is ideal for dual purification strategies designed to isolate full-length fusion proteins (2). Unique restriction sites are shown on the circle map.

- Skerra, A., Schmidt, T.G.M. (2000) *Meth. Enzymol.* 326, 271–304.
- Fiedler, M., Horn, C., Bandtlow, C., Schwab, M.E. (2002) *Protein Eng.* 15, 931–941.



pEx/Bac-3 cloning/expression regions

pIEx/Bac™-3 Restriction Sites

TB470 0807

Enzyme	# Sites	Locations	Enzyme	# Sites	Locations			
AatII	1	3251	EcoRV	1	1435			
Acc65I	1	2462	FspI	4	808 1850 5581 6754			
AccI	4	2829 5585 5958 6721	KasI	2	3800 4067			
AcuI	2	5014 6026	KpnI	1	2466			
AfeI	2	2928 4111	MfeI	4	2959 3304 3865 4230			
AgeI	1	3818	MluI	2	548 2068			
AhdI	2	648 5359	NaeI	2	932 6405			
AleI	2	2429 4049	NarI	2	3801 4068			
AlfI	2	2789 2823	NcoI	1	2393			
AlfI'	2	2789 2823	NgoMIV	2	930 6403			
Alol	1	6579	NheI	1	1802			
AlwNI	1	4882	NotI	1	2489			
ApaLI	3	2432 4780 6026	NsiI	1	333			
AscI	7	21 129 2027 2351 2595	PacI	2	2355 2595			
		3455 5531	PciI	1	4466			
BaeI	1	2446	PfoI	1	2279			
BaeI'	1	2479	PmeI	1	4346			
BamHI	1	2469	Ppil	3	5149 6023 6579			
BbeI	2	3804 4071	Ppil'	3	5181 5991 6547			
BbsI	2	647 3429	PsiI	6	469 558 1605 2995 4138			
BbvCI	1	565			6636			
BcgI	3	2679 3785 5864	PvuI	3	1403 5729 6775			
BcgI'	3	2713 3751 5898	PvuII	1	1972			
BciVI	5	731 1729 2064 4669 6196	SacI	1	2501			
BclI	3	1157 1167 1974	SalI	3	107 392 2482			
Bdal	4	2547 2581 4537 4571	SanDI	1	2455			
Bdal'	4	2547 2581 4537 4571	Scal	1	5839			
BglI	3	2191 5479 6747	SfoI	2	3802 4069			
BglII	1	1162	SgrAI	1	4058			
BlpI	1	2883	SmaI	1	2460			
BmgBI	1	3474	SnaBI	1	2970			
Bmri	2	2857 5399	SphI	2	1810 2893			
BmtI	1	1806	SspI	9	156 574 2169 2756 3064			
BpII	2	4236 4268			3261 4134 6163 6716			
BpII'	2	4236 4268	Swal	1	3060			
Bpml	1	5429	TaqII	4	4368 5707 5892 6606			
Bpu10I	1	565	TaqII'	3	308 6045 6062			
BpuEI	5	4238 4557 4855 5096 5964	TstI	2	2426 2701			
BsaI	2	434 5420	TstI'	2	2458 2733			
BsaXI	3	290 2696 6577	Tth111I	1	2819			
BsaXI'	3	320 2726 6547	XmaI	1	2458			
BsgI	3	3615 3909 4063	XmnI	2	3511 5958			
BsiWI	1	1980	ZraI	1	3249			
BsmBI	2	1098 3820	Enzymes that do not cut pIEx/Bac-3:					
BsmI	2	86 1971	AarI	AfII	AjuI	AjuI'	ApaI	AscI
BspHI	2	5186 6194	AsiSI	AvrII	BsaBI	BseRI	BspEI	BspMI
BsrBI	6	2059 3923 4218 4399 6200	BstAPI	BstEII	BstZ17I	FalI	FalI'	FseI
		6364	HindIII	HpaI	MscI	NdeI	NruI	
BsrDI	3	228 5420 5594	PasI	PfiMI	PmlI	PshAI	PspOMI	PspXI
BsrGI	4	198 917 2475 3300	PstI	RsrII	SacII	SapI	SbfI	SexAI
BssHII	2	1698 1959	SfiI	SpeI	SrfI	StuI	XbaI	XcmI
BssSI	3	1989 4639 6023	XhoI					
BstBI	1	2546						
BstXI	1	316						
Bsu36I	1	2566						
BtgZI	3	2842 3603 6503						
BtsI	3	1923 5759 5779						
Clal	3	1403 3942 4265						
CspCI	1	3192						
CspCI'	1	3157						
DraI	7	46 578 3060 4346 5225						
		5244 5936						
DraIII	1	6511						
DrdI	2	4574 6555						
EagI	2	2038 2489						
EarI	4	200 696 6154 6792						
Ecil	8	2112 3711 3717 3723 4528						
		4674 5502 6802						
EcoICRI	1	2499						
EcoNI	1	2449						
EcoP15I	3	4904 4962 5113						
EcoRI	5	1279 1381 1488 1561 1668						