
User Guide

Millex®-GV/HV Filter Unit (Sterile) with Durapore® Membrane

REF SLGV013SL*
SLHV013SL

Available languages

English	3	Svenska (Swedish)	27
Français (French)	6	Latviski (Latvian)	30
Italiano (Italian)	9	Magyar (Hungarian)	33
Deutsch (German)	12	Český (Czech)	36
Español (Spanish)	15	Polski (Polish)	39
Português (Portuguese)	18	Русский (Russian)	42
Nederlands (Dutch)	21	Norsk (Norwegian)	45
Dansk (Danish)	24		

*Brazil Product Registration
see page 2.

Definition of Symbols Used

Symbol	Definition	Symbol	Definition
	Sterilized using ethylene oxide		Use-by date
	Non-pyrogenic		Date of manufacture
	Consult instructions for use		Manufacturer
	Do not use if package is damaged and consult instructions for use		Caution
	Do not re-sterilize		CE conformity marking
	Do not re-use		Federal (USA) law restricts this device to sale by or on the order of a licensed healthcare professional.
	Catalogue number		Swiss authorised representative
	Batch code		Unique device identifier
	Temperature limit		

Made in Ireland
Merck Millipore Ltd.
Tullagreen, Carrigtwohill,
Co. Cork, Ireland
T45 KD29
an affiliate of Merck KGaA,
Darmstadt, Germany

MDSS CH GmbH
Laurenzenvorstadt 61
5000 Aarau
Switzerland

Brazilian Registration Holder (BRH)
Merck S/A
Rua Torre Eiffel, nº 100
Lote 3ª, Gleba A
Parque Rincão
Cotia/SP

*Brazil Product Registration
Number 81167709008

Millex®-GV/HV Filter Unit (Sterile) with Durapore® Membrane

Millex®-GV Filter **SLGV013SL** Qty: 100 pack

Millex®-HV Filter **SLHV013SL** Qty: 100 pack

- 13 mm
- Sterile
- Single use only
- Non-pyrogenic

Indications for Use/Purpose

The Millex®-GV/HV Filter Unit (Sterile) are medical devices intended for use as syringe filters to sterilize and/or clarify low volume solutions in direct patient care and pharmacy admixture applications. Devices are intended for medical and/or laboratory professional use only.

Introduction

These sterile 13 mm filters are syringe filters for the sterilization or clarification of aqueous solutions (column eluates, tissue culture additives, HPLC samples, protein pharmaceuticals, etc.) in volumes of 0.5 to 10 mL. These single-use, solvent-resistant filters are non-pyrogenic and non-toxic. They incorporate the low-protein binding Durapore® polyvinylidene fluoride (PVDF) membrane. These filters are suitable for use as syringe filters to sterile filter (GV), and/or clarify (GV/HV) low volume solutions in direct patient care and pharmacy admixture applications. Clinical benefits include the removal (GV) or reduction (HV) of microorganisms, particles, precipitates, and undissolved powders larger than the membrane's rated pore size from solutions for clinical use.

Chemical Compatibility


These filters are compatible with most aqueous solutions. Compatibility of solution should be tested before use. Visit [SigmaAldrich.com/MedMillex](https://www.sigmaaldrich.com/MedMillex) for more information.

WARNINGS

- Filter solutions with a temperature limit of 45 °C (113 °F).
- Syringe filters are for manual use only; do not use on automated systems.
- Not suitable for filtering high viscosity solutions or blood.
- Do not filter solutions containing 5 mg or less of active drug materials unless binding studies have been performed.
- To ensure sterility, do not use this product if the package is damaged.
- Use aseptic technique to maintain sterility removing product from primary package.
- Do not use with syringes smaller than 10 mL because pressures in excess of the maximum pressure rating may be reached, potentially causing damage to the filter and/or personal injury.
- Do not use for direct neuraxial applications which may come in contact with cerebrospinal fluid.
- Do not use the same syringe filter to filter solutions in both directions.
- Sudden loss of pressure could indicate failure of the filter.
- Single use only; do not re-use or resterilize.
- Make sure to wet the filter membrane thoroughly before injecting the solution; improperly wetted filters can become airlocked.
- Not suitable for use with medical pressure infusion equipment or gravity fed intravenous (IV) administration sets.
- Discard appropriately after single use.

Directions for Use

1. Wearing gloves using aseptic technique: fill syringe with solution to be filtered.

2. Remove cover from package.


3. Attach suitable syringe to filter inlet and remove assembly from package.

4. Hold syringe with filter pointing up and top off by pushing a few drops through to wet filter.

△ Excess solution may be hazardous and should be disposed with care.

5. Depress syringe plunger to filter solution into suitable container or compatible attachment on outlet.

Specifications and Performance Characteristics

Materials

Membrane	Hydrophilic Durapore® polyvinylidene fluoride (PVDF)
Pore size	Millex®-GV filter: 0.22 µm Millex®-HV filter: 0.45 µm
Housing	High density polyethylene (HDPE)

Dimensions

Inlet to outlet	21.1 mm (0.83 in.)
Diameter	14.7 mm (0.58 in.)
Filtration area	0.65 cm ² (0.10 in ²)

Temperature limit 45 °C (113 °F)

Pressure limit at 21 °C 5.2 bar (75 psi) inlet and differential

Filtration volume 0.5–10 mL

Hold-up volume ≤ 50 µL after air purge

Sterilization method Ethylene oxide gas

Connections Female Luer-Lok™ inlet; male Luer-slip outlet

Flow rate at 30 psi, 21 °C (2.1 bar, 70 °F) Millex®-GV filter: ≥ 7 mL/min

Flow rate at 10 psi, 21 °C (0.7 bar, 70 °F) Millex®-HV filter: ≥ 12 mL/min

Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

Contact Information

For the location of the office nearest you, go to SigmaAldrich.com/offices.

Disposal

Follow precautions for disposal of items contaminated with potentially infectious or hazardous biological material according to all applicable international, federal, state, and local regulations.

Technical Assistance

Visit the tech service page on our web site at SigmaAldrich.com/techservice.

Any serious incident of this device should be reported to manufacturer and competent authority of country where user is established.

Standard Warranty

The applicable warranty for the products listed in this publication may be found at SigmaAldrich.com/terms.

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Merck, Millipore, Millex, Durapore and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2023 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

The life science business of Merck operates
as MilliporeSigma in the U.S. and Canada.