

Technical Bulletin

Trypsin Inhibitor Assay Kit

Catalogue Number MAK620

Product Description

Trypsin (EC 3.4.21.4) is a digestive, serine protease that hydrolyzes dietary proteins in many eukaryotic and prokaryotic organisms. Trypsin predominantly cleaves peptide chains at the carboxyl side of lysine and arginine amino acids, but not before proline.

The Trypsin Inhibitor Assay Kit uses a fluorescein isothiocyanate (FITC)-labeled synthetic substrate. The fluorescein label is highly quenched. Upon digestion by trypsin present in the sample, the substrate is cleaved into smaller peptides, which abolishes the quenching of the fluorescence label. The fluorescence or fluorescence polarization (FP) of the FITC-labeled fragments is measured at $\lambda_{ex}/\lambda_{em} = 485/530$ nm. Inhibition is determined by the decrease in fluorescence.

Components

The kit is sufficient for 100 fluorometric assays in 96-well plates.

- Assay Buffer 10 mL
Catalogue Number MAK620A
- Substrate 100 μ L
Catalogue Number MAK620B

Reagents and Equipment Required but Not Provided

- Pipetting devices and accessories (e.g., multichannel pipettor)
- Fluorescent multiwell plate reader capable of $\lambda_{ex}/\lambda_{em} = 485/530$.
- Black plates with clear bottoms for fluorescence assays (Catalogue number

CLS3631 or equivalent) Cell culture or tissue culture treated plates are not recommended.

- 1.5 mL microcentrifuge tubes
- Purified trypsin (Catalogue number T8003)
- Control trypsin inhibitor
- Phenylmethanesulfonyl fluoride (PMSF) (Catalogue number P7626)

Precautions and Disclaimer

For R&D use only. Not for drug, household, or other uses. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

Storage/Stability

The product is shipped at ambient temperature. Store the kit at -20°C upon receipt.

Preparation Instructions

Prior to the assay, equilibrate all components to room temperature. This assay is based on a kinetic reaction. To ensure identical incubation time, addition of the Working Reagent to samples should be quick and mixing should be brief but thorough. Use of a multi-channel pipettor is recommended. Assays can be executed at room temperature.

Note: This kit does not include enzyme or inhibitor. Customers should provide enzyme or inhibitor of their choice.

Sample Preparation

To prevent autodigestion, pure trypsin solutions can be stored long term at -80°C

under acidic conditions (1 mM HCl, pH 3.0). Dilute purified trypsin in Assay Buffer. Trypsin should be kept on ice and used as soon as possible.

Dissolve test compounds in a solvent of choice (e.g. DMSO). It is prudent to first test the tolerance of the solvent by the enzyme of choice.

For porcine trypsin, the final DMSO concentration in the assay should be < 2%.

Procedure

Inhibitor Screening Reaction

1. Transfer 30 μL of trypsin into separate wells of a black, flat-bottom 96- well plate. Also, prepare 30 μL of trypsin ("Control") and 30 μL of Assay Buffer ("Blank") wells, respectively.
2. To the Control and Blank wells, add 10 μL of the solvent in buffer that the test compounds are dissolved in. For example, if the test compounds are dissolved in buffer containing 0.1% DMSO, add 10 μL of this solution to these wells.
3. To the remainder of the wells containing enzyme, add 10 μL of the test compounds. Tap plate to mix and incubate for 15 min at room temperature to allow the inhibitor to block Enzyme activity.
4. Prepare enough Working Reagent for all wells by mixing 1 μL Substrate and 99 μL Assay Buffer for each well. Add 80 μL of the Working Reagent to all wells. Briefly tap to mix. Incubate for 30 min at room temperature, protected from light.
5. Read the fluorescence or FP at $\lambda_{\text{ex/em}} = 485/530$ nm.

Results

Calculations

Trypsin inhibition for a test compound is calculated as follows:

$$\text{Enzyme Activity \%} = \frac{(F_{\text{Compound}} - F_{\text{Blank}})}{(F_{\text{Control}} - F_{\text{Blank}})} \times 100\%$$

Where:

Where F_{Compound} , F_{Control} , and F_{Blank} are the fluorescence or FP values at $\lambda_{\text{ex/em}} = 485/530$ nm of the test compound, no inhibitor "Control", and the no enzyme "Blank" wells, respectively.

Figure 1.

Trypsin was incubated with various concentrations of trypsin soybean inhibitor to generate a fluorescent IC50 curve.





Figure 2.

A second IC50 curve was generated for FP measurements.



References

1. Manea, M., Mezo, G., Hudecz, F., & Przybylski, M. (2007). Mass spectrometric identification of the trypsin cleavage pathway in lysylproline containing oligotuftsin peptides. *Journal of Peptide Science*. 13(4): 227–236.
2. López-Otín, C., Bond, J. S. (2008). Proteases: multifunctional enzymes in life and disease. *Journal of Biological Chemistry*, 283(45): 30433– 30437.
3. Shah, D., Mital, K. (2018). The Role of Trypsin: Chymotrypsin in Tissue Repair. *Advances in Therapy*, 35(1): 31–42.

Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

Technical Assistance

Visit the tech service page at SigmaAldrich.com/techservice.

Terms and Conditions of Sale

Warranty, use restrictions, and other conditions of sale may be found at SigmaAldrich.com/terms.

Contact Information

For the location of the office nearest you, go to SigmaAldrich.com/offices.

The life science business of Merck KGaA, Darmstadt, Germany
operates as MilliporeSigma in the U.S. and Canada.

MilliporeSigma, and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2026 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

mak620pis Rev 02/26

The Merck logo is displayed in a bold, red, sans-serif font.