

IRON(III)meso-TETRA-(N-METHYL-4-

ProductInformation

PYRIDINIUM)PORPHYRIN (Fe(III)T(4-N-MePy)P)

Product No. **I1648** Lot 75H8872

Store at 2-8°C CAS # [61943-73-5]

PRODUCT SUMMARY

Extinction coefficient (422nm, H_2O): 105 mM⁻¹cm⁻¹ Solubility: Fe(III)T(4-N-MePy)P is soluble at 10 mg/ml in H_2O .

Structure: See Figure 1.

Functional Assay

In the presence of O_2 and a reducing agent or terminal oxidants, Fe(III)T(4-N-MePy)P cleaves DNA at $(AT)_3^{-1}$.

References

- 1. Ward,B., Skorobogarty,A. and Dabrowiak,J.C. (1986) *Biochemistry*, **25**, 6875-6883
- Ward,B., Skorobogarty,A. and Dabrowiak,J.C. (1986) *Biochemistry*, 25, 7827-7833
- 3. Hui,X., Gresh,N. and Pullman,B. (1990) Nucl. Acids Res., **18**, 1109-1114
- Pasternack, R.F., Gibbs, E.J. and Villafranka, J.J. (1983) *Biochemistry*, 22, 5409-5417
- 5. Pasternack,R.F., Gibbs,E.J. and Villafranka,J.J. (1983) *Biochemistry*, **22**, 2406-2414
- Pasternack, R.F., Lee, H., Malek, P., Spencer, C., (1977) *J. Inorg. Nucl. Chem.* 39, 1865-1870
- 7. Byrnes,R.W., Fiel,R.J. and Datta-Gupta,N. (1988) *Chem.-Biol. Interact.* **67**, 225-242
- 8. Bernadou, J., Pratviel, G., Bennis, F. Girardet, M. and Meunier, B. (1989) *Biochemistry*, **28**, 7268-7275
- 9. Ding,L., Bernadou,J. and Meunier,B. (1991) *Bioconj. Chem.*, **2**, 201-206

BACKGROUND INFORMATION

Tetra-4-N-methylpyridinium porphyrins bind to DNA via groove binding and intercalative mechanisms¹⁻⁵. Intercalation is prohibited by axial ligation of the central metal ion by solvent molecules and/or counterions^{2,4,5}. These porphyrins are therefore restricted to binding to DNA via an electrostatic binding mechanism¹⁻⁵. Iron(III)T(4-N-MePy)P is an axially ligated ⁶ metallo tetracationic porphyrin which binds² and oxidatively cleaves¹ DNA from the minor groove of its binding site, minimally (AT)₃. Strand scission can be initiated by the addition of reducing agents and molecular oxygen or terminal oxidants to the porphyrin-DNA complex. The redox reagents include: ascorbate/O₂¹, DTT/O₂⁷, superoxide¹, iodosobenzene¹ and potassium peroxymonosulfate (oxone)^{8,9}.

FUNCTIONAL ASSAY

1 μ g of pBR322 DNA in 45 μ l of 10 mM Tris-HCl, pH 7.6, 1 mM EDTA containing from 1.0 mM to 0.1 nM Fe(III)T(4-N-MePy)P in 10 fold increments of dilution were incubated at 37 °C for 30 minutes. To initiate strand scission, 5 μ l of 20 mM oxone (Sigma Product No. P 1429) was added. The cleavage reaction was allowed to proceed for 30 minutes. Gel loading solution (Sigma Product No. G-2526, 40% (w/v) sucrose, 0.1 M EDTA, 0.5% (w/v) SDS and 0.05% (w/v) bromphenol blue) was then added to the reactions which were then immediately analyzed by agarose gel electrophoresis.

Lot specific results: At \geq 100 nM (Fe(III)T(4-N-MePy)P), the supercoiled band was completely degraded; at 10 nM (Fe(III)T(4-N-MePy)P), the supercoiled band was converted to >90% open circular DNA; at \leq 1 nM (Fe(III)T(4-N-MePy)P), the supercoiled band was comparable to the control which did not contain (Fe(III)T(4-N-MePy)P). In the absence of oxone no DNA degradation was observed.




Figure 1. Structure of (Fe(III)T(4-N-MePy)P)

Pcs2/00

Sigma brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.