## EnglishMillex®-SV Filter Unit (for low-binding applications)Cat. No. SLSV025LS (50/box)

- Single use only
- Sterile
- Non-pyrogenic

### Introduction

The bidirectionally supported Millex®-SV device is a low-binding, syringe-operated filter unit for the filtration of aqueous, proteinaceous solutions where minimal sample loss is desired. It will remove particles, precipitates, and undissolved powders larger than 5.0 microns (µm). The sterile Millex®-SV filter unit is non-pyrogenic and non-toxic. This single-use device consists of a 5.0 µm Durapore® membrane filter sealed in a plastic housing. Typical applications include particle or aggregate removal from protein solutions such as monoclonal antibodies, extracted cellular or tissue-based solutions, hormones, and drugs.

### Chemical Compatibility

The Millex®-SV filter unit is compatible with most aqueous solutions. It may be used to filter the agents listed below. This guide has been developed from technical publications, materials suppliers, and laboratory tests, and is believed to be reliable. However, because of variability in temperature, concentrations, duration of exposure, and other factors outside of our control which may affect the use of the unit, no warranty is given or is to be implied with respect to such information. Agents not listed below should be tested with the Millex®-SV filter unit prior to use.

|                          |                           |                        |                      |
|--------------------------|---------------------------|------------------------|----------------------|
| Acetic acid (glacial)    | Butyl alcohol             | Hydrofluoric acid      | Paraldehyde          |
| Acetic acid (5%)         | Carbon tetrachloride      | Hydrogen (gas)         | Pentane              |
| Aliphatic ethers         | Ethyl alcohol             | Hydrogen peroxide (3%) | Petroleum based oils |
| Ammonium hydroxide (6 N) | Ethylene glycol           | Hypo (photo)           | Petroleum ether      |
| Amyl acetate             | Formaldehyde              | Isobutyl alcohol       | Phenol (0.5%)        |
| Amyl alcohol             | Freon® solvent, TF or PCA | Isopropyl alcohol      | Pyridine             |
| Benzyl alcohol (1%)      | Glycerine (glycerol)      | Methyl alcohol         | Silicone oils        |
| Boric acid               | Helium                    | Mineral spirits        |                      |
| Brine (sea water)        | Hexane                    | Nitrogen               |                      |


### Directions for Use

### WARNINGS:

- Do not ensure sterility, do not use this product if the package is damaged.
- Do not use this product as an in-line filter for intravenous administration; it was not designed for long-term continuous use.
- Do not use with syringes smaller than 10 mL because pressures in excess of the maximum pressure rating may be reached, potentially causing damage to the filter unit and/or personal injury.

### CAUTIONS:

- Do not resterilize or reuse this filter unit, as Merck Millipore Ltd. cannot assure the sterility, integrity, and performance beyond a single use.
- Do not use the Millex®-SV unit at temperatures above 45 °C (113 °F).
- Do not use the same Millex®-SV unit to filter solutions in both directions.
- Do not use the Millex®-SV unit to filter emulsions or suspensions.
- Do not use the Millex®-SV unit to filter 5 mg or less of active drug materials unless binding studies have been performed.

### Specifications

|                                             |                                                                                 |
|---------------------------------------------|---------------------------------------------------------------------------------|
| <b>Materials</b>                            |                                                                                 |
| Membrane                                    | Hydrophilic Durapore® PVDF (polyvinylidene fluoride), type SV, pore size 5.0 µm |
|                                             | PVC (polyvinyl chloride)                                                        |
| <b>Housing</b>                              |                                                                                 |
| Inlet to outlet                             | 25 mm (0.98 in.)                                                                |
| Diameter                                    | 29 mm (1.14 in.)                                                                |
| Filtration area                             | 4 cm² (0.62 in²)                                                                |
| <b>Temperature limit</b>                    | 45 °C (113 °F)                                                                  |
| <b>Pressure limit at 25 °C</b>              | 5.2 bar (75 psi) inlet and differential                                         |
| <b>Filtration volume</b>                    | 1–100 mL                                                                        |
| <b>Hold-up volume</b>                       | ≤ 0.1 mL after air purge                                                        |
| <b>Sterilization method</b>                 | Ethylene oxide gas                                                              |
| <b>Connections</b>                          | Female Luer-Lok™ inlet, male Luer-slip outlet                                   |
| <b>Flow rate at 2.1 bar (30 psi), 25 °C</b> | ≥ 600 mL/min                                                                    |

### Notice

The information in this document is subject to change without notice and should not be construed as a commitment by Merck Millipore Ltd. ("Millipore") or an affiliate. Neither Merck Millipore Ltd. nor any of its affiliates assumes responsibility for any errors that may appear in this document.

### Technical Assistance

For more information, contact the office nearest you. In the U.S., call 1-800-MILLIPORE (1-800-645-5476). Outside the U.S., go to our web site at [www.millipore.com/offices](http://www.millipore.com/offices) for up-to-date worldwide contact information. You can also visit the tech service page on our web site at [www.millipore.com/techservice](http://www.millipore.com/techservice).

### Standard Product Warranty

The applicable warranty for the products listed in this publication may be found at [www.millipore.com/terms](http://www.millipore.com/terms) (within the "Terms and Conditions of Sale" applicable to your purchase transaction).

© 2012 EMD Millipore Corporation. All rights reserved.

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

Millipore, Millex®, and Durapore are registered trademarks of Merck KGaA.

Freon is a registered trademark of E. I. du Pont de Nemours and Co.

Luer-Lok is a trademark of Becton, Dickinson and Co.

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