


Product Information



71732 Sodium ionophore I

(ETH 227; *N,N',N''*-Triheptyl-*N,N',N''*-trimethyl-4,4',4''-propylidynetris(3-oxabutamide))
Selectophore®, function tested

Electrochemical Transduction

Ion-Selective Electrodes

Application 1 and Sensor Type¹⁻⁴

Assay of Na⁺ activity in whole blood, plasma, serum, undiluted urine and aqueous solutions with solvent polymeric membrane electrodes based on Sodium Ionophore I.

Recommended Membrane Composition

- 1.00 wt% Sodium Ionophore I ([71732](#))
- 66.00 wt% Bis(1-butylpentyl) decane-1,10-diyl diglutarate (ETH 469) ([30585](#))^{a)}
- 33.00 wt% Poly(vinyl chloride) high molecular weight ([81392](#))

^{a)} The use of bis(1-butylpentyl) adipate (BBPA) or bis(2-ethylhexyl)sebacate (DOS) leads to membrane electrodes of similar performance.

Recommended Cell Assembly

Reference || sample solution || liquid membrane | 0.1 M NaCl | AgCl, Ag

Electrode Characteristics and Function

Selectivity coefficients $\log K_{Na,M}^{Pot}$ determined by the mixed solution method.

	Required ^{b)} for blood	Required ^{b)} for urine	Found
$\log K_{Na,H}^{Pot}$	<-4.4	<-1.4	0.1
$\log K_{Na,K}^{Pot}$	<-0.6	<-2.1	-1.5
$\log K_{Na,Mg}^{Pot}$	<-1.2	<-0.6	-3.2
$\log K_{Na,Ca}^{Pot}$	<-1.3	<-0.6	-1.8
Stability:			
Drift [mV h ⁻¹]			0.03
Standard deviation [mV]	<-0.46		0.12
Reproducibility [mV]			0.23
Lifetime:			
log P _{TLC} ^{c)} ionophore	>8.4	>2.3	7.8
log P _{TLC} ^{c)} plasticizer	>12.8	>4.1	10.8

^{b)} 1% interference, worst case^{5,6}

^{c)} lipophilicity, determined by thin-layer chromatography⁷



Ion-Selective Field Effect Transistors

Application 1 and Sensor Type⁸⁻¹¹

Determination of sodium activity with Na⁺-selective field-effect transistor (ISFET). The gate membrane is based on „Urushi“ matrix with good durability.

Cocktail Composition

5.00 wt%	Sodium Ionophore I (71732)
44.77 wt%	Bis(2-ethylhexyl)phthalate (80030)
0.23 wt%	Potassium tetrakis(4-chlorophenyl)borate (60591)
50.0 wt%	Urushi (polymer from lacquer tree)

Electrode Characteristics and Function

Selectivity coefficients $\log K_{Mg,M}^{Pot}$.

$\log K_{Na,K}^{Pot}$	-1.5	$\log K_{Na,Ca}^{Pot}$	0.0
$\log K_{Na,NH_4}^{Pot}$	-1.6	$\log K_{Na,Mg}^{Pot}$	-2.4

Slope of linear regression: 53 mV/dec ($3 \cdot 10^{-4}$ to 10^{-1} M Na⁺)

¹ D. Ammann, P. Anker, E. Metzger, U. Oesch, W. Simon, in: Ion Measurements in Physiology and Medicine, Eds. M. Kessler, D. K. Harrison, J. Höper, Springer-Verlag, Berlin, Heidelberg 120 (1985).

² Neutral carrier based ion-selective electrode for the determination of Na⁺ in urine. H.-B. Jenny, D. Ammann, R. Dörig, B. Magyar, R. Asper, W. Simon, Mikrochim. Acta 74, 125 (1980).

³ Potentiometry of Na in undiluted serum and urine with use of an improved neutral carrier-based solvent polymeric membrane-electrode. P. Anker, H.-B. Jenny, U. Wuthier, R. Asper, D. Ammann, W. Simon Clin. Chem. 29, 1508 (1983).

⁴ D. Ammann, H.-B. Jenny, P. Anker, U. Oesch, W. Simon, in: Progress in Enzyme and Ion-Selective Electrodes, Eds. D. W. Lübbers, H. Acker, R. P. Buck, G. Eisenman, M. Kessler, W. Simon, Springer-Verlag, Berlin, Heidelberg, New York 21 (1981)

⁵ Ion selective electrodes in clinical chemistry. A. Lewenstam, Anal. Proc. 28, 106 (1991).

⁶ U. Oesch, P. Anker, D. Ammann, W. Simon, in: Ion-Selective Electrodes, ed. E. Pungor, I. Buzás, Akadémiai Kiadó, Budapest (1985).

⁷ Lifetime of neutral-carrier-based liquid membranes in aqueous samples and blood and the lipophilicity of membrane components. O. Dinten, U. E. Spichiger, N. Chaniotakis, P. Gehrig, B. Rusterholz, W. E. Morf, W. Simon, Anal. Chem. 63, 596 (1991).

⁸ Ion Sensitive Field Effect Transistors. J. Janata, R. Huber, Ion-Selective Electrode Rev. 1, 31 (1979).

⁹ The Operation of an ISFET as an Electronic device. P. Bergveld, Sensors and Actuators 1, 17 (1981).

¹⁰ Urushi matrix sodium, potassium, calcium and chloride-selective field-effect transistors. S. Wakida, M. Yamane, K. Higashi, K., Hiroy, Y. Ujihira, Sensors and Actuators 1, 412 (1990).

¹¹ S. Wakida, M. Yamane, K. Hiroy, Technical Digest of the 7th Sensor Symposium, 131 (1988).

