

Micropure®-EZ Enzyme Removers

■ Remove restriction and other enzymes from DNA in 30 seconds



Features

Fast

- 30-second enzyme removal at 14,000 x g

Easy to use

- 1. Load; 2. Spin; 3. Done.
- Negligible hands-on time
- One-step DNA purification
- Concentration or desalting with addition of Microcon

Safe

- Avoid phenol

High DNA recovery

- Over 85%

Convenient

- Ideal for multiple sample processing
- Microcentrifuge-compatible
- 14,000 x g rated

Compare with Alternatives

Although enzyme removal by heat inactivation after a reaction is simple and inexpensive, many enzymes cannot be completely inactivated by this method. Also, heat degradation of DNA may occur. Precipitation techniques and “bind and elute” protocols, requiring a series of steps, are labor-intensive. Extraction followed by precipitation techniques requires the use of organic solvents, which, if carried over, can inhibit subsequent enzyme reactions. By comparison, Micropure-EZ is simple, fast and assures high recovery. Just load the enzymatic reaction directly into Micropure-EZ and spin for 30 seconds. dsDNA, free of enzymes, is efficiently recovered in the filtrate (see Figure 1).

Micropure-EZ provides an easy, rapid means of removing restriction and other enzymes from double-stranded (ds) DNA. The membrane in the device has a high affinity for protein but not for dsDNA. In a single 30-second spin, enzymes from the reaction mix are selectively adsorbed and DNA is recovered in a vial — undiluted and enzyme-free. If concentration or desalting of the sample is required, Microcon® centrifugal filter devices may be substituted for the

filtrate vial, providing for DNA purification, concentration and desalting in a single unit.

Unlike precipitation or “bind and elute” purification, Micropure-EZ needs little hands-on time. This is especially useful when processing multiple samples. Also, DNA recovery is not affected by the binding capacity of an adsorption matrix. Micropure-EZ will deliver minimum 85% recovery of dsDNA, higher with a quick 5-10 µl rinse.

MILLIPORE

Concentrating and Desalting

When DNA concentration or buffer exchange are desired in addition to purification, all three steps can be performed in a single device by adding a Microcon centrifugal filter devices (see Figure 2). The combination of Micropure-EZ and Microcon allows DNA to be simultaneously separated from enzymes and concentrated. By adding buffer to the concentrated DNA, a second "diafiltration" spin can desalt the sample or exchange buffer. Result: DNA recoveries are superior to traditional techniques. No organic solvents or lengthy protocols are needed.

Specifications

Maximum initial sample volume
250 µL

Holdup volume, membrane, support
< 5 µL

Maximum relative centrifugal force
14,000 x g
(500 x g for Microcon YM-100)

Dimensions

Diameter
10.7 mm

Length
Micropure-EZ plus vial 40.9 mm
with Microcon 49.2 mm

Selection Guide

Intended use	Recommended device
Enzyme removal only	Micropure-EZ
Enzyme removal and DNA concentration/desalting	Micropure-EZ filter unit only and Microcon YM-30 or YM-50




Figure 1. Enzyme removal with Micropure-EZ

Figure 2. Simultaneous enzyme removal and concentration, buffer exchange or desalting, with Micropure-EZ and Microcon




Figure 3. Micropure-EZ removes enzyme activity without binding DNA or affecting subsequent DNA manipulation.

Lane 1: 1 kb DNA ladder

Lane 2: Uncut, untreated pBR322

Lane 3: pBR322 spun through Micropure-EZ (No detectable changes in DNA: No nick, cuts, etc.)

Lane 4: pBR322 incubated with 10 units of Pst I for 1 hr at 37 °C

Lane 5: pBR322 in the restriction buffer spun through Micropure-EZ, then digested with 20 units of Pst I for 1 hr at 37 °C (No enzyme inhibitors leaching from the membrane)

Lane 6: Reaction mixture containing 100 units of Pst I passed through Micropure-EZ and incubated with pBR322 for 18 hrs at 37 °C (No detectable enzyme activity).

Enzymes Removed by Micropure-EZ

In tests, all enzymes were removed from solutions containing DNA or RNA (in the case of RNase A). Five µg of bovine serum albumin were also present during removal of restriction enzymes. Removal was indicated by undetectable or inconsequential enzyme activity in filtrate (< 0.08% residual in the case of T4 polynucleotide kinase).

Enzyme	Challenge	Enzyme	Challenge
AMV reverse transcriptase	50 U	<i>EcoR</i> I	100 U
Calf intestinal alkaline phosphatase	10 U	<i>Hae</i> III	100 U
DNase I (bovine pancreas)	10 U	<i>Hinc</i> II	50 U
Exonuclease III (<i>E. coli</i>)	100 U	<i>Hind</i> III	100 U
MMLV reverse transcriptase	600 U	<i>Hpa</i> I	25 U
Mung bean nuclease	50 U	<i>Kpn</i> I	50 U
Proteinase K (Amresco)	5 µg	<i>Mbo</i> I	25 U
T4 DNA ligase	2,000 U*	<i>Mlu</i> I	50 U
T4 DNA polymerase	15 U	<i>Nco</i> I	50 U
T4 polynucleotide kinase**	50 U	<i>Nde</i> I	100 U
<i>Taq</i> DNA polymerase	5 U	<i>NgoM</i> I	50 U
Terminal deoxynucleotidyl transferase	45 U	<i>Nhe</i> I	25 U
<i>Acc</i> I	50 U	<i>Not</i> I	50 U
<i>Apa</i> I	100 U	<i>Nru</i> I	50 U
<i>Bam</i> H I	100 U	<i>Pst</i> I	100 U
<i>Bcl</i> I	50 U	<i>Sac</i> I	100 U
<i>Bgl</i> I	50 U	<i>Sac</i> II	100 U
<i>Bsi</i> W I	70 U	<i>Sal</i> I	100 U
<i>Bss</i> H II	20 U	<i>Sca</i> I	50 U
<i>Bst</i> N I	50 U	<i>Sph</i> I	25 U
<i>Dpn</i> I	100 U	<i>Sst</i> I	50 U
		<i>Xho</i> I	100 U

* New England Biolabs unit definition

** T4 polynucleotide kinase is not recommended with oligonucleotides (dsDNA only). The results suggest that this kinase mediates the binding of oligos to the membrane in Micropure-EZ device, causing oligo loss.

Enzymes Not Removed by Micropure-EZ

In tests, Micropure-EZ did not remove the indicated number of units of the listed enzymes. It may be effective in removing a lower number of units.

Enzyme	Challenge	Enzyme	Challenge
Bacterial alkaline phosphatase	0.6 U	<i>Eae</i> I	15 U
DNA polymerase I (Klenow)	20 U	<i>EcoR</i> V	50 U
Exonuclease I	50 U	<i>Hinf</i> I	50 U
<i>Pfu</i> DNA polymerase	2.5 U	<i>Msp</i> I	100 U
Ven TM DNA polymerase	4 U	<i>Pvu</i> I	25 U
Shrimp alkaline phosphatase	1 U	<i>Pvu</i> II	50 U
RNase A (bovine pancreas)	1 µg	<i>Sau</i> 3A I	20 U
<i>Apa</i> I	10 U	<i>Sfi</i> I	50 U
<i>Bgl</i> II	50 U	<i>Sma</i> I	25 U
<i>Bso</i> B I	50 U	<i>Xba</i> I	50 U
<i>Cla</i> I	25 U		

T4 RNA ligase is not recommended. Results suggest this ligase mediates binding of nucleic acids to the membrane in Micropure-EZ causing sample loss.

Ordering Information

Description	Membrane MW Cut-Off	Nucleotide Cut-off Double-stranded	Qty/Pk	Catalogue No.
Micropure-EZ Enzyme Removers			8	42528
			24	42549
			100	42530
Micropure-EZ Enzyme Removers Filter units only (for use with Microcon)			8	42531
			24	42532
			100	42533
Microcon Centrifugal UF Filter Devices (includes 2 vials with each device)				
Microcon YM-3	3,000	10	8	42420
			24	42403
			100	42404
Microcon YM-10	10,000	20	8	42421
			24	42406
			100	42407
Microcon YM-30	30,000	50	8	42422
			24	42409
			100	42410
Microcon YM-50	50,000	100	8	42423
			24	42415
			100	42416
Microcon YM-100	100,000	125	8	42424
			24	42412
			100	42413

To Place an Order or Receive Technical Assistance

For additional information call your nearest Millipore office:

In the U.S. and Canada, call toll-free **1-800-MILLIPORE (1-800-645-5476)**

In the U.S., Canada and Puerto Rico, fax orders to **1-800-MILLIFX (1-800-645-5439)**

On the Internet

<http://www.millipore.com>

E-mail: tech_service@millipore.com



MILLIPORE

Millipore, Amicon, Microcon and Micropure are registered trademarks of Millipore Corporation or an affiliated company.

Vent is a trademark of New England Biolabs, Inc.
Lit. No. PF174 Printed in U.S.A. 6/98 98-107

Copyright © 1998 Millipore Corporation or an affiliated company, Bedford, MA

Printed with soy inks on recycled paper.

Millipore Worldwide

AUSTRALIA

Tel. 1 800 222 111
or (02) 9888 8999
Fax (02) 9878 0788

AUSTRIA, CENTRAL EUROPE, C.I.S., AFRICA, MIDDLE EAST AND GULF
Tel. (43) 1 877-8926
Fax (43) 1 877-1654

BALTIC COUNTRIES
Tel. 358-9-804 5110
Fax 358-9-859 6616

BELGIUM AND LUXEMBOURG
Tel. (02) 726 88 40
Fax (02) 726 98 84

BRAZIL
Tel. (011) 548-7011
Fax (011) 548-7923

CANADA
Tel. 1-800-645-5476
Fax 1-800-645-5439

CHINA, PEOPLE'S REPUBLIC OF
Beijing:
Tel. (8610) 6500-8063
Fax (8610) 6500-7372
Guangzhou:
Tel. (8620) 8669-5293
Fax (8620) 8669-5302
Hong Kong:
Tel. (852) 2803-9111
Fax (852) 2513-0313
Shanghai:
Tel. (8621) 6373-7256
Fax (8621) 6320-0236

CZECH REPUBLIC
Tel. 02-205 138 41
02-205 138 42
Fax 02-205 14294

DENMARK
Tel. 46 59 00 23
Fax 46 59 13 14

FINLAND
Tel. (09) 804 5110
Fax (09) 859 6616

FRANCE
Tel. (01) 30.12.70.00
Fax (01) 30.12.71.80

GERMANY
Tel. (06196) 494-0
Fax (06196) 482237

HUNGARY
Tel. 1-205 9784
Fax 1-205 9792

INDIA
Tel. (91) 80-839 46 57
Fax (91) 80-839 63 45

ITALY
Vimodrone (Milano):
Tel. (02) 250781
Fax (02) 2650324
Roma:
Tel. (06) 5203600
Fax (06) 5295735

JAPAN
Tel. (03) 5442-9711
Fax (03) 5442-9736 Analytical
9737 BioProcess
9734 Lab Water
9738 Microelectr.

KOREA
Tel. (822) 551-0990
Fax (822) 551-0228 Analytical
0228 BioProcess
0227 Microelectr.

MALAYSIA
Tel. 603-757 1322
Fax 603-757 1711

MEXICO
Tel. (525) 576-9688
Fax (525) 576-8706

THE NETHERLANDS
Tel. 076-5022000
Fax 076-5022436

NORWAY
Tel. 22 67 82 53
Fax 22 66 04 60

POLAND
Tel. 22-669 12 25
22-663 70 31
Fax 22-663 70 33

PUERTO RICO
Tel. (787) 273-8495
Fax (787) 747-6553

SINGAPORE
Tel. (65) 842 1822
Fax (65) 842 4988

SPAIN AND PORTUGAL
Madrid:
Tel. 91 729 03 00
Fax 91 729 29 09
Barcelona:
Tel. 93 451 70 00
Fax 93 451 60 48

SWEDEN
Tel. 08-628 69 60
Fax 08-628 64 57

SWITZERLAND
Volkswil:
Tel. (01) 908-30-60
Fax (01) 908-30-80
Lausanne:
Tel. (021) 648 51 18
Fax (021) 648 51 37

TAIWAN
Taipei:
Tel. (886-2) 700-1742
Fax (886-2) 755-3267
Hsin Chu City:
Tel. (886-3) 571-0178
Fax (886-3) 572-9520

U.K. AND IRELAND
Tel. (01923) 816375
Fax (01923) 818297

U.S.A.
Tel. (781) 533-6000
Fax (781) 533-3110

IN ALL OTHER COUNTRIES
Millipore Intertech (U.S.A.)
Tel. +1 (781) 533-8622
Fax +1 (781) 533-8630