
# **MILLIPORE**



upstate
now part of Millipore

#### MAPK: Antibodies, Proteins, Kits, siRNA and Reagents

#### MAPK Signaling

- A Brief Review
- Pathways of MAPK

#### Classical MAPK

• FRK

#### Stress Induced MAPK

- JNK
- P38

#### MAPK Signaling Product Listing

• EGFR

#### Phosphorylation

 $\bullet$  Phoshoserine, threonine, and tyrosine, including  $4\text{G}10^{\text{\tiny{19}}}$  clone

# Upstate®, Chemicon® and Linco® are now part of Millipore

The goal of this combined company is to provide more innovative tools, services and application expertise that will improve your productivity. Our first priority is to support your work as a valued partner. Together we now offer:

- A broader range of products and services for markets including drug discovery, protein identification, antibody purification, molecular biology, stem cell research, water purification and general filtration.
- An increased number of applications and protocols that facilitate your development process.
- Extended global expertise, leadership and presence to serve you more effectively.
- Larger support and research development teams to deliver innovative products that cost-effectively impact your business.



### Antibodies, Proteins, Kits, siRNA and Reagents

The Mitogen-Activated Protein Kinase (MAPK) pathways transmit extracellular signals, through receptor tyrosine kinases (RTKs) or integrin stimulation, to both nuclear and cytoplasmic targets. The MAPK signaling cascade transmits and amplifies the signals of extracellular mitogens (growth factors, etc) through the subsequent activation of RTKs. This action creates a signaling cascade of kinases that results in the activation of numerous transcription factors that can permit a cell's entry into the cell cycle.

The MAPKs are comprised of four distinct families that transduce signals in response to numerous stimuli. They include the classical extracellular signal-related kinase 1 and 2 (ERK 1/2) pathway. The stress-related signaling pathways that encompass both the c-Jun N-terminal Kinase/Stress-Activated protein Kinase (JNK/SAPK) and the p38 MAPK pathways. Finally, the fourth MAPK pathway is the ERK5 (BMK1).

MAPKs are activated in a characteristic and highly specific phospho-relay system composed of three kinases (MEKK, MEK, and MAPK) which phosphorylate/activate each other in the cascading fashion described in the first paragraph. Different stimuli activate different MEKKs, which selectively activate their downstream target (MEK), which then ultimately phosphorylates and activates its MAPK.

MAPK is the most well-characterized signaling pathway in research, yet much is still unknown and ongoing studies for its implications in cancer and neurology continue to proliferate. Millipore is at the forefront of this ongoing research, providing complete solutions for active research in both academic and bio/pharma settings.

#### The Pathways of MAPK


The ERK pathway is activated by various growth factors (i.e. EGF, FGF, PDGF, VEGF, etc.) and integrin stimulation. The JNK and p38 pathways, also referred to as the stress-related MAPK, are regulated by stress inducing signals such as osmotic

stress and UV radiation as well as proinflammatory cytokines such as IL-1 and TNF $\alpha$ . It should be noted that although there is selective stimulation of each of these distinct pathways, there also is a great deal of cross-talk between them. For instance, p38-dependent activation of ERK has been demonstrated.

All of the MAPKs (ERK, JNK, p38MAPK, and ERK5) are activated through phosphorylation on both threonine and tyrosine residues at the Thr-X-Tyr dual phosphorylation motif, where X is Glu, Pro, and Gly for ERK, JNK, and p38, respectively. Once phosphorylated and active, they then go on to phosphorylate downstream targets in both the cytoplasm and the nucleus. One of their main targets of MAPK are transcription factors that are activated via this phosphorylation.

Many malignant cancers are characterized by the deregulation of MAPK signaling cascades. Cancerous cells do not respond to cell signaling that would normally result in cell cycle arrest or apoptosis (programmed cell death). In fact, constitutive MAPK signaling contributes to the evolution of some of the most lethal forms of cancer. In some cancers, MAPK is upregulated and results in the migration and invasion of the cancerous cells. Inhibition of the MAPK expression reduces migration and invasion. Many tumors show an increase, sustained activation of this pathway. This is one reason why the MAPK signaling pathway is a highly sought anti-cancer drug target. Among the cellular activities that MAPKs modulate are cell division, proliferation, survival, differentiation, apoptosis, motility and metabolism.

Millipore is your complete source for the study of MAPK signaling, particularly in these critical areas of disease research and drug discovery. We supply the highest quality cell signaling solutions and compounds available with our industry leading, fully guaranteed Upstate cell signaling products. As your single source provider, Millipore can help boost your productivity leaving more time for the science that is at the heart of your expertise.



#### Classical MAPK

The classical, canonical MAPK signaling pathway is an evolutionatily conserved pathway that controls the growth and survival of a large group of human tumors where it is shown to have increased and sustainable activation. Classical MAPK is involved in the control of many fundamental cellular processes that include cell proliferation, survival, differentiation, apoptosis, motility and metabolism. It accomplishes this through the promotion of cell proliferation and the prevention of apoptosis. In a normal cell, the signaling pathway is initiated by the binding of a growth factor, such as EGFR/ERBB, VEGF, PDGF, HGF, and FGF, to its respective receptor. This ligand binding activates RTK that then binds to the growth factor receptor bound protein 2 (Grb2) that subsequently binds to specific phosphorylated residues in the intracellular tail portion of the receptor. SOS (son of sevenless), a guanine nucleotide exchange factor, binds to Grb2 to activate the small G-protein Ras by catalyzing the replacement of GDP to GTP. It is in this GTP-bound active state that Ras activates the kinase activity of the Ser/Thr kinase Raf (MAP Kinase Kinase Kinase). Raf then binds to and activates MEK (MAP Kinase Kinase) by phosphorylating it on the two residue motif. Both of these two proteins - Ras and Raf - are mutated in a large number of human cancers. MEK, a dual kinase (a Ser/Thr and a Tyr kinase), phosphorylates and activates MAPK/ERK on the TxY motif.

The regulation of a large number of cellular processes are dependent upon the activation state of ERK, so controlling this pathway could have profound effects on various diseases. As noted above, many tumors show an increased, sustained activation of this pathway and Raf (BRAF and c-Raf/Raf1), the top kinase in this pathway, is an important modulator of this activity. In addition to stimulating the MEK/ERK pathway, It is widely held that Raf directly antagonizes the death promoting activity of apoptosis signaling-regulating kinase 1. Raf is thought to control apoptosis through two distinct mechanisms in both a MEK dependent and independent manner. In the MEK dependent manner, Raf activates ERK via MEK, which then phosphorylates RSK that leads to the phosphorylation and inactivation of the pro-apoptotic protein BAD and the activation of CREB.

Ras is a heavily studied target of both academic

and pharmaceutical research because of its implications in various pathways and diseases, as well as being mutated in a large number of human cancers. In its oncogenic mutated state, Ras is unable to hydrolyze GTP to GDP, thus staying in an active state and activating numerous pathways including the MAPK pathway through its activation of Raf, as well as others such as PI3 Kinase and RalGDS. For this reason it has become a popular drug target.

One path that the pharmaceutical industry has taken is to control Ras and its activity by finding its weakest feature. For its activation, Ras must localize to the plasma membrane, but interestingly, it lacks a transmembrane domain. To activate, Ras must first undergo a post-translational modification (PTM) known as prenylation at its C-terminal CAAX motif where C is a cysteine, the As are alipathic amino acids, and X can be any amino acid. Ras modification occurs in a controlled three step process. The first step in the process is the prenylation of the C in the CaaX motif that is initiated by the covalent attachment of farnesyl groups to the cysteine that is catalyzed by the heterodimer enzymes farnesyl transferases a and b. After this modification, the -aaX of the motif is proteolytically removed via Rce1 (Ras Converting Enzyme 1), a membrane associated endoprotease, by a mechanism that is still not fully understood. Finally, the C-terminal prenylcysteine is now methlylated by ICMT (Isoprenylcysteine Carboxymethyl Transferase). Many Ras targeting drugs have yet to pass clinical trials and their success in treating tumors associated with Ras activation remains unproven.

...20 years of experience backed by the most highly regarded scientific expertise.

Millipore's antibodies and kinases are widely used in drug discovery and our small G-protein kits have become one of the best selling in the market. With more than 20 years of experience backed by the most highly regarded scientific expertise, Millipore is the leader in complete cell signaling solutions.

The Raf family (A-Raf, B-Raf, and Raf1) is another heavily studied drug target because it

influences is so many downstream MAPK pathways. Raf also inhibits apoptosis by binding to and suppressing the activity of the pro-apoptotic kinase ASK1. Raf activates MEK by phosphorylating it on two serine residues. Some of the negative regulators of Raf include Akt and the serum/glucocorticoid inducible kinase. It both negatively and positively regulates large numbers of signaling targets. This development theory reasons that by selecting Raf as a drug target, the drug could inhibit the activity of a number of growth factor stimulations as opposed to a drug that inhibits only one of the growth factor receptors such as EGFR or VEGFR. It is thought that B-Raf might be the predominant activator of MEK and that Raf-1 has a role in protection against apoptosis; a process that does not require its kinase activity or its activation of MEK. This theory was further fueled by the finding of the larger amount of oncogenic BRAF mutations in human cancers.

MEK has also been implicated in a large number of human cancers. The two highly

homologous MEK proteins, MEK1 and MEK2, are dual-specific kinases that contain two consensus kinase motifs with one involved in the phosphorylation of Serine/Threonine residues and the other in the phosphorylation of tyrosine residues. MEK's only known substrate is ERK which it subsequently phosphorylates and activates in its TEY motif, phosphorylation sites Tyr185 followed by Thr183, that lay in the ERK activation loop.

The Raf/MEK/ERK pathway is also involved in focal adhesion and migration through its interactions with Paxillin. Paxillin is constitutively associated with MEK, and in response to hepatocyte growth factor (HGF), also recruits and binds activated Raf-1 and inactive ERK. ERK binds to the Src Paxillin phosphorylation site Tyr 1 18. This activates the MAPK signaling pathway specifically in newly forming focal adhesions that results in the ERK phosphorylation of Paxillin Ser83, a site that promotes its binding to FAK and to the downstream activation of Rac.

# Stress-Related MAPK

# Signaling Pathways

Millipore also provides signaling products for the study of stress-related MAPK pathways such as JNK and p38. Stress related pathways are activated by external stressors such as ultraviolet light and osmotic stress.

#### JNK

c-Jun N-terminal Kinase (JNK) was originally found to be stress-activated. When activated, JNK binds to and phosphorylates a number of transcription factors, most notable Jun, and are thus transcriptional regulators. C-Jun, a member of the AP-1 family of transcription factors, regulates cytokine gene expression. Along these lines, JNK inhibitors may be effective in the control of rheumatoid arthritis (RH), an autoimmune disease involving increased production of inflammatory cytokines. Inhibition of JNK signaling also enhances chemotherapy-induced inhibition of tumor cell growth, suggesting the JNKs may provide

a molecular target for the treatment of cancer.

Both JNK and p38 MAPK are often simultaneously activated in response to environmental changes. These two stress kinase signaling pathways have evolved to relay increasingly complex ranges of environmental stimuli. Just as ERK is activated by MEK, JNK is activated by MKK4 and MKK7. Once activated, it is known to phosphorylate a number of downstream targets including the transcription factors c-Jun, ATF-2, and ELK-1.

#### p38

The p38 MAPK family of Ser/Thr kinase proteins is made up of 4 isoforms, p38 $\alpha$  (MAPK14, SAPK2a), p38 $\beta$  (MAPK11, SAPK2b), p38 $\delta$  (MAPK13, SAPK4), and p38 $\gamma$  (MAPK12, SAPK3). p38 $\alpha$ , the first isoform discovered, was identified in LPS stimulated mouse macrophages and was found to have high sequence homology to HOG1 kinase.

The other isoforms share 74 %, 57 %, and 60 % sequence homology with the  $\alpha$  isoform, respectively. All four isoforms possess the TGY motif and are activated by the dual phosphorylation of the threonine and tyrosine residues in this motif, which are Thr180/Tyr182 in human p38 $\alpha$ . p38, especially the  $\alpha$  and  $\beta$ isoforms, are able to phosphorylate a wide array of targets including numerous kinases such as MAPKAP-K2, MAPKAP-K3, PRAK, MSK, and MNK and many transcription factors including STAT1, p53, NFATp, ATF-2, Elk-1 (although weakly), and MEF2a/c. They have traditionally been associated with stress, immune response, and regulation of cell proliferation and apoptosis. They also have a role in regulating the production of proinflammatory cytokines and have been associated with diseases such as asthma, rheumatoid arthritis,

diabetes, and inflammatory bowel disease. Just as ERK is phosphorylated and activated by MEK, p38 $\alpha$  MAP Kinase is activated predominantly by its upstream regulatory kinases MKK6 and MKK3 as well as MKK4 which also activates JNK, but is not a substrate for MKK7 as is JNK.

#### **References:**

Kolch, W. 2005. Nature Reviews Molecular Cell Biology. **6**:827-837

Wada, T. and Penninger, J.M. 2004. *Oncogene*. **23**:2838-2849. Seboli-Leopold, J. and Herrera, R. 2004. *Nature Reviews Cancer.* **4**:937-947.

Johnson, G. L. and Lapadat, R. 2002. *Science*. **298(5600)**: 1911-1912.


Gaestel, M. and Shi, Y. 2002. *Biol. Chem.*, Vol. **383**: 1519-1536.

Kolch, W. 2000. Biochem. J. 351:289-305.

Pomerance, M. et al. 2000. J. Biol. Chem. **275(51)**: 40539-40546.

Ishitani, T. et al. 1999. Nature. 399: 799-802.





# MAPK Products

**Note:** Additional species and applications may apply. For more information, call Tech Support at 800 437 7500. For international technical service, please look for your local listing in the back of this brochure.

#### **Tested Applications**

| Abbr. | Description                          | Abbr. | Description                     | Abbr.  | Description                                     |
|-------|--------------------------------------|-------|---------------------------------|--------|---|
| ABA   | Affinity Binding Assay               | HDAC  | Histone Deacetylase Assay       | NEUT   | Neutralizing                                    |
| ACT   | Activity Assay                       | HI    | Hemagglutination Inhibition     | NT     | Nitration                                       |
| ADH   | Stimulates ECM Adhesion              | HMT   | Histone Methyltransferase Assay | NUEX   | Nuclear Extraction                              |
| Al    | Agonist or Inhibitor                 | IAP   | Immunoaffinity Purification     | PA     | Phosphatase Assay                               |
| AMP   | DNA Amplification                    | IC    | Immunocytochemistry (Cells)     | PC     | Positive Control                                |
| APA   | Affinity Precipitation Assay         | ID    | Immunodiffusion                 | PCU    | Protein Clean-up                                |
| BD    | Beadlyte® Assay                      | IEP   | Immunoelectrophoresis           | PD     | Protein Determination                           |
| CA    | Caspase Assay                        | IF    | Immunofluorescence              | PIA    | Peptide Inhibition Assay                        |
| CC    | Culture Confirmation                 | IFIX  | Immunofixation                  | RIA    | Radioimmunoassay                                |
| ChIP  | Chromatin Immunoprecipitation        | IH    | Immunohistochemistry (Tissue)   | RNAi   | RNAi/siRNA/Gene Knockdown                       |
| CULT  | Cell Culture                         | IH(P) | Immunohistochemistry (Paraffin) | RPA    | Ribonuclease Protection Assay                   |
| DB    | Dot Blot                             | IND   | Induces Function                | RT-PCR | Reverse Transcriptase Polymerase Chain Reaction |
| EA    | Enzyme Assay                         | INHIB | Inhibits Activity/Function      | SW     | Software Needed                                 |
| ELISA | Enzyme Immunoassay (ELISA)           | IP    | Immunoprecipitation             | TFX    | Transfection                                    |
| EM    | Electron Microscopy                  | IPK   | IP-Kinase Assay                 | UC     | Uncharacterized Antiserum                       |
| EMSA  | Electrophoretic Mobility Shift Assay | IPX   | Immunoperoxidase Staining       | VVB    | Immunoblotting (Western)                        |
| FC    | Flow Cytometry (FACS)                | IRMA  | Immuno Radio-Metric Assay       | Web*   | Important additional product reactivity         |
| FUNC  | Affects Function                     | IT    | Immunotoxin                     |        | information available on datasheet              |
| GPA   | G-Protein Assay                      | KA    | Kinase Assay                    |        |   |
| HA    | Hemagluttination                     | LFA   | Lateral Flow Assay              |        |   |
| HAT   | Histone Acetyltransferase Assay      | NB    | Northern Blot                   |        |   |

#### **Tested Species Reactivity**

| Abbr. | Description          | Abbr. | Description       | Abbr. | Description                           |
|-------|----------------------|-------|-------------------|-------|---------------------------------------|
| A     | All Species          | Gr    | Gerbil            | R     | Rat                                   |
| Am    | Amphibian            | Gs    | Ground Squirrel   | Rb    | Rabbit                                |
| As    | Aspergillus          | Gt    | Goat              | Rc    | Raccoon                               |
| Av    | Avian                | Н     | Human             | rH    | Recombinant Human Protein             |
| В     | Bovine               | H-sp  | Human Only        | Rp    | Reptilian                             |
| Bab   | Baboon               | Ht '  | Hamster           | Sal   | Salamander                            |
| Bact  | Bacterial            | In    | Insect            | Seal  | Seal                                  |
| Bat   | Bat                  | Inv   | Invertebrates     | Sh    | Sheep                                 |
| Ca    | Canine (Dog)         | Kn    | Kangaroo          | Shk   | Shark                                 |
| Ch    | Chicken              | Lg    | Ligia             | SHm   | Syrian Hamster                        |
| Chp   | Chimpanzee           | Lz    | Lizard            | Shp   | Shrimp                                |
| Crb   | Crab                 | M     | Mouse             | Sį    | Schistosoma japonicum                 |
| Crf   | Crawfish             | Ma    | Mammals           | Sn    | Snail                                 |
| Di    | Dictyostelium        | Md    | Mule Deer         | Snk   | Snake                                 |
| Dr    | Drosophila           | Mi    | Mink              | Spd   | Spider                                |
| Ec    | E. coli Bacteria     | Mk    | Monkey            | Sqd   | Squid                                 |
| Ech   | Echinoderms          | M     | Mollusk           | Su    | Sea Urchin                            |
| Ecl   | Enterobacter cloacae | Nem   | Nematode          | T     | Tetrahymena                           |
| Elk   | Elk                  | Nr    | Neurospora crassa | Vo    | Vole                                  |
| Eq    | Equine (Horse)       | Ор    | Opposum           | Vrt   | Vertebrates                           |
| Eu    | Eukaryote            | Ox    | Ox                | ₩R.   | Most common vertebrate species tested |
| F     | Fish                 | Pl    | Green Plants      | Xn    | Xenopus                               |
| Fe    | Feline (Cat)         | Pm    | Primate           | Υ     | Yeast (S. cerevisiae)                 |
| Fg    | Frog                 | Pn    | Penicillium       | Zf    | Zebra Fish                            |
| Ft    | Ferret               | Po    | Porcine (Pig)     |       |                                       |
| Gp    | Guinea Pig           | Qu    | Quail             |       |                                       |

#### 14-3-3

The 14-3-3 proteins are a family (at least seven members in mammals) of linkers and adaptors that comprise as much as 1% of total cellular protein. The proteins serve many functions, binding to other cytoplasmic and nuclear proteins through phosphorylated serine or threonine residues in virtually all eukaryotic cells. These binding interactions frequently have regulatory significance in many key cellular processes such as cell cycle progression and apoptosis, cell adhesion, and cell morphology. Aberrations in expression of 14-3-3 isoform expression has been linked to several genetic disorders and to cancer.

| Description                     | Clone     | Species | <b>Applications</b> | Format | Host                | Quantity | Cat. No. |
|---------------------------------|-----------|---------|---------------------|--------|---------------------|----------|----------|
| Monoclonal Antibody             |           |         |                     |        |                     |          |          |
| Anti-14-3-3 γ                   | CG31-2B6  | HMR     | WB                  | Pur    | M lgG <sub>1</sub>  | 100 µg   | 05-639   |
| 14-3-3 Proteins                 | 8C3       | R Sh    | WB                  | Asc    | M lgG <sub>2a</sub> | 500 μL   | MAB3053  |
| Anti-14-3-3 s                   | CS112-2A8 | HMR     | WB IP               | Pur    | M lgG <sub>1κ</sub> | 200 µg   | 05-632   |
| Polyclonal Antibody             |           |         |                     |        |                     |          |          |
| Anti-14-3-3 β/ζ                 |           | Ma In   | WB                  | Pur    | Rb IgG              | 200 µg   | 06-351   |
| Anti-14-3-3 $\zeta/\gamma/\eta$ |           | VVR     | WB                  |        | Rb IgG              | 150 µL   | 06-408   |
| 14-3-3 Proteins                 |           | R       | WB                  | Serum  | Rb                  | 100 µL   | AB1671   |

#### AP-1 (Jun and Fos)

#### Fos

Fos is a component of the AP-1 transcription factor (heterodimer with Jun, or Jun homodimer). Whereas Jun is expressed constitutively but activated following phosphorylation, the Fos promoter contains a serum-response element, and Fos is expressed following mitogenic stimulation, with immediate-early kinetics. AP-1 activity thus differs based upon whether or not the cells are activated to express Fos.

|        | Description   | Clone | Species | Applications | Format | Host   | Quantity | Cat. No. |
|--------|---|-------|---------|--------------|--------|--------|----------|----------|
| Antibo | ody   |       |         |              |        |        |          |          |
|        | Anti-c-fos  | CF2   | Н       | IH not IH(P) | Sup    | M IgM  | 100 μL   | MAB1283  |
|        | Anti-Fos  |       | НМ      | IP WB IC     | APur   | Rb IgG | 100 µg   | 06-341   |
|        | Anti-c-fos  |       |         |              | APur   | Ch     | 100 µg   | AB9116   |
|        | Anti-c-fos  |       | R       | WB IH        | Pur    | Sh     | 500 µg   | AB1584   |
|        |   |       |         | IH(P) Web*   |        |        |          |          |
|        | Anti-c-fos  |       | H R     | WB IC        | Pur    | Sh     | 100 µg   | CBL440   |
| siRNA  |   |       |         |              |        |        |          |          |
|        | Fos SMARTpool® siRNA reagent                        |       | Н       | RNAi         |        |        | 5 nmol   | M-003265 |
|        | Fos siRNA/siAb™ Assay Kit                           |       | Н       | WB RNAi      |        |        | 1 kit    | 60-082   |
| Assay  |   |       |         |              |        |        |          |          |
|        | c-fos Transcription Factor Assay                    |       | H M R   | ACT          |        |        | 1 plate  | 70-545   |
|        | Jun/Fos Transcription Factor Assay                  |       | H M R   | ACT          |        |        | 1 plate  | 70-546   |
|        | AP-1 Family Transcription Factor Assay Colorimetric |       | H M R   | ACT          |        |        | 2 plates | 70-550   |

#### Jun

Jun is a component of the transcription factor AP-1, and can constitute AP-1 activity either as a homodimer or as a heterodimer with Fos. By contrast with Fos, Jun is expressed constitutively, and can be activated by phosphorylation of serines 63 and/or 73. These sites are targets for the SAPK1/JNK family of protein kinases, and their phosphorylation states can be monitored using antibodies specific for the phosphorepitopes. In addition, Jun is phosphorylated by GSK-3 on sites which inhibit Jun activity. Since GSK-3 is inactivated by Akt, the Akt pathway results in Jun activation.

| Description                   | Clo                          | ne  | Species | <b>Applications</b> | Format    | Host                | Quantity | Cat. No. |
|-------------------------------|------------------------------|-----|---------|---------------------|-----------|---------------------|----------|----------|
| Antibody                      |                              |     |         |                     |           |                     |          |          |
| Anti-c-jun                    | 6Ac                          | 6.2 | Н       | WB                  | Pur       | M IgG <sub>2a</sub> | 100 µg   | MAB3732  |
| Anti-phospho-c-Jun (Serz      | 73)                          |     | ΗМ      | WB IC               | Pur       | Rb IgG              | 100 μL   | 05-913   |
| Anti-c-Jun                    |                              |     | Av H M  | IP WB IC            | antiserum | Rb IgG              | 200 μL   | 06-225   |
| Anti-phospho-c-Jun (Ser       | 53)                          |     | HMR     | WB                  | Pur       | Rb IgG              | 100 μL   | 06-828   |
| Anti-phospho-c-Jun (Serz      | 73)                          |     | НМ      | IP WB IC            | APur      | Rb IgG              | 100 μL   | 06-659   |
| Anti-phospho-c-Jun (Thr9      | 1)                           |     | ΗМ      | WB                  | APur      | Sh IgG              | 100 μL   | 07-568   |
| Anti-phospho-c-Jun (Thr9      | 1/Thr93)                     |     | НМ      | WB                  | APur      | Sh IgG              | 100 μL   | 07-570   |
| Protein                       |                              |     |         |                     |           |                     |          |          |
| c-Jun (1-169)-GST, solu       | ole                          |     |         | KA                  |           |                     | 100 µg   | 14-195   |
| siRNA                         |                              |     |         |                     |           |                     |          |          |
| c-Jun siRNA/siAb Assa         | / Kit                        |     | Н       | WB RNAi             |           |                     | 1 kit    | 60-084   |
| c-Jun <i>SMART</i> pool siRN/ | A reagent                    |     | Н       | RNAi                |           |                     | 5 nmol   | m-003268 |
| Assay                         |                              |     |         |                     |           |                     |          |          |
| c-jun Transcription Facto     | or Assay                     |     | H M R   | ACT                 |           |                     | 1 plate  | sgt540   |
| Jun/Fos Transcription Fo      | actor Assay                  |     | HMR     | ACT                 |           |                     | 1 plate  | sgt546   |
| AP-1 Family Transcription     | on Factor Assay Colorimetric |     | H M R   | ACT                 |           |                     | 2 plates | sgt550   |

#### ASK1 (Apoptosis Signaling-regulating Kinase-1)

Apoptosis Signal-Regulating Kinase 1 (ASK1) is a Ser/Thr protein kinase and a member of the MAPKKK family. ASK1 is a component of stress response signaling and activates the JNK and SAPK signaling pathways.

|   | Description            | Clone | Species | <b>Applications</b> | Format | Host   | Quantity | Cat. No. |  |
|---|------------------------|-------|---------|---------------------|--------|--------|----------|----------|--|
| A | Antibody               |       |         |                     |        |        |          |          |  |
|   | Anti-ASK1              |       | Н       | WB                  | Pur    | Rb IgG | 200 µg   | 07-302   |  |
|   | Anti-ASK-1, C-terminus |       | Н       | WB                  | Pur    | Rb IgG | 100 µg   | AB16505  |  |

#### ATF-2

The ATF/CREB family of DNA binding proteins consists of transcription factors whose activity is regulated by specific signaling pathways, including the stress/cytokine/radiation inducible pathways that function through SAPK1/JNK activation. Two threonine residues at positions 69 and 71 in the N-terminal region of ATF-2 are essential in mediating adenovirus E1A-inducible transcriptional activation. These residues are efficiently phosphorylated *in vivo* and are excellent substrates for SAPK1/JNK *in vitro*.

| Description                  | Clone | Species | Applications | Format | Host               | Quantity | Cat. No. |
|------------------------------|-------|---------|--------------|--------|--------------------|----------|----------|
| Antibody                     |       |         |              |        |                    |          |          |
| Anti-phospho-ATF2 (Thr69/71) | AW65  | НМ      | IP WB        | Pur    | M IgG <sub>1</sub> | 200 µg   | 05-891   |
| Anti-ATF2                    |       | Н       | WB           | Pur    | Sh IgG             | 200 µg   | 06-326   |
| Anti-phospho-ATF2 (Thr69)    |       | Н       | KA           | Pur    | Sh IgG             | 200 µg   | 36-016   |
| Anti-phospho-ATF2 (Ser90)    |       | H M R   | WB           | Pur    | Sh IgG             | 200 µg   | 36-015   |
| Protein                      |       |         |              |        |                    |          |          |
| ATF2 (aa 19-96)              |       |         | KA           |        |                    | 1 mg     | 12-367   |
| ATF2, biotin conj.           |       |         | KA           |        |                    | 500 µg   | 12-432   |
| siRNA                        |       |         |              |        |                    |          |          |
| siRNA plasmid, pKD-ATF2-v3   |       |         | RNAi         |        |                    | 5 µg     | 62-172   |

#### BMK/ERK5

Big MAP kinase (BMK1), also known as Erk5, is a member of the MAP kinase family that is activated (phosphorylated on the Thr218 and Tyr220 residues) in cells in response to oxidative stress, hyperosmolarity and treatment with serum. EGF is a potent activator of BMK1, but in contrast to Erk1/2, the EGF-mediated activation of BMK1 occurs independently of Ras and requires Mek5. BMK1 is part of a distinct MAP-kinase pathway involved in EGF-induced cell proliferation and progression through the cell cycle.

|         | Description                            | Clone | Species                    | <b>Applications</b> | Format    | Host   | Quantity | Cat. No. |
|---------|--|-------|----------------------------|---------------------|-----------|--------|----------|----------|
| Antiboo | dy                                     |       |                            |                     |           |        |          |          |
|         | Anti-BMK1/Erk5                         |       | ΗМ                         | IP IPK WB           | antiserum | Rb IgG | 200 μL   | 07-039   |
|         | Anti-MEK5                              |       | Ca H M<br>Mk R Rb<br>Sh Xn | ELISA IP<br>VVB IH  | aPur      | Rb     | 50 µg    | AB3184   |
|         | Anti-phospho-BMK1/Erk5 (Thr218/Tyr220) |       | H M R                      | WB                  | APur      | Rb IgG | 200 μL   | 07-507   |

#### cAMP

Cyclic AMP (cAMP, adenosine 3', 5'-cyclic monophosphate) is a nucleotide that acts as a key second messenger in multiple signal transduction pathways. It is synthesized from ATP by the action of adenylate cyclase, and is inactivated by hydrolysis to 5'-AMP by the actions of phosphodiesterases. cAMP effects are mediated primarily by cAMP-dependent protein kinase (PKA), result in cAMP being responsible for the regulation of many physiological processes such as metabolism, cell growth and differentiation, gene transcription, ion transport and ion channel function.

|        | Description                                 | Clone | Species | Applications | Format | Host | Quantity     | Cat. No. |
|--------|---|-------|---------|--------------|--------|------|--------------|----------|
| Polycl | onal Antibody                               |       |         |              |        |      |              |          |
|        | Anti-cyclic AMP (cAMP)                      |       | А       | ELISA RIA IH | Serum  | Rb   | 100 µL       | 09-176   |
|        | Anti-cyclic AMP (cAMP)                      |       | Α       | ELISA RIA IH | Serum  | Rb   | 50 μL        | AB306    |
|        | Anti-cAMP                                   |       | А       | RIA          | Serum  | Rb   | 1,000 assays | AB505    |
| Assay  |   |       |         |              |        |      |              |          |
|        | cAMP HTS Immunoassay Kit (Chemiluminescent) |       |         | ELISA        |        |      | 192 assays   | 17-418   |

#### Cdc42

cdc42 is a member of the Rho family of small GTP-binding proteins, and plays roles in cytoskeletal Actin organization as well as transformation. Rac is activated downstream of cdc42, or independently by Ras signaling, and in turn Rac stimulates Rho. Activation of cdc42 stimulates filopodia formation, whereas Rac results in the formation of lamellipodia, and Rho stimulates formation of actin stress fibers. Effectors of the cdc42/Rac pathway include NADPH Oxidase (Phox) and the SAPK1/JNK family of kinases. The latter pathway is activated by the PAK kinases, which are direct targets of cdc42.

|       | Description                          | Clone | Species | Applications | Format | Host    | Quantity | Cat. No. |
|-------|--------------------------------------|-------|---------|--------------|--------|---------|----------|----------|
| Antib | ody                                  |       |         |              |        |         |          |          |
|       | Cdc42                                |       |         |              | Pur    | Μ       | 100 µg   | MAB3707  |
|       | Anti-Rac / Cdc42                     |       | BHMR    | IP WB        | Pur    | Rb      | 50 µg    | AB3302   |
|       | Anti-Cdc42                           |       | H M R   | WB           | Pur    | Rb      | 100 µg   | AB4201   |
| Recor | nbinant Protein                      |       |         |              |        |         |          |          |
|       | Cdc42, recombinant human full length |       |         |              | Pur    | E. coli | 20 µg    | SGT211   |

| Descr   | iption                                  | Clone | Species | Applications | Format | Host | Quantity  | Cat. No. |
|---------|---|-------|---------|--------------|--------|------|-----------|----------|
| Assay   |   |       |         |              |        |      |           |          |
| cdc42   | 2 Activation Assay Kit                  |       |         | ABA          |        |      | 1 kit     | 17-286   |
| Anti-co | dc42 Immunoblotting Kit                 |       |         | WB           |        |      | 1 kit     | 17-299   |
| Cdc4    | 2/Rac Activation Assay                  |       |         |              |        |      | 30 assays | SGT445   |
| cDNA    |   |       |         |              |        |      |           |          |
| cdc42   | 2 cDNA (wt) in pUSEamp                  |       |         | TFX          |        |      | 5 µg      | 21-191   |
| cdc42   | 2 cDNA (dominant negative) in pUSEamp   |       |         | TFX          |        |      | 5 µg      | 21-192   |
| cdc42   | 2 cDNA (activated) in pUSEamp           |       |         | TFX          |        |      | 5 µg      | 21-197   |
| Protein |   |       |         |              |        |      |           |          |
| Rac/a   | cdc42 Assay Reagent (PAK1 PBD, agarose) |       |         | ABA          |        |      | 300 µg    | 14-325   |

#### Crk

CRKL belongs to a family of adaptor proteins containing SH2/SH3 domains that couple to effector proteins including p130CAS, SOS, Abl, and Arg. CRKL plays a role in the mediation of Integrin- and Rac-dependent cell migration and in Ras-dependent activation of JNK. The protein is also involved in T-cell activation.

| Description         | Clone | Species | Applications | Format | Host   | Quantity | Cat. No. |
|---------------------|-------|---------|--------------|--------|--------|----------|----------|
| Monoclonal Antibody |       |         |              |        |        |          |          |
| Anti-CrkL           | 5-6   | Н       | IP WB IH     | Pur    | M lgG  | 200 µg   | 05-414   |
| Polyclonal Antibody |       |         |              |        |        |          |          |
| Anti-CrkL           |       | Н       | WB           | Pur    | Rb IgG | 200 µg   | 07-620   |
| Anti-CrkL, NT       |       | Н       | WB           | Pur    | Rb IgG | 200 µg   | 07-621   |
| Protein             |       |         |              |        |        |          |          |
| Crk (120-225)-GST   |       |         | KA           |        |        | 500 µg   | 14-468   |

#### **CREB**

CREB was identified as the transcription factor which binds to cyclic AMP response elements (CRE). Its activity is inducible by phosphorylation, which can be effected by PKA, CaM K IV, and perhaps other protein kinases. Upon activation, CREB dimerizes, binds to DNA, and recruits histone acetyl transferases such as CBP/p300 to chromatin. Phosphorylation state-specific antibodies specific for phosphoserine 133 can provide a surrogate readout of CREB activation.

|        |                                 |       |          |                     |        | -                   |          |          |
|--------|---------------------------------|-------|----------|---------------------|--------|---------------------|----------|----------|
|        | Description                     | Clone | Species  | <b>Applications</b> | Format | Host                | Quantity | Cat. No. |
| Mono   | clonal Antibody                 |       |          |                     |        |                     |          |          |
|        | Anti-CREB                       | NL904 | H M R    | IP WB               | Pur    | Rb IgG              | 200 μL   | 05-767   |
|        | Anti-phospho-CREB (Ser 133)     | 10E9  | H M R    | WB                  | Pur    | M lgG <sub>1κ</sub> | 100 µg   | 05-667   |
|        | Anti-phospho-CREB (Ser 133      | 634-2 | M R      | IF WB IC            | Pur    | M lgG               | 200 µg   | 05-807   |
|        | Anti-CREB                       |       | H M R    | ELISA WB            | Pur    | M lgG <sub>1</sub>  | 100 µg   | MAB5432  |
| Polycl | onal Antibody                   |       |          |                     |        |                     |          |          |
|        | Anti-CREB                       |       | НМ       | IP WB               | Pur    | Rb IgG              | 200 µg   | 06-863   |
|        | Anti-phospho-CREB (Ser 133)     |       | H Ht M R | EMSA IP             | Pur    | Rb IgG              | 200 μL   | 06-519   |
|        |                                 |       |          | WB IH               |        |                     |          |          |
|        | Anti-CREB                       |       | H M R    | ELISA EMSA          | Pur    | Rb                  | 100 µL   | AB3006   |
|        |                                 |       |          | IP WB               |        |                     |          |          |
|        | Anti-CREB, pSer133              |       | H M Mk R | WB                  | Serum  | Rb                  | 100 µL   | AB3442   |
| Peptid | les                             |       |          |                     |        |                     |          |          |
|        | CREB (123-136)                  |       |          | KA                  |        |                     | 50 µg    | 12-490   |
|        | CREB Immunizing Peptide         |       |          | WB                  |        |                     | 50 µg    | 12-377   |
|        | phospho CREB Immunizing Peptide |       |          | WB                  |        |                     | 50 µg    | 12-378   |
| siRNA  | 4                               |       |          |                     |        |                     |          |          |
|        | CREB SMARTpool siRNA reagent    |       | Н        | RNAi                |        |                     | 5 nmol   | M-003619 |
|        | CREB siRNA/siAb Assay Kit       |       | Н        | WB RNAi             |        |                     | 1 kit    | 60-012   |
|        | siRNA plasmid, pKD-CREB-v2      |       | Н        | RNAi                |        |                     | 5 µg     | 62-011   |
| Assay  | ,                               |       |          |                     |        |                     |          |          |
|        | CREB Transcription Factor Assay |       |          | ACT                 |        |                     | 1 plate  | 70-575   |
|        |                                 |       |          |                     |        |                     |          |          |

#### Elk

Elk-1 is a member of the ternary complex factor (TCF) family of transcription factors. It is phosphorylated and activated by the MAP kinase pathway in the C-terminal transcription factor domain on Serine 383. This phosphorylation is critical to it activation.


| Description                           | Clone | Species | Applications | Format | Host | Quantity | Cat. No. |
|---------------------------------------|-------|---------|--------------|--------|------|----------|----------|
| Polyclonal Antibody                   |       |         |              |        |      |          |          |
| Anti-ELK-1, phospho-specific (Ser383) |       | Н       | WB           | APur   | Rb   | 100 µL   | AB3809   |

### **FGFR**

Epidermal growth factor receptor (EGFR) is one of four members of the EGFR family of receptor tyrosine kinases. Signaling through these receptors is induced by binding of EGF to the extracellular domain, which results in receptor homo- and hetero-dimerization, trans- and autophosphorylation of tyrosine residues within the cytoplasmic domain and ultimately to the activation of the classical MAPK signaling pathway. EGFR plays a key role in the regulation of essential normal cellular processes and in the pathophysiology of hyperproliferative diseases such as cancer. It is known to be essential for mediation of both proliferative and survival signals to cells. Activation of the EGFR

signaling pathway has been linked with increased cell proliferation, angiogenesis, metastasis and decreased apoptosis. EGFR is found in most solid tumors and its autophosphorylation activates the intrinsic kinase activity toward heterologous substrates, as well as creating docking sites for adapter proteins for multiple adapter proteins to bind and nucleate signaling complexes that activate the Ras, PI3 Kinase, and PLC pathways. Many EGFR mutations are directly related to various diseases and thus are the targets of numerous drugs.

Millipore offers numerous products for the study of EGFR, particularly in the areas of overexpression in malignancies such as lung cancer. For instance, the amino acid substitution L858R is one of several




heterozygous mutations that have been identified in Non-Small-Cell Lung Cancer (NSCLC) patients who have clinical responses to the EGFR inhibitor Iressa® preparation. There is some evidence that these mutations result in elevated activity and enhanced sensitivity to Iressa preparation. In patients with tumors bearing Iressa-sensitive mutations, resistant subclones containing an additional EGFR mutation, T790M, emerge in the presence of the drug. It has been shown experimentally that the T790M mutation leads to high-level functional resistance to Iressa preparation. In patients with tumors bearing Iressa-sensitive mutations (eg. L858R, L861Q), resistant subclones containing the T790M mutation emerge in the presence of the drug.

| Description                         | Quantity | Cat. No. | Description                                 | Quantity   | Cat. No. |
|-------------------------------------|----------|----------|---|------------|----------|
| Anti-EGFR, neutralizing             | 250 µg   | 05-101   | Anti-phospho-EGFR (Tyr 1 0 68)              | 100 µg     | 04-285   |
| Anti-EGFR                           | 250 µg   | 05-104   | Anti-phospho-EGFR (Tyr 1 086)               | 100 pL     | 04-340   |
| Anti-EGFR, Rb Mab                   | 100 μL   | 04-338   | Anti-phospho-EGFR (Tyr 1 086)               | 100 µg     | 04-286   |
| Anti-EGFR, Rb Mab                   | 100 μL   | 04-337   | Anti-phospho-EGFR (Tyr 1 173), Rb Mab       | 100 μL     | 04-341   |
| Anti-EGFR (aa 1140-1160)            | 100 µg   | 04-287   | Anti-phospho-EGFR (Tyr 1 173)               | 50 µg      | 05-483   |
| Anti-EGFR (extracellular domain)    | 100 µg   | 04-289   | Anti-EGFR (non-phospho-Tyr 1 173)           | 50 µg      | 05-484   |
| Anti-EGFR (C-Terminus)              | 100 µg   | 04-290   | Anti-Epidermal Growth Factor                | 200 μL     | MAB126   |
| Anti-phospho-EGFR (Thr654)          | 100 µg   | 04-282   | Anti-EGF Receptor, intracellular, activated | 100 µg     | MAB3052  |
| Anti-phospho-EGFR (Thr669)          | 100 µg   | 04-281   | Anti-EGFR, extracellular, EGF binding site  | 50 µg      | MAB88910 |
| Anti-phospho-EGFR (Tyr845)          | 100 µg   | 04-283   | Anti-EGFR, N-terminal                       | 100 µg     | CBL416   |
| Anti-phospho-EGFR (Tyr 1 045)       | 100 µg   | 04-284   | Anti-EGFR, N-terminal                       | 100 assays | CBL416F  |
| Anti-phospho-EGFR (Ser 1047)        | 100 µg   | 04-280   | Anti-EGFR, Cytosolic                        | 100 µg     | CBL417   |
| Anti-phospho-EGFR (Tyr 1068), RbMab | 100 μL   | 04-339   | Anti-EGFR, Cytosolic                        | 100 assays | CBL417F  |

#### EGF/EGFR (erbB-1) (Epidermal Growth Factor Receptor)


EGF exerts its signals through the EGF Receptor and its hetero-oligomerization partners, members of the ErbB family. The Receptors are transmembrane tyrosine kinases which are activated as a consequence of ligand-induced dimerization. Every EGF molecule has two receptor-binding sites, so the activated ligand-receptor complex consists of two molecules of EGF and two receptor molecules. Receptor autophosphorylation activates the intrinsic kinase activity towards heterologous substrates, as well as creating docking sites for adapter proteins to bind and nucleate signaling complexes that activate the Ras, P13 Kinase, and PLC pathways. Overexpression of the EGF Receptor is typical in most solid tumors.

|   | Description             | Clone | Species | Applications | Format | Host               | Quantity | Cat. No. |
|---|-------------------------|-------|---------|--------------|--------|--------------------|----------|----------|
| I | Monoclonal Antibody     |       |         |              |        |                    |          |          |
|   | Anti-EGFR, neutralizing | LA1   | Н       | IP WB        | Pur    | M lgG <sub>1</sub> | 250 µg   | 05-101   |
|   |                         |       |         | IC NEUT      |        |                    |          |          |



#### Immunofluorescence Analysis


Left: A431 cells were stained with 2  $\mu$ g/mL of anti-EGF Receptor (neutralizing), clone LA1, Alexa Fluor 488 conjugate (green) and DAPI (blue).



#### **Immunoblot Analysis**


Right: Non-reduced A431 cell lysate probed with anti-EGF Receptor (1 µg/mL)

| Anti-EGFR                             | LA22 | Н  | IP WB IC Pur   | M IgG <sub>2a</sub> | 250 µg | 05-104 |
|---------------------------------------|------|----|----------------|---------------------|--------|--------|
| Anti-EGFR, Rb Mab                     |      | Н  | WB IP IC FC    | Rb IgG              | 100 pL | 04-338 |
| Anti-EGFR, Rb Mab                     |      | Н  | WB ELISA       | Rb IgG              | 100 µL | 04-337 |
| Anti-EGFR (aa 1140-1160)              |      | Н  | WB ELISA       | M lgG <sub>1κ</sub> | 100 µg | 04-287 |
| Anti-EGFR (extracellular domain)      |      | Н  | WB ELISA       | M lgG <sub>1κ</sub> | 100 µg | 04-289 |
| Anti-EGFR (C-Terminus)                |      | Н  | WB ELISA       | M lgG <sub>1κ</sub> | 100 µg | 04-290 |
| Anti-phospho-EGFR (Thr654)            |      | Н  | WB ELISA       | M lgG <sub>1κ</sub> | 100 µg | 04-282 |
| Anti-phospho-EGFR (Thr669)            |      | Н  | WB ELISA       | M lgG <sub>1κ</sub> | 100 µg | 04-281 |
| Anti-phospho-EGFR (Tyr845)            |      | Н  | WB ELISA       | M lgG <sub>1κ</sub> | 100 µg | 04-283 |
| Anti-phospho-EGFR (Tyr 1 045)         |      | Н  | WB ELISA       | M lgG <sub>1κ</sub> | 100 µg | 04-284 |
| Anti-phospho-EGFR (Ser 1 047)         |      | Н  | WB ELISA       | M lgG <sub>1κ</sub> | 100 µg | 04-280 |
| Anti-phospho-EGFR (Tyr 1 068), Rb/Mab |      | Н  | WB ELISA       | Rb IgG              | 100 pL | 04-339 |
| Anti-phospho-EGFR (Tyr 1 068)         |      | Н  | WB ELISA       | M lgG <sub>1κ</sub> | 100 µg | 04-285 |
| Anti-phospho-EGFR (Tyr 1 086)         |      | Н  | WB ELISA       | Rb IgG              | 100 pL | 04-340 |
| Anti-phospho-EGFR (Tyr 1 086)         |      | Н  | WB ELISA       | M lgG1k             | 100 µg | 04-286 |
| Anti-phospho-EGFR (Tyr 1 173), Rb Mab |      | Н  | WB IP IF IH(P) | Rb IgG              | 100 pL | 04-341 |
| Anti-phospho-EGFR (Tyr 1 173)         | 9H2  | ΗМ | IP WB FC IF    | M lgG1k             | 50 µg  | 05-483 |



#### Immunofluorescence Analysis

Left: A431 cells were untreated (Top) and treated (Bottom) with EGF then stained with 2  $\mu$ g/mL of anti-phospho-EGFR (Tyr1173), clone 9H2, Alexa Fluor 555 conjugate (red) and DAPI (blue).



#### Immunoblot Analysis

Right: 20 micrograms of EGF-stimulated (lanes 1 and 3) and unstimulated (lanes 2 and 4) HFF cell lysate was probed with Anti-EGF Receptor (Cat. No. 06-129, lanes 1 and 2) or Anti-phospho-EGFR (Tyr1173), clone 9H2 (lanes 3 and 4) at 0.5 µg/mL.

| Anti-EGFR (non-phospho-Tyr1 173)            | 20G3   | ΗМ | IP WB                     |      | $\text{M IgG}_{1\kappa}$ | 50 µg      | 05-484   |  |
|---|--------|----|---------------------------|------|--------------------------|------------|----------|--|
| Anti-Epidermal Growth Factor                | 144-8  | ΗМ | IP IH                     | SPur | M lgG <sub>1</sub>       | 200 μL     | MAB126   |  |
| Anti-EGF Receptor, intracellular, activated | 74     | Н  | ELISA VVB<br>IC IH        | Pur  | M lgG₁                   | 100 µg     | MAB3052  |  |
| Anti-EGFR, extracellular, EGF binding site  | EGFR-1 | Н  | ELISA FC IP<br>RIA IH INF |      | M lgG₁                   | 50 µg      | MAB88910 |  |
| Anti-EGFR, N-terminal                       | EGFR1  | Н  | FC IP IC                  | Pur  | M lgG <sub>2b</sub>      | 100 µg     | CBL416   |  |
| Anti-EGFR, N-terminal                       | EGFR1  | Н  | FC IF                     | FITC | M lgG <sub>2b</sub>      | 100 assays | CBL416F  |  |

|         | Description  | Clone | Species        | Applications     | Format | Host                     | Quantity   | Cat. No. |
|---------|--|-------|----------------|------------------|--------|--------------------------|------------|----------|
|         | Anti-EGFR, Cytosolic                                     | F4    | Eq H           | ELISA IP WB      | Pur    | M lgG <sub>1</sub>       | 100 µg     | CBL417   |
|         |  |       |                | IC IH IH(P)      |        |                          |            |          |
|         | Anti-EGFR, Cytosolic                                     | F4    | Eq H           | IC IH IH(P)      | FITC   | M lgG <sub>1</sub>       | 100 assays | CBL417F  |
|         | Anti-c-erbB-1, a.a. 580-591                              | C11   | Н              | WB IH(P)         | Pur    | M lgG <sub>1</sub>       | 100 µg     | CBL754   |
| Polyclo | onal Antibody  |       |                |                  |        |                          |            |          |
|         | Anti-EGF, neutralizing                                   |       | Μ              | WB NEUT          |        | Rb IgG                   | 1 mg       | 06-102   |
|         | Anti-EGFR  |       | H M R<br>Ca Ch | IP WB<br>IC IH   |        | Sh IgG                   | 200 µg     | 06-129   |
|         | Anti-EGFR  |       | НМ             | IP WB            |        | Rb IgG                   | 200 µg     | 06-847   |
|         | Anti-Epidermal Growth Factor                             |       | Н              | RIA IH           | Serum  | Rb                       | 100 µL     | AB1910   |
|         | Anti-phospho-EGFR (Tyr845)                               |       | Н              | WB IF IHC,       | APur   | Rb IgG                   | 100 µL     | 07-820   |
|         | Anti-phospho-EGFR (Tyr992)                               |       | Н              | WB               | APur   | Rb IgG                   | 100 µL     | 07-821   |
|         | Anti-phospho-EGFR (Tyr 1069)                             |       | Н              | WB IP            | APur   | Rb IgG                   | 100 µL     | 07-715   |
|         | Anti-phospho-EGFR (Tyr1086)                              |       | Н              | WB               | APur   | Rb IgG                   | 100 µL     | 07-818   |
|         | Anti-phospho-EGFR (Tyr1 148)                             |       | Н              | WB               | APur   | Rb IgG                   | 100 µL     | 07-819   |
| Antibo  | dy Conjugates  |       |                |                  |        |                          |            |          |
|         | Anti-EGFR, neutralizing, Alexa Fluor® 488 Conjugate      | LA1   | Н              | FC WB IC NEUT    | Pur    | M lgG₁                   | 125 µg     | 16-246   |
|         | Anti-EGFR, neutralizing, Alexa Fluor 555 Conjugate       | LA1   | Н              | FC WB<br>IC NEUT | Pur    | M lgG <sub>1</sub>       | 125 µg     | 16-247   |
|         | Anti-phospho-EGFR (Tyr 1 173), Alexa Fluor 488 Conjugate | 9H2   | НМ             | FC IF            |        | M lgG <sub>1κ</sub>      | 25 µg      | 16-244   |
|         | Anti-phospho-EGFR (Tyr1173), Alexa Fluor 555 Conjugate   | 9H2   | НМ             | FC IF            |        | $\text{M IgG}_{1\kappa}$ | 25 µg      | 16-245   |
| Kinase  | , 9  |       |                |                  |        |                          |            |          |
|         | EGFR, active   |       |                |                  |        |                          | 10 µg      | 14-531   |
|         | EGFR (L858R), active                                     |       |                |                  |        |                          | 10 µg      | 14-626   |
|         | EGFR (L861Q), active                                     |       |                |                  |        |                          | 10 µg      | 14-627   |
|         | EGFR (T790M), active                                     |       |                |                  |        |                          | 10 µg      | 14-725   |
|         | EGFR (T790M, L858R), active                              |       |                |                  |        |                          | 10 µg      | 14-721   |
| Growt   | h Factor   |       |                |                  |        |                          |            |          |
|         | EGF, culture grade                                       |       |                | CULT             |        |                          | 100 µg     | 01-101   |
|         | EGF, receptor grade                                      |       |                | CULT             |        |                          | 100 µg     | 01-102   |
|         | EGF, human recombinant                                   |       |                | CULT             |        |                          | 100 µg     | 01-107   |
|         | EGF  |       |                | CULT             |        |                          | 500 μg     | 01-407   |
|         | Epidermal Growth Factor, mouse tissue culture grade      |       |                |                  | SPur   |                          | 100 µg     | EA140    |
|         | Epidermal Growth Factor, recombinant human               |       |                |                  | Pur    |                          | 500 µg     | GF001    |
|         | Epidermal Growth Factor, recombinant mouse               |       |                |                  | Pur    |                          | 500 µg     | GF123    |
| cDNA    |  |       |                |                  |        |                          |            |          |
|         | EGFR cDNA (wt) in pUSEamp                                |       |                | TFX              |        |                          | 5 μg       | 21-176   |
| siRNA   |  |       |                |                  |        |                          |            |          |
|         | EGFR siRNA/siAb Assay Kit                                |       | Н              | WB RNAi          |        |                          | 1 kit      | 60-015   |
|         | EGFR SMARTpool siRNA reagent                             |       | Н              | RNAi             |        |                          | 5 nmol     | M-003114 |

#### erbB-2/HER2

v-ErbB was identified as the viral counterpart of the EGF Receptor, and its relatives have adopted this nomenclature (or the HER nomenclature for human proteins). These proteins are activated by Heregulin, Neuregulin, and Amphiregulin. All the family members are capable of heterodimerization, and together activate downstream signaling pathways. erbB-2/HER-2/neu is amplified in highly aggressive breast cancers. ErbB-3 has a natural "substitution" in an invariant Lys residue in the ATP-binding site. It is therefore inactive as a kinase, and can only become phosphorylated by heterodimerization with another erbB-family member. All members of this family, when Tyrphosphorylated, can nucleate the formation of signaling complexes through adapter proteins.

| Description                                  | Clone | Species  | Applications Form | nat Host | Quantity | Cat. No. |
|--|-------|----------|-------------------|----------|----------|----------|
| Antibody                                     |       |          |                   |          |          |          |
| Anti-erbB2 (intracellular domain aa 860-880) |       | Н        | WB ELISA          | M lgG    | 100 µg   | 04-291   |
| Anti-erbB2 (phospho-Ser 1113)                |       | Н        | WB ELISA          | M lgG    | 100 µg   | 04-292   |
| Anti-erbB2 (phospho-Thr 686)                 |       | Н        | WB ELISA          | M lgG    | 100 µg   | 04-293   |
| Anti-erbB2 (phospho-Tyr 1112)                |       | Н        | VVB ELISA         | M lgG    | 100 µg   | 04-294   |
| Anti-erbB-2/HER-2                            |       | H M R MK | WB IP IH Pur      | Rb IgG   | 200 µg   | 06-562   |

| Description                             | Clone | Species | Applications | Format | Host   | Quantity | Cat. No. |
|---|-------|---------|--------------|--------|--------|----------|----------|
| Anti-phopsho-erbB-2 / HFR-2 (Tyr 1 248) |       | H R     | WB FC IC     | APur   | Rb IaG | 100 ul   | 06-229   |



#### Immunofluorescence Analysis

Left: 2  $\mu$ g/mL of this lot showed positive immunostaining for phospho-erbB2/HER-2 in EGF treated A431 cells.



#### Immunoblot Analysis

Right: Lysates from non-stimulated (lane 1) and heregulin-stimulated (lane 2) MCF-7 cells were probed with Anti-phospho-Erb B2/HER-2 (Tyr1248) (1:2000 dilution).

|                | nopsho-erbB-2/HER-2 (Tyr1248),   |      | H R | FC IC   | APur | Rb IgG | 100 µg | 16-235   |
|----------------|----------------------------------|------|-----|---------|------|--------|--------|----------|
|                | a Fluor 488 Conjugate            |      |     |         |      |        |        |          |
| Anti-ph        | nopsho-erbB-2/HER-2 (Tyr 1 248), |      | H R | FC IC   | APur | Rb IgG | 50 µg  | 16-236   |
| Alex           | a Fluor 555 Conjugate            |      |     |         |      |        |        |          |
| Anti-c-        | erbB-2                           | CB11 | Н   | IH(P)   | Sup  | M lgG1 | 200 μL | MAB1282  |
| Anti-c-        | erbB-2                           |      | H R | WB IH   | Pur  | Sh     | 100 µg | AB1368   |
| erbB-2         | ? siRNA/siAb Assay Kit           |      | Н   | WB RNAi |      |        | 1 kit  | 60-016   |
| erbB-2         | ? SMARTpool siRNA reagent        |      | Н   | RNAi    |      |        | 5 nmol | M-003126 |
| erbB-3/HER3    |                                  |      |     |         |      |        |        |          |
| Monoclonal A   |                                  |      |     |         |      |        |        |          |
| IVIONOCIONAL A | ANTIDOGV                         |      |     |         |      |        |        |          |

| Monoclonal Antibody                   |           |          |                                 |        |                             |        |         |
|---------------------------------------|-----------|----------|---------------------------------|--------|-----------------------------|--------|---------|
| Anti-erbB-3/HER-3                     | 2F12      | BHMR     | IP WB IH                        |        | M lgG <sub>2α<b>κ</b></sub> | 200 µg | 05-390  |
| Anti-erbB-3/HER-3                     | H3.105.5  | НМ       | IP NEUT                         |        | M lgG <sub>1</sub>          | 100 µg | 05-471  |
| Anti-c-erbB-3                         | RTJI/2E11 | Н        | IP WB IH(P)                     | Pur    | M IgM                       | 100 µg | MAB4021 |
| Anti-c-erbB-3, cytoplasmic domain     | RTJ2      | Н        | ELISA IP<br>WB IH               | Pur    | M lgG <sub>1</sub>          | 100 µg | MAB4023 |
| Anti-c-erbB-2, a.a. 1238-1255         | N3/D10    | Н        | FC WB<br>IC IH(P)               | Pur    | M lgG <sub>1</sub>          | 100 µg | CBL755  |
| Anti-c-erbB-2, extracellular domain   | 300G9     | Н        | IP WB IH IH(                    | P) Pur | M lgG1                      | 100 µg | CBL772  |
| Anti-c-erbB-3                         | SGP. 1    | Н        | ELISA FC IP IH not WB not IH(P) | Pur    | M lgG₁                      | 100 µg | MAB406  |
| Polyclonal Antibody                   |           |          |                                 |        |                             |        |         |
| Anti-phospho-erbB-2/HER-2 (Tyr 1 248) |           | Н        | WB                              |        | Rb IgG                      | 100 μL | 06-229  |
| Anti-erbB-2/HER-2                     |           | H M Mk R | IP WB IH                        | Rb IgG | 200 µg                      | 06-562 |         |
| Anti-erbB-4/HER-4                     |           | Н        | IP WB                           |        | Rb IgG                      | 200 µg | 06-572  |

#### erbB-4/HER4

v-ErbB was identified as the viral counterpart of the EGF Receptor, and its relatives have adopted this nomenclature (or the HER nomenclature for human proteins). These proteins are activated by Heregulin, Neuregulin, and Amphiregulin. All the family members are capable of heterodimerization, and together activate downstream signaling pathways. erbB-2/HER-2/neu is amplified in highly aggressive breast cancers. ErbB-3 has a natural "substitution" in an invariant Lys residue in the ATP-binding site. It is therefore inactive as a kinase, and can only become phosphorylated by heterodimerization with another erbB-family member. All members of this family, when Tyrphosphorylated, can nucleate the formation of signaling complexes through adapter proteins.

| Description                    | Clone    | Species | Applications Format | Host                | Quantity | Cat. No. |
|--------------------------------|----------|---------|---------------------|---------------------|----------|----------|
| Antibody                       |          |         |                     |                     |          |          |
| Anti-erbB-4/HER-4              | H4.72.8  | Н       | FC IP NEUT          | M lgG <sub>2a</sub> | 100 µg   | 05-478   |
| erbB-4 siRNA/siAb Assay Kit    |          | Н       | WB RNAi             |                     | 1 kit    | 60-080   |
| erbB-4 SMARTpool siRNA reagent |          | Н       | RNAi                |                     | 5 nmol   | M-003128 |
| Anti-c-erbB-4                  | HFR1/2G4 | Н       | IP WB IH(P) Pur     | M laGal             | 100 ua   | MAB4025  |

#### Erk/MAP Kingse


MAP Kinases, or Erks, are Ser/Thr kinases which phosphorylate PXS/TP motifs in many different proteins, and are activated principally in response to mitogenic stimulation. MAP Kinases comprise MAP kinase 1/Erk1, or p44 MAP kinase, and MAP kinase 2/Erk2, or p42 MAP kinase. Both are activated by MEK1 or MEK2, by dual phosphorylation of a threonine and tyrosine residue in the activation loop (TEY motif). Either phosphorylation alone can induce an electrophoretic mobility shift, but both are required for activation of the kinase. This dual phosphorylation is efficiently detected by phosphorylation state-specific antibody directed to the pTEPY motif. No known mutations exist which can lead to constitutive activation of MAP kinases. Once activated, MAP kinases phosphorylate a broad spectrum of substrates, including cytoskeletal proteins, translation regulators, transcription factors, and the Rsk family of protein kinases.

|   | Description            | Clone | Species  | <b>Applications</b> | Format | Host                | Quantity | Cat. No.       |
|---|------------------------|-------|----------|---------------------|--------|---------------------|----------|----------------|
| I | Monoclonal Antibody    |       |          |                     |        |                     |          |                |
|   | Anti-MAP Kinase 2/Erk2 | 1B3B9 | Av H M R | IP VVB              | Pur    | M lgG <sub>2a</sub> | 200 µg   | 05-1 <i>57</i> |

|        | Description                           | Clone | Species   | Applications | Format | Host               | Quantity | Cat. No. |
|--------|---------------------------------------|-------|-----------|--------------|--------|--------------------|----------|----------|
|        | Anti-phospho-MAP Kinase 1/2 (Erk 1/2) | 12D4  | Av H M R  | IP WB        | Pur    | M lgG $_{1\kappa}$ | 50 µg    | 05-481   |
|        | Anti-phospho-MAP Kinase 1/2 (Erk1/2)  | AW39  | H R       | BD WB        | Pur    | Rb IgG             | 100 μL   | 05-797   |
|        | (Thr 185/Tyr 187)                     |       |           |              |        |                    |          |          |
|        | Anti-MAP Kinase 1/Erk 1, CT           |       | M R       | IPK WB       | Pur    | Rb IgG             | 100 μL   | 05-957   |
| Polycl | onal Antibody                         |       |           |              |        |                    |          |          |
|        | Anti-MAP Kinase 1/2 (Erk1/2), CT      |       | Av H M R  | IP WB IC     | APur   | Rb IgG             | 100 µg   | 06-182   |
|        |                                       |       | Sh Xn Ech |              |        |                    |          |          |



# 



|                     | Immunoblot Analysis     |
|---------------------|-------------------------|
| Right: NIH 3T3 cell | lysate probed with      |
| anti-MAP Kinase     | $1/2$ (0.5 $\mu$ g/mL). |

| Anti-MAP Kinase 2/Erk2                        | H M R     | IP WB       | Pur  | Rb IgG | 200 µg  | 06-333          |
|---|-----------|-------------|------|--------|---------|-----------------|
| Anti-phospho-MAP Kinase 1/2 (Erk 1/2)         | H M R     | IP WB       | APur | Sh IgG | 100 µg  | 06-642          |
| (Tyr180), polyclonal                          |           |             |      |        |         |                 |
| Anti-phospho-MAP Kinase 1/2 (Erk 1/2)         | H R WR    | WB          | Pur  | Rb IgG | 200 µg  | 07-467          |
| Anti-ERK1/2                                   | HMRB      | ELISA IP    | Pur  | Rb     | 50 µg   | AB3053          |
|   | Ch Dr Sh  | WB IC       |      |        |         |                 |
|   | Xn Ech Ml |             |      |        |         |                 |
| Anti-ERK1                                     | BHM       | ELISA IP WB | APur | Rb     | 100 µg  | AB3189          |
|   | Po R Rb   |             |      |        |         |                 |
| Anti-ERK1/2, phospho-specific (Thr202/Tyr204) | H M R Xn  | WB          | APur | Rb     | 100 pL  | AB3826          |
| Antibody Conjugate                            |           |             |      |        |         |                 |
| Anti-MAP Kinase 1/2 (Erk1/2), agarose         |           | IAP IP      |      |        | 50 µg   | 16-111          |
| Assay   |           |             |      |        |         |                 |
| MAP Kinase/Erk Assay Kit                      |           | KA          |      |        | 1 kit   | 17-133          |
| MAP Kinase/Erk Sampler Pack                   |           | KA          |      |        | 1 kit   | 17-171          |
| MAP Kinase/Erk Assay Kit, non-radioactive     |           | KA          |      |        | 1 kit   | 1 <i>7</i> -191 |
| MAP Kinase/Erk Immunoprecipitation Kinase     |           | KA          |      |        | 1 kit   | 1 <i>7</i> -192 |
| Assay Kit, non-radioactive                    |           |             |      |        |         |                 |
| MAP Kinase (ERK1/2) Activity Assay            |           | ACT         |      |        | 1 plate | SGT415          |
| MAP Kinase/Erk Substrate Cocktail I           |           | KA          |      |        | 1 mL    | 20-115          |
| MAP Kinase/Erk Substrate Cocktail II          |           | KA          |      |        | 1 mL    | 20-166          |
| MAP Kinase 2/Erk2 cDNA (wt) in pUSEamp        |           | TFX         |      |        | 50 µg   | 21-112          |
| MAP Kinase 1/Erk1 KinEASE™ FP                 |           | KA          |      |        | 1 kit   | 32-044          |
| Fluorescein Green Assay                       |           |             |      |        |         |                 |
| MAP Kinase 2/Erk2 KinEASE FP                  |           | KA          |      |        | 1 kit   | 32-045          |
| Fluorescein Green Assay                       |           |             |      |        |         |                 |
| MAP Kinase 1/Erk1 KinEASE                     |           | KA          |      |        | 1 kit   | 32-124          |
| FP-645nm FarRed Assay                         |           |             |      |        |         |                 |
| MAP Kinase 2/Erk2 KinEASE                     |           | KA          |      |        | 1 kit   | 32-125          |
| FP-645nm FarRed Assay                         |           |             |      |        |         |                 |
| rotein  |           |             |      |        |         |                 |
| MAP Kinase 2/Erk 2 (K52R), inactive           |           | KA          |      |        | 100 µg  | 14-696          |
| MAP Kinase 1/Erk1 agarose, unactive           |           | KA          |      |        | 250 µg  | 14-121          |
| MAP Kinase 2/Erk2, active                     |           | KA          |      |        | 10 µg   | 14-550          |
| MAP Kinase 1/Erk1, active                     |           | KA          | Pur  |        | 10 µg   | 14-439          |
| MAP Kinase 2/Erk2, unactive                   |           | KA          | Pur  |        | 50 µg   | 14-198          |
|   |           |             |      |        |         |                 |

|      | Description                               | Clone | Species | Applications Format | Host | Quantity | Cat. No. |
|------|---|-------|---------|---------------------|------|----------|----------|
| siRN | A   |       |         |                     |      |          |          |
|      | MAP Kinase 1/Erk2 siRNA/siAb Assay Kit    |       | Н       | WB RNAi             |      | 1 kit    | 60-029   |
|      | siRNA plasmid, pKD-MAP Kinase 1/Erk2-v1   |       | HM      | RNAi                |      | 5 µg     | 62-046   |
|      | siRNA plasmid, pKD-MAP Kinase/Erk1-v4     |       | Н       | RNAi                |      | 5 µg     | 62-149   |
|      | siRNA plasmid, pKD-MAP Kinase 1/Erk2-v5   |       |         | RNAi                |      | 5 µg     | 62-191   |
|      | siRNA plasmid, pKD-MAP Kinase/Erk1-v3     |       |         | RNAi                |      | 5 µg     | 62-192   |
|      | MAP Kinase 1/Erk2 SMARTpool siRNA reagent |       | Н       | RNAi                |      | 5 nmol   | M-003555 |

#### **FGF**

Fibroblast Growth Factors (FGFs) are mitogens for a wide variety of fibroblast and epithelial cells. They signal through receptor tyrosine kinases of the FGF Receptor (FGFR) family. The extracellular portions of these receptors have lg-like loops, and mutations occurring in these loops have been associated with developmental disorders. Disruption of disulfide bonds in the Ig-loops allows for adjacent receptor molecules to form inter-molecular disulfide bonds, dimerizing the receptor and activating it in the absence of ligand. FGF-2 (basic FGF) is a potent inducer of angiogenesis, and acts by stimulating VEGF expression in endothelial cells.

| •      |  |                |         | •            |        |                     |          |          |
|--------|--|----------------|---------|--------------|--------|---------------------|----------|----------|
|        | Description  | Clone          | Species | Applications | Format | Host                | Quantity | Cat. No. |
| Mond   | oclonal Antibody                                   |                |         |              |        |                     |          |          |
|        | Anti-FGF-2/basic FGF                               | bFM-1          | BHMR    | ria Neut     |        | M $\lg G_{1\kappa}$ | 500 µg   | 05-117   |
|        |  | (neutralizing) |         |              |        |                     |          |          |
|        | Anti-FGF-2/basic FGF                               | bFM-2          | BHMR    | RIA WB IH    |        | M lgG $_{1\kappa}$  | 500 µg   | 05-118   |
|        | Anti-FGFR1   | 19B2           | H R     | IP WB        |        | M IgG <sub>1</sub>  | 250 µg   | 05-149   |
|        | Anti-Fibroblast Growth Factor basic                |                | Н       | IF VVB       | Pur    | M IgG <sub>1</sub>  | 150 µg   | MAB120   |
|        | Anti-Fibroblast Growth Factor Receptor             | VBS1           | HMR     | ELISA IP WB  | Pur    | M IgM               | 100 µg   | MAB125   |
|        |  |                | B Ch    | IH(P) INHIB  |        |                     |          |          |
| Polyc  | lonal Antibody                                     |                |         |              |        |                     |          |          |
|        | Anti-Fibroblast Growth Factor basic                |                | Н       | ELISA WB     | Pur    | Rb                  | 100 µg   | AB1435   |
|        | Anti-Fibroblast Growth Factor basic                |                | Ma      | WB IC IH(P)  | Serum  | Rb                  | 100 μL   | AB1458   |
|        | Anti-Fibroblast Growth Factor basic                |                | Ma      | WB IC IH(P)  | Serum  | Rb                  | 100 μL   | AB1459   |
|        | Anti-Fibroblast Growth Factor basic                |                | HMR     | IH           | Serum  | Sh                  | 100 μL   | AB5396   |
|        | Anti-Fibroblast Growth Factor basic                |                | HMR     | IH           | APur   | Sh                  | 50 µg    | AB5396P  |
| Protei | n  |                |         |              |        |                     |          |          |
|        | FGF-2 / basic FGF, human recombinant               |                |         | CULT         |        |                     | 25 µg    | 01-106   |
|        | FGF-2/basic FGF, carrier-free                      |                |         | CULT         |        |                     | 25 µg    | 01-114   |
|        | FGF-1/acidic FGF, human recombinant                |                |         | CULT         |        |                     | 25 µg    | 01-116   |
|        | FGF-7/KGF, human recombinant                       |                |         | CULT         |        |                     | 10 µg    | 01-118   |
|        | Fibroblast Growth Factor basic peptide,            |                |         |              | Pur    |                     | 1 mg     | FA011    |
|        | Synthetic, brain derived (1-24)                    |                |         |              |        |                     |          |          |
|        | Fibroblast Growth Factor acidic, recombinant human |                |         |              | Pur    |                     | 50 µg    | GF002    |
|        | Fibroblast Growth Factor basic, recombinant human  |                |         |              | Pur    |                     | 50 µg    | GF003    |
|        | Fibroblast Growth Factor-4, recombinant human      |                |         |              | Pur    |                     | 25 µg    | GF098    |
|        |  |                |         |              |        |                     |          |          |

#### FNTA/FNTB (Farnesyl Transferase a and b subunits)

Many membrane-associated proteins must undergo a complicated post-translational processing in the endoplasmic reticulum in order to become active. This processing occurs on a c-terminal CAAX motif in members of the Ras/Rho family of G-Proteins, the  $\gamma$  subunit of heterotrimeric G-Proteins, and nuclear Lamin. The first step in processing is prenylation (farnesylation by Farnesyl Transferase) of the Cys residue. Prenylated proteins are then cleaved (by Rce 1 protease) to remove the amino acids following the prenylated cysteine. The last step in processing is methylation of the neo-c-terminus of the prenylated protein by lcmt (Isoprenylcysteine carboxylmethyltransferase). The enzymes involved in this process are thought to be targets for development of new anti-cancer drugs.

|   | Description                                       | Clone | Species | <b>Applications</b> | Format | Host               | Quantity | Cat. No. |  |
|---|---|-------|---------|---------------------|--------|--------------------|----------|----------|--|
| A | Antibody  |       |         |                     |        |                    |          |          |  |
|   | Anti-Farnesyl Transferase $\alpha$ subunit (FNTA) | 2418  | Н       | WB                  | Pur    | M IgG <sub>1</sub> | 100 mL   | 04-470   |  |
|   | Anti-Farnesyl Transferase β subunit (FNTB)        |       | Н       | WB                  | APur   | Rb IgG             | 100 mL   | 09-121   |  |

#### Fos

see AP-1 above

#### Fyn

Fyn is a Src-family non-receptor tyrosine kinase. It is involved in T-cell and B-cell activation as well as in keratinocyte differentiation. In T-cells, Fyn associates with the T-cell Receptor and Thy-1. Lck and Fyn localize to different cellular regions, suggesting that these two members of the Src family regulate different aspects of T-cell activation.

| Description         | Clone | Species | Applications Format | Host               | Quantity | Cat. No. |  |
|---------------------|-------|---------|---------------------|--------------------|----------|----------|--|
| Monoclonal Antibody |       |         |                     |                    |          |          |  |
| Anti-Fyn            | 18    | H M R   | WB IP IPK IH Pur    | M lgG $_{1\kappa}$ | 100 µg   | 05-436   |  |
| Anti-Fyn, Rb Mab    |       | H M R   | IF WB IH(P) Serum   | Rb IgG             | 100 pL   | 04-353   |  |

| Description                            | Clone | Species  | Applications | Format    | Host   | Quantity | Cat. No. |
|--|-------|----------|--------------|-----------|--------|----------|----------|
| Polyclonal Antibody                    |       |          |              |           |        |          |          |
| Anti-Fyn                               |       | Av H M R | IP WB        | antiserum | Rb IgG | 250 µL   | 06-133   |
| Anti-Fyn                               |       | НМ       | IP WB        | Serum     | Rb     | 50 µL    | AB1378   |
| Protein                                |       |          |              |           |        |          |          |
| Fyn, active                            |       |          | KA           |           |        | 10 µg    | 14-441   |
| Fyn, partially Pured                   |       |          | KA           |           |        | 300 mU   | 14-107   |
| siRNA                                  |       |          |              |           |        |          |          |
| Fyn SMARTpool siRNA reagent            |       | Н        | RNAi         |           |        | 5 nmol   | M-003140 |
| Fyn siRNA/siAb Assay Kit               |       | Н        | WB RNAi      |           |        | 1 kit    | 60-050   |
| siRNA plasmid, pKD-Fyn-v3              |       | Н        | RNAi         |           |        | 5 µg     | 62-087   |
| siRNA plasmid, pKD-Fyn-v6              |       | Н        | RNAi         |           |        | 5 µg     | 62-088   |
| Assay                                  |       |          |              |           |        |          |          |
| Fyn KinEASE FP Fluorescein Green Assay |       |          | KA           |           |        | 1 kit    | 32-069   |
| Fyn KinEASE FP-645nm FarRed Assay      |       |          | KA           |           |        | 1 kit    | 32-149   |
| General Reagents                       |       |          |              |           |        |          |          |
| Assay Dilution Buffer I (ADBI)         |       |          | KA           |           |        | 1 mL     | 20-108   |
| Assay Dilution Buffer II (ADBII)       |       |          | KA           |           |        | 1 mL     | 20-111   |
| Magnesium/ATP Cocktail                 |       |          | KA           |           |        | 1 mL     | 20-113   |
| Phosphate-Citrate Buffer, 5X           |       |          |              |           |        | 8 mL     | 20-143   |
| Assay Dilution Buffer, 5X              |       |          | KA           |           |        | 1 mL     | 20-145   |
| Mg2+ Lysis/Wash Buffer, 5X             |       |          | GPA          |           |        | 18 mL    | 20-168   |
| Protein Phosphatase Dilution Buffer    |       |          | PA           |           |        | 1 mL     | 20-169   |
| 100X GTPγS, 10mM                       |       |          | GPA          |           |        | 50 µL    | 20-176   |
| 100X GDP, 100mM                        |       |          | GPA          |           |        | 50 µL    | 20-177   |
| pNPP Ser/Thr Assay Buffer              |       |          | PA           |           |        | 20 mL    | 20-179   |
| pNPP Tyr Assay Buffer                  |       |          | PA           |           |        | 20 mL    | 20-180   |
| Tris Assay Dilution Buffer, 10X        |       |          | KA           |           |        | 1 mL     | 20-181   |
| Kinase Assay Blocking Buffer           |       |          | KA           |           |        | 10 mL    | 20-189   |
| TBS, 20X                               |       |          |              |           |        | 50 mL    | 20-190   |
| 10% BSA in TBS                         |       |          | KA           |           |        | 10 mL    | 20-191   |
| Kinase Assay Clearing Buffer           |       |          | KA           |           |        | 10 mL    | 20-193   |
| Dithiothreitol, 1M                     |       |          |              |           |        | 450 µL   | 20-265   |
| 1M MgCl <sub>2</sub>                   |       |          |              |           |        | 500 μL   | 20-303   |
| 50mM AMP                               |       |          |              |           |        | 50 µL    | 20-304   |
| 1 Om/M ATP                             |       |          |              |           |        | 300 µL   | 20-306   |
| 0.5M EDTA, pH 7.2                      |       |          |              |           |        | 2 mL     | 20-307   |
| 1 OX Detection Buffer                  |       |          |              |           |        | 3 mL     | 20-308   |

#### Grb2

Growth Factor Receptor Binding proteins (GRB) associate with activated receptor tyrosine kinases through their SH2 domains. GRB1 is more commonly known as PI3 Kinase p85. GRB2 is an adapter protein which recruits other proteins such as the Ras GEF SOS. Other GRB proteins include GRB7, and GRB10.

| Description         | Clone | Species    | <b>Applications</b> | Format | Host               | Quantity | Cat. No. |  |
|---------------------|-------|------------|---------------------|--------|--------------------|----------|----------|--|
| Antibody            |       |            |                     |        |                    |          |          |  |
| Anti-GRB2           | 3F2   | B H Ht M R | WB IH               | Pur    | M lgG $_{1\kappa}$ | 100 µg   | 05-372   |  |
| Anti-GRB2 (C-term)  |       | Н          | WB                  |        | Rb IgG             | 100 µL   | 04-358   |  |
| Anti-Grb2           | 2GB04 | B H M R Xn | IF IH(P)            | Pur    | M lgG <sub>1</sub> | 100 µg   | MAB1127  |  |
| GRB2, agarose conj. |       |            | IP WB               |        |                    | 100 µg   | 14-128   |  |

#### Heregulin

Ī

Heregulin  $\alpha$  and  $\beta$  are members of the Neuregulin family of ligands for Erb B or HER family receptor tyrosine kinases. Heregulin- $\beta$  is a specific ligand for HER-3 and HER-4, though HER-2 also becomes activated as a result of heregulin stimulation, likely by hetero-oligomerization with HER-4. *In vivo*, Heregulin promotes differentiation of mammary epithelial cells into secretory lobuloalveoli, a process involving Rab3A, which is activated by heregulin stimulation.

|        | Description                           | Clone | Species | <b>Applications</b> | Format | Host   | Quantity | Cat. No. |
|--------|---------------------------------------|-------|---------|---------------------|--------|--------|----------|----------|
| Antibo | ody                                   |       |         |                     |        |        |          |          |
|        | Heregulin-β, EGF Domain               |       |         | CULT                |        |        | 100 µg   | 01-201   |
|        | Anti-Heregulin-1 precursor            |       | ΗМ      | WB                  | Pur    | Rb IgG | 200 µg   | 07-494   |
|        | Anti-Herstatin                        | CW001 | Н       | WB IH               |        | M lgG  | 200 µg   | 05-582   |
|        | Herstatin Detection Kit (ELISA based) |       |         | ELISA               |        |        | 1 kit    | 17-354   |

#### CMT

Icmt (Isoprenylcysteine carboxylmethyltransferase) methylates the cleaved form of Ras/Rho small G-protein family members and nuclear Lamin. Methylation is part of the required post-translational modification (prenylation by FT, cleavage by Rce 1, methylation by Icmt) required for these proteins to become active.

|       | Description  | Clone | Species | Applications | Format | Host   | Quantity | Cat. No. |
|-------|--|-------|---------|--------------|--------|--------|----------|----------|
| Antib | oody   |       |         |              |        |        |          |          |
|       | Anti-Isoprenylcysteine Carboxymethyl Transferase (IC | MT)   | Н       | WB           | APur   | Rb IgG | 100 mL   | 09-119   |

#### ILK

Integrin-Linked Kinase (ILK) is a 59 kDa Ser/Thr kinase with four Ankyrin-like repeats. It was cloned from a human placental cDNA library using the yeast two-hybrid system and the cytoplasmic domain of  $\beta$ 1 Integrin as bait. The association of ILK and  $\beta$ 1 Integrin has been reported *in vivo*. The kinase activity appears to be reduced following Integrin activation, and overexpression of p59ILK inhibits cellular adhesion to Integrin substrates.

| Description                           | Clone  | Species | <b>Applications</b> | Format | Host    | Quantity | Cat. No. |
|---------------------------------------|--------|---------|---------------------|--------|---------|----------|----------|
| Antibody                              |        |         |                     |        |         |          |          |
| ILK                                   |        | H M R   | WB IP IPK IF        | APur   | Rb IgG  | 100 µg   | 06-592   |
| ILK cDNA (hyperactive) in pUSEamp     |        |         | TFX                 |        |         | 5 µg     | 21-185   |
| ILK cDNA (inactive, R211A) in pUSEamp |        |         | TFX                 |        |         | 5 µg     | 21-182   |
| ILK cDNA (inactive, S343A) in pUSEamp |        |         | TFX                 |        |         | 5 µg     | 21-183   |
| ILK cDNA (wt) in pUSEamp              |        |         | TFX                 |        |         | 5 µg     | 21-184   |
| ILK siRNA/siAb Assay Kit              |        | Н       | WB RNAi             |        |         | 1 kit    | 60-116   |
| ILK SMARTpool siRNA reagent           |        |         | RNAi                |        |         | 5 nmol   | M-004499 |
| ILK                                   | 65.1.9 | H M R   | WB IP IC            | Pur    | M laGab | 200 µg   | 05-575   |

#### **IGF** and **IGF** Receptor

IGF-I and -II signal through the IGF-I Receptor, which is homologous to the Insulin Receptor. The high-affinity IGF-II Receptor does not play a direct role in signaling, but regulates the concentration of free IGF-II. The IGFs are involved in skeletal growth, and are essential for prevention of apoptosis. Serum levels of free IGFs are kept low by the action of IGF binding proteins (IGFBPs), which sequester the IGFs. Overexpression of IGFBPs may induce apoptosis, presumably by reduction of free IGF; IGFBP levels are also altered in some cancers. The IGF-I Receptor is not as mitogenic as some other growth factor receptors, but its ability to activate the PI3 Kinase pathway, through the Insulin Receptor Substrate (IRS) proteins, is critical for mediating cell survival.

|        | Description  | Clone  | Species             | Applications        | Format    | Host                     | Quantity | Cat. No. |
|--------|--|--------|---------------------|---------------------|-----------|--------------------------|----------|----------|
| Grow   | th Factor Protein  |        |                     |                     |           |                          |          |          |
|        | IGF-II   |        |                     | CULT                |           |                          | 25 µg    | 01-142   |
|        | IGF-I (resistant to IGFBPs)  |        |                     | CULT                |           |                          | 25 µg    | 01-189   |
|        | IGF-I  |        |                     | CULT                |           |                          | 25 µg    | 01-208   |
|        | IGF-I, biotin conj.  |        |                     | WB                  | HPLC      |                          | 2 µg     | 01-212   |
|        | Insulin-like Growth Factor-I, recombinant human                        |        |                     |                     | Pur       |                          | 100 µg   | GF006    |
|        | Insulin-like Growth Factor-II, recombinant human                       |        |                     |                     | Pur       |                          | 50 µg    | GF007    |
|        | Insulin-like Growth Factor-I, recombinant mouse                        |        |                     |                     | Pur       |                          | 50 µg    | GF121    |
| Mono   | clonal Antibody  |        |                     |                     |           |                          |          |          |
|        | Anti-IGF-II  | S1F2   | H R                 | ELISA WB<br>IC NEUT |           | M lgG <sub>1</sub>       | 500 µg   | 05-166   |
|        | Anti-IGF-I   | Sm1.2  | Av H M R            | IP WB IH<br>NEUT    |           | $\text{M IgG}_{1\kappa}$ | 200 µg   | 05-172   |
|        | Anti-IGF-IR  | JBW902 | ΗМ                  | IP WB IH            | Pur       | M lgG <sub>2a</sub>      | 200 µg   | 05-656   |
|        | Anti-phospho-IGF-IR (Tyr1131)/InsR (Tyr1158)                           | JY202  | ΗМ                  | IP WB               | Pur       | M lgG <sub>2bκ</sub>     | 200 µg   | 05-694   |
|        | Anti-Insulin-like Growth Factor-I Receptor, $lpha$ -Subunit            | 24-31  | Н                   | IP IH(P)            | Pur       | M lgG <sub>1</sub>       | 100 µg   | MAB1120  |
|        | Anti-Insulin-like Growth Factor-I Receptor, $lpha$ -Subunit exon $7/8$ | 24-57  | H Rb                | INHIB               | Pur       | M lgG <sub>1</sub>       | 100 µg   | MAB1122  |
|        | Anti-Insulin-like Growth Factor-l Receptor, β-Subunit                  | 1-2    | H R                 | WB                  | Pur       | M lgG <sub>2b</sub>      | 100 µg   | MAB1123  |
|        | Anti-Insulin-like Growth Factor-I                                      | M23    | Н                   | IH(P)               | Pur       | M lgG <sub>1</sub>       | 100 µg   | CBL52    |
|        | Anti-Insulin-like Growth Factor-II                                     | W2H1   | Н                   | ELISA RIA IH        | Pur       | M lgG <sub>1</sub>       | 100 µg   | CBL82    |
| Polycl | onal Antibody  |        |                     |                     |           |                          |          |          |
|        | Anti-IGFBP-1   |        | Н                   | IP WB IC            | antiserum | Rb IgG                   | 200 μL   | 06-106   |
|        | Anti-IGFBP-2   |        | H M R B<br>Eq Gt Po | IP WB IC            | antiserum | Rb IgG                   | 200 μΙ   | 06-107   |
|        | Anti-IGFBP-3   |        | Н                   | IP WB               | antiserum | Rb IgG                   | 200 μL   | 06-108   |
|        | Anti-IGFBP-4   |        | ВНМРо               | IP WB               | antiserum | Rb IgG                   | 200 μL   | 06-109   |
|        | Anti-IGFBP-5   |        | H R                 | IP WB               | antiserum | Rb IgG                   | 200 μL   | 06-110   |
|        | Anti-IGF-IRα subunit   |        | ΗМ                  | IP WB IC            |           | Ch IgY                   | 250 µg   | 06-429   |
|        | Anti-IGFBP-1   |        | MR                  | IP WB               | Pur       | Rb IgG                   | 200 µg   | 06-580   |
|        | Anti-Insulin-like Growth Factor-I                                      |        | Н                   | ELISA WB            | Pur       | Rb                       | 100 µg   | AB1437   |
|        | Anti-Insulin-like Growth Factor-II                                     |        | Н                   | ELISA WB            | Pur       | Rb                       | 100 µg   | AB1438   |
|        |  |        |                     |                     |           |                          |          |          |

|        | Description   | Clone | Species | Applications | Format | Host | Quantity  | Cat. No. |
|--------|---|-------|---------|--------------|--------|------|-----------|----------|
|        | Anti-Insulin-like Growth Factor-I                           |       | Μ       | ELISA WB     | APur   | Gt   | 50 µg     | AB2131P  |
|        |   |       |         | NEUT         |        |      |           |          |
|        | Anti-Insulin-like Growth Factor Binding Protein 6           |       |         |              | APur   | Ch   | 100 µg    | AB9148   |
|        | Anti-Insulin-like Growth Factor-1 Receptor, $lpha$ -Subunit |       | Н       | IP WB        | Pur    | Rb   | 100 µg    | CBL257   |
|        | IGFBP-3   |       |         | ABA          |        |      | 5 µg      | 12-131   |
| Assay  |   |       |         |              |        |      |           |          |
|        | CpG WIZ® H19-IGF2 Amplification Kit                         |       | Н       |              |        |      | 25 assays | S7842    |
|        | IGF-IR KinEASE FP Fluorescein Green Assay                   |       |         | KA           |        |      | 1 kit     | 32-070   |
|        | IGF-IR KinEASE FP-645nm FarRed Assay                        |       |         | KA           |        |      | 1 kit     | 32-150   |
| iRNA   |   |       |         |              |        |      |           |          |
|        | siRNA plasmid, pKD-IGF-IR-v2                                |       | Н       | RNAi         |        |      | 5 µg      | 62-114   |
|        | siRNA plasmid, pKD-IGF-IR-v1                                |       | Н       | RNAi         |        |      | 5 µg      | 62-115   |
|        | IGF-1R SMARTpool siRNA reagent                              |       |         | RNAi         |        |      | 5 nmol    | M-003012 |
| Protei | n   |       |         |              |        |      |           |          |
|        | IGF-IR (δ1-958), active                                     |       |         | KA           |        |      | 20 µg     | 14-465   |
|        | IGF-IRtide  |       |         |              |        |      | 1 mg      | 12-527   |
|        | IGFBP-3   |       |         | ABA          |        |      | 5 µg      | 12-131   |
|        |   |       |         |              |        |      |           |          |

#### Insulin and Insulin Receptor

The Insulin Receptor (IR) is synthesized as a single polypeptide, which is subsequently cleaved to generate an extracellular  $\alpha$ -chain and a transmembrane and intracellular  $\beta$ -chain, tethered together by disulfide bonds. The  $\beta$ -chain has multiple tyrosine phosphorylation sites, including three autophosphorylation sites at its activation loop. The overall structure of the IR is highly homologous to the IGF-I Receptor, except in their c-termini, where the two proteins diverge somewhat. The IR signals primarily by phosphorylating the Insulin Receptor Substrate (IRS) family of proteins, which creates docking sites for SH2-domain containing proteins. Insulin signaling is highly dependent on the PI3 Kinase pathway and Akt, which appear to mediate the functions of insulin.

|         | Description   | Clone    | Species   | Applications | Format | Host               | Quantity | Cat. No. |
|---------|---|----------|-----------|--------------|--------|--------------------|----------|----------|
| Recom   | binant Growth Factor                                |          |           |              |        |                    |          |          |
|         | Insulin (Arg-Insulin)                               |          |           | CULT         |        |                    | 10 mg    | 01-207   |
|         | Incelligent™ SG, recombinant human                  |          |           |              | Pur    |                    | 1 g      | 4502-01  |
|         | Incelligent AF, recombinant human                   |          |           |              | Pur    |                    | 1 g      | 4506-01  |
| Monoc   | lonal Antibody                                      |          |           |              |        |                    |          |          |
|         | Anti-Insulin  | MAB 1    | Н         | ELISA RIA    | Pur    | M IgG <sub>1</sub> | 1 mg     | CBL71    |
|         | Anti-Pro-Insulin C Peptide                          | C-PEP-01 | Н         | IH           | Pur    | M IgG <sub>1</sub> | 1 mg     | CBL94    |
|         | Anti-Insulin Receptor, $\alpha$ subunit             | 47-9     | B H Po Sh | INHIB not IP | Pur    | M IgG <sub>1</sub> | 100 µg   | MAB1137  |
|         |   |          | not M not | R not WB     |        |                    |          |          |
|         | Anti-Insulin Receptor, $\alpha$ subunit, azide free | 47-9     | B H Po Sh | INHIB not IP | Pur*   | M IgG <sub>1</sub> | 100 µg   | MAB1137Z |
|         |   |          | not M not | R not WB     |        |                    |          |          |
|         | Anti-Insulin Receptor, $\alpha$ subunit             | 83-7     | ВНРо      | ELISA FC KA  | Pur    | M IgG <sub>1</sub> | 100 µg   | MAB1138  |
|         |   |          | Rb Sh     | IH(P) not WE | 3      |                    |          |          |
|         | Anti-Insulin Receptor, $\beta$ subunit, C-terminus  | CT-3     | HMMkR     | ELISA KA     | Pur    | M IgG <sub>1</sub> | 100 µg   | MAB1139  |
|         |   |          |           | WB IH(P)     |        |                    |          |          |
| Polyclo | onal Antibody                                       |          |           |              |        |                    |          |          |
|         | Anti-Insulin Receptor, β subunit                    |          | НМ        | IP WB        | APur   | Rb IgG             | 100 µL   | 07-724   |
|         | Anti-phospho-IR (Tyr972)                            |          | НМ        | WB           | APur   | Rb IgG             | 100 µL   | 07-838   |
|         | Anti-phospho-IR/IGF1R (Tyr1158)                     |          | Н         | WB IC IH     | APur   | Rb IgG             | 100 µL   | 07-839   |
|         | Anti-phospho-IR/IGF1R (TyrpY1162/1163)              |          | Н         | WB IC IH     | APur   | Rb IgG             | 100 µL   | 07-840   |
|         | Anti-phospho-IR/IGF1R (Tyr1158/Tyr1162/Tyr1163      | )        | НМ        | WB           | APur   | Rb IgG             | 100 µL   | 07-841   |
|         | Anti-Insulin  |          | ВНРо      | WB IH        | Serum  | Gp                 | 1 mL     | AB3440   |
|         | Anti-Insulin  |          | Н         | ELISA WB     | APur   | Ch                 | 100 µg   | AB3455   |
|         | Anti-Insulin Receptor B, exon 11                    |          | Н         | RIA          | Pur    | Rb                 | 100 µg   | CBL245   |
|         | Anti-Insulin Receptor, C-terminus                   |          | Н         | IP WB        | Pur    | Rb                 | 100 µg   | CBL258   |
| Assay   |   |          |           |              |        |                    |          |          |
|         | Insulin R KinEASE FP-645nm FarRed Assay             |          |           | KA           |        |                    | 1 kit    | 32-151   |
|         | Insulin R KinEASE FP Fluorescein Green Assay        |          |           | KA           |        |                    | 1 kit    | 32-071   |
| Protein |   |          |           |              |        |                    |          |          |
|         | Insulin Receptor, active                            |          |           | KA           |        |                    | 10 µg    | 14-466   |
|         |   |          |           |              |        |                    |          |          |

#### IRSp53

IRSp53 is a recently identified substrate for the Insulin Receptor tyrosine kinase which is predominantly expressed in mammalian brain and is involved in guidance of the neuronal growth cone. At least five isoforms resulting from alternative splicing have been found; splicing produces proteins with alternative c-terminal sequences that modulate subcellular localization. IRSp53 contains a CRIB motif that allows it to bind and activate the small G-Protein cdc42, an SH3 domain that couples it to actin. when overexpressed, IRSp53 causes extensive formation of filopodia.

|   | Description | Clone | Species | Applications | Format | Host  | Quantity | Cat. No. |
|---|-------------|-------|---------|--------------|--------|-------|----------|----------|
| Δ | ntibody     |       |         |              |        |       |          |          |
|   | Anti-IRSp53 |       | H M R   | WB           | Serum  | R IgG | 100 µL   | 07-786   |

#### JIP (JNK Interacting Protein)

JIP, JNK Interacting protein, is a scaffolding protein that binds to JNK. It was originally identified as cytoplasmic inhibitor of JNK. It is also known to bind to MAP Kinase Kinase 7, MKP7, and members of the MLK family. It is believed to be essential for JNK activation in certain cell types including neurones.

| Description        | Clone | Species | Applications Forma | t Host | Quantity | Cat. No. |
|--------------------|-------|---------|--------------------|--------|----------|----------|
| Antibody           |       |         |                    |        |          |          |
| Anti-JIP1          |       |         |                    | M lgG  | 100 µg   | 04-481   |
| Anti-phospho-JIP 1 |       |         |                    |        |          |          |
| JIP2               |       |         |                    |        |          | 04-482   |

#### JNK/SAPK

The JNK/SAPK1 kinases, like the other MAPK-like kinases, are thought to phosphorylate multiple substrates and regulate many processes, including proliferation (in some cell types) and apoptosis. The SAPK2/3 family is widely referred to as the p38 family. These kinases are activated by stresses, most notably inflammatory cytokines, irradiation, and certain toxins such as anisomycin, and arsenite. The activating kinases of SAPK2/3 are SKK2/MEK3 for SAPK2a and 2b, and MKK6 for SAPK3. The targets of the SAPK2/3 family include the MAPKAP kinases 2 and 3/3pK. In addition, SKK4 is related to this family, exhibiting 60% identity, and is activated by MKK6.

| Description                                | Clone    | Species | Applications | Format | Host   | Quantity | Cat. No. |
|--|----------|---------|--------------|--------|--------|----------|----------|
| noclonal Antibody                          |          |         |              |        |        |          |          |
| Anti-JNK2                                  |          | H M R   | WB ELISA     | APur   | M lgG  | 100 µg   | 05-986   |
| Anti-JNK3/SAPK1b                           | C05T     | HR      | WB           | Pur    | M lgG  | 100 μL   | 05-893   |
| yclonal Antibody                           |          |         |              |        |        |          |          |
| Anti-JNK/SAPK1                             |          | H M R   | IP WB        | Pur    | Rb IgG | 200 µg   | 06-748   |
| Anti-phospho-JNK (Thr 183/Tyr 185, Thr 221 | /Tyr223) | HMR     | IP WB        | Pur    | Rb IgG | 200 µg   | 07-175   |
| Anti-JNK1                                  |          | MR      | IP KA        | Pur    | Sh     | 250 µg   | AB4081   |
| Anti-JNK2                                  |          | HMR     | ELISA WB     | Serum  | Rb     | 50 μL    | AB8910   |
|  |          |         | not IP       |        |        |          |          |
| NA   |          |         |              |        |        |          |          |
| JNK2/SAPK1a siRNA/siAb Assay Kit           |          | Н       | WB RNAi      |        |        | 1 kit    | 60-099   |
| siRNA plasmid, pKD-JNK2a2/SAPK1a-v1        |          | Н       | RNAi         |        |        | 5 µg     | 62-097   |
| siRNA plasmid, pKD-JNK2a2/SAPK1a-v5        |          | ΗМ      | RNAi         |        |        | 5 µg     | 62-098   |
| JNK2 SMARTpool siRNA                       |          | Н       | RNAi         |        |        | 5 nmol   | M-003515 |
| tein                                       |          |         |              |        |        |          |          |
| JNK1a1/SAPK1c, active                      |          |         | KA           |        |        | 10 µg    | 14-327   |
| JNK1a1/SAPK1c, unactive                    |          |         | KA           |        |        | 50 µg    | 14-328   |
| JNK2a2/SAPK1a, active                      |          |         | KA           |        |        | 10 µg    | 14-329   |
| JNK3/SAPK1b, active                        |          |         | KA           |        |        | 10 µg    | 14-501   |
| JNK3/SAPK1b, unactive                      |          |         | KA           |        |        | 50 µg    | 14-523   |
| JNK3/SAPK1b (K55R), inactive,              |          |         | KA           |        |        | 200 µg   | 13-126   |
| Streptavidin Binding Peptide               |          |         |              |        |        |          |          |
| JNK3tide                                   |          |         |              |        |        | 1 mg     | 12-528   |
| MKK6/SKK3, active                          |          |         | KA           |        |        | 10 µg    | 14-303   |

#### Jun

see "AP-1" above

#### LIMK

The Actin cytoskeleton undergoes extensive changes during cell morphogenesis and motility. Rho, the small GTPase, is involved with regulating these processes. LIM Kinase 1 phosphorylates the Actin-depolymerizing factor, Cofilin, thus regulating Actin cytoskeletal reorganization. It appears that LIM Kinase 1 is activated by Rho and its downstream protein kinase ROK (Rho-associated Kinase).

| Description                           | Clone | Species | Applications | Format | Host   | Quantity | Cat. No. |
|---------------------------------------|-------|---------|--------------|--------|--------|----------|----------|
| Antibody                              |       |         |              |        |        |          |          |
| Anti-LIM Kinase 1                     |       | НМ      | WB           | Pur    | Rb     | 100 µg   | ab3814   |
| Anti-phospho-LIMK 1/2 (Tyr507/Thr508) |       | Н       | WB           | APur   | Rb IgG | 100 µL   | 07-850   |
| Protein                               |       |         |              |        |        |          |          |
| LIM Kinase 1, active                  |       |         | KA           |        |        | 10 µg    | 14-656   |

|       | Description                        | Clone | Species | Applications Format Host | Quantity | Cat. No.     |
|-------|------------------------------------|-------|---------|--------------------------|----------|--------------|
|       | LIM Kinase 1, active, mouse        |       |         | KA                       | 15 µg    | 14-457       |
|       | LIM Kinase 1, unactive             |       |         | KA                       | 50 µg    | 14-659       |
| siRNA | A                                  |       |         |                          |          |              |
|       | siRNA plasmid, pKD-LIM Kinase 1-v2 |       |         | RNAi                     | 5 μg     | 62-256       |
|       | siRNA plasmid, pKD-LIM Kinase 1-v3 |       |         | RNAi                     | 5 μg     | 62-257       |
|       | AdenoSilence™ RNAi Kit LIMK1       |       | Н       | RNAi                     | 3 vials  | GAL10103     |
|       | AdenoSilence RNAi Virus LIMK1-v1   |       | Н       | RNAi                     | 500 μL   | GAL10103-V1  |
|       | AdenoSilence RNAi Virus LIMK1-v2   |       | Н       | RNAi                     | 500 μL   | GAL10103-V2  |
|       | AdenoSilence RNAi Virus LIMK1-v3   |       | Н       | RNAi                     | 500 μL   | GAL10103-V3  |
|       | AdenoSilence RNAi Kit LIMK2        |       | Н       | RNAi                     | 3 vials  | GAL10104     |
|       | AdenoSilenc RNAi Virus LIMK2-v1    |       | Н       | RNAi                     | 500 μL   | GAL10104-V1  |
|       | AdenoSilence RNAi Virus LIMK2-v1 1 |       | Н       | RNAi                     | 500 μL   | GAL10104-V11 |
|       | AdenoSilence RNAi Virus LIMK2-v2   |       | Н       | RNAi                     | 500 μL   | GAL10104-V2  |

#### **MAPKAP Kingse**

MAPKAP Kinases are Ser/Thr kinases activated by phosphorylation by the SAPK2/3, or p38 family of stress-activated kinases. MAPKAP Kinase 2 and 3 phosphorylate substrates with the sequence H-X-R-X-X-S where H is any hydrophobic residue, including CREB (Ser 133) and HSP27. MAPKAP kinases are similar to Rsk kinases, which are activated by mitogen-signaling pathways, rather than stress-signaling pathways.

| Description                     | Clone | Species | <b>Applications Format</b> | Host | Quantity | Cat. No. |  |
|---------------------------------|-------|---------|----------------------------|------|----------|----------|--|
| Antibody                        |       |         |                            |      |          |          |  |
| MAPKAP Kinase Substrate Peptide |       |         | KA                         |      | 1 mg     | 12-240   |  |
| Rsk3, active                    |       |         | KA                         |      | 10 µg    | 14-462   |  |
| Rsk4, active                    |       |         | KA                         |      | 10 µg    | 14-702   |  |

#### Rsk1 (MAPKAP Kinase 1a)

Rsk1, Rsk2, and Rsk3 are Ser/Thr kinases which are activated by MAP Kinases, and are alternatively referred to as MAPKAP kinases 1a, 1b, and 1c. Following activation, Rsk isoforms translocate to the nucleus, where Rsk2 phosphorylates Ser10 of Histone H3, a critical step in the stimulation of immediate-early gene expression. Rsk1 is able to phosphorylate Ser133 of CREB, activating it as a transcription factor. Rsk1 activation can be monitored by phosphorylation state specific antibodies directed against sites involved in its activation.

| Description                              | Clone  | Species | Applications | Format | Host   | Quantity | Cat. No. |
|--|--------|---------|--------------|--------|--------|----------|----------|
| Antibody                                 |        |         |              |        |        |          |          |
| RSK1 (N-term)                            |        |         |              |        |        |          | 04-417   |
| Rsk1 Phospho (pT359/pS363)               |        |         |              |        |        |          | 04-419   |
| Rsk1 Phospho (pS380)                     |        |         |              |        |        |          | 04-418   |
| Anti-Rsk1/MAPKAP Kinase 1a               |        | H M R   | IP WB        |        | Rb IgG | 100 µg   | 06-668   |
| Anti-phospho-Rsk1/MAPKAP Kinase 1a (Thr3 | 59)    |         | WB           |        | Rb lgG | 100 µL   | 07-526   |
| Anti-phospho-Rsk1/MAPKAP Kinase 1a (Ser3 | 363)   | H M R   | WB           |        | Sh IgG | 100 µg   | 06-824   |
| Anti-Rsk 1                               |        | H M R   | ELISA IP WB  | APur   | Rb     | 100 µg   | AB3190   |
| siRNA plasmid, pKD-Rsk1/MAPKAP Kinase    | l a-v5 | Н       | RNAi         |        |        | 5 µg     | 62-106   |
| Rsk1/MAPKAP Kinase 1a siRNA/siAb Assay   | / Kit  | Н       | WB RNAi      |        |        | 1 kit    | 60-052   |
| Rsk1/MAPKAP Kinase 1a                    |        | Н       | RNAi         |        |        | 5 nmol   | M-003025 |
| SMARTpool siRNA reagent                  |        |         |              |        |        |          |          |
| protein                                  |        |         |              |        |        |          |          |
| S6 Kinase/Rsk Substrate Peptide 1        |        |         | KA           |        |        | 2 mg     | 12-124   |
| Rsk1/MAPKAP Kinase 1a, active            |        |         | KA           |        |        | 10 µg    | 14-479   |
| Rsk1/MAPKAP Kinase 1a, active            |        |         | KA           |        |        | 10 µg    | 14-509   |
| Rsk2/MAPKAP Kinase 1b, active            |        |         | KA           |        |        | 10 µg    | 14-480   |
| siRNA                                    |        |         |              |        |        |          |          |
| Rsk1/MAPKAP Kinase 1a SMARTpool          |        | Н       | RNAi         |        |        | 5 nmol   | M-003025 |
| siRNA reagent                            |        |         |              |        |        |          |          |

#### Rsk2 (MAPKAP Kinase 1b)

| Antib | ody   |     |         |        |        |          |
|-------|---|-----|---------|--------|--------|----------|
|       | Rsk2 (C-term)                               |     |         |        |        | 04-420   |
|       | Rsk2 (N-term)                               |     |         |        |        | 04-421   |
|       | Anti-Rsk2/MAPKAP Kinase 1b                  | HMR | IP WB   | Rb IgG | 100 µg | 06-918   |
| siRN  | 4   |     |         |        |        |          |
|       | siRNA plasmid, pKD-Rsk2/MAPKAP Kinase 1b-v4 | Н   | RNAi    |        | 5 µg   | 62-109   |
|       | siRNA plasmid, pKD-Rsk2/MAPKAP Kinase 1b-v5 | Н   | RNAi    |        | 5 µg   | 62-110   |
|       | Rsk2/MAPKAP Kinase 1b siRNA/siAb Assay Kit  | Н   | WB RNAi |        | 1 kit  | 60-100   |
|       | Rsk2/MAPKAP Kinase 1b SMARTpool             | Н   | RNAi    |        | 5 nmol | M-003026 |
|       | siRNA reagent                               |     |         |        |        |          |

23

| Description                           | Clone | Species   | Applications | Format | Host   | Quantity | Cat. No. |
|---------------------------------------|-------|-----------|--------------|--------|--------|----------|----------|
| MAPKKAP-K2                            |       |           |              |        |        |          |          |
| Antibody                              |       |           |              |        |        |          |          |
| Anti-MAPKAPK2                         | 7H4.2 | Н         | IP WB        | Pur    | M lgG  | 100 µg   | 04-539   |
| Anti-phospho-MAPKAPK2                 |       | Н         | WB           | Pur    | M lgG  | 100 µg   | 05-1016  |
| Anti-phospho-MAPKAP Kinase 2 (Thr222) |       | H Ht M Rb | KA WB        |        | Rb IgG | 200 µg   | 07-155   |
| Anti-MAPKAP Kinase 2, long form       |       | H M       | WB           |        | Rb IgG | 100 µg   | 07-241   |
| Anti-MAPKAP Kinase 2, agarose conj.   |       |           | ELISA WB     |        | Rb IgG | 200 μL   | 16-203   |
| MAPKAP Kinase 2 Immunoprecipitation   |       |           | KA           |        |        | 1 kit    | 17-297   |
| Kinase Assay Kit, non-radioactive     |       |           |              |        |        |          |          |
| siRNA                                 |       |           |              |        |        |          |          |
| siRNA plasmid, pKD-MAPKAP Kinase 2-v1 |       | HMR       | RNAi         |        |        | 5 µg     | 62-143   |
| siRNA plasmid, pKD-MAPKAP Kinase 2-v4 |       | Н         | RNAi         |        |        | 5 µg     | 62-144   |
| Rsk2/MAPKAP Kinase 1b SMARTpool       |       | Н         | RNAi         |        |        | 5 nmol   | M-003026 |
| siRNA reagent                         |       |           |              |        |        |          |          |
| МАРККАР-КЗ                            |       |           |              |        |        |          |          |
| Antibody                              |       |           |              |        |        |          |          |
| Anti-MAPKAP Kinase 3 (3pK)            |       | Н         | IP           |        | Sh IgG | 100 µg   | 06-601   |

#### **MBP**

siRNA plasmid, pKD-Rsk1/MAPKAP Kinase 1a-v3

Myelin Basic Protein is a widely used substrate for the assay of several protein kinases, including MAP kinase and PKC. MAP kinase phosphorylates a single site on MBP (Thr98), and the extent of phosphorylation of this site can be monitored by use of a phosphorylation state-specific antibody, providing the basis of a non-radioactive assay for MAP kinase activity.

RNAi

62-125

5 µg

Н

| Description                                 | Clone | Species | Applications | Format | Host                | Quantity | Cat. No.        |  |  |  |  |
|---|-------|---------|--------------|--------|---------------------|----------|-----------------|--|--|--|--|
| Aonoclonal Antibody                         |       |         |              |        |                     |          |                 |  |  |  |  |
| Anti-phospho-MBP                            | P12   | WR      | KA WB        | Pur    | M lgG               | 100 µg   | 05-429          |  |  |  |  |
| Anti-MBP                                    | SKB3  | H M R   | WB           | Pur    | M lgG <sub>1κ</sub> | 50 µg    | 05-675          |  |  |  |  |
| Anti-phospho-MBP                            | BK403 | WR      | WB           | Pur    | M $\lg G_{1\kappa}$ | 200 µg   | 05-705          |  |  |  |  |
| Anti-phospho-MBP                            | BK403 |         | WB           | Pur    | M $\lg G_{1\kappa}$ | 1 mg     | 05-705MG        |  |  |  |  |
| Antibody Conjugate                          |       |         |              |        |                     |          |                 |  |  |  |  |
| Anti-phospho-MBP, HRP conj.                 |       |         | ELISA KA WE  | }      |                     | 50 μL    | 16-198          |  |  |  |  |
| Anti-phospho-MBP, HRP conj.                 |       |         |              |        |                     | 50 μL    | 16-206          |  |  |  |  |
| Protein                                     |       |         |              |        |                     |          |                 |  |  |  |  |
| MBP (myelin basic protein), bovine, Pured   |       |         | KA           |        |                     | 10 mg    | 13-104          |  |  |  |  |
| MBP, Dephosphorylated                       |       |         | KA           |        |                     | 5 mg     | 13-110          |  |  |  |  |
| MBP, biotin conj.                           |       |         | KA           |        |                     | 1 mg     | 13-111          |  |  |  |  |
| Assay                                       |       |         |              |        |                     |          |                 |  |  |  |  |
| MBP Assay Kit, 96-well,                     |       |         | KA           |        |                     | 1 kit    | 1 <i>7</i> -353 |  |  |  |  |
| Chemiluminescence Detection                 |       |         |              |        |                     |          |                 |  |  |  |  |
| MBP Coated Microplate, 96 Well              |       |         | KA           |        |                     | 2 plates | 30-011          |  |  |  |  |
| KinEASE FP Fluorescein Green Assay - Module | 2 -   |         | KA           |        |                     | 1 kit    | 32-002          |  |  |  |  |
| MBP (Thr98) Kinases                         |       |         |              |        |                     |          |                 |  |  |  |  |
| KinEASE FP-645nm FarRed Assay Module 2 -    |       |         | KA           |        |                     | 1 assay  | 32-005          |  |  |  |  |
| MBP (Thr98) Kinases                         |       |         |              |        |                     |          |                 |  |  |  |  |

#### MEK

MAP Kinase/Erk Kinase (MEK), alternatively known as MKK, is a true dual-specificity kinase, in that it phosphorylates the MAP kinases on both the Thr and Tyr of the activation motif TEY. In vitro, the tyrosine phosphorylation is favored, whereas in vivo both phosphorylation events appear to occur simultaneously. This suggests that an additional factor is present in cells to facilitate the reaction. MEK1 and MEK2 are activated by phosphorylation of two serine residues (Ser218/222 in MEK1 and Ser222/226 in MEK2), which are substrates for the Raf family of kinases. Mutation of the phosphorylation sites from Ser to Asp creates a protein with constitutive kinase activity, which when expressed in cells is able to cause transformation.

|    | Description   | Clone | Species           | Applications          | Format | Host   | Quantity | Cat. No. |
|----|---|-------|-------------------|-----------------------|--------|--------|----------|----------|
| Мо | noclonal Antibody                                   |       |                   |                       |        |        |          |          |
|    | Anti-phospho-MEK1 (Ser218/222)/MEK2<br>(Ser222/226) |       | H M R             | WB                    | Pur    | Rb IgG | 100 µL   | 05-747   |
|    | Anti-MEK1, NT                                       | C12T  | Ca H M<br>Mk R Rb | IP WB                 | Pur    | Rb IgG | 100 μL   | 05-925   |
|    | Anti-MEK1 CT, Rb Mab                                |       | H R               | FC IP VVB<br>IC IH(P) | Serum  | Rb IgG | 100 µL   | 04-376   |
|    | Anti-MEK2 NT, Rb Mab                                |       | H R               | FC IP WB              | Serum  | Rb IgG | 100 µL   | 04-377   |

| Description   | Clone | Species                    | Applications       | Format | Host   | Quantity | Cat. No.    |
|---|-------|----------------------------|--------------------|--------|--------|----------|-------------|
| Polyclonal Antibody                                     |       |                            |                    |        |        |          |             |
| Anti-phospho-MEK1 (Ser298)                              |       | H M R                      | WB                 | APur   | Rb IgG | 200 μL   | 07-339      |
| Anti-phospho-MEK1 (Ser218/222)/MEK2<br>(Ser222/226)     |       | H M R                      | WB                 | Pur    | Rb IgG | 200 µg   | 07-461      |
| Anti-MEK1   |       | H M R                      | IP IPK VVB         | Pur    | Rb IgG | 200 µg   | 07-641      |
| Anti-phospho-MEK1 (Thr292)                              |       | НМ                         | WB                 | APur   | Rb IgG | 100 µL   | 07-852      |
| Anti-phospho-MEK1 (Thr386)                              |       | H M R                      | WB                 | APur   | Rb IgG | 100 µL   | 07-853      |
| Anti- phospho-MEK2 (Thr394)                             |       | Н                          | WB                 | APur   | Rb IgG | 100 µL   | 07-854      |
| Anti-phospho-MEK2 (Thr394)                              |       | M R                        | WB                 | APur   | Rb IgG | 100 µL   | 07-855      |
| Anti-phospho-MKK3/6 (Ser189/Thr193)/<br>(Ser207/Thr211) |       | Н                          | WB                 | APur   | Rb lgG | 100 µL   | 07-856      |
| Anti-phospho-MKK4 (Ser257/Thr261)                       |       | Н                          | WB                 | APur   | Rb IgG | 100 µL   | 07-857      |
| Anti-phospho-MKK7 (Ser271/Thr275)                       |       | ΗМ                         | WB                 | APur   | Rb IgG | 100 µL   | 07-858      |
| Anti-MEK1, N-terminus                                   |       | H M R Xn                   | IP WB IH           | Pur    | Rb     | 50 µg    | AB3183      |
| Anti-MEK5   |       | Ca H M<br>Mk R Rb<br>Sh Xn | ELISA IP<br>VVB IH | APur   | Rb     | 50 µg    | AB3184      |
| Anti-MKK6   |       | H M R Xn                   | IP WB              | Pur    | Rb     | 50 µg    | AB3185      |
| Anti-MEK1/2, phospho-specific (Ser218/222)              |       | Н                          | WB                 | APur   | Rb     | 100 µL   | AB3810      |
| Anti-MEK1/2, phospho-specific (Ser218/222)              |       | Н                          | WB                 | APur   | Rb     | 25 µL    | AB3810-25UL |
| Anti-MEK1, phospho-specific (Thr386)                    |       | ΗМ                         | WB                 | APur   | Rb     | 100 µL   | AB4209      |
| Anti-MEK1, phospho-specific (Thr292)                    |       | ΗМ                         | WB                 | APur   | Rb     | 100 µL   | AB4210      |
| Assay and cDNA  |       |                            |                    |        |        |          |             |
| MEK1 Assay Kit  |       |                            | KA                 |        |        | 1 kit    | 17-157      |
| MEK1 cDNA (wt) in pUSEamp                               |       |                            | TFX                |        |        | 5 µg     | 21-106      |
| MEK1 cDNA (wt) in pUSEamp                               |       |                            | TFX                |        |        | 5 µg     | 21-107      |
| MEK1 cDNA (activated) in pUSEamp                        |       |                            | TFX                |        |        | 5 µg     | 21-119      |
| MEK1 cDNA (activated) in pUSEamp                        |       |                            | TFX                |        |        | 5 µg     | 21-120      |
| siRNA   |       |                            |                    |        |        |          |             |
| MEK1 siRNA/siAb Assay Kit                               |       | Н                          | WB RNAi            |        |        | 1 kit    | 60-096      |
| siRNA plasmid, pKD-MEK1-v2                              |       | H M R                      | RNAi               |        |        | 5 µg     | 62-152      |
| MEK1 SMARTpool siRNA reagent                            |       | Н                          | RNAi               |        |        | 5 nmol   | M-003571    |
| Protein   |       |                            |                    |        |        |          |             |
| MEK1 (6His), unactive                                   |       |                            | KA                 |        |        | 50 µg    | 14-706      |
| MEK1, active  |       |                            | KA                 | Pur    |        | 10 µg    | 14-429      |
| MEK1, active, rabbit                                    |       |                            | KA                 |        |        | 10 µg    | 14-206      |
| MEK1 (K97R), inactive                                   |       |                            | EA                 |        |        | 50 µg    | 14-737      |
| MEK1, unactive  |       |                            | KA                 |        |        | 50 µg    | 14-420      |
| MEK1, unactive, rabbit                                  |       |                            | IPK KA             | Pur    |        | 50 µg    | 14-205      |
| MEK2, active  |       |                            | KA                 |        |        | 10 µg    | 14-528      |
| MEK2 (K101M), inactive                                  |       |                            | KA                 |        |        | 50 µg    | 14-709      |
| MEK2, unactive  |       |                            | KA                 |        |        | 25 µg    | 14-541      |
| MEK2, unactive  |       |                            | KA                 |        |        | 50 µg    | 14-532      |

#### **MEKK**

MEKK1 (Mitogen-Activated Protein Kinase Kinase Kinase 1 or MAP3K1) mediates the responses to a number of metabolic and mitogenic signals and is involved in JNK-mediated apoptosis (MEKK1 seems to have an anti-apoptotic role). MEKK1 also exhibits E3 Ubiquitin Ligase activity towards Erk2. MEKK3 (MAP3K3) regulates the activity of MEK5 and BMK1/Erk5 during growth factor-induced cellular stimulation. BMK1 (Big Mitogen-Activated Protein Kinase, aka Erk5) is a member of the MAP Kinase family whose cellular activity is elevated in response to growth factors, oxidative stress, and hyperosmolar conditions. MEKK3 regulates BMK1 activity through its regulation of MEK5 activity. MEKK3 may also regulate NFκB-dependent gene transcription.

#### Antibody

| Anti-MEKK1                  | H M R   | WB not IP | APur | Rb     | 50 µg  | AB3186 |  |
|-----------------------------|---------|-----------|------|--------|--------|--------|--|
| Anti-MEKK2 CT, Rb Mab       | H not M | FC IP WB  | Pur  | Rb IgG | 100 µL | 04-378 |  |
|                             | not R   | IC IH(P)  |      |        |        |        |  |
| Anti-MEKK3                  | H M R   | WB        |      | Rb IgG | 100 μL | 07-691 |  |
| Anti-MEKK3 (N-term)         | Н       | WB        |      | Rb IgG | 100 μL | 04-379 |  |
| Protein                     |         |           |      |        |        |        |  |
| MEKK1, active               |         | KA        | Pur  |        | 10 µg  | 14-196 |  |
| siRNA                       |         |           |      |        |        |        |  |
| siRNA plasmid, pKD-MEKK2-v1 |         | RNAi      |      |        | 5 µg   | 62-260 |  |
|                             |         |           |      |        |        |        |  |

#### MKK

Members of the p38 Mitogen-Activated Protein Kinase family are activated by treatment of cells with cytokines and by exposure to environmental stress. The effects of these stimuli on p38 MAP Kinase are mediated by the MAP Kinase Kinases (MKKs) MKK3, MKK4, and MKK6. Studies have shown that knocking out MKK3 protein expression causes defects in the response of fibroblasts to the proinflammatory cytokine TNF. MKK3 appears to be a critical component of a TNF-stimulated signaling pathway that causes increased expression of inflammatory cytokines.

| causes increased expression of inflammatory cytokines.     |       |         |                   |        |        |                |          |
|--|-------|---------|-------------------|--------|--------|----------------|----------|
| Description  | Clone | Species | Applications      | Format | Host   | Quantity       | Cat. No. |
| NKK3   |       |         |                   |        |        |                |          |
| Antibody   |       |         |                   |        |        |                |          |
| Anti-MKK3 NT, Rb Mab                                       |       | Н       | FC IP WB          | Pur    | Rb IgG | 100 µL         | 04-381   |
| siRNA plasmid, pKD-MKK3-v3                                 |       |         | RNAi              |        | . 0    | 5 µg           | 62-261   |
| siRNA plasmid, pKD-MKK3-v4                                 |       |         | RNAi              |        |        | 5 μg           | 62-262   |
| siRNA plasmid, pKD-MKK3-v4                                 |       |         | RNAi              |        |        | 5 µg           | 62-262   |
| siRNA plasmid, pKD-MKK3-v3                                 |       |         | RNAi              |        |        | 5 μg           | 62-261   |
|  |       |         |                   |        |        | 10             |          |
| AKK4   |       |         |                   |        |        |                |          |
| Antibody   |       |         |                   |        |        |                |          |
| Anti-MKK4 CT, Rb Mab                                       |       | H M R   | FC IP WB<br>IH(P) | Pur    | Rb IgG | 100 μL         | 04-382   |
| Anti-phospho-MKK4 (Ser257/Thr261)                          |       | Н       | WB                |        | Rb IgG | 100 pL         | 07-857   |
| rotein   |       |         |                   |        |        |                |          |
| MKK4/SKK1, active  |       |         | KA                | Pur    |        | 10 µg          | 14-377   |
| MKK4/SKK1, unactive  |       |         | KA                |        |        | 50 µg          | 14-378   |
| irna   |       |         |                   |        |        |                |          |
| siRNA plasmid, pKD-MKK4/SKK1-v2                            |       |         | RNAi              |        |        | 5 µg           | 62-263   |
| siRNA plasmid, pKD-MKK4/SKK1-v4                            |       |         | RNAi              |        |        | 5 µg           | 62-264   |
| siRNA plasmid, pKD-MKK4/SKK1-v4                            |       |         | RNAi              |        |        | 5 µg           | 62-264   |
| siRNA plasmid, pKD-MKK4/SKK1-v2                            |       |         | RNAi              |        |        | 5 μg           | 62-263   |
| NKK6   |       |         |                   |        |        |                |          |
| Antibody   |       |         |                   |        |        |                |          |
| Anti-MKK6/SKK3   |       | H M R   | WB                | Pur    | Rb IgG | 200 µg         | 07-417   |
|  |       | H M R   | IP VVB            |        | _      | . 0            |          |
| Anti-MKK6 NT, Rb Mab<br>rotein                             |       | □ /V\ K | IF VVD            | Pur    | Rb IgG | 100 pL         | 04-383   |
| MKK6/SKK3, active  |       |         | KA                |        |        | 10 µg          | 14-303   |
| MKK6/SKK3, unactive  |       |         | KA                |        |        | 10 pg<br>50 µg | 14-304   |
| MKK6/SKK3 (S599D, T603D), active                           |       |         | KA                |        |        | 30 рд<br>10 рд | 14-537   |
| iRNA   |       |         | IV-A              |        |        | то ру          | 14307    |
| MKK6 SMARTpool siRNA reagent                               |       | Н       | RNAi              |        |        | 5 nmol         | M-003967 |
| MKK6 siRNA/siAb Assay Kit                                  |       | Н       | WB RNAi           |        |        | 1 kit          | 60-094   |
| siRNA plasmid, pKD-MKK6/SKK3-v1                            |       | ''      | RNAi              |        |        | 5 µg           | 62-200   |
| siRNA plasmid, pKD-MKK6/SKK3-v2                            |       |         | RNAi              |        |        | 5 µg           | 62-201   |
| siRNA plasmid, pKD-MKK6/SKK3-v2                            |       |         | RNAi              |        |        | 5 µg           | 62-201   |
| sint of plasma, profitting, sinter ve                      |       |         | KI W II           |        |        | 3 pg           | 02 201   |
| NKK7   |       |         |                   |        |        |                |          |
| Antibody   |       |         |                   |        |        |                |          |
| Anti-MKK7/SKK4   |       | Н       | WB                |        | Rb IgG | 100 µg         | 07-399   |
| Anti-phospho-MKK7/SKK4 (Thr275/Ser277)                     |       | Н       | WB                |        | Rb IgG | 100 pL         | 36-010   |
| Anti-phospho-MKK7/SKK4 (Thr275)                            |       | Н       | WB                | Pur    | Rb IgG | 200 µg         | 36-013   |
| rotein   |       |         |                   |        |        |                |          |
| MKK $7lpha$ 1, active                                      |       |         | KA                | APur   |        | 25 µg          | 14-301   |
| MKK7 $\beta$ 1, active                                     |       |         | KA                |        |        | 10 µg          | 14-543   |
| MKK7 $eta$ 1, unactive                                     |       |         | KA                |        |        | 50 µg          | 14-542   |
| irna   |       |         |                   |        |        |                |          |
| MKK7 siRNA/siA Assay Kit                                   |       | Н       | WB RNAi           |        |        | 1 kit          | 60-095   |
| siRNA plasmid, pKD-MKK7-v2                                 |       | ΗМ      | RNAi              |        |        | 5 µg           | 62-154   |
| siRNA plasmid, pKD-MKK7-v5                                 |       | ΗМ      | RNAi              |        |        | 5 µg           | 62-155   |
| MKK7 SMARTpool siRNA reagent                               |       | Н       | RNAi              |        |        | 5 nmol         | M-004016 |
| NLCK (Myosin Light Chain Kinase)                           |       |         |                   |        |        |                |          |
| Antibody   |       | וח ח    | \                 | Λ D    | DL I C | 100            | 07.040   |
| Anti-phospho-Myosin Light Chain Kinase<br>(MLCK) (Ser1760) |       | B Rb    | VVB               | APur   | Rb IgG | 100 μL         | 07-860   |
| MLCK, active   |       |         | KA                |        |        | 10 µg          | 14-638   |

#### MKP-1

MKP-1 belongs to a protein-tyrosine phosphatase family that dephosphorylats MAP Kinase 2 / ERK2 on both Thr183 and Tyr185 sites. MKP-1 is induced by oxidative stress and heat shock-inducible gene and is thought to play a regulatory role in the MAPK and stress-activated protein kinase signaling pathways.

|   | Description | Clone | Species | Applications | Format | Host   | Quantity | Cat. No. |
|---|-------------|-------|---------|--------------|--------|--------|----------|----------|
| A | Antibody    |       |         |              |        |        |          |          |
|   | Anti-MKP1   |       | НМ      | WB           | Pur    | Rb IgG | 100 µg   | 07-535   |

#### MSK-1

MSK1 (mitogen-and-stress-activated protein kinase 1) contains two protein kinase domains in a single polypeptide. It is activated *in vitro* by MAPK2/Erk2 or SAPK2/p38. Endogenous MSK1 is activated in 293 cells by either growth factor/phorbol ester stimulation, or by exposure to UV radiation and oxidative and chemical stress. MSK1 may mediate the growth factor and stress-induced activation of CREB.

| Description                | Clone | Species | Applications Format | Host   | Quantity | Cat. No. |
|----------------------------|-------|---------|---------------------|--------|----------|----------|
| Antibody                   |       |         |                     |        |          |          |
| Anti-phospho-MSK1 (Ser376) |       | Н       | WB                  | Rb IgG | 100 µL   | 04-384   |

#### p38 ( $\alpha$ , $\beta$ , $\delta$ , $\gamma$ )

The p38 family refers to the SAPK2/3 Kinases, a subfamily of JNK/SAPK family. p38/SAPK Kinases are also activated by stresses, most notably inflammatory cytokines, irradiation, and certain toxins such as anisomycin, and arsenite. The activating kinases of SAPK2/3 are SKK2/MEK3 for SAPK2a and 2b, and MKK6 for SAPK3. The targets of the SAPK2/3 family include the MAPKAP Kinases 2 and 3/3pK. In addition, SKK4 is related to this family, exhibiting 60% identity, and is activated by MKK6. SAPKs have been implicated in many pathological conditions, including inflammation, cancer, and neurodegenerative diseases.

|       | Description   | Clone | Species | Applications | Format    | Host   | Quantity  | Cat. No. |
|-------|---|-------|---------|--------------|-----------|--------|-----------|----------|
| lono  | clonal Antibody   |       |         |              |           |        |           |          |
|       | Anti-p38/SAPK2  | 2F11  | H M R   | IP WB        | Pur       | M lgG  | 200 µg    | 05-454   |
| olycl | onal Antibody   |       |         |              |           |        |           |          |
|       | Anti-phospho-p38α (Thr180/Tyr182)   |       | H M R   | WB IF IP IPK | APur      | Rb IgG | 100 μL    | 09-470   |
|       | Anti-p38a   |       | HMRB    | WB IP IPK KA | APur      | Rb IgG | 100 µg    | 09-272   |
|       | Anti-p388/SAPK4   |       | H R     | WB           | Pur       | Rb IgG | 200 µg    | 07-603   |
|       | Anti-p388/SAPK4   |       | H R Rb  | IP WB        | APur      | Sh IgG | 100 µg    | 06-652   |
|       | Anti-p38γ/SAPK3   |       | M R     | WB           | Pur       | Rb IgG | 200 µg    | 07-139   |
|       | Anti-p38γ/SAPK3   |       | Н       | WB           | Pur       | Rb IgG | 200 µg    | 07-474   |
|       | Anti-p38γ/SAPK3   |       |         | WB           | antiserum | Rb IgG | 200 μL    | 07-508   |
| oteii | n   |       |         |              |           |        |           |          |
|       | p38 $lpha$ /SAPK2a, active  |       |         | KA           | Pur       |        | 10 µg     | 14-210   |
|       | p38 $lpha$ /SAPK2a, active  |       |         | KA           |           |        | 10 µg     | 14-251   |
|       | p38 $lpha$ /SAPK2a, unactive  |       |         | KA           |           |        | 50 µg     | 14-252   |
|       | p38α/SAPK2a (6His-tag), active  |       |         | KA           |           |        | 10 µg     | 14-587   |
|       | p38α/SAPK2a (6His-tag), unactive  |       |         | KA           |           |        | 50 µg     | 14-588   |
|       | p38 $\alpha$ /SAPK2a (T106M), active  |       |         | KA           | Pur       |        | 10 µg     | 14-687   |
|       | p38β2/SAPK2b2, active   |       |         | KA           |           |        | 10 µg     | 14-253   |
|       | p38γ/SAPK3, active  |       |         | KA           |           |        | 10 µg     | 14-246   |
|       | p38δ/SAPK4, active  |       |         | KA           |           |        | 10 µg     | 14-249   |
| RNA   |   |       |         |              |           |        |           |          |
|       | siRNA plasmid, pKD-p38α/SAPK2a-v1   |       | Н       | RNAi         |           |        | 5 µg      | 62-294   |
|       | siRNA plasmid, pKD-p38α/SAPK2a-v3   |       | H M R   | RNAi         |           |        | 5 µg      | 62-295   |
|       | siRNA plasmid, pKD-p38β2/SAPK2b2-v1   |       |         | RNAi         |           |        | 5 µg      | 62-208   |
|       | siRNA plasmid, pKD-p38β2/SAPK2b2-v4   |       |         | RNAi         |           |        | 5 µg      | 62-209   |
|       | siRNA plasmid, pKD-Sapk2b-v1  |       | Н       | RNAi         |           |        | 5 µg      | 62-328   |
|       | siRNA plasmid, pKD-Sapk2b-v2  |       | Н       | RNAi         |           |        | 5 µg      | 62-329   |
|       | p38γ SMARTpool siRNA reagent  |       | Н       | RNAi         |           |        | 5 nmol    | M-003590 |
|       | siRNA plasmid, pKD-p38γ/SAPK3-v3  |       | Н       | RNAi         |           |        | 5 µg      | 62-118   |
|       | siRNA plasmid, pKD-p38&/SAPK4-v3  |       |         | RNAi         |           |        | 5 µg      | 62-220   |
|       | siRNA plasmid, pKD-p38&/SAPK4-v4  |       |         | RNAi         |           |        | 5 µg      | 62-221   |
|       | p38/SAPK2 SMARTpool siRNA reagent   |       |         | RNAi         |           |        | 5 nmol    | M-003512 |
|       | p38/SAPK2 siRNA/siAb Assay Kit  |       | Н       | WB RNAi      |           |        | 1 kit     | 60-032   |
| ssay  |   |       |         |              |           |        |           |          |
|       | Non-radioactive p38 Kinase Assay Kit  |       | H M R   | Activity     |           |        | 40 assays | SGT455   |
|       | p38α/SAPK2a Assay Kit, 50 Assays  |       |         | KA           |           |        | 1 kit     | 17-169   |
|       | p $38lpha/SAPK2a$ KinEASE FP Fluorescein Green Assa   | ly    |         | KA           |           |        | 1 kit     | 32-046   |
|       | p38α/SAPK2a KinEASE FP-645nm FarRed Assay   |       |         | KA           |           |        | 1 kit     | 32-126   |
|       | p38β2/SAPK2b2 KinEASE FP Fluorescein Green A  | ssay  |         | KA           |           |        | 1 kit     | 32-047   |
|       | Professional Control of the Control |       |         |              |           |        |           |          |
|       | p38β2/SAPK2b2 KinEASE FP-645nm FarRed Assay   | у     |         | KA           |           |        | 1 kit     | 32-127   |

Lysate fr 1 c trai (Thr 1)

| Description                                 | Clone | Species | Applications Format | Host | Quantity | Cat. No. |
|---|-------|---------|---------------------|------|----------|----------|
| p388/SAPK4 KinEASE FP-645nm FarRed Assay    |       |         | KA                  |      | 1 kit    | 32-129   |
| p38γ/SAPK3 KinEASE FP Fluorescein Green Ass | ay    |         | KA                  |      | 1 kit    | 32-048   |
| p38γ/SAPK3 KinEASE FP-645nm FarRed Assay    |       |         | KA                  |      | 1 kit    | 32-128   |
| Inhibitors                                  |       |         |                     |      |          |          |
| p38/SAPK2 Inhibitor (SB 202190)             |       |         | KA                  |      | 1 mg     | 19-134   |
| p38/SAPK2 Inhibitor (SB 203580)             |       |         | KA                  |      | 1 mg     | 19-135   |

#### PAK

p21-activated kinases (PAK) are Ser/Thr kinases, which are activated by, and are effectors of Rho family of GTPases. As such, they are implicated in Actin polymerization, and activation of the stress-activated kinase cascades. PAK1 is phosphorylated by cdk5, resulting in its inhibition. PAK2 can be activated by Caspase-cleavage in addition to p21-binding, implicating it in cytoskeletal changes during apoptosis. PAK3 is expressed at high levels in post-mitotic neurons of the developing post-natal cerebral cortex and hippocampus, and is the mutant gene thought responsible for non-syndromic X-linked mental retardation.

| Description                            | Clone        | Species  | Applications   | Format | Host   | Quantity | Cat. No. |
|--|--------------|----------|----------------|--------|--------|----------|----------|
| Antibody                               |              |          |                |        |        |          |          |
| Anti-PAK1, Rb Mab                      |              | H M R    | IF IP WB IH(P  |        |        | 100 μL   | 04-394   |
| Anti-PAK3, Rb Mab                      |              | H M R    | IF IP VVB IH(P |        | Rb IgG | 100 μL   | 04-395   |
| Anti-PAK3, NT                          |              | M R      | IP WB          | Pur    | Rb IgG | 200 µg   | 06-902   |
| Anti-PAK1 (Ser144)/PAK2 (Ser141), phos | pho-specific | Н        | WB             | APur   | Rb     | 100 μL   | AB3833   |
| Anti-PAK1 (Ser199/204)/PAK2 (Ser192,   | /197),       | Н        | WB             | APur   | Rb     | 100 μL   | AB3834   |
| phospho-specific                       |              |          |                |        |        |          |          |
| Anti-PAK2, phospho-specific (Ser20)    |              | Н        | WB IH(P)       | APur   | Rb     | 100 μL   | AB3836   |
| Anti-PAK 1                             |              | H M Mk R | IP WB IH(P)    | APur   | Rb     | 100 μL   | AB3844   |
| Anti-PAK 1/2/3                         |              | H M Mk R | WB             | APur   | Rb     | 100 μL   | AB3845   |
| PAK-PBD GST Protein, Rac/Cdc42 bindin  | g domain     |          |                | Pur    |        |          | SGT224   |
| Protein                                |              |          |                |        |        |          |          |
| PAK2, active                           |              |          | KA             |        |        | 10 µg    | 14-481   |
| PAK3, active                           |              |          | KA             |        |        | 10 µg    | 14-683   |
| PAK4, active                           |              |          | KA             |        |        | 10 µg    | 14-584   |
| PAK5, active                           |              |          | KA             |        |        | 10 µg    | 14-699   |
| PAK6, active                           |              |          | KA             |        |        | 10 µg    | 14-633   |
| Assay                                  |              |          |                |        |        |          |          |
| PAK2 KinEASE FP Fluorescein Green Assa | У            |          | KA             |        |        | 1 kit    | 32-018   |
| PAK2 KinEASE FP-645nm FarRed Assay     |              |          | KA             |        |        | 1 kit    | 32-098   |

#### Paxillin

Paxillin is a 68 kDa focal adhesion protein that is phosphorylated on tyrosine residues in response to Integrin stimulation, growth factor stimulation, and oncogenic transformation (v-Src, v-Crk and BCR/ABL). Paxillin binds other focal adhesion proteins such as Vinculin, Talin, Tensin and FAK.

|        | Description                          | Clone | Species  | Applications | Format | Host               | Quantity | Cat. No. |
|--------|--------------------------------------|-------|----------|--------------|--------|--------------------|----------|----------|
| Mono   | clonal Antibody                      |       |          |              |        |                    |          |          |
|        | Anti-Paxillin                        | 5H11  | HMRBAv   | WB IP IC     | Pur    | M lgG <sub>1</sub> | 250 µg   | 05-417   |
|        | Paxillin                             | 349   | H M R    | WB ELISA     | Pur    | M lgG <sub>1</sub> | 100 µg   | MAB3060  |
|        | Anti-phospho-Paxillin (Tyr 1 1 8)    | 30    | НМ       | WB           | Pur    | M lgG <sub>1</sub> | 100 µg   | MAB1145  |
|        | Paxillin, phospho-specific [Tyr31]   | M102  | Н        | WB ELISA     |        | M lgG <sub>1</sub> | 100 μL   | MAB1146  |
| Polycl | onal Antibody                        |       |          |              |        |                    |          |          |
|        | Paxillin                             |       | H M R Rb | WB IP        | Pur    | Rb                 | 50 µL    | AB3794   |
|        |                                      |       | Ca Ch    | IC IH        |        |                    |          |          |
|        | Paxillin (Tyr118)                    |       | Μ        | WB           | APur   | Rb IgG             | 100 μL   | 07-733   |
|        | Paxillin [pS126]                     |       |          |              | APur   |                    |          | AB3793   |
|        | Paxillin, phospho-specific [Ser 178] |       | Н        | WB ELISA     |        | Rb IgG             | 100 pL   | AB1962   |
|        | Phospho-Paxillin                     |       | H M R    | WB IH(P)     | APur   | Rb                 | 100 pL   | AB3837   |
|        |                                      |       |          |              |        |                    |          |          |

#### PI3 Kinase

P13 Kinase (Phosphatidylinositol 3-Kinase) is responsible for phosphorylation of the 3 position of the inositol ring of PI[4,5]P2, to generate PI[3,4,5]P3, a potent second messenger required for survival signaling, and insulin action. P13 Kinase is a heterodimeric complex composed of an 85 kDa regulatory subunit and a 110 kDa catalytic subunit. Tyrosine phosphorylation of growth factor receptors creates docking sites for binding of p85 (through its SH2 domains) on the receptors; p85 brings with it p110, which is then proximal to its phospholipid substrate on the membrane. P13 Kinase is also activated by Ras, and by the  $\beta$ : $\gamma$  subunits of heterotrimeric G-proteins. P13 Kinase is inhibitable by wortmannin, a useful tool for the study of the P13 Kinase signaling pathway.

|   | Description                 | Clone  | Species | Applications | Format | Host                | Quantity | Cat. No. |
|---|-----------------------------|--------|---------|--------------|--------|---------------------|----------|----------|
| 1 | Monoclonal Antibody         |        |         |              |        |                     |          |          |
|   | Anti-Pl3 Kinase, p85, N-SH3 | AB6    | ΗМ      | IP WB IC     | Pur    | M lgG <sub>1κ</sub> | 100 µg   | 05-212   |
|   | Anti-Pl3 Kinase, p85, N-SH2 | UB93-3 | H M R   | IP WB        | Sup    | M lgG               | 200 μL   | 05-217   |
|   | Anti-Pl3 Kinase, p110δ      | AW103  | Н       | WB           | Pur    | M lgG <sub>1κ</sub> | 200 µg   | 05-703   |

# Phosphorylation and Phosphotyrosine, clone 4G10®

#### **Phosphorylation**


Many aspects of cell activity are controlled by reversible phosphorylation of target proteins that are mediated by protein kinases and their subsequent dephosphorylation by phosphatases. In many cases, mutations or misregulation of members of these two families of enzymes, kinase and phosphatases, are indicated in various disease states. Protein kinases play crucial roles in the regulation of many cellular events such as signal transduction, cell cycle progression, targeted proteolysis, protein trafficking, cytoskeletal organization, and gene expression. A protein's activity may be enabled or inhibited by the phosphorylation of specific serine, threonine, or tyrosine amino acids by protein kinases resulting in altered conformations.

#### Anti-Phosphotyrosine, clone 4G10


### The Gold Standard in Anti-Phosphotyrosine Antibodies

The development of the anti-Phosphotyrosine clone 4G10 in 1989 was a monumental discovery for researchers as it facilitated the study of tyrosine kinase activity by removing the need for radioactivity and converting the method to a simple blotting procedure. Millipore's clone 4G10 was the first and is still the most highly cited and validated phosphotyrosine antibody to consistently produce high quality, reproducible data. Over the past 17 years, clone 4G10 has been validated by thousands of scientific and medical researchers in virtually every application and tyrosine target including western blot (WB), flow cytometry (FC), fluorescence polarization (FP), ELISA, Homogeneous Time Resolved Fluorescence (HTRF®), immunocytochemistry (ICC), immunohistochemistry (IHC) and immunoprecipitation (IP). Clone 4G10 is now a critical tool for researchers in areas as diverse as serine/threonine phospho-proteins (including p38, Akt, and ERK 1/2), G-proteins, and neurological targets. For your full range of research needs, clone 4G10 is available as a number of conjugates including agarose, biotin, FITC, and HRP as well as the newly released

Catch and Release® Phosphotyrosine, clone 4G10 for easier, cleaner immunoprecipitation of your tyrosine phosphorylated material. In addition to its unparalleled characterization and applicability, Clone 4G10 is well known as one of the most sensitive phosphotyrosine antibodies, detecting at least twice as many phosphoproteins in a side-by-side comparison with PY20 and PT-66 by independent investigators. This sensitivity is also observed in other applications, including FP, where clone 4G10 and recombinant clone 4G10 from Millipore were shown to be more sensitive that other anti-phosphotyrosine antibodies tested.



EGF treated A431 cells that are triple stained with anti-Phosphotyrosine, clone 4G10 (green), anti-actin Alexa Fluor® 555, (red), and Dapi nuclei, (blue).



Clone 4G10 IHC Stain: Paraffin embedded poorly-differentiated adenocarcinoma, stained with anti-phosphotyrosine, clone 4G10 (Cat. No. 05-321) at 10 µg/mL.

#### **Antibodies**

| Description                             | Quantity | Cat. No. |
|---|----------|----------|
| Phosphotyrosine, clone 4G10             | 100 µg   | 05-321   |
| Phosphotyrosine, clone 4G10             | 1 mg     | 05-321MG |
| Phosphotyrosine, clone 4G10             | 50 µg    | 05-321X  |
| Phosphotyrosine, recombinant clone 4G10 | 1 mg     | 05-777   |

#### Conjugates

| Description   | Quantity  | Cat. No. |
|---|-----------|----------|
| Phosphotyrosine, clone 4G10, agarose conjugate          | 1 mg      | 16-101   |
| Phosphotyrosine, clone 4G10, HRP conjugate              | 10 blots  | 16-105   |
| Phosphotyrosine, clone 4G10, FITC conjugate             | 100 µg    | 16-104   |
| Phosphotyrosine, clone 4G10, biotin conjugate           | 100 µg    | 16-103   |
| Phosphotyrosine, recombinant clone 4G10, HRP conjugate  | 10 blots  | 16-184   |
| Phosphotyrosine, recombinant clone 4G10,                | 1 mg      | 16-199   |
| agarose conjugate                                       |           |          |
| Phosphotyrosine, recombinant clone 4G10, FITC conjugate | 100 µg    | 16-205   |
| Catch and Release Phosphotyrosine, clone 4G10           | 50 assays | 17-502   |

For a complete listing of products, please visit our website at **www.millipore.com**.

| D        | Description                             | Clone   | Species  | Applications          | Format    | Host                | Quantity | Cat. No. |
|----------|---|---------|----------|-----------------------|-----------|---------------------|----------|----------|
| A        | Anti-Pl3 Kinase, p110α, Rb Mab          |         | Н        | IF IP WB              |           | Rb IgG              | 100 µL   | 04-399   |
| A        | Anti-Pl3 Kinase, p110β, Rb Mab          |         | Н        | FC IP WB              |           | Rb IgG              | 100 µL   | 04-400   |
| A        | Anti-Pl3 Kinase, p110δ;, Rb Mab         |         | Н        | IF IP WB              |           | Rb IgG              | 100 µL   | 04-401   |
| A        | Anti-Pl3 Kinase, p110γ, Rb Mab          |         | Н        | FC IP WB              |           | Rb IgG              | 100 μL   | 04-402   |
| A        | Anti-Pl3 Kinase, p85α, Rb Mab           |         | H M R    | FC IP VVB<br>IC IH(P) |           | Rb IgG              | 100 µL   | 04-403   |
| A        | Anti-Pl 3-Kinase p85α                   | 8-2D-4D | H M R    | FC IF IP              | Pur       | M lgG <sub>1</sub>  | 100 µg   | MAB1143  |
| Polyclon | nal Antibody                            |         |          |                       |           |                     |          |          |
| A        | Anti-Pl3 Kinase, p85                    |         | H M Mk R | IP WB                 | antiserum | Rb IgG              | 125 µL   | 06-195   |
| A        | Anti-P13 Kinase, p85, N-SH2 domain      |         | H M Mk R | IP WB                 | Pur       | Rb IgG              | 250 µg   | 06-496   |
| A        | Anti-Pl3 Kinase, p85                    |         | H M Mk R | IP WB                 | Pur       | Rb IgG              | 250 µg   | 06-497   |
| A        | Anti-Pl3 Kinase, p110b                  |         | B H M R  | IP                    | Pur       | Rb IgG              | 200 µg   | 06-568   |
| A        | Anti-Pl3 Kinase, p101                   |         | H R      | WB                    | Pur       | Rb IgG <sub>1</sub> | 200 µg   | 07-281   |
| A        | Anti-Pl3 Kinase, p110a                  |         | H WR     | IP                    | Pur       | Rb IgG              | 200 µg   | 07-658   |
| A        | Anti-PI4-Kinase b                       |         | H M R    | IP WB IC              | Pur       | Rb IgG              | 200 µg   | 06-578   |
| siRNA    |   |         |          |                       |           |                     |          |          |
| Р        | P13 Kinase p85a siRNA/siAb Assay Kit    |         | Н        | WB RNAi               |           |                     | 1 kit    | 60-048   |
| S        | siRNA plasmid, pKD-PI3 Kinase, p85-v3   |         |          | RNAi                  |           |                     | 5 μg     | 62-222   |
| Р        | P13 Kinase p85a SMARTpool siRNA reagent |         | Н        | RNAi                  |           |                     | 5 nmol   | M-003020 |
| Protein  |   |         |          |                       |           |                     |          |          |
| Р        | PI3 Kinase (p110β/p85a), active         |         |          | KA                    |           |                     | 20 µg    | 14-603   |
| Р        | P13 Kinase (p1108/p85a), active         |         |          | KA                    |           |                     | 10 µg    | 14-604   |
| Р        | Pl3 Kinase (p110α/p85a), active         |         |          | KA                    |           |                     | 20 µg    | 14-602   |
| Р        | P13 Kinase (p120γ)                      |         |          | KA                    |           |                     | 20 µg    | 14-558   |
| $\vee$   | Wortmannin, (PI3 Kinase inhibitor)      |         |          | KA                    |           |                     | 1 mg     | 12-338   |
| Assay    |   |         |          |                       |           |                     |          |          |
| Р        | PI 3-Kinase HTRF Assay                  |         |          | KA                    |           |                     | 1 plate  | 33-016   |
| P        | PI 3-Kinase HTRF Assay                  |         |          | KA                    |           |                     | 5 plates | 33-017   |

#### PKA

Protein kinase A is perhaps the prototype Ser/Thr kinase, which is activated by the second messenger cyclic AMP (cAMP). The PKA holoenzyme consists of two catalytic subunits complexed with two regulatory subunits. In addition, the regulatory subunits bind various AKAP proteins, which localize the holoenzyme to specific locations within the cell. Binding of cAMP to the regulatory subunits results in release of catalytic subunits, which phosphorylate an extremely broad range of substrates, including the transcription factor CREB. PKA activity in a lysate can be neutralized with the pseudosubstrate PKI.

|        | Description  | Clone | Species          | Applications | Format | Host   | Quantity | Cat. No.        |
|--------|--|-------|------------------|--------------|--------|--------|----------|-----------------|
| Antibo | ody  |       |                  |              |        |        |          |                 |
|        | Anti-PKA, RII Subunits                             |       | H R              | IP WB IC     | Pur    | Gt IgG | 500 µg   | 06-411          |
|        | Anti-phospho-PKA, RII (Ser96)                      |       | M R              | WB           | Pur    | Rb IgG | 100 µg   | 06-704          |
|        | Anti-PKA, NT                                       |       | B H Ht<br>M Po R | WB           | Pur    | Rb IgG | 200 µg   | 06-903          |
|        | Anti-phospho-PKA catalytic subunits a/b (Thr197)   |       | M                | WB           | APur   | Rb IgG | 100 µL   | 07-867          |
|        | Anti-phospho-PKA Catalytic b subunit (Ser338)      |       | Μ                | WB           | APur   | Rb IgG | 100 µL   | 07-868          |
|        | Anti-phospho-PKA, Regulatory subunit IIb (Ser1 14) |       | Μ                | WB           | APur   | Rb IgG | 100 µL   | 07-869          |
| Assay  |  |       |                  |              |        |        |          |                 |
|        | PKA Assay Kit                                      |       |                  | KA           |        |        | 1 kit    | 1 <i>7</i> -134 |
|        | PKA KinEASE FP Fluorescein Green Assay             |       |                  | KA           |        |        | 1 kit    | 32-020          |
|        | PKA KinEASE FP-645nm FarRed Assay                  |       |                  | KA           |        |        | 1 kit    | 32-100          |
| Reage  | ents   |       |                  |              |        |        |          |                 |
|        | PKA Inhibitor Cocktail                             |       |                  | KA           |        |        | 1 mL     | 20-114          |
|        | PKA Inhibitor Peptide                              |       |                  | KA           |        |        | 1 mL     | 20-120          |
|        | PKA/PKC Inhibitor Cocktail                         |       |                  | KA           |        |        | 1 mL     | 20-129          |
|        | PKA/CaMK Inhibitor Cocktail                        |       |                  | KA           |        |        | 1 mL     | 20-132          |
| siRNA  | N.   |       |                  |              |        |        |          |                 |
|        | PKA siRNA/siAb Assay Kit                           |       | Н                | WB RNAi      |        |        | 1 kit    | 60-115          |
|        | siRNA plasmid, pKD-PKAa-v2                         |       | Н                | RNAi         |        |        | 5 µg     | 62-069          |
|        | siRNA plasmid, pKD-PKAa-v3                         |       | Н                | RNAi         |        |        | 5 µg     | 62-070          |
|        | Anti-cAMP-Dependent Protein Kinase,                |       | H M Mk           | WB IH        | APur   | Rb     | 100 µg   | AB1612          |

|        | Description                              | Clone | Species | Applications | Format | Host           | Quantity | Cat. No. |
|--------|--|-------|---------|--------------|--------|----------------|----------|----------|
|        | Regulatory Subunit I- $\beta$ , internal |       |         |              |        |                |          |          |
|        | Anti-cAMP-Dependent Protein Kinase,      |       | HMMk    | ELISA FC IP  | APur   | Rb             | 100 µg   | AB1613   |
|        | Regulatory Subunit II- $lpha$            |       |         | WB IH        |        |                |          |          |
|        | Anti-cAMP-Dependent Protein Kinase,      |       | НМ      | WB IH        | APur   | Rb             | 100 µg   | AB1614   |
|        | Regulatory Subunit II-β, internal        |       | Mk Ech  |              |        |                |          |          |
|        | PKA SMARTpool siRNA reagent              |       |         | RNAi         |        |                | 5 nmol   | M-004649 |
|        | Anti-cAMP-Dependent Protein Kinase II,   | RS30  | ВНРо    | ELISA IF WB  | Pur    | $M \lg G_{2a}$ | 100 µg   | MAB1697  |
|        | Regulatory Subunit $lpha/eta$            |       |         |              |        |                |          |          |
| Protei | n  |       |         |              |        |                |          |          |
|        | PKA, catalytic subunit                   |       |         | KA           |        |                | 10 µg    | 14-114   |
|        | PKA, catalytic subunit, recombinant      |       |         | KA           |        |                | 10 µg    | 14-440   |
|        | PKA Substrate Peptide                    |       |         | KA           |        |                | 2 mg     | 12-152   |
|        | PKA Substrate Peptide, biotin conj.      |       |         | KA           |        |                | 500 µg   | 12-394   |
|        | PKA Inhibitor peptide                    |       |         | KA           |        |                | 2 mg     | 12-151   |

#### PKC

Members of the PKC family fall into three categories: conventional (PKC  $\alpha$ ,  $\beta$ I,  $\beta$ II,  $\gamma$ ), Novel ( $\delta$ ,  $\epsilon$ ,  $\eta$ ,  $\mu$ ), and atypical ( $\lambda$ ,  $\xi$ ). Conventional PKC isoforms are dependent on Ca<sup>2+</sup> and diacylglycerol (DAG) for their activity, whereas novel PKC isoforms are Ca<sup>2+</sup>-independent, and atypical PKC isoforms are Ca<sup>2+</sup>- and DAG-independent. All PKC isoforms are dependent on membrane phospholipid, although some may become activated in a soluble form. PKC family members are broad-specificity kinases, which phosphorylate many proteins, the best known of which are MARCKS and nuclear lamins. The role of phorbol ester activators of PKC as tumor promoters clearly implicates PKC as activators of signaling pathways, including the MAP kinase pathway.

| Description                 | Clone | Species | Applications Format | Host | Quantity | Cat. No. |
|-----------------------------|-------|---------|---------------------|------|----------|----------|
| Antibody                    |       |         |                     |      |          |          |
| PKC, active                 |       |         | KA                  |      | 1 µg     | 14-115   |
| Assays                      |       |         |                     |      |          |          |
| PKC Assay Kit               |       |         | KA                  |      | 1 kit    | 17-139   |
| PKC/CaMK Inhibitor Cocktail |       |         | KA                  |      | 1 mL     | 20-119   |
| PKC Substrate Cocktail      |       |         | KA                  |      | 1 mL     | 20-131   |
| PKC Lipid Activator         |       |         | KA                  |      | 1 mL     | 20-133   |

#### $PKC\alpha$

#### Antibody

| Anti-PKC $lpha$                        | M4    | $B \; H \; M \; R \; Rb$ | IP WB NEUT | Pur  | M IgG <sub>1</sub> | 100 µg | 05-154   |
|--|-------|--------------------------|------------|------|--------------------|--------|----------|
| Anti-Protein Kinase $Clpha$            | 1F3.2 | Н                        | ELISA WB   | Asc  | M IgM              | 100 μL | MAB3074  |
| Anti-phospho-PKC $\alpha$ (Ser657)     |       | $B \; H \; M \; R \; Rb$ | WB         | Pur  | Rb IgG             | 200 µg | 06-822   |
| Anti-PKC $\alpha$ , $\beta$ , $\gamma$ |       | $B \; H \; M \; R \; Rb$ | WB         | Pur  | Rb IgG             | 200 µg | 06-870   |
| Anti-phospho-PKC $\alpha$ (Thr638)     |       | НМ                       | WB         | APur | Rb IgG             | 100 μL | 07-871   |
| PKC $lpha$ siRNA/siAb Assay Kit        |       | Н                        | WB RNAi    |      |                    | 1 kit  | 60-110   |
| siRNA plasmid, pKD-PKCα-v4             |       | Н                        | RNAi       |      |                    | 5 µg   | 62-104   |
| siRNA plasmid, pKD-PKCα-v6             |       | ΗМ                       | RNAi       |      |                    | 5 µg   | 62-105   |
| PKC $lpha$ SMARTpool siRNA reagent     |       |                          | RNAi       |      |                    | 5 nmol | M-003523 |
| PKC $\alpha$ , active                  |       |                          | KA         |      |                    | 10 µg  | 14-484   |

#### РКСβ

#### Antibody

| Anti-PKC $\alpha$ , $\beta$ , $\gamma$ | BHMRRb WB | Pur Rb Ig(                   | G 200 μg 06-87 | 0          |
|--|-----------|------------------------------|----------------|------------|
| Anti-PKC ζ                             | H M R WB  | antiserum Rb lg <sup>0</sup> | G 200 μL 07-26 | 04         |
| Anti-phospho-PKCβ I (Thr642)           | H WB      | APur Rb Ig <sup>Q</sup>      | G 100 μL 07-87 | 72         |
| Anti-phospho-PKCβ II (Thr641)          | H WB      | APur Rb Igo                  | G 100 μL 07-87 | '3         |
| Anti-phospho-PKCβ I&II (Thr500)        | H WB      | APur Rb Ig <sup>Q</sup>      | G 100 μL 07-87 | <b>'</b> O |
| PKCβI, active                          | KA        |                              | 10 µg 14-50    | )3         |
| PKCβII, active                         | KA        |                              | 10 µg 14-49    | 06         |

#### РКСδ

#### Antibody

| Anti-PKCδ                           | H M R | WB       | Pur  | Rb IgG | 200 µg | 06-990 |
|-------------------------------------|-------|----------|------|--------|--------|--------|
| Anti-phospho-PKCδ (Ser645)          | Н     | WB       | APur | Rb IgG | 100 µL | 07-874 |
| Anti-phospho-PKC& (Ser664)          | Н     | WB       | APur | Rb IgG | 100 µL | 07-875 |
| Anti-PKC δ                          | ΗМ    | IP WB IH | Pur  | Sh     | 100 µg | AB1685 |
| siRNA plasmid, pKD-PKC <b>δ</b> -v3 | Н     | RNAi     |      |        | 5 µg   | 62-071 |

|                        | Description   | Clone                | Species           | Applications    | Format        | Host        | Quantity         | Cat. No.           |
|------------------------|---|----------------------|-------------------|-----------------|---------------|-------------|------------------|--------------------|
|                        | siRNA plasmid, pKD-PKC <b>&amp;</b> -v6   |                      | H M R             | RNAi            |               |             | 5 µg             | 62-072             |
|                        | PKC $\delta$ , active   |                      |                   | KA              |               |             | 10 µg            | 14-504             |
| ίζγ                    |   |                      |                   |                 |               |             |                  |                    |
| ntibo                  | odv   |                      |                   |                 |               |             |                  |                    |
|                        | Anti-PKCα, β, γ   |                      | B H M R Rb        | WB              | Pur           | Rb IgG      | 200 µg           | 06-870             |
|                        | Anti-phospho-PKCγ (Thr514)  |                      | Н                 | WB              | APur          | Rb IgG      | 100 µL           | 07-878             |
|                        | Anti-phospho-PKCγ (Thr655)  |                      | Н                 | WB              | APur          | Rb IgG      | 100 µL           | 07-879             |
|                        | Anti-phospho-PKCγ (Thr674)  |                      | Н                 | WB              | APur          | Rb IgG      | 100 µL           | 07-880             |
|                        | PKCγ, active  |                      |                   | KA              |               | , ,         | 10 µg            | 14-483             |
| ΚCε                    |   |                      |                   |                 |               |             |                  |                    |
| <b>ιιτ</b> ε<br>Intibo | -4.   |                      |                   |                 |               |             |                  |                    |
| MIIDO                  | pay<br>Anti-phospho-PKCε (Ser729)   |                      | H M R Rb          | WB              | Pur           | Rb IgG      | 200 µg           | 06-821             |
|                        | Anti-PKC <sub>E</sub>   |                      | H M R             | IP WB IH        | Pur           | Rb IgG      | 200 µg           | 06-991             |
|                        | PKCE siRNA/siAb Assay Kit   |                      | Н                 | WB RNAi         | FUI           | kb igG      | 200 pg<br>1 kit  | 60-111             |
|                        | siRNA plasmid, pKD-PKCE-v1  |                      | Н                 | RNAi            |               |             | 5 µg             | 62-121             |
|                        | siRNA plasmid, pKD-PKCε-v5  |                      | Н                 | RNAi            |               |             | 5 µg             | 62-121             |
|                        | PKCe SMARTpool siRNA reagent  |                      | 11                | RNAi            |               |             | 5 pg<br>5 nmol   | M-004653           |
|                        | PKCε, active  |                      |                   | KA              |               |             | 10 µg            | 14-518             |
|                        | Treet, delive   |                      |                   | 10.1            |               |             | 10 49            | 14010              |
| KCη                    |   |                      |                   |                 |               |             |                  |                    |
| Antibo                 | ody   |                      |                   |                 |               |             |                  |                    |
|                        | PKC η, active   |                      |                   | KA              |               |             | 10 µg            | 14-497             |
| KCı                    |   |                      |                   |                 |               |             |                  |                    |
| Antibo                 | ody   |                      |                   |                 |               |             |                  |                    |
|                        | Anti-phospho-PKCι (Thr555)/PKCλ (Thr563)  |                      | НМ                | WB              | APur          | Rb IgG      | 100 µL           | 07-881             |
|                        | PKC ι, active   |                      |                   | KA              |               | Ü           | 10 µg            | 14-505             |
| u.c                    |   |                      |                   |                 |               |             |                  |                    |
| КСμ                    |   |                      |                   |                 |               |             |                  |                    |
| Antibo                 |   |                      | 11.44             | ) A /D          | 4.0           | D O         | 100 1            | 07.000             |
|                        | Anti-phospho-PKCµ/PKD (Ser742)  |                      | НМ                | WB              | APur          | Rb IgG      | 100 µL           | 07-882             |
|                        | PKCµ, active  |                      |                   | KA              |               |             | 10 µg            | 14-508             |
| КСτ                    |   |                      |                   |                 |               |             |                  |                    |
| ntibo                  | ody   |                      |                   |                 |               |             |                  |                    |
|                        | Anti-phospho-PKCτ (Ser676)  |                      | Н                 | WB              | APur          | Rb IgG      | 100 µL           | 07-883             |
|                        | Anti-phospho-PKCτ (Ser695)  |                      | Н                 | WB              | APur          | Rb IgG      | 100 µL           | 07-884             |
|                        | Anti-phospho-PKCτ (Thr538)  |                      | Н                 | WB              | APur          | Rb IgG      | 100 µL           | 07-885             |
|                        | PKC $\tau$ , active   |                      |                   | KA              |               |             | 10 µg            | 14-444             |
| ΚCζ                    |   |                      |                   |                 |               |             |                  |                    |
| ncs<br>Antibo          | ody   |                      |                   |                 |               |             |                  |                    |
| VIIII DC               | pay<br>PKCζ siRNA/siAb Assay Kit  |                      | Н                 | WB RNAi         |               |             | 1 kit            | 60-112             |
|                        | siRNA plasmid, pKD-PKC5-v1  |                      | Н                 | RNAi            |               |             | 5 µg             | 62-073             |
|                        | siRNA plasmid, pKD-PKC5-v3  |                      | Н                 | RNAi            |               |             | 5 µg             | 62-073             |
|                        | PKC\$ SMARTpool siRNA reagent   |                      |                   | RNAi            |               |             | 5 pg<br>5 nmol   | M-003526           |
|                        | PKC ζ, active   |                      |                   | KA              |               |             | 10 µg            | 14-525             |
|                        | 5, 40110  |                      |                   |                 |               |             | 10 kg            | 14020              |
|                        | protein Ser/Thr phosphatases (PP1) have a high<br>we to Inhibitor-2 and microcystin LR. | specific activity ag | jainst phosphoryl | ase, and are re | esistant to o | kadaic acid | up to concentrat | ions of 50 nM, but |
|                        | Description   | Clone                | Species           | Applications    | Format        | Host        | Quantity         | Cat. No.           |
| olycl                  | onal Antibody   |                      | 2424102           |                 |               |             | ~~~~~~ <u>~</u>  |                    |
| / 51                   | Anti-PP1 a  |                      | B H M R Rb        | \/\/R           | Pur           | Rh IaG      | 200 ua           | 06-221             |

|      | Description                       | Clone | Species    | Applications | Format | Host   | Quantity | Cat. No. |
|------|-----------------------------------|-------|------------|--------------|--------|--------|----------|----------|
| Poly | clonal Antibody                   |       |            |              |        |        |          |          |
|      | Anti-PP l $lpha$                  |       | B H M R Rb | WB           | Pur    | Rb IgG | 200 µg   | 06-221   |
|      | Anti-PP1 $\delta$                 |       | BHMR       | IP WB        | APur   | Rb IgG | 200 μL   | 07-270   |
|      | Anti-PP l $lpha$                  |       | ΗМ         | WB           | APur   | Rb IgG | 200 μL   | 07-273   |
|      | Anti-Protein Phosphatase 1 $lpha$ |       | BHMR       | IP WB        | APur   | Rb IgG | 100 µg   | AB4082   |
|      | Anti-Protein Phosphatase 1 β      |       | BHMR       | IP WB        | Pur    | Rb     | 100 µg   | AB4083   |
|      | Anti-Protein Phosphatase 1 γ1     |       | BHMR       | IP WB        | Pur    | Rb     | 100 µg   | AB4084   |
|      | Anti-Protein Phosphatase 1 γ2     |       | BHMR       | IP WB        | Pur    | Sh     | 100 µg   | AB4085   |
|      | Anti-Protein Phosphatase 1 γ1     |       | BHMR       | IP VVB       | Pur    | Rb     | 100 µg   | AB4084   |

| Description |        | Clone | Species | Applications | Format | Host | Quantity | Cat. No. |
|-------------|--------|-------|---------|--------------|--------|------|----------|----------|
| Assay       |        |       |         |              |        |      |          |          |
| PP1/PP2A T  | oolbox |       |         | PA           |        |      | 1 kit    | 17-301   |
| Protein     |        |       |         |              |        |      |          |          |
| PP1α        |        |       |         | KA           | Pur    |      | 10 µg    | 14-595   |
| PP1, Pured  |        |       |         | PA           |        |      | 10 units | 14-110   |
| Inhibitor 2 |        |       |         | PA           |        |      | 100 µg   | 14-162   |

#### Prenylation

see FNTA, FNTB, ICMT, and Rce-1

#### PTEN

The PTEN (MWAC1) tumor suppressor is a lipid phosphatase that removes the 3' phosphate from PI(3,4,5)P3. Overexpression of PTEN inhibits cell migration, and expression of antisense PTEN mRNA enhances migration. Integrin-mediated cell spreading is down-regulated by wild-type PTEN.

|   | 0     | '       | 0                        | ,      | / 1                |          |                 |
|---|-------|---------|--------------------------|--------|--------------------|----------|-----------------|
| Description   | Clone | Species | <b>Applications</b>      | Format | Host               | Quantity | Cat. No.        |
| Monoclonal Antibody                                 |       |         |                          |        |                    |          |                 |
| Anti-PTEN   | 6H2.1 | H M R   | WB IP IHC IC             | 0      | M lgG              | 100 µg   | 04-035          |
| Anti-PTEN, Rb Mab                                   |       | H M R   | WB IH(P)<br>IF FC IP     |        | Rb IgG             | 100 µL   | 04-409          |
| Polyclonal Antibody                                 |       |         |                          |        |                    |          |                 |
| Anti-phospho-PTEN (Ser370)                          |       | Μ       | WB                       | APur   | Rb IgG             | 100 µL   | 07-889          |
| Anti-phospho-PTEN (Ser385)                          |       | Н       | WB                       | APur   | Rb IgG             | 100 µL   | 07-890          |
| Anti-phospho-PTEN (Ser380/Thr382/<br>Thr383/Ser385) |       | ΗМ      | WB                       | APur   | Rb IgG             | 100 µL   | 07-891          |
| Assay   |       |         |                          |        |                    |          |                 |
| PTEN Malachite Green Assay Kit                      |       |         | PA                       |        |                    | 1 kit    | 1 <i>7-</i> 351 |
| PTEN Enzyme Assay Buffer, 5X<br>iRNA                |       |         | PA                       |        |                    | 1 mL     | 20-165          |
| PTEN SMARTpool siRNA reagent                        |       | Н       | RNAi                     |        |                    | 5 nmol   | M-003023        |
| Anti-PTEN, C-terminus                               | A2b1  | H M R   | IP WB IC<br>IH not IH(P) | Pur    | M lgG <sub>1</sub> | 100 µg   | MAB4037         |
| Protein   |       |         |                          |        |                    |          |                 |
| PTEN, active  |       |         | PA                       |        |                    | 10 µg    | 14-488          |
| PTP   |       |         |                          |        |                    |          |                 |
| Polyclonal Antibody                                 |       |         |                          |        |                    |          |                 |
| Anti-PTP-PEST                                       | AG25  | ΗМ      | WB IP IC                 | Asc    | M IgG <sub>1</sub> | 100 pL   | MAB3739         |
| PTPMEG-2  |       |         | KA                       | Pur    |                    | 10 µg    | 14-592          |
| LMPTP-B   |       |         | PA                       |        |                    | 10 µg    | 14-620          |
| LMPTP-A   |       |         | PA                       |        |                    | 10 µg    | 14-619          |
| PTPMEG-1  |       |         | PA                       |        |                    | 10 µg    | 14-642          |
| HePTP   |       |         | KA                       | Pur    |                    | 10 µg    | 14-593          |
| TCPTP   |       |         | PA                       |        |                    | 10 µg    | 14-646          |
| TCPTP   |       |         | PA                       |        |                    | 10 µg    | 14-646          |

#### PTP-1D

see SHP-2 below

#### Rac

Rac is a member of the Rho family of small GTP-binding proteins, and plays roles in cytoskeletal actin organization as well as transformation. Rac is activated downstream of cdc42, or independently by Ras signaling, and in turn Rac stimulates Rho. Activation of Rac results in the formation of lamellipodia, whereas cdc42 stimulates filopodia formation, and Rho stimulates formation of actin stress fibers. Effectors of the Rac pathway include NADPH Oxidase and the SAPK1/JNK family of kinases. The latter pathway is activated by the PAK kinases, which are direct targets of Rac.

| Description                    | Clone            | Species | Applications | Format | Host                | Quantity | Cat. No. |
|--------------------------------|------------------|---------|--------------|--------|---------------------|----------|----------|
| Monoclonal Antibody            |                  |         |              |        |                     |          |          |
| Anti-Rac 1                     | 23A8             | HMR     | IP WB IH     | Pur    | M IgG <sub>2b</sub> | 250 µg   | 05-389   |
| Polyclonal Antibody            |                  |         |              |        |                     |          |          |
| Anti-Rac2                      |                  | НМ      | IP WB APA    | Serum  | Rb IgG              | 50 μL    | 07-604   |
| Anti-phospho-Rac 1/cdc42 (Ser. | 71)              | Н       | WB           | APur   | Rb IgG              | 100 μL   | 07-896   |
| Anti-Rac1/Cdc42, phospho-spe   | ecific (Ser71)   | Н       | WB           | APur   | Rb                  | 100 μL   | AB3838   |
| Anti-Rac 1                     |                  | HMR     | WB           | Pur    | Rb                  | 100 µg   | AB4202   |
| Recombinnat Protein            |                  |         |              |        |                     |          |          |
| Rac/cdc42 Assay Reagent (PAI   | <1 PBD, agarose) |         | ABA          |        |                     | 300 µg   | 14-325   |

|       | Description                               | Clone | Species  | <b>Applications</b> | Format | Host                | Quantity  | Cat. No.        |
|-------|---|-------|----------|---------------------|--------|---------------------|-----------|-----------------|
| Assay |   |       |          |                     |        |                     |           |                 |
|       | Rac1 Activation Assay Kit                 |       |          | ABA                 |        |                     | 1 kit     | 1 <i>7</i> -283 |
|       | Rac2 Activation Assay Kit                 |       |          | ABA                 |        |                     | 30 assays | 1 <i>7</i> -369 |
|       | Rac1 Activation Assay (96 well)           |       |          | ABA                 |        |                     | 96 assays | 17-450          |
|       | Rac1/Cdc42 Activation Assay (96 well)     |       |          | ABA                 |        |                     | 96 assays | 1 <i>7-</i> 452 |
| cDNA  |   |       |          |                     |        |                     |           |                 |
|       | Rac 1 cDNA (activated) in pUSEamp         |       |          | TFX                 |        |                     | 5 μg      | 21-193          |
|       | Rac 1 cDNA (dominant negative) in pUSEamp |       |          | TFX                 |        |                     | 5 μg      | 21-199          |
|       | Rac 1 cDNA Allelic Pack                   |       |          | TFX                 |        |                     | 1 kit     | 17-309          |
|       | Rac 1 cDNA (wt) in pUSEamp                |       |          | TFX                 |        |                     | 5 μg      | 21-200          |
| siRNA | A   |       |          |                     |        |                     |           |                 |
|       | Rac1 siRNA/siAb Assay Kit                 |       | Н        | WB RNAi             |        |                     | 1 kit     | 60-037          |
|       | Rac1 SMARTpool siRNA reagent              |       | Н        | RNAi                |        |                     | 5 nmol    | M-003560        |
|       | Anti-Rac 1                                | 102   | H M R    | IF WB IH            | Pur    | M IgG <sub>2b</sub> | 100 µg    | MAB3735         |
|       |   |       | Ca Ch Dr |                     |        |                     |           |                 |

#### Raf

The Raf proteins (Raf-1, A-Raf, B-Raf) are Ser/Thr kinases with homology to the PKC family, containing an n-terminal regulatory domain and a c-terminal catalytic domain. Members of the Raf family bind to activated Ras GTPase, which results in Raf translocation to the plasma membrane and activation. Activated Raf proteins phosphorylate MEKs, and are therefore the principle transducers of signals from Ras to MAP kinase. In addition to its role in mitogenesis, Raf-1 may play a role in regulation of apoptosis and cell cycle progression. Activation of Raf-1 involves phosphorylation of Ser338/339 and Tyr340/341. Activating mutations of B-Raf that disrupt its auto-inhibition loop have been implicated in a number of cancers, including melanoma and colon cancer.

| Description       | on   | Clone | Species | Applications   | Format | Host               | Quantity       | Cat. No.         |
|-------------------|--|-------|---------|----------------|--------|--------------------|----------------|------------------|
| Monoclonal Antil  | body                                       |       |         |                |        |                    |                |                  |
| Anti-Raf-1        |  | AM223 |         | IP WB          | Pur    | Rb IgG             | 100 µg         | 05-739           |
| Anti-Raf-1,       | , Rb Mab                                   |       |         | WB IP IF IH(P) |        | Rb IgG             | 100 μL         | 04-412           |
| Anti-phosp        | pho-Raf-1 (Ser259), Rb Mab                 |       |         | WB IF IH(P)    |        | Rb IgG             | 100 μL         | 04-411           |
| Anti-phosp        | pho-Raf-1 (Ser338)                         |       | Μ       | WB             | Pur    | R IgG <sub>1</sub> | 100 µg         | 05-534           |
| Anti-phosp        | pho-Raf-1 (Ser338)                         |       | Μ       | WB             | Pur    | M lgG              | 200 µg         | 05-538           |
| Anti-phosp        | pho-Raf-1 (Ser338)                         |       |         | WB             | Pur    | M lgG              | 1 mg           | 05-538MG         |
| Anti-C-raf        | (phospho-Ser 621)                          |       | НМ      | WB             | Pur    |                    |                | 04-306           |
| Polyclonal Antibo | ody  |       |         |                |        |                    |                |                  |
| Anti-Raf-1        |  |       | Н       | IPK WB         | Pur    | Rb IgG             | 200 µg         | 07-396           |
| Anti-B-Raf        |  |       | H M R   | IPK WB         | Pur    | Rb IgG             | 200 µg         | 07-453           |
| Anti-B-Raf,       | NT   |       | H M R   | KA WB          | APur   | Rb IgG             | 200 μL         | 07-583           |
| Anti-phosp        | pho-Raf1 (Ser259)                          |       | Н       | WB             | APur   | Rb IgG             | 100 µL         | 07-811           |
| Anti-phosp        | pho-Raf1 (Ser43)                           |       | Н       | WB             | APur   | Rb IgG             | 100 μL         | 07-812           |
| Anti-phosp        | pho-Raf-1 (Ser621)                         |       | Н       | WB             | APur   | Rb lgG             | 100 μL         | 07-813           |
| Anti-phosp        | pho-Raf-1 (Ser338/Tyr340)                  |       | Н       | WB             | APur   | Rb IgG             | 100 μL         | 07-814           |
| Anti-phosp        | oho-Raf-1 (Tyr340/Tyr341)                  |       | Н       | WB             | APur   | Rb lgG             | 100 μL         | 07-815           |
| Assay             |  |       |         |                |        |                    |                |                  |
| Raf-1 Kina        | ase Cascade Assay Kit                      |       |         | IPK KA         |        |                    | 1 kit          | 17-357           |
| B-Raf Kina        | ase Cascade Assay Kit                      |       |         | IPK KA         |        |                    | 1 kit          | 17-358           |
| B-Raf Kina        | ase Assay Kit, Chemiluminescence Detectio  | n     |         | KA WB          |        |                    | 1 kit          | 17-359           |
| Raf-1 Kina        | ase Assay Kit, Chemiluminescence Detection | n     |         | KA             |        |                    | 1 kit          | 17-360           |
| Protein           |  |       |         |                |        |                    |                |                  |
| Raf-1-RBD         | GST Protein, ras binding domain            |       |         |                | Pur    |                    | 300 µg         | SGT223           |
| cDNA              |  |       |         |                |        |                    |                |                  |
| Raf-1 cDN         | NA (wt) in pUSEamp                         |       |         | TFX            |        |                    | 5 µg           | 21-111           |
| siRNA             |  |       |         |                |        |                    |                |                  |
| Raf-1 siRN        | NA/siAb Assay Kit                          |       | Н       | WB RNAi        |        |                    | 1 kit          | 60-109           |
| siRNA pla         | asmid, pKD-Raf-1-v4                        |       | H M R   | RNAi           |        |                    | 5 µg           | 62-085           |
| siRNA pla         | asmid, pKD-Raf-1-v6                        |       | Н       | RNAi           |        |                    | 5 µg           | 62-086           |
| Raf-1 SMA         | ARTpool siRNA reagent                      |       |         | RNAi           |        |                    | 5 nmol         | M-003601         |
|                   |  |       |         |                |        |                    |                |                  |
| Protein           |  |       |         |                |        |                    |                |                  |
|                   | 99E), active                               |       |         | KA             |        |                    | 10 µg          | 14-557           |
| B-Raf (V59        | 99E), active<br>415), active               |       |         | KA<br>KA       |        |                    | 10 µg<br>10 µg | 14-557<br>14-530 |

#### Ral

RalA and RalB are Rasfamily small G-Proteins that can be activated through a Ras-dependent mechanism. Mutants of Ras that are defective in activating all known Ras effectors except Ral GEF are still able to induce invasive phenotypes in transfected cells, suggesting a role for Ral in these processes.


| •  |       |         |                     |               |        |          |          |
|--|-------|---------|---------------------|---------------|--------|----------|----------|
| Description                              | Clone | Species | <b>Applications</b> | Format        | Host   | Quantity | Cat. No. |
| Antibody                                 |       |         |                     |               |        |          |          |
| Anti-RalB                                | 25    | H M R   | WB                  | Pur           | M lgG  | 100 µL   | 309      |
| Protein                                  |       |         |                     |               |        |          |          |
| Ral Assay Reagent (Ral BP1, agarose)     |       |         | ABA                 |               |        | 300 µg   | 14-415   |
| Assays                                   |       |         |                     |               |        |          |          |
| RalA Activation Assay Kit                |       |         | ABA                 |               |        | 1 kit    | 17-300   |
| RalB Activation Assay Kit                |       |         | ABA                 |               |        | 1 kit    | 17-439   |
| Ral Activation Assay Buffer, 5X          |       |         | ABA                 |               |        | 18 mL    | 20-196   |
| cDNA                                     |       |         |                     |               |        |          |          |
| RalA cDNA (activated) in pUSEamp         |       |         | TFX                 |               |        | 5 µg     | 21-189   |
| RalA cDNA (dominant negative) in pUSEam  | р     |         | TFX                 |               |        | 5 µg     | 21-187   |
| RalA cDNA (wt) in pUSEamp                |       |         | TFX                 |               |        | 5 µg     | 21-190   |
| Ran                                      |       |         |                     |               |        |          |          |
| Antibody                                 |       |         |                     |               |        |          |          |
| Ran                                      |       | H Xn    | WB IP IC            | antiserum     | Rb IgG | 200 µL   | 07-517   |
| RanGEF (RCC1)                            |       | 11741   | WB IP IC            | antiserum     | Rb IgG | 200 pt   | 07-519   |
|  |       |         | ,,,,,,,,            | GIIIIOOI OIII | 90     | 200 pt   | 0, 0.,   |
| RanBP17                                  |       |         |                     |               |        |          |          |
| Antibody                                 |       |         |                     |               |        |          |          |
| RanBP17                                  |       | Н       | WB ELISA            | APur          | Rb     | 100 µg   | AB3488   |
| Rap                                      |       |         |                     |               |        |          |          |
| Antibody                                 |       |         |                     |               |        |          |          |
| Rap1 Activation Assay Kit                |       |         | ABA                 |               |        | 1 kit    | 17-321   |
| Rap 1 Assay Reagent (Ral GDS-RBD, agaros | e)    |         | ABA                 |               |        | 650 µg   | 14-455   |
|  | -1    |         |                     |               |        | 222 FB   | 00       |
| Rap1GAP                                  |       |         |                     |               |        |          |          |
| Antibody                                 |       |         |                     |               |        |          |          |
| Rap1GAP (N-term)                         |       |         |                     |               |        |          | 04-413   |
|  |       |         |                     |               |        |          |          |

#### Ras

Ras proteins are small GTP-binding proteins which contain all GTPase and effector functions within a single polypeptide. Three isoforms of Ras exist, Ki-Ras, Ha-Ras, and N-Ras, with distinct expression patterns but similar signaling activity. Ras is palmitoylated and farnesylated at the c-terminus, anchoring it in the membrane. In resting cells, Ras is loaded with GDP; activation by growth factor stimulation of receptors (for example) recruits Guanine nucleotide Exchange Factors (GEFs) which induce exchange of GDP with GTP. In its GTP-bound form, Ras binds and activates effector proteins, such as Raf, RalGDS, or Pl3 Kinase. GTP hydrolysis returns Ras to its inactive state; hydrolysis is accelerated by GTPase Activating Proteins (GAPs), such as RasGAP or NF-1.

|       | Description                               | Clone | Species    | Applications | Format | Host                    | Quantity | Cat. No. |
|-------|---|-------|------------|--------------|--------|-------------------------|----------|----------|
| Mon   | oclonal Antibody                          |       |            |              |        |                         |          |          |
|       | Anti-Ras                                  | RAS10 | H M R      | ELISA FC IP  | Pur    | M lgG $_{2a\mathbf{k}}$ | 100 µg   | 05-516   |
|       |   |       |            | WB IC IH     |        |                         |          |          |
|       | Anti-Ha-Ras                               | MC57  | Н          | WB           | Pur    | Rb IgG                  | 100 pL   | 05-775   |
|       | Anti-RasGAP                               | B4F8  | B H Ht M R | IP WB IC     | Pur    | M lgG                   | 200 µg   | 05-178   |
| Poly  | clonal Antibody                           |       |            |              |        |                         |          |          |
|       | Anti-RasGAP                               |       | Av H M R   | IP WB        | Pur    | Rb IgG                  | 200 µg   | 06-157   |
| Prote | ein                                       |       |            |              |        |                         |          |          |
|       | Ras-GST, agarose conj.                    |       |            | ABA          |        |                         | 100 µg   | 14-139   |
|       | Ras Assay Reagent (Raf-1 RBD, agarose)    |       |            | ABA          |        |                         | 600 µg   | 14-278   |
|       | Ras, recombinant human full length        |       |            |              | Pur    |                         | 100 µg   | SGT213   |
| cDN   | A   |       |            |              |        |                         |          |          |
|       | H-Ras cDNA (wt) in pUSEamp                |       |            | TFX          |        |                         | 5 µg     | 21-102   |
|       | H-Ras cDNA (activated) in pUSEamp         |       |            | TFX          |        |                         | 5 µg     | 21-103   |
|       | H-Ras cDNA (dominant negative) in pUSEamp |       |            | TFX          |        |                         | 5 µg     | 21-104   |
|       |   |       |            |              |        |                         |          |          |

| Description                             | Clone | Species | Applications Format | Host | Quantity  | Cat. No. |  |
|---|-------|---------|---------------------|------|-----------|----------|--|
| Assay                                   |       |         |                     |      |           |          |  |
| Ras GTPase Activation ELISA Kit-96 well |       | H M R   | ELISA GPA           |      | 96 assays | 17-424   |  |



#### **Ras Activation Titration**

Increasing amounts of EGF stimulated HeLa whole cell extracts were titrated in a 2 fold serial dilutions down to 1.25 mg of total cellular protein using the Ras GTPase Activation EUSA Kit.

| Ras Activation Assay Kit    |       |       | ABA      |     |                     | 1 kit  | 17-218   |
|-----------------------------|-------|-------|----------|-----|---------------------|--------|----------|
| H-Ras cDNA Allelic Pack     |       |       | TFX      |     |                     | 1 kit  | 17-267   |
| siRNA                       |       |       |          |     |                     |        |          |
| Ras siRNA/siAb Assay Kit    |       | Н     | WB RNAi  |     |                     | 1 kit  | 60-093   |
| siRNA plasmid, pKD-Ras-v1   |       |       | RNAi     |     |                     | 5 μg   | 62-214   |
| siRNA plasmid, pKD-Ras-v2   |       |       | RNAi     |     |                     | 5 μg   | 62-215   |
| Ras SMARTpool siRNA reagent |       | Н     | RNAi     |     |                     | 5 nmol | M-004142 |
| Anti-H-Ras                  | 7D7.2 | H M R | ELISA WB | Pur | M IgG <sub>2b</sub> | 100 µg | MAB3291  |

#### Prenylation

see, FNTA, FNTB, ICMT, and Rce-1

#### Rce1 (Ras converting enzyme 1)

Rce1 is a metalloproleinase that cleaves prenylated members of the Ras/Rho family of small G-proteins and nuclear Lamin. Cleavage of prenylated proteins is part of the required post-translational modification (prenylation by FT, cleavage by Rce1, methylation by lcmt) required for these proteins to become active.

|    | Description | Clone | Species | Applications | Format | Host   | Quantity | Cat. No. |
|----|-------------|-------|---------|--------------|--------|--------|----------|----------|
| An | tibody      |       |         |              |        |        |          |          |
|    | Anti-Rce-1  |       | Н       | WB           | APur   | Rb IgG | 100 mL   | 09-120   |

#### Rheb

Rheb is a member of the Ras family of GTPases, primarily expressed in the brain. It is a target of the TSC2 (Tuberous Sclerosis Complex 2) GAP (GTPase Activating Protein), and is thus integral in mTOR signaling. Two Rheb proteins have been identified, Rheb2 (aka Rheb, hRheb1 or Rheb2) and Rheb11 (aka hRheb2).

|      | Description | Clone | Species | <b>Applications</b> | Format | Host   | Quantity | Cat. No. |
|------|-------------|-------|---------|---------------------|--------|--------|----------|----------|
| Anti | body        |       |         |                     |        |        |          |          |
|      | Anti-RHEB 1 |       | Н       | WB                  | Pur    | Rb IgG | 100 µL   | 09-247   |

#### Rho

The Rho-family of small GTP-binding proteins includes Rac and cdc42, as well as several Rho isoforms. Collectively, this family is involved in the regulation of cytoskeletal changes that accompany signal transduction, and is also involved in the activation of the stress-activated MAP kinase cascades.

| Description                               | Clone | Species | Applications | Format  | Host               | Quantity | Cat. No. |
|---|-------|---------|--------------|---------|--------------------|----------|----------|
| Monoclonal Antibody                       |       |         |              |         |                    |          |          |
| Rho (-A, -B, -C)                          | 3L74  | H R     | WB           | Pur     | Rb IgG             | 100 µL   | 05-822   |
| Rho (-A, -B, -C)                          | 55    | HMR     | WB IC        | Pur     | M lgG <sub>1</sub> | 200 µg   | 05-778   |
| RhoE/Rnd3                                 | 4     | ΗМ      | WB IC        | ascites | M lgG              | 200 μL   | 05-723   |
| RhoG                                      |       | HM      | WB           |         | M lgG              | 100 µg   | 04-486   |
| RhoGAP p190                               | D2D6  | HMRMk   | WB IP        | Pur     | M lgG              | 200 µg   | 05-378   |
| Polyclonal Antibody                       |       |         |              |         |                    |          |          |
| Rho                                       |       | Н       | WB           | Pur     | Rb                 | 100 µg   | AB3884   |
| RhoGDI                                    |       | НМВ     | WB IP        | Pur     | Rb IgG             | 200 µg   | 06-730   |
| Protein                                   |       |         |              |         |                    |          |          |
| Rho Assay Reagent (Rhotekin RBD, agarose) |       |         | ABA          |         |                    | 650 µg   | 14-383   |
| Rhotekin, GST fusion protein              |       |         |              |         |                    | 500 µg   | 14-662   |
| Rho-GST                                   |       |         |              | Pur     |                    | 20 µg    | SGT212   |

|       | Description                              | Clone | Species | Applications | Format | Host | Quantity | Cat. No. |
|-------|--|-------|---------|--------------|--------|------|----------|----------|
| Assay | /  |       |         |              |        |      |          |          |
|       | Rho Activation Assay Kit                 |       |         | ABA          |        |      | 1 kit    | 17-294   |
| cDN/  | A  |       |         |              |        |      |          |          |
|       | RhoA cDNA (activated) in pUSEamp         |       |         | TFX          |        |      | 5 µg     | 21-195   |
|       | RhoA cDNA (dominant negative) in pUSEamp |       |         | TFX          |        |      | 5 µg     | 21-196   |
|       | RhoA cDNA (wt) in pUSEamp                |       |         | TFX          |        |      | 5 µg     | 21-194   |
| DI O  | A.D.                                     |       |         |              |        |      |          |          |

#### RhoGAP

#### Antibody

| Anti-RhoGAP p190 | D2D6 | H M Mk R | IP WB Pui | M IgG               | 200 µg | 05-378  |
|------------------|------|----------|-----------|---------------------|--------|---------|
| Anti-Rin         | 14G7 | Н        | WB Pui    | M IgG <sub>1</sub>  | 100 µg | MAB3744 |
| Anti-Rit         | 27G2 | НМ       | WB Pui    | M IgG <sub>2b</sub> | 100 µg | MAB3743 |

#### **RKIP**

RKIP, or Raf Kinase Inhibitory Protein is also known as PEBP (Phosphatidylethanolamine-Binding Protein). RKIP mediates the interaction between TAK1 and MEK3. RKIP is a cytoplasmic serine protease that inhibits Raf activity, and is also thought to inhibit Thrombin.

| Description                           | Clone | Species | Applications | Format | Host   | Quantity | Cat. No. |
|---------------------------------------|-------|---------|--------------|--------|--------|----------|----------|
| Antibody                              |       |         |              |        |        |          |          |
| Anti-RKIP                             |       | H M R   | IP WB        |        | Rb IgG | 200 µg   | 07-137   |
| ROK/ROCK                              |       |         |              |        |        |          |          |
| Antibody                              |       |         |              |        |        |          |          |
| ROCK (1113) cleavage-specific product |       | Н       | VVB          | APur   | Rb IgG | 100 pL   | 07-903   |
| ROCK-1                                |       | ΗМ      | WB           | APur   | Rb     | 100 µg   | AB3885   |
| ROKa/ROCK-II                          |       | R       | WB IP        | Pur    | Rb IgG | 200 µg   | 07-443   |
| ROKa/ROCK-II                          | A9W4  | ΗR      | WB IP        | Pur    | Rb IgG | 100 µL   | 05-841   |
| siRNA plasmid, pKD-ROKa/ROCK-II-v4    |       | H M R   | RNAi         |        |        | 5 µg     | 62-075   |
| siRNA plasmid, pKD-ROKa/ROCK-II-v6    |       | HMR     | RNAi         |        |        | 5 µg     | 62-076   |
| ROKa/ROCK-II, active                  |       |         | KA           |        |        | 10 µg    | 14-451   |
| ROKb/ROCK-I, active                   |       |         | KA           |        |        | 10 µg    | 14-601   |

#### SHC

SHC is expressed as three alternatively spliced adapter proteins which share an n-terminal PTB domain, a central glycine-rich sequence, and a c-terminal SH2 domain. SHC proteins are tyrosine phosphorylated by receptor tyrosine kinases and bind to phosphotyrosine residues on those receptors. SHC also binds GRB2, and in the case of many receptor tyrosine kinases, GRB2 requires SHC as an intermediate to bind the receptor.

| Description                         | Clone          | Species  | Applications | Format    | Host              | Quantity | Cat. No. |
|-------------------------------------|----------------|----------|--------------|-----------|-------------------|----------|----------|
| Antibody                            |                |          |              |           |                   |          |          |
| phospho-SHC (Tyr317)                | 1 <i>5</i> E11 | ΗМ       | WB           |           | M $lgG_{1\kappa}$ | 100 µg   | 05-668   |
| Anti-shc (phospho-Tyr 239/pTyr 240) |                |          |              |           | M lgG             | 100 µg   | 04-312   |
| Anti-shc/p66 (N-Terminus)           |                |          |              |           | M lgG             | 100 µg   | 04-313   |
| Anti-shc/p66 (phospho-Ser 36)       |                |          |              |           | M lgG             | 100 µg   | 04-314   |
| phospho-SHC (Tyr239)                |                | Н        | WB ELISA     | APur      | Rb IgG            | 200 μL   | 07-209   |
| phospho-SHC (Tyr317)                |                | ΗМ       | WB           | APur      | Rb IgG            | 200 μL   | 07-206   |
| SHC                                 |                | H M R Po | WB IP IH     | Pur       | Rb IgG            | 250 µg   | 06-203   |
| SHC                                 |                | H M R    | WB IP        | antiserum | Rb IgG            | 200 μL   | 07-150   |
| SHC [p66; SH2 domain Protein C1]    |                |          |              | Pur       |                   | 100 μL   | AB3824   |

#### SOS

The mammalian homolog of Drosophila Son-of-Sevenless is a GDP/GTP exchange factor for Ras proteins, which has an SH3-binding motif that associates with GRB2. Together with GRB2, it associates with activated receptor tyrosine kinases in the plane of the membrane, and activates Ras. Expression of SOS 1 that is modified to constitutively reside at the membrane is transforming in most cells.

|   | Description | Clone | Species | Applications | Format | Host   | Quantity | Cat. No. |
|---|-------------|-------|---------|--------------|--------|--------|----------|----------|
| 1 | Antibody    |       |         |              |        |        |          |          |
|   | SOS1, CT    |       | НМ      | WB IP        | Pur    | Rb IgG | 200 µg   | 07-337   |

Description Clone Species Applications Format Host Quantity Cat. No

#### **STAT**

STATs (Signal-Transducing Activators of Transcription) are the principle substrates of JAK kinases, and mediators of cytokine signaling. STAT proteins have SH2 domains, and once phosphorylated by JAK kinases, STATs dimerize in a head-to-tail fashion, through their SH2 domains. STAT1 is a critical mediator of IFN signaling and STAT3 is involved in multiple signaling pathways, and mice deficient in it are not viable. In addition to tyrosine phosphorylation, many STAT proteins are targets of Ser/Thr kinases, and those phosphorylation events potentiate STAT-induced transcription. The phosphorylation state of STATs can be monitored by phosphorylation state-specific antibodies. Please see our website or Transcription brochure for a complete listing of STAT products.

| STAT1      | Description  | Clone | Species  | Applications        | Format | Host                | Quantity | Cat. No. |
|------------|--|-------|----------|---------------------|--------|---------------------|----------|----------|
| Antibo     | •  | Cione | Species  | Applications        | romiui | 11031               | Quality  | cui. No. |
| , ,,,,,,,, | Anti-STAT 1 , CT                                   |       | НМ       | EMSA IP WB          | Pur    | Rb IgG              | 250 µg   | 06-501   |
|            | Anti-STAT1   | 79    | H M R    | WB                  | Pur    | M lgG <sub>1</sub>  | 100 µg   | 05-987   |
|            | Anti-STAT1a  |       | НМ       | IP WB               | Pur    | Rb                  | 100 µg   | ab16951  |
|            | Anti-phospho-STAT1 (Tyr701)                        |       | Ch H M R | WB                  | APur   | Rb IgG              | 100 μL   | 07-307   |
|            | Anti-STAT1, phospho-specific (Tyr701)              |       | Н        | WB                  | APur   | Rb                  | 100 μL   | ab3892   |
|            | Anti-phospho-STAT1 (Ser727)                        |       | НМ       | WB                  | Serum  | Rb                  | 200 μL   | 07-714   |
|            | STATPAK, Anti-STAT (1, 2, 3, 5A, 5B) Miniature Set |       |          | EMSA IP<br>WB IC    |        |                     | 1 kit    | 17-176   |
| Assay      |  |       |          | 115.10              |        |                     |          |          |
|            | STAT 1 a EZ-TFA                                    |       | H M R    | Activity            |        |                     | 1 plate  | 70-535   |
| STAT3      |  |       |          |                     |        |                     |          |          |
| Antibo     | ody  |       |          |                     |        |                     |          |          |
|            | Anti-phospho-STAT3 (Tyr705)                        | 9E12  | НМ       | IP WB               |        | M lgG <sub>1κ</sub> | 50 µg    | 05-485   |
|            | Anti-phospho-STAT3 (Ser727)                        |       | Μ        | WB                  | APur   | Rb IgG              | 100 μL   | 07-703   |
|            | Anti-STAT3   |       | H M R    | ChIP EMSA  IP WB IC | Pur    | Rb IgG              | 200 µg   | 06-596   |
|            | Anti-STAT3, C-terminus                             |       | H M R    | EMSA IP WB          | Pur    | Rb                  | 100 µg   | AB3162   |
| Assay      |  |       |          |                     |        |                     |          |          |
|            | Anti-STAT3, phospho-specific (Ser727)              | 6E4   | НМ       | IP WB IC IH(P       | ')     | M lgG <sub>2a</sub> | 100 μL   | MAB3705  |
|            | STAT3 Transcription Factor Assay                   |       | H M R    |                     |        |                     | 1 plate  | SGT530   |

#### TAK

TAK1 (TGF- $\beta$  activated kinase) is a member of the MAPKKK family and its kinase activity is stimulated in response to TGF- $\beta$ , bone morphogenic protein (BMP) and ceramide. Activation of TAK1 is dependent on TAB1 (TAK1 binding protein 1), and overexpression of both together results in NIK-independent activation of NF- $\kappa$ B. TAK1 is an activator of SKK1/MEK4, but also phosphorylates SKK2/MEK3 and SKK3/MEK6 *in vitro*. Preliminary evidence suggests that mutations of TAK1 may play a role in lung cancer development.


|   | Description                                    | Clone | Species | Applications | Format | Host  | Quantity | Cat. No. |
|---|--|-------|---------|--------------|--------|-------|----------|----------|
| 1 | Antibody                                       |       |         |              |        |       |          |          |
|   | Anti-C-TAK1 (cdc25C associated protein kinase) |       | НМ      | WB           |        | M lgG | 200 µg   | 05-680   |

#### **VEGF and VEGF Receptor**

VEGF is a dimeric ligand, and is among the most potent angiogenic mitogens. VEGF is secreted by tumor cells and other cells exposed to hypoxia. Expression of VEGF is stimulated by FGF-2, and it activates FIt-1 and KDR, the two high-affinity receptors for VEGF. The VEGF Receptors stimulate the Ras-MAPK pathway, suggesting that it signals as a conventional receptor tyrosine kinase.

| Description                     | Clone | Species  | <b>Applications</b> | Format | Host                | Quantity | Cat. No. |
|---------------------------------|-------|----------|---------------------|--------|---------------------|----------|----------|
| Monoclonal Antibody             |       |          |                     |        |                     |          |          |
| Anti-VEGF                       | JH121 | H M R Rb | IP IH NEUT          |        | M IgG <sub>1</sub>  | 100 µg   | 05-443   |
| Anti-VEGF, recognizes 121, 165, | VG1   | Н        | ELISA WB            |        | M $\lg G_{1\kappa}$ | 100 µg   | MAB3734  |
| and 189 a.a. isoforms           |       |          | IH(P)               |        |                     |          |          |
| Anti-VEGF                       | CH-10 | HR       | ELISA WB            | Pur    | M IgG <sub>1</sub>  | 100 µg   | MAB1665  |
| Anti-VEGF Receptor-1            |       | ΗМ       | ELISA WB            | Pur    | M IgG <sub>1</sub>  | 50 μL    | MAB1664  |
| VEGFR-1 (N-term)                |       |          |                     |        |                     |          | 04-431   |

| Description  | Clone        | Species | Applications         | Format    | Host                | Quantity | Cat. No. |
|--|--------------|---------|----------------------|-----------|---------------------|----------|----------|
| Anti-KDR/Flk-1/VEGFR2  | CH-11        | НМ      | ELISA WB<br>IC IH    |           | M lgG <sub>1</sub>  | 100 µg   | 05-554   |
| Anti-phospho-KDR/Flk-1/VEGFR2 (Tyr                                   | 1054) D1W    | ΗМ      | WB                   | Pur       | Rb IgG              | 100 µL   | 05-894   |
| Anti-Flt-1   | BK302        | Н       | WB                   | Pur       | M IgG <sub>1κ</sub> | 200 µg   | 05-696   |
| Anti-VEGF Receptor-2   | 4H3B6H9      | Μ       | ELISA FC<br>IP WB    | Pur       | R IgG <sub>2b</sub> | 100 µg   | MAB1147  |
| Anti-VEGF Receptor-2   | 89B3A5       | Μ       | FC IP WB             | Pur       | R IgG <sub>2a</sub> | 100 µg   | MAB1669  |
| Anti-VEGF Receptor-3, extracellular do                               | omain 9D9F9  |         | FC IF IP<br>WB IH(P) | Asc       | M lgG <sub>1</sub>  | 100 μL   | MAB3757  |
| Polyclonal Antibody  |              |         |                      |           |                     |          |          |
| Anti-VEGF, N-terminus  |              | Н       | ELISA IP<br>RIA WB   | Pur       | Rb                  | 500 µg   | CBL42    |
| Anti-VEGF  |              | НМ      | ELISA WB<br>NEUT     | Pur       | Rb                  | 100 µg   | AB1442   |
| Anti-VEGF  |              | Μ       | ELISA WB             | Pur       | Rb                  | 100 µg   | AB1876   |
| Anti-VEGF  |              | Μ       | ELISA WB<br>NEUT     | APur      | Rb                  | 50 µg    | AB2142P  |
| Anti-phospho-VEGFR2/FLK-1/KDR (Se                                    | er 1 188)    | Н       | WB                   | APur      | Rb                  | 100 µg   | 09-021   |
| Anti-Flt-1, CT   |              | Н       | WB                   |           | Ch IgY              | 200 µg   | 06-670   |
| Anti-KDR/Flk-1/VEGFR2  |              | ΗМ      | IP WB IH             | antiserum | Rb IgG              | 200 μL   | 07-158   |
| Anti-phospho-KDR/Flk-1/VEGFR2 (Tyr                                   | 1214)        | ΗМ      | WB                   | Pur       | Rb IgG              | 200 µg   | 07-374   |
| Anti-VEGF Receptor-3   |              | Н       | WB IC IH(P)          | Pur       | Rb                  | 100 µg   | AB1875   |
| Anti-VEGF Receptor-3   |              | Μ       | ELISA WB             | APur      | Rb                  | 50 µg    | AB3127   |
| Anti-VEGF Receptor-2, phospho-specif                                 | fic (Tyr996) | Н       | IP WB                | APur      | Rb                  | 100 μL   | AB3847   |
| Anti-Flt-3, extracellular region                                     |              | ΗМ      | FC IP WB IC          |           | Rb IgG              | 200 µg   | 06-646   |
| Anti-Flt-3, cytoplasmic domain                                       |              | ΗМ      | IP WB IC             |           | Rb IgG              | 200 µg   | 06-647   |
| Protein  |              |         |                      |           |                     |          |          |
| VEGF, recombinant  |              |         | CULT                 |           |                     | 10 µg    | 01-185   |
| Vascular Endothelial Growth Factor, recombinant human, 165aa isoform | 1            |         |                      | Pur       |                     | 10 µg    | GF025    |
| Vascular Endothelial Growth Factor 1<br>recombinant mouse            | 64,          |         |                      | Pur       |                     | 10 µg    | GF060    |
| Flt-3 Ligand, recombinant human                                      |              | Н       |                      | Pur       |                     | 10 µg    | GF038    |
| VEGF Receptor-3, control peptide for                                 | AB3127       |         | PIA                  | Pur       |                     | 100 µg   | AG659    |
| siRNA  |              |         |                      |           |                     |          |          |
| VEGF SMARTpool siRNA reagent   |              | Н       | RNAi                 |           |                     | 5 nmol   | M-003550 |
| KDR/Flk-1/VEGFR2 siRNA/siAb Asso                                     | ay Kit       | Н       | WB RNAi              |           |                     | 1 kit    | 60-104   |
| KDR/Flk-1/VEGFR2 SMARTpool siRN                                      | JA reagent   |         | RNAi                 |           |                     | 5 nmol   | M-003148 |
| Flt-1 SMARTpool siRNA reagent  |              | Н       | RNAi                 |           |                     | 5 nmol   | M-003136 |
| Flt-3 SMARTpool siRNA reagent  |              |         | RNAi                 |           |                     | 5 nmol   | M-003137 |
| Flt-3 siRNA/siAb Assay Kit   |              | Н       | WB RNAi              |           |                     | 1 kit    | 60-106   |
| Assay  |              |         |                      |           |                     |          |          |
| ChemiKine™ Vascular Endothelial Gro<br>Competitive ELISA             | wth Factor,  | H R     | ELISA                |           |                     | 1 kit    | CYT132   |
| ChemiKine Vascular Endothelial Grow<br>Competitive ELISA             | vth Factor,  | Μ       | ELISA                |           |                     | 1 kit    | CYT133   |
| ChemiKine Vascular Endothelial Grow<br>Sandwich EUSA                 | vth Factor,  | Н       | ELISA                |           |                     | 1 kit    | CYT214   |



Millipore provides innovative tools, services and application expertise that improve the scientist's productivity in the pharmaceutical, academic, life science, applied science, clinical and government research sectors. Our product range spans water purification, protein research, analytical sample preparation, cell-based assays and high throughput screening. Please visit www.millipore.com/bioscience for more details.

#### Millipore Worldwide

**UNITED STATES AUSTRALIA FRANCE GERMANY SWEDEN** UK TO PLACE AN ORDER: TO PLACE AN ORDER OR TO PLACE AN ORDER: TO PLACE AN ORDER: TO PLACE AN ORDER: TO PLACE AN ORDER: TEL 040 22 42 00 TEL 0805 300 576 TEL 0800 759 8965 TEL 800 437 7500 FOR TECH SUPPORT: TEL 0800 0190 333 FAX 0800 0190 444 FAX 951 676 9209 TEL 1800 252 265 FAX 0805 300 579 FAX 0800 759 8967 03 9839 2000 FOR TECH SUPPORT: FOR TECH SUPPORT: FOR TECH SUPPORT: FOR TECH SUPPORT: FAX 03 9887 391 TEL 0805 300 578 TEL 0800 759 8966 TEL 0800 0190 555 TEL 800 437 7500

For additional contact information, visit us online at www.millipore.com.

### **MILLIPORE**

Millipore, Upstate, Chemicon, Linco, Beadlyte, Catch and Release, CpG WIZ and 4G10 are registered trademarks of Millipore Corporation.

KinEASE, AdenoSilence, siRNA/siAb, Incelligent, ChemiKine are trademarks of Millipore Corporation.

SMARTpool is a registered trademark of Dharmacon, Inc.

Iressa is a registered trademark of Astra Zeneca UK LTD.

Alexa Fluor is a registered trademark of Molecular Probes, Inc.

HTRF is a registered trademark of Cisbio International.

Lit. No. PB1726EN00 Printed in U.S.A. 1779-0906

© 2007 Millipore Corporation, Billerica, MA 01821 All rights reserved.

