

3050 Spruce Street
Saint Louis, Missouri 63103 USA
Telephone 800-325-5832 • (314) 771-5765
Fax (314) 286-7828
email: techserv@sial.com
sigma-aldrich.com

# **ProductInformation**

Lectin from *Phaseolus vulgaris* Phytohemagglutinin PHA-P Cell Culture Tested

Product Number **L 9132** Storage Temperature 2-8 °C

# **Product Description**

PHA-P is a mixture of PHA-E (MW = 128 kDa) and PHA-L (MW = 126 kDa).

Lectins are proteins or glycoproteins of non-immune origin that agglutinate cells and/or precipitate complex carbohydrates. Lectins are capable of binding glycoproteins even in presence of various detergents. The agglutination activity of these highly specific carbohydrate-binding molecules is usually inhibited by a simple monosaccharide, but for some lectins, di, tri, and even polysaccharides are required.

Lectin PHA-P is not inhibited easily by monosaccharides, but may be inhibited by oligosaccharides.

Lectins are isolated from a wide variety of natural sources, including seeds, plant roots and bark, fungi, bacteria, seaweed and sponges, mollusks, fish eggs, body fluids of invertebrates and lower vertebrates, and from mammalian cell membranes. The precise physiological role of lectins in nature is still unknown, but they have proved to be very valuable in a wide variety of applications *in vitro*, including:

- blood grouping and erythrocyte polyagglutination studies.
- 2. mitogenic stimulation of lymphocytes.
- 3. lymphocyte subpopulation studies.
- 4. fractionation of cells and other particles.
- histochemical studies of normal and pathological conditions.

Sigma offers a range of lectins suitable for the above applications. Most Sigma lectins are highly purified by affinity chromatography, but some are offered as purified or partially purified lectins, suitable for specific applications.

Many of the lectins are available conjugated to (conjugation does not alter the specificity of the lectin):

- 1. fluorochromes (for detection by fluorimetry).
- enzymes (for enzyme-linked assays).
- 3. insoluble matrices (for use as affinity media).

Please refer to the table for general information on the most common lectins.

This lectin has been tested with cell lines to verify the product is not cytotoxic.

#### **Procedure**

A general agglutination procedure using this lectin with 96 well plates is as follows:

- Prepare a lectin solution of 1 mg/ml in PBS buffer, pH 6.8.
- 2. Pipette 50  $\mu$ l of fresh PBS into each well and add 50  $\mu$ l of the lectin solution into the first well.
- 3. Serial dilutions are made by pipetting 50 μl from each successive well into the next well.
- 4. Blood type A with a 2% hematocrit is used as the substrate.
- Pipette 50 μl of blood into each well.
- 6. Visually determine agglutination.

# **Precautions and Disclaimer**

For Laboratory Use Only. Not for drug, household or other uses.

#### **Preparation Instructions**

This lectin is soluble in phosphate buffered saline, pH 7.2 (1 mg/ml).

## Storage/Stability

Aggregation is thought to occur in the presence of high concentrations of 2-mercaptoethanol.

Lastin	MAA/ /L-D - \	O l	Mitogenic		
Lectin	MW (kDa)	Subunits	Blood Group	Sugar	Activity
Abrus precatorius	134	4	_	gol	+
Agglutinin Abrin A (toxin)	60	4 2		gal	
	63.8			gal	
Abrin B (toxin)	58.5	2(αβ)		gal	
Agarius bisporus		_	_	β-gal(1→3)galNAc	
Anguilla anguilla	40	2	H	α-L-Fuc	
Arachis hypogaea	120	4	T	β-gal(1→3)galNAc	
Artocarpus integrifolia Bandeiraea simplicifolia	42	4	Т	α-gal→OMe	+
BS-I	114	4	A, B	$\alpha$ -gal, $\alpha$ -galNAc	
BS-I-A <sub>4</sub>	114	4	Α	α-galNAc	
BS-I-B <sub>4</sub>	114	4	В	α-gal	
BS-II	113	4	acq, B, Tk, T	glcNAc	
Bauhinia purpurea	195	4	_	β-gal(1→3)galNAc	+
Caragana arborescens	60; 120 <sup>a</sup>	2/4	_	galNAc	
Cicer arietinum	44	2	_	fetuin	
Codium fragile	60	4	_	galNAc	
Concanavalin A	102	4	_	$\alpha$ -man, $\alpha$ -glc	+_
Succinyl-Concanavalin A	51	2	_	$\alpha$ -man, $\alpha$ -glc	+ <sup>b</sup>
Cytisus scoparius	_	_	_	galNAc, gal	
Datura stramonium	86	$2(\alpha\beta)$	_	(glcNAc) <sub>2</sub>	
Dolichos biflorus	140	4	$A_1$	α-galNAc	
Erythrina corallodendron	60	2	_	β-gal(1→4)glcNAc	+
Erythrina cristagalli	56.8	2(αβ)	_	β-gal(1→4)glcNAc	
Euonymus europaeus	166	$4(\alpha\beta)$	B, H	α-gal(1→3)gal	+
Galanthus nivalis	52	4	(h)	non-reduc. α-man	
Glycine max	110	4	_	galNAc	+ <sup>c</sup>
Helix aspersa	79	_	Α	galNAc	
Helix pomatia	79	6	Α	galNAc	
Lathyrus odoratus	40-43	4(αβ)	_	$\alpha$ -man	+
Lens culinaris	49	2	_	α-man	+
Limulus polyphemus	400	18	_	NeuNAc	
Bacterial agglutinin	_	_	_	galNAc, glcNAc	
Lycopersicon esculentum	71	_	_	(glcNAc) <sub>3</sub>	
Maackia amurensis	130	2(αβ)	0	sialic acid	+
Maclura pomifera	40-43	$2(\alpha\beta)$	_	$\alpha$ -gal, $\alpha$ -galNAc	
Momordica charantia	115-129	$4(\alpha\beta)$	_	gal, galNAc	
Naja mocambique mocambique		— ( - 1 · 7 ·	_	_	
Naja naja kaouthia	_	_	_	_	
Narcissus pseudonarcissus	26	2	(h)	α-D-man	
Perseau americana	_	_	· · · / -	_	
Phaseolus coccineus	112	4	_	_	
Phaseolus limensis	247(II)	8	Α	galNAc	+
	124(III)	4	, .	gao	
Phaseolus vulgaris	_ · (···)	-			
PHA-E	128	4	_	oligosaccharide	+
PHA-L	128	4	_	oligosaccharide	+
PHA-P	.20	·		Singoodooridindo	•
PHA-M					

			Mitogenic		
<u>Lectin</u>	MW (kDa)	Subunits	Blood Group	Sugar	Activity
Phytolacca americana	32	_	_	(glcNAc) <sub>3</sub>	+
Pisum sativum	49	4(αβ)	_	α-man	+
Pseudomonas aeruginosa PA-I	13-13.7	_	_	gal	+ <sup>c</sup>
Psophocarpus tetragonolobus	35	1	_	galNAc, gal	
Ptilota plumosa	65; 170	_	В	$\alpha$ -gal	
Ricinus communis					
Toxin, RCA <sub>60</sub>	60	2	_	galNAc, β-gal	
Toxin, RCA <sub>120</sub>	120	4	_	β-gal	
Sambucus nigra	140	4(αβ)	_	αNeuNAC(2→6)gal	+ <sup>c</sup>
-				galNAc	
Solanum tuberosum	50; 100 <sup>a</sup>	1, 2	_	(glcNAc)₃	
Sophora japonica	133	4	A, B	β-galNAc	
Tetragonolobus purpureas	120(A)	4	Н	$\alpha$ -L-fuc	
	58(BA)	2	Н	$\alpha$ -L-fuc	
	117(C)	4	Н	$\alpha$ -L-fuc	
Triticum vulgaris	36	2	_	(glcNAc) <sub>2</sub> , NeuNAc	+
Ulex europaeus					
UEA I	68	_	Н	$\alpha$ -L-fuc	
UEA II	68	_	_	(glcNAc) <sub>2</sub>	
Vicia faba	50	$4(\alpha\beta)$	_	man, glc	+
Vicia sativa	40	$4(\alpha\beta)$	_	glc, man	+
Vicia villosa	139	4	$A_{1+}T_n$	galNAc	
$A_4$	134	4	$A_1$	galNAc	
$B_4$	143	4	$T_n$	galNAc	
Vigna radiata	160	4	_	α-gal	
Viscum album	115	$4(\alpha\beta)$	_	β-gal	
Wisteria floribunda	68	2	_	galNAc	

<sup>&</sup>lt;sup>a</sup> Concentration-dependent molecular weight

### References

 Rueben, L., et al., Activities of lectins and their immobilized derivatives in detergent solutions. Implications on the use of lectin affinity chromatography for the purification of membrane glycoproteins. Biochemistry, 16, 1787-1794 (1977).

IRB/MWM/JRC/NSB/SAG 11/02

b Non-agglutinating and mitogenic

<sup>&</sup>lt;sup>c</sup> Mitogenic for neuraminidase-treated lymphocytes