

Product Information

00254 *Fluoroselect*™ Glycerol Kit

Product description

Glycerol (glycerin or glycerine, C₃H₅(OH)₃) is widely used in foods, beverages and pharmaceutical formulations. It is also a main by-product of biodiesel production. Simple, direct and automation ready procedures for measuring glycerol concentrations find wide applications.

Sigma-Aldrich's glycerol assay uses a single working reagent that combines glycerol kinase, glycerol phosphate oxidase and color reactions in one step. The color intensity of the reaction product at 570 nm, or fluorescence intensity at $\lambda_{ex/em} = 530/590$ nm is directly proportional to glycerol concentration in the sample.

Detection ranges and limits

Linear detection range: 2 -100 μ M, 18-920 μ g/dL or 0.18-9.20 ppm.

Sensitive and accurate. Use as little as 10 μ L samples. Simple and convenient. The procedure involves addition of a single working reagent and incubation for 20 min at room temperature.

Improved reagent stability. The optimized formulation has greatly enhanced the reagent and signal stability.

Equipment required but not included

[Z805491-1EA](#) *FluoroSELECT*™ Single channel fluorometer λ_{ex} 530 nm; λ_{em} 590 nm

[Z805823-100EA](#) Glass vials for *FluoroSELECT*™ fluorometer

Components

1. 24 mL Assay Buffer
2. 500 μ L Enzyme Mix
3. 250 μ L ATP
4. 220 μ L Dye Reagent
5. 100 μ L Glycerol Standard

The kit is sufficient for approximately 200 assays.

Storage conditions

Store at -20°C

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

Procedure

1. Standard. Prepare 1 mM standard by mixing 5 μ L of the provided standard with 495 μ L distilled H₂O.

Then mix 50 μ L of the 1mM standard with 450 μ L H₂O to obtain 100 μ M glycerol standard. In separate mini-glass tubes add 10 μ L H₂O ("Blank"), 10 μ L 100 μ M glycerol standard ("Std"), and 10 μ L sample.

2. Prepare enough working reagent for all assay tubes, by mixing per tube: 100 μ L assay buffer, 2 μ L enzyme mix, 1 μ L ATP and 1 μ L dye reagent in a clean Eppendorf tube. Then add 100 μ L working reagent to each tube and mix. Incubate for 20 min in the dark.

3. Switch on the reader. To calibrate the reader, place the "Blank" tube into the sample holder. Press "Calibrate", "Assay 1", then "Blank". Reader starts measuring. Press left arrow on "<-Std ->", until the window shows "100.000". Place the "Std" tube into the sample holder. Press "Measure". The reader shows "Calibration Finished". Press "Return".

4. Measure. Place the sample tube into the sample holder. Press "Measure", "Assay 1", "Measure". The glycerol concentration (μ M) will be displayed in the window. Record the data, or press "Save" to save the data for later retrieval. Press "Return" and then press "Measure" for the next sample.