



## Product Information

### N-ACETYL-ILE-GLU-THR-ASP-7-AMIDO-4-METHYLCOUMARIN

Product Number **A 4188**

Storage Temperature  $-20^{\circ}\text{C}$

(Ac-IETD-AMC)

#### Product Description

Appearance: White Powder

Formula Weight: 675.7

Molecular Formula  $\text{C}_{31}\text{H}_{41}\text{N}_5\text{O}_{12}$

Purity:  $\approx 97\%$  by HPLC.

Fluorescent substrate for caspase-8 and granzyme B.

- Fluorometric detection when AMC is cleaved from peptide (excitation wavelength 380 nm emission wavelength 460 nm)
- Molar Extinction coefficient:  $1.78 \times 10^4$  liter/mole cm in EtOH.
- Absorption maximum: 354 nm in EtOH
- AMC is soluble in DMF or absolute ethanol
- Sensitivity of enzyme assay relevant to fluorescence detection is greatly increased over 4-Methoxy-2-Naphthylamide
- AMC by the Ames Test has been shown to be a non-mutagenic chemical

#### Preparation Instructions

Soluble in DMF or absolute ethanol at 20 mM.

#### Storage/Stability

Store tightly sealed and desiccated at  $-20^{\circ}\text{C}$ . Allow powder to reach room temperature before opening vial. May be stored desiccated in solid form at room temperature for one year. Store DMSO/DMF solutions at  $-20^{\circ}\text{C}$  for up to 6 months.

#### Procedure

##### Fluorometric Enzyme Assay

- Buffer: 100 mM HEPES, pH 7.5, 20 % (v/v) glycerol, 5 mM DTT, 0.5 mM EDTA
  - Substrate: 20 mM stock solution of Ac-IETD-AMC in DMF
  - Enzyme: Cell lysate or purified enzyme solution ( $\sim 15$  ng enzyme)
  - Fluorescence Standard: 80  $\mu\text{M}$  free AMC (Product Number A 9891) in DMF
1. Add 10  $\mu\text{l}$  of enzyme to 470  $\mu\text{l}$  buffer. Mix. Incubate at  $30^{\circ}\text{C}$  for 30 minutes.
  2. With fluorometer adjusted to 380 nm excitation and 460 nm emission, add 20  $\mu\text{l}$  of substrate to enzyme solution.
  3. Record increase in fluorescence (FLU) per minute from  $T_0$  to  $T_{\text{end}}$  where the fluorescence generated at  $T_{\text{end}}$  is significantly different from that of  $T_0$ .
  4. Calculate the ? FLU/min. from the linear portion of the curve.
  5. Record fluorescence units (FLU) generated by 10  $\mu\text{l}$ , 20  $\mu\text{l}$ , and 30  $\mu\text{l}$  free AMC and 490  $\mu\text{l}$  (1.6  $\mu\text{M}$ ), 480  $\mu\text{l}$  (3.2  $\mu\text{M}$ ), and 470  $\mu\text{l}$  (4.8  $\mu\text{M}$ ) buffer solution, respectively. These solutions contain 0.8, 1.6 and 2.4 nanomoles, respectively, of free AMC product in 0.5 ml of solution.
  6. Graph the fluorescence units (FLU) vs.  $\mu\text{M}$  the amount of free AFC (nanomoles). The standard curve is the best line connecting the data points. Determine the value of fluorescent units per nanomole (FLU/nmole) of free AMC from the standard curve.
  7. Calculate activity as follows:

$$1 \text{ unit of activity (nmole/min/ml)} = \frac{(\text{?FLU/min}) \times (\text{dilution factor})}{(\text{FLU/nmole}) \times (\text{Vol.})}$$

**?FLU/min** = value determined for enzyme assay in step 4

**Dilution factor** = any dilution of original protein sample prior to addition to reaction.

**FLU/ nmole** = value determined from standard curve in step 6

**Vol.** = volume in ml of enzyme solution in the reaction

#### References

1. Thornberry, N.A., et al., A combinatorial approach defines specificities of members of the caspase family and granzyme B. Functional relationships established for key mediators of apoptosis. *J. Biol. Chem.*, **272**, 17907-17911 (1997).
2. Han, Z., et al., A sequential two-step mechanism for the production of the mature p17:p12 form of caspase 3 in vitro. *J. Biol. Chem.*, **272**, 13432-13436 (1997).

lpg 11/00