# **Product Information**

#### sigma-aldrich.com

3050 Spruce Street, Saint Louis, MO 63103 USA Tel: (800) 521-8956 (314) 771-5765 Fax: (800) 325-5052 (314) 771-5757 email: techservice@sial.com sigma-aldrich.com

#### Apyrase, High Activity from potato (S. tuberosum)

Catalog Number **A2230** Storage Temperature –20 °C

CAS RN 9000-95-7 EC 3.6.1.5

# **Product Description**

Apyrase has adenosine 5'-triphosphatase and adenosine 5'-diphosphatase activities. There are at least two isoenzymes present in different varieties of *S. tuberosum*:<sup>1,2</sup> one with a high ATPase/ADPase ratio (~10) and another with a low ratio (~1). The two isozymes are glycoproteins with approximately the same molecular mass, 49 kDa (gel filtration)<sup>2</sup> and 45 kDa ( $\gamma$ -ray inactivation).<sup>3</sup>

$$ATP \rightarrow ADP + P_i \rightarrow AMP + 2 P_i$$

Isolectric point:2

pl = 8.74 (Pimpernel isoform) pl = 6.69 (Desiree isoform)

Divalent metal ions are required for activity and best activity is observed with calcium ion at 5 mM.

# Optimal pH:<sup>2</sup> pH 6.0, (50% activity at pH 4.6 and 8.4) For hydrolysis of organic di- and triphosphates, the optimal pH is 6, and for inorganic substrates, the optimal pH is 5.1.

This product is derived from red potatoes and supplied as a lyophilized powder containing  $\geq$ 30% protein with the balance potassium succinate buffer salts. It contains predominately the Desiree isoform.

Specific Activity:

ATPase ≥600 units/mg protein ADPase >50% of base activity Acid phosphatase ≤2% of base activity

Unit Definition: One unit will liberate 1.0  $\mu$ mole of inorganic phosphate from ATP or ADP per minute at pH 6.5 at 30 °C.

#### Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

#### **Preparation Instructions**

This product is soluble in water (1 mg/ml).

To prepare solutions of <1 mg/ml, dissolve the product in HEPES buffer, pH 7.5, containing 1 mM MgC1<sub>2</sub>, 1 mM DTT, 1 mM EDTA, and 1 mg/ml bovine serum albumin.

## Storage/Stability

Store the product at -20 °C. When stored at -20 °C, the enzyme retains activity for at least one year.

An enzyme solution ( $\geq$ 1 mg/ml in water) retains activity stored frozen in aliquots. This solution retains activity between pH 5–7. Outside this pH range, activity is lost rapidly.

Repeated freeze-thaw cycles and room temperature exposure for several hours will result in loss of activity. A solution of apyrase stored at 2–8 °C will gradually form a black insoluble precipitate with nearly the same activity as the soluble form.

## References

- 1. Molnar, J., and Lorand, L., Studies on Apyrases. Arch. Biochem. Biophys., **93**, 353 (1961).
- Kettlun, A., *et al.*, Properties of Two Apyrases from Solanum tuberosum, Phytochemistry, **21**, 551 (1982).
- 3. Traverso-Cori, A., *et al.*, Arch. Biochem. Biophys., **109**, 173 (1965).

RBG,MAM 01/08-1

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.

