

Product Information

β -Galactosidase from *E. coli*, recombinant overexpressed in *E. coli*


Product Number **G 3153**
Storage Temperature $-20\text{ }^{\circ}\text{C}$

CAS[#] 9031-11-2
EC 3.2.1.23

Synonyms: Lactase; β -D-Galactoside galactohydrolase;
 β -D-Galactopyranosidase; β -Lactosidase

Product Description

The enzyme cleaves terminal galactose residues, which are β -1 \rightarrow 4 linked to a monosaccharide, glycopeptide, or oligosaccharide.



Molecular weight: \sim 540 kDa (tetramer)

pH Optimum:¹ 6.8–7.4

The enzyme requires Mg^{2+} for optimal activity. The enzymatic assay performed by Sigma uses 1 mM MgCl_2 in the reaction mixture.

Inhibitors:

2-mercaptoethanol¹
citrate¹
 Zn^{2+} ions²
ethanolamine¹

The enzyme is a lyophilized powder containing phosphate buffer salts and sucrose. It does not contain any substances which may interfere with the derivatization of $-\text{NH}_2$ or $-\text{SH}$ groups, for example, ammonium salts, primary amines, 2-mercaptoethanol, etc. Thus, an aqueous solution of the lyophilized material can be directly used for derivatization, such as an enzyme label for IgG, without prior dialysis or gel filtration.³ The protein contains 12–20 moles of free thiol per mole of enzyme.

Specific activity: \geq 500 units/mg protein

Unit Definition: One unit will hydrolyze 1 μ mole of 2-nitrophenyl β -D-galactopyranoside to 2-nitrophenol and D-galactose per minute at the pH 7.3 and $37\text{ }^{\circ}\text{C}$.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

Storage/Stability

It is recommended to store the product at $-20\text{ }^{\circ}\text{C}$, under nitrogen as originally sealed. If stored properly, no significant decrease in either specific activity or in the number of free thiol groups will occur within 12 months. Storage at ambient temperature during shipment does not have any detrimental effect on the product.

References

1. Wallenfels, K., and Weil, R., The Enzymes, Third Edition (Boyer, P.D., ed.), **7**, 617–663 (1972).
2. Wallenfels, K., and Malhotra, O.M., The Enzymes, Second Edition (Boyer, P.D., *et al.*, eds.), **4**, 409–430 (1960).
3. O'Sullivan, M.J., *et al.*, Ann. Clin. Biochem., **16**, 221-240 (1979).

AE,MAM 02/05-1

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.