



3050 Spruce Street
Saint Louis, Missouri 63103 USA
Telephone (800) 325-5832 (314) 771-5765
Fax (314) 286-7828
email: techserv@sial.com
sigma-aldrich.com

ProductInformation

pcDNA 3.1 (-) MycHisC Expression Vector

Product Number **E0648**

Store at -20°C

Product Summary:

Package Size: 20 mg

Lyophilized in 10 mM Tris-HCl and
1 mM EDTA pH 7.5

Plasmid Size: 5.5 kb

Antibiotic Resistance: Ampicillin and Neomycin

Resuspend the expression vector in 0.2 µm filtered
water.

Description:

pcDNA3.1 (-)/myc-his C is an expression vector that
contains the strong CMV enhancer-promoter for high
level expression of recombinant proteins.

pcDNA3.1(-)/myc-his C expression vector gene for the
expression of a myc-his tag for detection of the
recombinant protein.

The pcDNA3.1(-)/myc-his C expression vector
contains a multiple cloning site for easy cloning.

Restriction endonuclease digestion of supercoiled
pcDNA 3.1(-)/myc-his C yields the following
fragments: Apa I (5500 bp), Xba I 5,500) and SnaB I
(5,088, 412 bp).

pcDNA3.1(-) /myc-his C expression vector has a
multiple cloning site in an opposite orientation from the
pcDNA3.1(+). The vector is supplied in one of three
reading frames to facilitate in frame cloning with a C-
terminal tag containing a polyhistidine metal-binding
region and the myc (c-myc) epitope. The human
cytomegalovirus immediate-early (CMV) promoter
provides high-level expression in many mammalian
cells. Also, the vector will replicate in cell lines which
are infected with SV40.

T7 Multiple Cloning Site

Nhe I
Pme I
Apa I
Xba I
Xho I
Not I
BstX I
EcoR V
EcoR I
BstX I
BamH I
Asp718 I
Kpn I
Hind III
Apa I
Xba I
Myc epitope
His₆
Term
Pme I

Additional Information:

To ensure proper expression of your recombinant
protein, you must clone your gene in frame with the
C-terminal peptide. The vector can be used with two
additional vectors to facilitate cloning of the gene in
frame.

The expression vectors are:

**E0398 pc DNA3.1(-) MycHis A
and E0523 pcDNA3.1(-) MycHisB**

Features of pcDNA3.1 (-)/Myc-His C	Benefits of pcDNA3.1(-)/Myc-His C
Human cytomegalovirus (CMV) immediate-early promoter/enhancer.	Permits efficient, high-level expression of recombinant protein.
T7 promoter/priming site	Allows for <i>in vitro</i> transcription in the sense orientation and sequencing through the insert.
Multiple cloning site	Allows insertion of your gene and facilitates cloning in frame with the polyhistidine C-terminal tag.
Myc epitope (c-myc)	Allows detection of the recombinant protein.
C-terminal polyhistidine tag	Permits purification of the recombinant protein on metal chelating resins.
pcDNA3.1/BGH reverse priming site	Permits sequencing through the insert.
Bovine growth hormone (BGH) polyadenylation signal	Efficient transcription termination and polyadenylation of mRNA.
F1 origin	Allows rescue of single-DNA.
Neomycin (G418) resistance gene	Selection of stable transfectants in mammalian cells
SV40 polyadenylation signal	Efficient transcription termination and polyadenylation of mRNA
ColE1 origin (pUC-derived)	High copy number replication and growth in <i>E. coli</i>
Ampicillin resistance gene (β -lactamase)	Selection of vector in <i>E. coli</i>
SV40 early promoter and origin	Allows efficient, high-level expression of neomycin resistance gene and episomal replication in cells expressing SV40 large T antigen.

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.