

User Manual

Ready-to-Assay™ D4 Dopamine Receptor Frozen Cells

HTS223RTA

FOR RESEARCH USE ONLY

Not for use in diagnostic procedures. Not for Human or Animal Consumption

Product Overview

Ready-to-Assay™ GPCR frozen cells are designed for simple, rapid calcium assays with no requirement for intensive cell culturing. The freezing conditions have been optimized to provide cells with high viability and functionality post-thaw. The user simply thaws the cells and resuspends them in media, dispenses cell suspension into assay plates and, following overnight recovery, assays for calcium response.

The neurotransmitter dopamine is involved in a wide variety of neuroendocrine, emotional, cognitive and locomotor functions. These activities of dopamine are mediated by a group of 5 G-protein-coupled receptors, 2 of which (D₁ and D₅) couple to G_s and 3 of which (D₂, D₃, D₄) couple to G_i (Missale et al., 1998). The D₄ receptor is highly polymorphic, particularly in the third cytoplasmic loop, which contains from 2 to 7 repeats of 16 amino acids. Individuals with 7 repeats have increased risk of developing ADHD (Thapar et al., 2005). Mice engineered to lack D₄ have reduced locomotor activity, increased sensitivity to drugs of abuse that elevate dopamine (cocaine, methamphetamine and ethanol), and reduced response to novelty (Rubinstein et al., 1997; Dulawa et al., 1999). Cloned human D₄-expressing cell line is made in the Chem-5 host, which supports high levels of recombinant D₄ expression on the cell surface and contains optimized levels of a recombinant promiscuous G protein to couple the receptor to the calcium signaling pathway. Thus, the cell line is an ideal tool for screening for agonists, antagonists, and modulators at D₄.

Materials Provided

- Pack contains 2 vials of mycoplasma-free cells, 1 mL per vial. Store in liquid N₂.
- 50 mL of Media Component. Store at 4 °C (-20°C for prolonged storage).

Use Restrictions

Please see User Agreement (Label License) for further details. One such restriction is that the contents of the supplied vial(s) are limited to a single use and shall not be propagated and/or re-frozen by licensee.

GMO

This product contains genetically modified organisms.
Este producto contiene organismos genéticamente modificados.
Questo prodotto contiene degli organismi geneticamente modificati.
Dieses Produkt enthält genetisch modifizierte Organismen.
Ce produit contient organismes génétiquement des modifiés.
Dit product bevat genetisch gewijzigde organismen.
Tämä tuote sisältää geneettisesti muutettuja organismeja.
Denna produkt innehåller genetiskt ändrade organismer.

Applications

Calcium Flux Assays

Applications Data




Figure 1. Representative data for activation of D4 receptor. Calcium flux in D4 -expressing Chem-5 cell line induced by Dopamine. D4 -expressing Chem-5 cells were loaded with a calcium dye, and calcium flux in response to the indicated ligand(s), 4-fold serial dilution with each concentration performed in duplicate, was determined on a Molecular Devices FLIPR Tetra® with ICCD camera. Maximal fluorescence signal obtained in this experiment was 11,600 RLU (Relative Light Units).

Table 1. EC₅₀ values of D4 -expressing Chem-5 cells

Ligand	Assay	Potency (nM)	Reference
Dopamine	Calcium Flux	10	Internal Data

Assay Setup

1. Immediately upon receipt, thaw cells or place cells in liquid nitrogen.
2. Thaw cells rapidly by removing from liquid nitrogen and immediately immersing in a 37 °C water bath. Immediately after ice has thawed, sterilize the exterior of the vial with 70% ethanol.
3. Add 1mL of pre-warmed Media Component to each vial of cells. Place contents from two vials into a 15 mL conical tube and bring the volume to 10 mL of Media Component.
4. Centrifuge the cell suspension at 190 x g for four minutes.
5. Remove supernatant and add 10.5 mL of pre-warmed Media Component to resuspend the cell pellet.
6. Seed cell suspension into appropriate assay microplate (100 µL/well for 96-well plate, 25 µL/well for 384-well plate).
7. When seeding is complete, place the assay plate at room temperature for 30 minutes.
8. Move assay plate to a humidified 37 °C 5% CO₂ incubator for 24 hours.
9. After 24-hour incubation, remove assay plate from the incubator and wash sufficiently with Hank's Balanced Salt Solution (HBSS) supplemented with 20 mM HEPES, 2.5 mM Probenecid at pH 7.4 to remove all trace of Media Component.
10. Prepare Fluo-8, AM (AAT Bioquest: 21080) Ca²⁺ dye by dissolving 1mg of Fluo-8 NW in 200 µL of DMSO. Once dissolved place 10 µL of Fluo-8 NW Ca²⁺ dye solution into 10 mL of HBSS 20 mM HEPES, 2.5 mM Probenecid pH 7.4 buffer and apply to assay microplate (Ca²⁺ dye at 10 µL /10 mL is sufficient for loading one (1) microplate).
11. Set-up FLIPR to dispense 3x ligand to appropriate wells in the assay plate. Set excitation wavelength at 470-495 nm (FLIPR Tetra®) or 485 nm (FLIPR1, FLIPR2, FLIPR3) and emission wavelength at 515-565 nm

(FLIPRTETRA) or emission filter for Ca²⁺ dyes (FLIPR1, FLIPR2, FLIPR3). Set pipet tip height to 5 µL below liquid level and dispense rate to 75 µL/sec (96-well format) or 50 µL/sec (384-well format). Set up plate layout and tip layout for each individual experiment. Set time course for 180 seconds, with ligand addition at 10 seconds.

12. Ligands are prepared in non-binding surface Corning plates (Corning 3605: 96-well or Corning 3574: 384-well).
13. After the run is complete, negative control correction is applied and data analyzed utilizing the maximum statistic.

Assay Materials

All items may be purchased at SigmaAldrich.com unless otherwise noted.

Description	Catalogue Number
HBSS	SH30268.02 (Hyclone)
HEPES 1M Stock	TMS-003-C
Probenicid	P8761
Quest Fluo-8™, AM	21080 (AAT Bioquest)
Dopamine ligand	H8502
Non-binding white plates (for ligand prep)	3605 (96-well)/3574 (384-well) (Corning)
Black (clear bottom) tissue-culture treated plates	3904 (96-well)/3712 (384-well) (Corning)

FLIPR Settings

Settings for FLIPR^{TETRA}® with ICCD camera option.

Read Mode	Fluorescence
Ex/Em	Ex470_495/Em515_575
Camera Gain	2000
Gate Open	6%
Exposure Time	0.53
Read Interval	1 second
Dispense Volume	50 µl (25 µl for 384-well)
Dispense Height	25 µl (50 µl for 384-well)
Dispense Speed	75 µl L/sec (50 µl for 384-well)
Expel Volume	0 µl
Analysis	Subtract Bias Sample 1

Host Cell

Chem-5, an adherent rat hematopoietic cell line expressing endogenous G-15 protein as well as an exogenous proprietary promiscuous Ga protein.

Exogenous Gene Expression

DRD4.4 cDNA (Accession Number: NM_000797; see CODING SEQUENCE below) expressed from a proprietary pHS plasmid.

Coding Sequence

```
ATGGGGAACCGCAGCACCCGCGACCGGACGGGCTGCTGGCTGGGCGCGGGCCGCGGGGCATCTGCGGGGCATCTGCGGGGCTG - 90
1 - M G N R S T A D A D G L L A G R G P A A G A S A G A S A G L - 30
91 - GCTGGGCAGGGCGCGGCGCTGGTGGGGCGTGTGCTCATCGGCGCGGTGCTCGCGGGAACTCGCTCGTGTGCGTGAGCGTGGCC - 180
31 - A G Q G A A A L V G G V L L I G A V L A G N S L V C V S V A - 60
181 - ACCGAGCGCGCCCTGCAGACGCCCACCAACTCCTTCATCGTGAGCCTGGCGGCCCGACCTCCTCCTCGCTCTCCTGGTGCTGCCGCTC - 270
61 - T E R A L Q T P T N S F I V S L A A A D L L L A L L V L P L - 90
271 - TTCGTCTACTCCGAGGTCCAGGGTGGCGCGTGGCTGTGAGCCCCGCTGTGCGACGCCCTCATGGCCATGGACGTATGCTGTGCACC - 360
91 - F V Y S E V Q G G A W L L S P R L C D A L M A M D V M L C T - 120
361 - GCCTCCATCTTCAACCTGTGCGCCATCAGCGTGGACAGGTTTCGTGGCCGTGGCCGTGCGCTACAAACCGGACGGTGGGAGCCGC - 450
121 - A S I F N L C A I S V D R F V A V A V P L R Y N R Q G G S R - 150
451 - CGGCAGCTGTGCTCATCGGCGCCACGTGGTGTGCTGCCGCGCGGTGGCGGCCCGCTACTGTGCGGCCTCAACGACGTGCGCGGCCGC - 540
151 - R Q L L L I G A T W L L S A A V A A P V L C G L N D V R G R - 180
541 - GACCCCGCGTGTGCGCCTGGAGGACCGCGACTACGTGGTCTACTCGTCCGTGTGCTCCTTCTTCTACCTGCCCCGCTCATGCTGTG - 630
181 - D P A V C R L E D R D Y V V Y S S V C S F F L P C P L M L L - 210
631 - CTCTACTGGGCCACGTTCCGCGGCCGTCAGCGCTGGGAGGTGGCACGTGCGGCCAAGCTGCACGGCCGCGCGCCCCGCGACCCAGCGGC - 720
211 - L Y W A T F R G L Q R W E V A R R A K L H G R A P R R P S G - 240
721 - CCTGGCCCCGCTTCCCCACGCCACCCGCGCCCCGCTCCCCAGGACCCCTGCGGCCCGACTGTGCGCCCCCGCGCCCGGCTTCCC - 810
241 - P G P P S P T P P A P R L P Q D P C G P D C A P P A P G L P - 270
811 - CGGGTCCCTGCGGCCCGACTGTGCGCCCGCGGCCAGCCTCCCCAGGACCCCTGTGGCCCCGACTGTGCGCCCCCGCGCCCGGC - 900
271 - R G P C G P D C A P A A P S L P Q D P C G P D C A P P A P G - 300
901 - CTCCCCCGGACCCCTGCGGCTCCAACGTGTGCTCCCCCGACGCCGTGAGAGCCGCGGCTCCACCCAGACTCCACCGCAGACCCGC - 990
301 - L P P D P C G S N C A P P D A V R A A A L P P Q T P P Q T R - 330
991 - AGGAGCGCGGTGCCAAGATCACCGCCGGGAGCGCAAGGCCATGAGGGTCTGCCGGTGGTGGTGGGGCCTTCTGCTGTGCTGGACG -1080
331 - R R R R A K I T G R E R K A M R V L P V V V G A F L L C W T - 360
1081 - CCCTTCTCGTGGTGACATCACGAGCGCTGTGTCTGCTGCTCCGTGCCCCCGGGCTGGTCAGCGCGTACCTGGCTGGGCTAC -1170
361 - P F F V V H I T Q A L C P A C S V P P R L V S A V T W L G Y - 390
1171 - GTCAACAGCGCCCTCAACCCGTCATCTACACTGTCTTCAACGCCGAGTTCGCCAACGTCTCCGCAAGGCCCTGCGTGCCTGTGCTGA -1260
391 - V N S A L N P V I Y T V F N A E F R N V F R K A L R A C C Stp - 420
```

Related Products

- HTSCHEM-1RTA Ready-to-Assay™ Chem-1 host frozen cells (control cells).
Note: Chem-5 cells are derived from Chem-1 cells.
- HTS223M ChemoScreen™ D4 Dopamine receptor membrane Preparation Recombinant Human D4 Dopamine Receptor.

References

1. Chemel BR et al. (2006) WAY-100635 is a potent dopamine D4 receptor agonist. *Psychopharmacology* 188: 244-251.
2. Dulawa SC et al. (1999) Dopamine D4 receptor-knock-out mice exhibit reduced exploration of novel stimuli. *J. Neurosci.* 19: 9550-9556.
3. Millan MJ et al. (1998) S 18126 ([2-[4-(2,3-dihydrobenzo [1,4] dioxin-6-yl) piperazin-1-yl methyl] indan-2-yl]), a potent, selective and competitive antagonist at dopamine D4 receptors: an in vitro and in vivo comparison with L 745,870 (3-(4-[4-chlorophenyl]piperazin-1-yl)methyl-1H-pyrrolo[2, 3b]pyridine) and raclopride. *J. Pharmacol Exp. Ther.* 287: 167-186.
4. Missale C et al. (1998) Dopamine receptors: From structure to function. *Physiol. Rev.* 78.

User Agreement (Label License)

In addition to the General Terms and Conditions section, these specific terms also apply for Ready-to-Assay™ D4 Dopamine Receptor Frozen Cells, Product No. HTS223RTA

BY USING THE THIS PRODUCT LICENSED TO YOU ("LICENSEE") HEREUNDER, YOU ARE HEREBY REPRESENTING THAT YOU HAVE THE RIGHT AND AUTHORITY TO LEGALLY BIND YOURSELF OR YOUR COMPANY, AS APPLICABLE, AND ARE CONSENTING TO BE LEGALLY BOUND BY ALL OF THE TERMS OF THIS USER AGREEMENT ("AGREEMENT"). IF YOU DO NOT AGREE TO ALL THESE TERMS, DO NOT USE THE PRODUCT, AND IMMEDIATELY RETURN SUCH PRODUCTS TO THE APPLICABLE SELLER FOR A REFUND. This is a legal agreement between Licensee and Millipore governing use of the Ready-to-Assay™ Cells products and/or any accompanying operating/use protocols (the "Product(s)") provided to Licensee.

LICENSEE shall obtain no ownership interest in the Product or use/culture protocols accompanying the Product other than through the perpetual limited license granted herein. If the Product is licensed through an authorized Millipore distributor, Licensee shall be obligated to disclose its identity to Eurofins Pharma Bioanalytics Services US Inc. to insure compliance with this User Agreement.

Limited License and Restrictions. Pursuant to the terms and conditions of this Agreement, Millipore conveys to Licensee the non-exclusive and non-transferable right to use the Licensed Product only for Research Purposes conducted by Licensee (whether Licensee is an academic user or a for-profit entity). "Research Purposes" means any biological research and development application or use, including without limitation, developing, demonstrating or validating biological assays, life sciences and/or pharmaceutical research. "Research Purposes" excludes applications outside biology (including but not limited to consumer electronics or materials sciences), and specifically excludes the following applications of whatever kind or nature: Clinical Diagnostics (any use of a product or service for clinical diagnosis where data from an individual's sample is given to such individual or used for the purpose of diagnosis or treatment of a medical condition in such individual, where that result may be used in the treatment of such individual), therapeutics, clinical imaging, environmental testing and cosmetics. Contents of the supplied vial(s) are limited to a single use and shall not be propagated and/or re-frozen by licensee. Licensee cannot sell or otherwise transfer (a) this Product or (b) materials made using this Product to a third party. Licensee may transfer information or materials made through use of this Product to a scientific collaborator, provided that such transfer is not for the commercial purposes, and that such collaborator agrees in writing: (a) not to transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for Research Purposes and not for commercial purposes. Commercial purposes means any activity by a user of the Product for consideration that may include, but is not limited to: (1) operating a service business that uses the Products to develop information or data which is resold for research and development applications; (2) use of the Product in manufacturing; (3) use of the Product for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the Product, whether or not such Product is resold for use in research. Licensee expressly represents and warrants to Millipore that Licensee will properly test and use any Product purchased from Millipore or its affiliated companies in accordance with the practices of a reasonable person who is an expert in the field and in strict compliance with all applicable laws and regulations, now and hereinafter enacted. Licensee agrees to comply with instructions, if any, furnished by Millipore relating to the use of the Product and to not misuse the Product in any manner. Licensee shall not reverse engineer, disassemble or modify the Product or create any derivative works of the written materials accompanying the Product, including but not limited to any material data sheets or similar materials with respect to the Products' specifications. Licensee acknowledges that Millipore or its affiliated companies retains ownership of all patents, copyrights, trademarks, trade secrets and other proprietary rights relating to or residing in the Product or any portion thereof.

Term and Termination. This Agreement commences upon Licensee's use of the Products and shall remain in effect in perpetuity unless terminated sooner as set forth hereunder. Millipore may terminate this Agreement immediately if Licensee breaches any provision herein. Upon any such termination, all rights granted to Licensee hereunder will immediately terminate, and Licensee shall immediately cease using the Product and, at Millipore's option, return or destroy all Products (certifying such destruction to Millipore in writing).

Assignment. Licensee shall not sublicense, assign (by operation of law or otherwise) or otherwise transfer this Agreement or any of the rights or licenses granted under this Agreement without the prior written consent of Millipore. Any attempted assignment, sublicense or transfer by Licensee without such consent shall be null and void.

Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

Technical Assistance

Visit the tech service page at SigmaAldrich.com/techservice.

Terms and Conditions of Sale

Warranty, use restrictions, and other conditions of sale may be found at SigmaAldrich.com/terms.

Contact Information

For the location of the office nearest you, go to SigmaAldrich.com/offices.

The life science business of Merck KGaA, Darmstadt, Germany
operates as MilliporeSigma in the U.S. and Canada.

MilliporeSigma, and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2023 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

Document Template 20769660 Ver 4.0

00072270 Ver 1.0, Rev 15MAR2023, SS

**MILLIPORE
SIGMA**