


Product Information

Anisomycin from *Streptomyces griseolus*

Product Number **A 9789**
Storage Temperature 2-8 °C

CAS RN: 22862-76-6

Synonym: Flagecidin



Molecular Formula: C₁₄H₁₉O₄N
Molecular Weight: 265.30

Melting Point: 140-141 °C¹

E_{mM} (224 nm) = 10.8¹
 E_{mM} (277nm) = 1.8¹
 E_{mM} (283nm) = 1.6¹

$[\alpha]_D^{23} = -30^\circ$ (methanol)¹

Product Description

Anisomycin is a pyrrolidine antibiotic isolated from culture filtrates of *Streptomyces griseolus*. It is primarily an antiprotozoal agent with little antibacterial or antifungal activity. It is a potent, structurally specific and reversible inhibitor of protein and DNA synthesis in eukaryotic organisms including HeLa cells.² In rabbit reticulocytes and in HeLa cells greater than 99% inhibition of protein synthesis was achieved by using 10⁻⁵ M anisomycin. RNA synthesis was unaffected by anisomycin at 10⁻⁶ M.² At minimum concentrations ranging from 1.5 µg/mL to 3.1 µg/mL, anisomycin inhibits bacteria-free cultures of *Trichomonas vaginalis*, *Trichomonas foetus* and *Endamoeba histolytica*.³ Anisomycin has been extensively used in studying the apoptosis process. It affects both pro- and anti-apoptotic mechanisms, depending on the

concentration used. It was found to be a stimulator of two anti-apoptotic proteins Akt and Bcl-2⁴ and an activator of JNKs^{5,6} and P38, which participate in stress and apoptotic responses.⁶ Injections of Anisomycin affect the acquisition of classically conditioned responses in intermediate cerebellum during learning⁷ and the long-term memory in a spatial memory task.⁸

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

Preparation Instructions

Anisomycin is soluble at 25 mg/mL in ethanol and 50 mg/mL in DMSO. Moderately soluble in water at 2 mg/ml. For higher solubility in aqueous solution (100 mg/ml), the pH should be lowered to 5.0. It is also soluble in lower alcohols, esters, ketones, chloroform; slightly soluble in benzene, toluene and hexane.¹ Sigma routinely tests the solubility at 20 mg/mL in methanol, yielding a clear, faint yellow solution.

Storage/Stability

Store the product desiccated at 2-8 °C. Under these conditions the product is stable for 4 years. Aqueous solutions of anisomycin are most stable at neutral pH. Slightly acidic solutions lose activity at relatively slow rates, while alkaline solutions lose activity more rapidly.³ DMSO solutions are stable for at least one month at 2-8 °C.

References

1. Merck Index, 13th ed., No. 673, 113 (2001).
2. Grollman, A.P., J. Biol. Chem., **242**, 3226-3233 (1967).
3. Tanner, F.W. et al., Antibiotics Annual, 809-812 (1954-1955).
4. Torocsik, B. and Szeberenyi, J., Biochem. Biophys. Res. Commun., **278**, 550-6 (2000)
5. Curtin, J.F. and Cotter T.G., Br. J. Cancer, **87**, 1188-94 (2002).

6. Torocsik, B. and Szeberenyi. J., Eur. J. Neurosci., **12**, 527-32 (2000).
7. Bracha, V., et al., Brain Res., **788**, 169-78 (1998).
8. Meiri, N, and Rosenblum, K., Brain Res., **789**, 48-55 (1998).

NDH/PHC 12/04

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.