

# INS-1 832/3 Rat Insulinoma Cell Line

Cancer Cell Line

Cat. # SCC208

Pack size:  $\geq 1 \times 10^6$

viable cells/vial

FOR RESEARCH USE ONLY.  
NOT FOR USE IN DIAGNOSTIC PROCEDURES.  
NOT FOR HUMAN OR ANIMAL CONSUMPTION.  
THIS PRODUCT CONTAINS GENETICALLY MODIFIED ORGANISMS.

Store in liquid nitrogen



Data Sheet

page 1 of 4

## Background:

Glucose-stimulated insulin secretion (GSIS) is potentiated by pancreatic beta cells and is critical to the physiological control of blood glucose levels. Insulin secretion is impaired in type 2 diabetes, and insight into the mechanisms and regulation of insulin secretion is fundamental to understanding the roles of beta cells in metabolic disease.

The rat insulinoma cell line INS-1<sup>1</sup> is a well-established model for studies of pancreatic islet beta-cell function; however, the GSIS response of INS-1 cells may decrease over time. The INS-1 832/3 cell line is a subclone of INS-1 that was selected for robust GSIS, producing and secreting both rat and human insulin. INS-1 832/3 harbors a human insulin expression cassette allowing for human insulin secretion to be maintained over extended passages with selection<sup>2</sup>. INS-1 832/3 cells are widely used to study the mechanisms of cellular insulin secretion, storage and synthesis. INS-1 832/3 cells may be characterized by granular staining for synaptotagmin, as described for the parental cell line<sup>3</sup>. The INS-1 832/3 cell line exhibits the unique feature of additional insulin secretion in response to natural incretin hormones such as glucagon-like peptide 1 (GLP1), pituitary adenylate cyclase-activating peptide (PACAP), and gastric inhibitory peptide (PIP)<sup>4</sup>, making INS-1 832/3 a valuable tool for physiologically relevant investigations of insulin regulation.

## Source

INS-1 832/3 is a derivative of INS-1 cells originally established from an x-ray induced insulinoma in rat<sup>1</sup>. The INS-1 832/3 cell line is a subclone of INS-1 that was stably transfected with a CMV promoter-human insulin expression plasmid carrying a geneticin (G418)-resistance marker for selection<sup>2</sup>.


## Storage and Handling

INS-1 832/3 Rat Insulinoma Cell Line should be stored in liquid nitrogen. The cells can be cultured for at least 10 passages after initial thawing without significantly affecting the cell marker expression and functionality.

## Quality Control Testing

- Each vial contains  $\geq 1 \times 10^6$  viable cells.
- Cells are tested negative for infectious diseases by a Mouse/Rat Comprehensive CLEAR panel by Charles River Animal Diagnostic Services.
- Cells are verified to be of rat origin and negative for inter-species contamination from mouse, chinese hamster, Golden Syrian hamster, human and non-human primate (NHP) as assessed by a Contamination CLEAR panel by Charles River Animal Diagnostic Services.
- Cells are negative for mycoplasma contamination.

## Representative Data



**Figure 1.** INS-1 832/3 cells one day after thawing in a T75 flask (A, 10X magnification), (B, 20X magnification). Cells express actin (Phalloidin, C) and Synaptotagmin (D, red).

## References

1. *Endocrinology* 1992; 130(1): 167-178.
2. *Diabetes* 2000; 49(3): 424-30.
3. *EMBO J* 1997; 16(19):5837-5846.
4. *J Biol Chem* 2008; 283(43): 28909-28917.

Please visit [www.millipore.com](http://www.millipore.com) for additional product information and references.

Submit your published journal article and credit toward future purchases. Visit [www.millipore.com/publicationrewards](http://www.millipore.com/publicationrewards) to learn more!

## Protocols

**Note:** Extensive passaging may cause cells to lose the expression of the human insulin gene and the population may become heterogeneous for expression of human insulin. However, the loss of human insulin expression does not affect cell function measured as GSIS. G418 at 0.3 mg/mL may be added to apply selective pressure when growing larger batches of cells.

### Thawing Cells

1. Do not thaw the cells until the recommended medium is on hand. Cells can grow on normal tissue cultureware surfaces without any additional coating.

**INS-1 832/3 Expansion Medium:** Cells are thawed and expanded in RPMI-1640 (Sigma Cat. No. R0883) supplemented with 2 mM L-Glutamine (Cat. No. TMS-002-C), 1 mM sodium pyruvate (Cat. No. TMS-005-B), 10 mM HEPES (Cat. No. TMS-003-C), 0.05 mM  $\beta$ -mercaptoethanol (Cat. No. ES-007-E) and 10% FBS (Cat. No. ES-009-B).

**Note:**  $\beta$ -mercaptoethanol is critical for the continued propagation of the cell line and should not be omitted from the culture medium.

2. Remove the vial of frozen INS-1 832/3 cells from liquid nitrogen and incubate in a 37°C water bath. Closely monitor until the cells are completely thawed. Maximum cell viability is dependent on the rapid and complete thawing of frozen cells.

**IMPORTANT: Do not vortex the cells.**

3. As soon as the cells are completely thawed, disinfect the outside of the vial with 70% ethanol. Proceed immediately to the next step.
4. In a laminar flow hood, use a 1 or 2 mL pipette to transfer the cells to a sterile 15 mL conical tube. Be careful not to introduce any bubbles during the transfer process.
5. Using a 10 mL pipette, slowly add dropwise 9 mL of INS-1 832/3 Expansion Medium (Step 1 above) to the 15 mL conical tube.

**IMPORTANT: Do not add the entire volume of media all at once to the cells. This may result in decreased cell viability due to osmotic shock.**

6. Gently mix the cell suspension by slowly pipetting up and down twice. Be careful not to introduce any bubbles.

**IMPORTANT: Do not vortex the cells.**

7. Centrifuge the tube at 300 x g for 2-3 minutes to pellet the cells.
8. Decant as much of the supernatant as possible. Steps 5-8 are necessary to remove residual cryopreservative (DMSO).
9. Resuspend the cells in 15 mL of INS-1 832/3 Expansion Medium.
10. Transfer the cell mixture to a T75 tissue culture flask.
11. Incubate the cells at 37°C in a humidified incubator with 5% CO<sub>2</sub>.

### Subculturing Cells

1. Carefully remove the medium from the T75 tissue culture flask containing the confluent layer of INS-1 832/3 cells.
2. Rinse the flask with 10 mL 1X PBS. Aspirate after the rinse.
3. Apply 5-7 mL of Accutase or trypsin-EDTA solution and incubate in a 37°C incubator for 3-5 minutes.
4. Inspect the flask and ensure the complete detachment of cells by gently tapping the side of the flask with the palm of your hand.
5. Add 5-7 mL of INS-1 832/3 Expansion Medium to the plate.
6. Gently rotate the flask to mix the cell suspension. Transfer the dissociated cells to a 15 mL conical tube.
7. Centrifuge the tube at 300 x g for 3-5 minutes to pellet the cells.
8. Discard the supernatant, then loosen the cell pellet by tapping the tip of the tube with a finger.
9. Apply 2-5 mL of INS-1 832/3 Expansion Medium to the conical tube and resuspend the cells thoroughly.

**IMPORTANT: Do not vortex the cells.**

10. Count the number of cells using a hemocytometer.
11. Plate the cells to the desired density. Typical split ratio is 1:6.

### Cryopreservation of Cells

INS-1 832/3 rat insulinoma cell line may be frozen in the expansion medium plus 10% DMSO using a Nalgene slow freeze Mr. Frosty container.

## INS-1 832/3 Rat Insulinoma Cell Line

Cat # SCC208

### **Glucose stimulated insulin secretion (GSIS):**

For each condition, triplicate wells should be plated.

1. Cells are seeded at a density of  $0.5 \times 10^6$ /well in 1 mL medium for a 24-well plate, or  $1 \times 10^6$ /well in 2 mL of medium for a 12-well plate. For best results in glucose stimulated insulin secretion, cells need to be confluent.
2. After 2 days, the medium should be changed, and the assay should be performed on day 3.
3. The glucose stimulated insulin secretion is performed in HBSS (Hepes balanced salt solution): 114 mmol/L NaCl, 4.7 mmol/L KCl, 1.2 mmol/L  $\text{KH}_2\text{PO}_4$ , 1.16 mmol/L  $\text{MgSO}_4$ , 20 mmol/L HEPES, 2.5 mmol/L  $\text{CaCl}_2$ , 25.5 mmol/L  $\text{NaHCO}_3$ , and 0.2% bovine serum albumin, pH 7.2. (Hohmeier et al. 2000, Diabetes 49:424-430).
4. For the assay, wash the cells twice with HBSS + 2.5 mM glucose. The first wash is just a quick rinse, for the second wash, leave the HBSS on for 1 hr.

Caution: The cells do not attach firmly to the plates and will wash off very easy, if solutions are added with too much force.

5. After 1 hrs, the secretagogues diluted in HBSS are added for 2 hr. For a 24-well plate, add 1mL/well and 1.5 or 2 mL/well to a 12-well plate.
6. After 2 hrs, remove the solution for insulin ELISA. Keep in mind that these cells secrete a mixture of rat and human insulin and an insulin ELISA which is cross reactive between these two species will give you higher sensitivity – especially if you want to adapt the assay for a 96-well format. We recommend the Chemi-rodent-insulin ELISA from ALPCO Cat. No. 80-INSMR-CH01.
7. Wash the cells twice with PBS and lyse cells with RIPA buffer to measure total cellular protein for normalization of secreted insulin to protein.
8. Conditions used for quality control:
  1. HBSS+2.5 mM glucose (basal insulin secretion)
  2. HBSS+15 mM glucose (stimulated insulin secretion)
  3. HBSS+15 mM glucose +50 nM GLP-1

■ antibodies ■ Multiplex products ■ biotools ■ cell culture ■ enzymes ■ kits ■ proteins/peptides ■ siRNA/cDNA products

**Please visit [www.millipore.com](http://www.millipore.com) for additional product information, test data and references**

EMD Millipore Corporation, 28820 Single Oak Drive, Temecula, CA 92590, USA 1-800-437-7500

Technical Support: T: 1-800-MILLIPORE (1-800-645-5476) • F: 1-800-437-7502

**FOR RESEARCH USE ONLY.** Not for use in diagnostic procedures. Not for human or animal consumption. Purchase of this Product does not include any right to resell or transfer, either as a stand-alone product or as a component of another product. Any use of this Product for purposes other than research is strictly prohibited.

EMD Millipore®, the M mark, Upstate®, Chemicon®, Linco® and all other registered trademarks, unless specifically identified above in the text as belonging to a third party, are owned by Merck KGaA, Darmstadt, Germany. Copyright ©2008-2018 Merck KGaA, Darmstadt, Germany. All rights reserved.



We Buy 100% Certified Renewable Energy

**ACADEMIC USE AGREEMENT**  
(subject to local law)

**THIS PRODUCT MAY ONLY BE USED BY INDIVIDUALS EMPLOYED BY AN ACADEMIC INSTITUTION AND IS INTENDED SOLELY TO BE USED FOR ACADEMIC RESEARCH, WHICH IS FURTHER DEFINED BELOW. BY OPENING THIS PRODUCT, YOU (“PURCHASER”) HEREBY REPRESENT THAT YOU HAVE THE RIGHT AND AUTHORITY TO LEGALLY BIND YOURSELF AND/OR YOUR EMPLOYER INSTITUTION, AS APPLICABLE, AND CONSENT TO BE LEGALLY BOUND BY THE TERMS OF THIS ACADEMIC USE AGREEMENT. IF YOU DO NOT AGREE TO COMPLY WITH THESE TERMS, YOU MAY NOT OPEN OR USE THE PRODUCT AND YOU MUST CALL MILLIPORESIGMA (“SELLER”) CUSTOMER SERVICE (1-800-645-5476) TO ARRANGE TO RETURN THE PRODUCT FOR A REFUND.**

“Product” means INS-1 832/3 Rat Insulinoma Cell Line (SCC208)

“Academic Research” means any internal *in vitro* research use by individuals employed by an academic institution. Academic Research specifically excludes the following uses of whatever kind or nature:

- Re-engineering or copying the Product
- Making derivatives, modifications, or functional equivalents of the Product
- Obtaining patents or other intellectual property rights claiming use of the Product
- Using the Product in the development, testing, or manufacture of a Commercial Product
- Using the Product as a component of a Commercial Product
- Reselling or licensing the Product
- Using the Product in clinical or therapeutic applications including producing materials for clinical trials
- Administering the Product to humans
- Using the Product in collaboration with a commercial or non-academic entity

“Commercial Product” means any product intended for: (i) current or future sale; (ii) use in a fee-for-service; or (iii) any diagnostic, clinical, or therapeutic use.

Access to the Product is limited solely to those officers, employees, and students of PURCHASER’s academic institution who need access to the Product to perform Academic Research. PURCHASER shall comply with all applicable laws in its use and handling of the Product and shall keep it under reasonably safe and secure conditions to prevent unauthorized use or access.

These use restrictions will remain in effect for as long as PURCHASER possesses the Product.

**COMMERCIAL OR NON-ACADEMIC ENTITIES INTERESTED IN PURCHASING OR USING THE PRODUCT MUST CONTACT [licensing@emdmillipore.com](mailto:licensing@emdmillipore.com) AND AGREE TO SEPARATE TERMS OF USE PRIOR TO USE OR PURCHASE.**

**GMO**

This product contains genetically modified organisms.  
Este producto contiene organismos genéticamente modificados.  
Questo prodotto contiene degli organismi geneticamente modificati.  
Dieses Produkt enthält genetisch modifizierte Organismen.  
Ce produit contient des organismes génétiquement modifiés.  
Dit product bevat genetisch gewijzigde organismen.  
Tämä tuote sisältää geneettisesti muutettuja organismeja.  
Denna produkt innehåller genetiskt ändrade organismer.

■ antibodies ■ Multiplex products ■ biotools ■ cell culture ■ enzymes ■ kits ■ proteins/peptides ■ siRNA/cDNA products

Please visit [www.millipore.com](http://www.millipore.com) for additional product information, test data and references

EMD Millipore Corporation, 28820 Single Oak Drive, Temecula, CA 92590, USA 1-800-437-7500

Technical Support: T: 1-800-MILLIPORE (1-800-645-5476) • F: 1-800-437-7502

**FOR RESEARCH USE ONLY.** Not for use in diagnostic procedures. Not for human or animal consumption. Purchase of this Product does not include any right to resell or transfer, either as a stand-alone product or as a component of another product. Any use of this Product for purposes other than research is strictly prohibited.

EMD Millipore®, the M mark, Upstate®, Chemicon®, Linco® and all other registered trademarks, unless specifically identified above in the text as belonging to a third party, are owned by Merck KGaA, Darmstadt, Germany. Copyright ©2008-2018 Merck KGaA, Darmstadt, Germany. All rights reserved.



We Buy 100% Certified  
Renewable Energy