

Data Sheet

# MOC1 Mouse Oral Squamous Cell Carcinoma (OSCC)

Cancer Cell Line

**SCC469****Pack Size:  $\geq 1 \times 10^6$  viable cells/vial****Store in liquid nitrogen.****FOR RESEARCH USE ONLY****Not for use in diagnostic procedures. Not for human or animal consumption.**

## Background

Oral squamous cell carcinoma (OSCC) comprises 95% of all forms of head & neck cancer. The five-year survival rate for OSCC patients is approximately 50%; this poor patient outcome has not changed significantly in the last few decades. More recently, it has been suggested that immune cells may play a role in tumor-host interactions in tumorigenicity. The lack of available syngeneic mouse models has, however, hindered progress in studying this role. Xenograft models exist, but require the use of immune-compromised mice, and the necessary components of adaptive immunity are therefore missing. The MOC (Mouse Oral Cancer) cell lines provide a syngeneic head and neck squamous cell carcinoma (HNSCC) model that can be transplanted into immune-competent C57BL/6 mice.<sup>1</sup> These cell lines are important models for facilitating studies on immune cell infiltration of tumors.

The MOC cell lines (MOC1, MOC2, MOC22) exhibit diverse *in vivo* tumor growth phenotypes. MOC1 and 22 demonstrate indolent *in vivo* growth phenotypes. MOC2 has a more aggressive growth phenotype – an injection of as few as 10,000 cells sufficient to form tumors in mouse models. The more aggressive MOC2 line generates higher CD4+ T cell levels. CD44, which has been implicated as a cancer stem cell marker that characterizes tumors with higher resistance to therapeutic treatment, is also expressed at a higher level in MOC2 cells as compared with the MOC1 cell line. The MOC1 cell line, which has a less aggressive growth pattern, generates a tumor microenvironment with higher CD8+ T cell levels compared to MOC2. This mirrors data that associates a higher CD8+ T-cell count with more favorable outcomes in HNSCC patients.<sup>1</sup>

The MOC cell lines have been impactful in elucidating CXCL14 downregulation in metastatic lymph nodes. Previously, CXCL14 expression had been associated with T cell-associated mechanisms for tumor suppression. Tumor Infiltrating Lymphocyte (TIL) studies using MOC1 and MOC2 have found a cell-specific association in Oral Squamous Cell Carcinomas (OSCCs) involving CXCL14 downregulation resulting in TIL increase which leads to tumor suppression.<sup>2</sup> The MOC lines have been utilized to compare and identify potential clinical scenarios that are applicable to human OSCC. The MOC group of cell lines contained common HNSCC mutations in addition to human related OSCC driver pathways.

## Source

Cell line was derived from a carcinogen-induced oral tumor in an immunocompetent C57BL/6 mouse.<sup>1</sup>

## Short Tandem Repeat

M18-3: 16	M4-2: 20.3	M6-7: 15	M19-2: 13	M1-2: 19	M7-1: 25.2	M1-1: 16
M3-2: 14	M8-1: 16	M2-1: 16	M15-3: 22.3	M6-4: 18	M11-2: 16	M17-2: 15
M12-1: 17	M5-5: 18	MX-1: 27	M13-1: 17			


## Quality Control Testing

- MOC1 cells are verified to be of mouse origin and negative for human, rat, Chinese hamster, Golden Syrian hamster, and nonhuman primate interspecies contamination, as assessed by a Contamination Clear panel by Charles River Animal Diagnostic Services.
- Cells tested negative for infectious diseases against a Mouse Essential CLEAR panel by Charles River Animal Diagnostic Services.
- Cells tested negative for mycoplasma.

## Storage and Handling

MOC1 cells should be stored in liquid nitrogen until use. The cells can be cultured for at least 10 passages after initial thawing without significantly affecting the cell marker expression and functionality.

## Representative Data



**Figure 1.** **A.** Bright-field image of MOC1 cells one day after thaw in a T75 flask. **B.** (MDCI10000) Cell confluency was assessed throughout the culture using the Millicell® Digital Cell Imager. **C.** (CPHCC360KIT) Cell counting was performed using the Scepter™ 3.0 handheld automated cell counter using 60  $\mu\text{m}$  senss. **D.** (49409) Cells express actin, **E.** (SAB5700076) CD44, **F.** (Fisher Scientific PA5-38014) CD133 and **G.** (Fisher Scientific PA5-32127) ALDH1A1.

**NOTE:** Product catalog numbers indicated in ( ) can be purchased at [SigmaAldrich.com](https://www.sigmaaldrich.com) unless otherwise stated.

## Protocols

### Thawing the Cells

MOC cell lines proliferate very rapidly. Upon thaw, it will take approximately 2-4 days for the cells to reach confluence. A passage ratio of 1:12 from a 80% confluent T175 flask) will take 2-4 days to reach 80% confluency.

1. Do not thaw the cells until the recommended medium is on hand. Cells can grow on standard tissue cultureware surfaces without any additional coating. Cells are thawed and expanded in MOC Expansion Medium comprising DMEM/F12 (DF-042-B) containing 10% FBS (ES-009-B), and 2 mM L-Glutamine (TMS-002-C).
2. Remove the vial of frozen MOC1 cells from liquid nitrogen and incubate in a 37 °C water bath. Closely monitor until the cells are completely thawed. Maximum cell viability is dependent on the rapid and complete thawing of frozen cells.

**IMPORTANT:** Do not vortex the cells.

3. As soon as the cells are completely thawed, disinfect the outside of the vial with 70% ethanol. Proceed immediately to the next step.
4. In a laminar flow hood, use a 1- or 2-mL pipette to transfer the cells to a sterile 15 mL conical tube. Be careful not to introduce any bubbles during the transfer process.
5. Using a 10 mL pipette, slowly add dropwise 9 mL of MOC Expansion Medium (Step 1 above) to the 15 mL conical tube.  
**IMPORTANT:** Do not add the entire volume of media all at once to the cells. This may result in decreased cell viability due to osmotic shock.
6. Gently mix the cell suspension by slowly pipetting up and down twice. Be careful not to introduce any bubbles.  
**IMPORTANT:** Do not vortex the cells.
7. Centrifuge the tube at 300 x g for 2-3 minutes to pellet the cells.
8. Decant as much of the supernatant as possible. Steps 5-8 are necessary to remove residual cryopreservative (DMSO).
9. Resuspend the cells in 25–30 mL of MOC Expansion Medium.
10. Transfer the cell mixture to a T175 tissue culture flask.
11. Incubate the cells at 37 °C in a humidified incubator with 5% CO<sub>2</sub>.

### Subculturing the Cells

MOC1 cells are highly adherent and will require higher Accutase® volumes with longer incubation times.

1. Do not allow the cells to grow to confluency. MOC1 cells should be passaged at ~80-85% confluency.
2. Carefully remove the medium from the T175 tissue culture flask containing the 80% confluent layer of MOC1 cells.
3. Rinse the flask with 15 mL 1X PBS. Aspirate after the rinse.
4. Apply 10 mL of Accutase® and incubate in a 37 °C incubator for 10-15 minutes.
5. Inspect the flask and ensure the complete detachment of cells by gently tapping the side of the flask with the palm of your hand.

**NOTE:** This cell line requires longer incubations and higher Accutase® volumes compared to most cancer cell lines.

6. Add 10 mL of MOC Expansion Medium to the plate.
7. Gently rotate the flask to mix the cell suspension. Transfer the dissociated cells to a 50 mL conical tube.
8. Centrifuge the tube at 300 x g for 3-5 minutes to pellet the cells.
9. Discard the supernatant, then loosen the cell pellet by tapping the tip of the tube with a finger.
10. Apply 2-5 mL of MOC Expansion Medium to the conical tube and resuspend the cells thoroughly. Large cell clumps may be broken up by gentle trituration.

**IMPORTANT:** Do not vortex the cells.

11. Count the number of cells using a hemocytometer or a Scepter™ 3.0 handheld automated cell counter.
12. Plate the cells to the desired density. Typical split ratio is 1:10–1:12.

## Cryopreservation of the Cells

MOC1 cells may be frozen in MOC Expansion Medium supplemented with 10% DMSO using a Nalgene® slow freeze Mr. Frosty™ container.

## References

1. Judd NP, Allen CT, Winkler AE, Uppaluri R. 2012. Comparative analysis of tumor-infiltrating lymphocytes in a syngeneic mouse model of oral cancer. *Otolaryngol Head Neck Surg.* 147(3): 493–500.
2. Parikh A, Shin J, Faquin W, Lin DT, Tirosh I, Sunwoo JB, Puram SV. 2020. Malignant cell-specific CXCL14 promotes tumor lymphocyte infiltration in oral cavity squamous cell carcinoma. *J Immunother Cancer.* 8(2): e001048.
3. Onken MD, Winkler AE, Kanchi K-L, Chalivendra V, Law JH, Rickert CG, Kallogjeri D, Judd NP, Dunn GP, Piccirillo JF, et al. 2014. A surprising cross-species conservation in the genomic landscape of mouse and human oral cancer identifies a transcriptional signature predicting metastatic disease. *Clin Cancer Res.* 20(11): 2873-2884.

## Academic Use Agreement

Subject to local law

THIS PRODUCT MAY ONLY BE USED BY INDIVIDUALS EMPLOYED BY AN ACADEMIC INSTITUTION AND IS INTENDED SOLELY TO BE USED FOR ACADEMIC RESEARCH, WHICH IS FURTHER DEFINED BELOW. BY OPENING THIS PRODUCT, YOU ("PURCHASER") HEREBY REPRESENT THAT YOU HAVE THE RIGHT AND AUTHORITY TO LEGALLY BIND YOURSELF AND/OR YOUR EMPLOYER INSTITUTION, AS APPLICABLE, AND CONSENT TO BE LEGALLY BOUND BY THE TERMS OF THIS ACADEMIC USE AGREEMENT. IF YOU DO NOT AGREE TO COMPLY WITH THESE TERMS, YOU MAY NOT OPEN OR USE THE PRODUCT AND YOU MUST CALL CUSTOMER SERVICE (1-800-645-5476) TO ARRANGE TO RETURN THE PRODUCT FOR A REFUND.

"Product" means MOC1 Mouse Oral Squamous Cell Carcinoma (OSCC) Cell Line (SCC469).

"Academic Research" means any internal in vitro research use by individuals employed by an academic institution. Academic Research specifically excludes the following uses of whatever kind or nature:

- Re-engineering or copying the Product
- Making derivatives, modifications, or functional equivalents of the Product
- Obtaining patents or other intellectual property rights claiming use of the Product
- Using the Product in the development, testing, or manufacture of a Commercial Product
- Using the Product as a component of a Commercial Product
- Reselling or licensing the Product
- Using the Product in clinical or therapeutic applications including producing materials for clinical trials
- Administering the Product to humans
- Using the Product in collaboration with a commercial or non-academic entity

"Commercial Product" means any product intended for: (i) current or future sale; (ii) use in a fee-for-service; or (iii) any diagnostic, clinical, or therapeutic use.

Access to the Product is limited solely to those officers, employees, and students of PURCHASER's academic institution who need access to the Product to perform Academic Research. PURCHASER shall comply with all applicable laws in its use and handling of the Product and shall keep it under reasonably safe and secure conditions to prevent unauthorized use or access.

These use restrictions will remain in effect for as long as PURCHASER possesses the Product.

COMMERCIAL OR NON-ACADEMIC ENTITIES INTERESTED IN PURCHASING OR USING THE PRODUCT MUST CONTACT [licensing@emdmillipore.com](mailto:licensing@emdmillipore.com) AND AGREE TO SEPARATE TERMS OF USE PRIOR TO USE OR PURCHASE.

---

## Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

## Technical Assistance

Visit the tech service page at [SigmaAldrich.com/techservice](https://SigmaAldrich.com/techservice).

## Terms and Conditions of Sale

Warranty, use restrictions, and other conditions of sale may be found at [SigmaAldrich.com/terms](https://SigmaAldrich.com/terms).

## Contact Information

For the location of the office nearest you, go to [SigmaAldrich.com/offices](https://SigmaAldrich.com/offices).

The life science business of Merck KGaA, Darmstadt, Germany  
operates as MilliporeSigma in the U.S. and Canada.

Merck Millicell, Scepter and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2023 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

Document Template 20306518 Ver 6.0

00149644 Ver 1.0, Rev 17NOV2023, RC, AB

