



Product Information

Concanavalin A from *Canavalia ensiformis* (Jack bean)

Product Number **C 7275**
Storage Temperature -0 °C

Product Description

CAS Number: 11028-71-0
Molecular Weight: The monomeric molecular weight is 25.5 kDa. At pH 5.5 Con A exists as a dimer and at pH>7 it exists as a tetramer.¹
Extinction Coefficient: $E^{1\%} = 11.4$ (280 nm, 100 mM NaCl)²
pI: 4.5, 4.7, 5.1, 5.5³
Synonyms: Con A, Concanavalin A from Jack bean

Concanavalin A isolated from Jack bean is a lectin, although it is not a glycoprotein like most other phytohemagglutinins.⁴ The protein does not contain any cysteine residues.⁵ Con A is a metalloprotein and requires a transition metal ions such as manganese and calcium for saccharide binding to occur. Each subunit has a binding site for one Mn^{2+} , one Ca^{2+} , and one oligo- or polysaccharide. Removal of these cations by dialysis under acidic conditions abolishes the carbohydrate binding activity. This protein will bind to the glucosyl and mannosyl residues at the nonreducing terminus of polysaccharides and glycoproteins. Unmodified hydroxyl groups at the C3, C4, and C6 positions of D-glucopyranosyl or D-mannopyranosyl rings are considered to be essential for binding. α -glucopyranosides and α -mannopyranosides bind more strongly than the β -forms.⁶

Concanavallin A is a stimulator of lymphocyte blast formation⁷ and possesses mitogenic activity.^{8,9} It also agglutinates red blood cells and complexes with blood group sugars.¹⁰ Concanavalin A, when immobilized to a chromatography resin, may also be used to purify glycoproteins. The bound glycoprotein may be eluted with 500 mM α -methyl-D-mannoside, which has a high affinity for Concanavalin A.¹¹

Precautions and Disclaimer

For Laboratory Use Only. Not for drug, household or other uses.

Preparation Instructions

This product is soluble in water (10 mg/ml), yielding a slightly hazy, faint yellow solution.

Storage/Stability

Solutions of Concanavalin A in 1 M NaCl are stable when stored at -20 °C.

References

1. Senear, D.F., and Teller, D.C., Thermodynamics of concanavalin A dimer-tetramer self-association: sedimentation equilibrium studies. *Biochemistry*, **20**, 3076-3083 (1981).
2. Agrawal, B.B., and Goldstein, I.J., Protein-carbohydrate interaction. VII. Physical and chemical studies on concanavalin A, the hemagglutinin of the jack bean. *Arch. Biochem. Biophys.*, **124**, 218-229 (1968).
3. Entlicher, G., et al., Studies on phytohemagglutinins. VIII. Isoelectric point and multiplicity of purified concanavalin A. *Biochim. Biophys. Acta*, **236**, 795-797 (1971).
4. Olson, M.O., and Liener, I.E., Some physical and chemical properties of concanavalin A, the phytohemagglutinin of the jack bean., *Biochemistry*, **6**, 105-111 (1967).
5. Cunningham, B.A., et al., The covalent and three dimensional structure of concanavalin A. II. Amino acid sequence of cyanogen bromide fragment F3. *J. Biol. Chem.*, **250**, 1503-1512 (1975).
6. Goldstein, I.J., et al., Protein carbohydrate interaction. II. Inhibition studies on the interaction of concanavalin A with polysaccharides. *Biochemistry*, **4**, 876-83 (1965).

7. Perlman, P., et al., Inhibition of cytotoxicity of lymphocytes by concanavalin A in vitro. *Science*, **168**, 1112-1115 (1970).
8. Kitao, T., and Hattori, K., Concanavalin A as a carrier of daunomycin. *Nature*, **265**, 81-82 (1977).
9. Douglas, S.D., et al., Human lymphocyte response to phytoimitogens in vitro: normal, agammaglobulinemic and paraproteinemic individuals. *J. Immunol.*, **103**, 1185-1195 (1969).
10. Clark, A.E., and Denborough, M., The interaction of concanavalin A with blood-group-substance glycoproteins from human secretions. *Biochem. J.*, **121**, 811-816 (1971).
11. Campbell, K.P., and MacLennan, D.H., Purification and characterization of the 53,000 dalton glycoproteins from the sarcoplasmic reticulum. *J. Biol. Chem.*, **256**, 4626-4632 (1981).

TMG/AJH 11/02

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.