Supelco®

www.sigmaaldrich.com

1.14729.0001

1.14729.0007

Spectroquant®

Phosphate Cell Test

P

for the determination of orthophosphate and total phosphorus
USEPA equivalent for drinking water and wastewater

1. Method

In sulfuric solution orthophosphate ions react with molybdate ions to form molybdophosphoric acid. Ascorbic acid reduces this to phosphomolybdenum blue (PMB) that is determined photometrically.

The method is analogous to EPA 365.2+3, APHA 4500-P E, and DIN EN ISO 6878.

2. Measuring range and number of determinations

Measuring range	Number of determinations
0.5 - 25.0 mg/l PO ₄ -P	25
1.5 - 76.7 mg/l PO ₄ ³⁻	
1.1 - 57.3 mg/l P ₂ O ₅	

For programming data for selected photometers / spectrophotometers see www.sigmaaldrich.com/photometry.

3. Applications

This test measures only orthophosphate. Samples must be decomposed by digestion before total phosphorus can be measured (see section 6).

Sample material:

Groundwater and surface water, seawater
Drinking water
Wastewater
Nutrient solutions for fertilization
Soils after appropriate sample pretreatment
Food after appropriate sample pretreatment

4. Influence of foreign substances

This was checked individually in solutions containing 15 and 0 mg/l PO₄-P. The determination is not yet interfered with up to the concentrations of foreign substances given in the table. Cumulative effects were not checked; such effects can, however, not be excluded.

Concentrations of foreign substances in mg/l or %					
Ag ⁺	1000	F ⁻	250	Pb ²⁺	100
AsO ₄ ³⁻	1.0	Fe ³⁺	1000	S ²⁻	10
Ca ²⁺	1000	Hg ²⁺	50	SiO ₃ ²⁻	1000
Cd ²⁺	1000	Mg ²⁺	1000	SO ₃ ²⁻	1000
CN ⁻	1000	Mn ²⁺	1000	Zn ²⁺	1000
Cr ³⁺	1000	NH ₄ ⁺	1000	Na-acetate	5 %
Cr ₂ O ₇ ²⁻	25	Ni ²⁺	1000	NaCl	20 %
Cu ²⁺	1000	NO ₂ ⁻	1000	NaNO ₃	10 %
				Na ₂ SO ₄	10 %

Reducing agents interfere with the determination.

¹⁾ tested with nonionic, cationic, and anionic surfactants

²⁾ **A higher COD may impair the efficacy of the digesting mixture in the determination of total phosphorus and thus result in false-low readings. Up to a maximum of 1500 mg/l COD, this can be avoided by adding 2 doses of reagent P-1K instead of 1.**

5. Reagents and auxiliaries

Please note the warnings on the packaging materials!

The test reagents are stable up to the date stated on the pack when stored closed at +15 to +25 °C.

Package contents:

1 bottle of reagent P-1K
1 bottle of reagent P-2K
1 bottle of reagent P-3K
25 reaction cells
1 green dose-metering cap
1 blue dose-metering cap
1 sheet of round stickers for numbering the cells

Other reagents and accessories:

MQuant® Phosphate Test, Cat. No. 1.10428,
measuring range 10 - 500 mg/l PO₄³⁻ (3.3 - 163 mg/l PO₄-P)
MQuant® Universal indicator strips pH 0 - 14, Cat. No. 1.09535
Sulfuric acid 0.5 mol/l Titripur®, Cat. No. 1.09072
Spectroquant® CombiCheck 20, Cat. No. 1.14675
Spectroquant® CombiCheck 80, Cat. No. 1.14738
Phosphorus (total) standard solution CRM, 4.00 mg/l PO₄-P,
Cat. No. 1.25047
Phosphorus (total) standard solution CRM, 15.0 mg/l PO₄-P,
Cat. No. 1.25048
Hydrochloric acid 25 % for analysis EMSURE®, Cat. No. 1.00316

Pipette for a pipetting volume of 1.0 ml
Thermoreactor

6. Preparation

- Use only phosphate-free detergents to rinse glassware. Otherwise fill with hydrochloric acid (approx. 10 %) and leave to stand for several hours.

At the first use **replace the screw caps of the reagent bottles P-1K and P-3K by the corresponding dose-metering caps:**

Reagent P-1K: green dose-metering cap
Reagent P-3K: blue dose-metering cap

Hold the respective reagent bottle **vertically** and, at each dosage, press the slide **all the way** into the dose-metering cap. **Before each dosage** ensure that the slide is **completely retracted**.

Reclose the reagent bottles with the corresponding screw caps at the end of the measurement series, since the function of the reagents is impaired by the absorption of atmospheric moisture.

- Analyze immediately after sampling.

- Digestion for the determination of total phosphorus (Wear eye protection!):**

Pretreated sample	1.0 ml	Pipette into a reaction cell.
Reagent P-1K	1 dose ¹⁾	Add, close the cell tightly , and mix.

Heat the cell at 120 °C in the preheated thermoreactor for 30 min. Allow the closed cell to cool to room temperature in a test-tube rack.

Do not cool with cold water!

¹⁾ in the case of high COD values: 2 doses

- Check the phosphate content with the MQuant® PhosphateTest. Samples containing more than 25.0 mg/l PO₄-P must be diluted with distilled water **prior to** digestion.
- The pH must be within the range 0 - 10.** Adjust, if necessary, with sulfuric acid.
- Filter turbid samples.

7. Procedure

Pretreated sample (10 - 35 °C)	1.0 ml	Pipette into a reaction cell and mix or - after digestion for total phosphorus - shake the tightly closed cell vigorously after cooling.
Reagent P-2K ¹⁾	5 drops ²⁾	Add, close the cell tightly, and mix.
Reagent P-3K ¹⁾	1 dose	Add, close the cell tightly, and shake vigorously until the reagent is completely dissolved .

Leave to stand for 5 min (reaction time), then measure the sample in the photometer.

¹⁾ In the case of high chloride contents, it is recommended to switch the sequence of the reagents P-2K and P-3K.

²⁾ **Hold the bottle vertically while adding the reagent!**

Notes on the measurement:

- For photometric measurement the cells must be clean. Wipe, if necessary, with a clean dry cloth.
- Measurement of turbid solutions yields false-high readings.
- The pH of the measurement solution must be within the range 0.80 - 0.95.
- The color of the measurement solution remains stable for at least 60 min after the end of the reaction time stated above.

8. Analytical quality assurance

recommended before each measurement series
To check the photometric measurement system (test reagents, measurement device, handling) and the mode of working the phosphorus (total) standard solutions CRM, (see section 5) - for the determination of total phosphorus - or Spectroquant® CombiCheck 20 and 80 can be used. Besides a **standard solution** with 8.0 mg/l PO₄-P (CombiCheck 20) or, respectively, 15.0 mg/l PO₄-P (CombiCheck 80) these articles also contain an **addition solution** for determining sample-dependent interferences (**matrix effects**).

Additional notes see under **www.qa-test-kits.com**.

For quality and batch certificates for Spectroquant® test kits see the website, where you will find all data in production control, that are determined in accordance with ISO 8466-1 and DIN 38402 A51.

9. Notes

- Reclose the reagent bottles immediately after use.
- Information on disposal can be obtained at www.disposal-test-kits.com.**

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© 2021 Merck KGaA, Darmstadt, Deutschland und/oder Tochterunternehmen. Alle Rechte vorbehalten. Das Bunte M ist eine Marke der Merck KGaA, Darmstadt, Deutschland. Alle anderen Marken sind Eigentum der jeweiligen Inhaber. Ausführliche Informationen zu Markennamen sind über öffentlich zugängliche Informationsquellen erhältlich.