



Product Information

D₁ DOPAMINE RECEPTOR, HUMAN Expressed in SF9 Cells

Product Number **D-178**

Product Description

Dopamine receptors were initially divided into two general categories on the basis of differences in receptor pharmacology and biochemical mechanisms of signal transduction. With the application of the techniques of molecular biology, two predominant dopamine receptors, D₁ and D₂, were cloned. Later other dopamine receptors with homology to either the D₁ or D₂ receptor were identified. Thus, at present, two families of vertebrate dopamine receptors (designated as D₁-like and D₂-like) are recognized. The D₁-like family consists of the D₁ and D₅ receptors while the D₂-like family consists of the D₂, D₃ and D₄ receptors.

The D₁ and D₂ receptors occur in sufficiently high concentrations that they can be studied *in situ*. The D₃, D₄ and D₅ receptors occur in such low concentrations that study of them *in situ* is difficult. Thus, the majority of study of these receptors has been accomplished using cell lines cloned to express these receptors.

Reagents

D₁ Dopamine Receptor, Human (Sf9) is provided as purified membranes suspended in 50 mM Tris-HCl at pH 7.4, containing 10% glycerol and 1% bovine serum albumin (BSA).

Storage/Stability

Store tightly sealed at -80 °C. When stored in its original packaging solution, the membranes retain their original specific activity for several months.

Procedure

Incubation Buffer

50 mM Tris-HCl at pH 7.4, containing 5.0 mM MgCl₂, 5.0 mM EDTA, 5.0 mM KCl and 1.5 mM CaCl₂

Typical Binding Protocol

1. Membranes

Dilute in incubation buffer (0.5 ml of membranes to 24.5 ml of incubation buffer).

2. Assay mixture

500 µl of diluted membranes

20 µl of [³H]-radioligand in buffer

20 µl of incubation buffer or unlabeled ligand in buffer

Radioligand: [³H]-SCH-23390 at a final concentration of 1.6 nM for competition studies.

Unlabeled ligand: (+)-Butaclamol (Cat. No. D-033) at a final concentration of 10 µM.

3. Incubation time

90 minutes at 27 °C.

4. Separation

Over GF/C filter (5 mm diameter, presoaked in 0.3% polyethylenimine, Cat. No. P-182) then wash 4x with 500 µl of ice cold 50 mM Tris-HCl at pH 7.4.

Results

Typical affinities using standard binding assay above. Results may vary from lot to lot.

Ligand	Affinity (K _i) (nM)
[³ H]-SCH-23390	1.3 (K _d)
(+)-Butaclamol (D-033)	1.6
SKF 83566 (S-110)	3.6
Haloperidol (H-100)	38
Spiperone (D-050)	261

References

1. Jarvie, K.R., et al., "Molecular cloning, stable expression and desensitization of the human dopamine D_{1b}/D₅ receptor." *J. Recept. Res.* **13**, 573-590 (1993).
2. Sunahara, R.K., et al., "Cloning of the gene for a human D₅ receptor with higher affinity for dopamine than D₁." *Nature* **350**, 614-619 (1991).
3. MacKenzie, R.G., et al., "A D₁/D₂ chimeric dopamine receptor mediates a D₁ response to a D₂-selective agonist." *FEBS Lett.* **323**, 59-62 (1993).

SMS 10/00

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.