Millipore®

SigmaAldrich.com

CE

0459

Rx only

English	Português	Lietuviškai
Français	Ελληνικά	Magyar
Italiano	Dansk	Český
Deutsch	Svenska	Русский
Español	Latviski	Norsk

Dutch (Nederlands) and Polish (Polski) available in User Guide online.
Go to the product page at [SigmaAldrich.com](#) to download.

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada.

PRO5200 Rev 08/21

Symbol	Definition	Symbol	Definition	
	Federal (U.S.A.) law restricts this device to sale by or on the order of a licensed healthcare professional.		CE conformity marking	
	Catalogue number		Date of manufacture	
	Do not reuse		Manufacturer	
	Use-by date		Do not use if package is damaged	
	Batch code		Sterilized using ethylene oxide	

Made in Ireland
Merck Millipore Ltd.
Tullagreen,
Carrigrohilly,
Co. Cork, IRL

The initial M, Millipore, Millex, Durapore, and Sigma Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2021 Merck KGaA, Darmstadt, Germany and/or its affiliates. All rights reserved.

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Hergestellt in Irland

English

Millex® (13 mm) Filter with Durapore® Membrane

- Single use only
- Sterile
- Non-pyrogenic

Indications for Use/Purpose
The Millex®-GV/HV filters are medical devices intended for use as a syringe filter to sterilize (GV only) and/or clarify low volume solutions in direct patient care and pharmacy admixture applications.

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a licensed healthcare professional.

Introduction
The sterile 13 millimeter (mm) Millex®-GV and HV filters are syringe-operated filters for the sterilization or clarification of aqueous solutions (column eluates, tissue culture additives, HPLC samples, protein pharmaceuticals, etc.) in volumes of 0.5 to 10 milliliters (mL). These solvent-resistant, single-use filters incorporate the low-protein binding Durapore® polyvinylidene fluoride membrane. The Millex® filters are suitable for use as syringe filters to sterile filter (GV), and/or clarify (GV/HV) low volume solutions in direct patient care and pharmacy admixture applications.

Chemical Compatibility
Millex® filters with Durapore® membrane are compatible with most aqueous solutions. Based on information from technical publications, materials suppliers, and laboratory tests, Merck Millipore Ltd. believes that the agents listed in the following table are safe to use with Millex® filters. However, because of the effects of variability in temperature, concentrations, duration of exposure, and other factors outside of our control, Merck Millipore Ltd. does not provide or imply a warranty with respect to this information. Agents that are not listed should be tested with the Millex® filters prior to use.


Acetic acid (5%)	Cyclohexanone	Hexane	Mineral spirits
Acetic acid, glacial	Developers (photo)	Hydrochloric acid	Pentane
Acetonitrile	Ethyl alcohol	Hydrogen (gas)	Petroleum based oils
Ammonium hydroxide (6 N)	Ethylene glycol	Hydrofluoric acid	
Amyl alcohol	Formaldehyde	Hypo (photo)	Silicone oils
Benzyl alcohol (1%)	Freon® solvent, TF or PCA	Isobutyl alcohol	Sodium hydroxide
Boric acid	Gasoline	Isopropyl alcohol	Trichloroacetic acid
Brine (sea water)	Glycerine (glycerol)	Kerosene	Trichloroethane
Butyl alcohol	Helium	Methyl alcohol	


Directions for Use
WARNINGS


- To ensure sterility, do not use this product if the package is damaged.
- Do not use this product as an in-line filter for intravenous administration; it was not designed for long-term continuous use.
- Do not use with syringes smaller than 10 mL because pressures in excess of the maximum pressure rating may be reached, potentially causing damage to the filter and/or personal injury.


CAUTIONS


- Do not resterilize or reuse the Millex® filter, as Merck Millipore Ltd. cannot assure the sterility, integrity, and performance beyond a single use.
- Do not use the Millex® filter at temperatures above 45 °C (113 °F).
- Do not use the same Millex® filter to filter solutions in both directions.
- Do not use the Millex® filter to filter emulsions or suspensions.
- Do not use the Millex® filter to filter 5 mg or less of active drug materials unless binding studies have been performed.

Specifications	
Materials	
Membran	Low binding hydrophilic polyvinylidene fluoride (PVDF)
Pore size	Millex®-GV filter: 0.22 µm Millex®-HV filter: 0.45 µm High density polyethylene (HDPE)
Housing	
Dimensions	
Inlet to outlet	21.1 mm (0.83 in.)
Diameter	14.7 mm (0.58 in.)
Filtration area	0.65 cm² (0.10 in²)
Temperature limit	45 °C (113 °F)
Pressure limit at 21 °C	5.2 bar (75 psi) inlet and differential
Filtration volume	0.5–10 mL
Hold-up volume	≤ 50 µL after air purge
Sterilization method	Ethylene oxide gas
Connections	Female Luer-Lok™ inlet; male Luer-slip outlet
Flow rate at 30 psi, 21 °C (2.1 bar, 70 °F)	Millex®-GV filter: ≥ 7 mL/min
Flow rate at 10 psi, 21 °C (0.7 bar, 70 °F)	Millex®-HV filter: ≥ 12 mL/min

Notice
We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.
The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

Contact Information
For the location of the office nearest you, go to [SigmaAldrich.com/offices](#).
Technical Assistance
Visit the tech service page on our web site at [SigmaAldrich.com/techservice](#).
Standard Warranty
The applicable warranty for the products listed in this publication may be found at [SigmaAldrich.com/terms](#).

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