

Product Information

ANTI- BCL-10 (MK-17)

Developed in Rabbit, Affinity Isolated Antibody

Product Number **B 0431**

Product Description

Anti-Bcl-10 is developed in rabbit using a synthetic peptide corresponding to amino acid residues 1-17 of human Bcl-10, conjugated to KLH with glutaraldehyde, as immunogen. The corresponding sequence in rat and mouse differs by 2 amino acids. The antibody is affinity-purified using the immunizing peptide immobilized on agarose.

Anti-Bcl-10 specifically recognizes Bcl-10 by immunoblotting and immunoprecipitation (approximately 32 kDa). By immunoblotting, staining of Bcl-10 is inhibited with the immunizing peptide. The antibody may also be used for detection of Bcl-10 by immunohistochemistry. The epitope(s) recognized by the antibody is compatible with routine formalin-fixation and paraffin embedding. The antibody reacts with Bcl-10 of human, rat and mouse origin.

Bcl-10, an N-terminal CARD (Caspase Recruitment Domain) containing protein, is also designated CIPER, mE10, cE10, CARMEN and CLAP.¹⁻⁶ Bcl-10 is a cellular homolog of the equine herpesvirus-2 protein E-10 (vCLAP). It was implicated in the regulation of apoptosis by interacting with caspase 9, enhancing procaspase 9 processing and triggering its activation when overexpressed in the cell.^{3,7}

Bcl-10 cellular overexpression induces JNK, p38 and NF- κ B activation. Deregulation of Bcl-10 expression was also demonstrated to be involved in cellular oncogenesis.^{1,5} In mice, Bcl-10 plays an important role in the immune system functioning and in the development of the central nervous system while its roles in the *in vivo* execution of cell death and oncogenesis are not clear.⁸

Mucosa-associated lymphoid tissue (MALT) B lymphomas with the t(1;14)(p22;q32) are associated with overexpression and constitutive activity of Bcl-10. Such

tumors contain a variety of mutations, most of which result in truncations either in the CARD domain or carboxy-terminal to it. Bcl-10 mutations are also found in cases of follicular lymphoma and diffuse large B cell lymphoma.⁹ Mutations of the Bcl-10 gene do not appear to play a major role in the pathogenesis of human solid neoplasms or leukemias.

In normal tissues Bcl-10 is detectable in lymphoid organs and in the cytoplasm of mammary gland cells. On the other hand, both nuclear and cytoplasmic expression are detected in MALT lymphomas especially those with the trans location t(1;14)(p22;q32).¹⁰ Bcl-10 protein was reported to bind itself, TRAF1, TRAF2, TRAF5 and CARD9.^{4,7,11} Overexpressed Bcl-10 protein was shown to be arranged in cytoplasmic filaments in cultured cells and reported to be essential for recruitment of several signal transducer molecules such as TRADD and RIP.¹²

Reagent

Anti-Bcl-10 is provided as affinity isolated antibody in 0.01 M phosphate buffered saline, pH 7.4, containing 1 % bovine serum albumin and 15 mM sodium azide.

Precautions and Disclaimer

Due to the sodium azide content a material safety sheet (MSDS) for this product has been sent to the attention of the safety officer of your institution. Consult the MSDS for information regarding hazards and safe handling practices.

Storage/Stability

For continuous use, store at 2 °C to 8 °C for up to one month. For extended storage, freeze in working aliquots. Repeated freezing and thawing is not recommended. Storage in "frost-free" freezers is not recommended. If slight turbidity occurs upon prolonged storage, clarify the solution by centrifugation before use. Working dilution samples should be discarded if not used within 12 hours.

Product Profile

A minimum working dilution of 1:500 is determined by immunoblotting using whole extracts of human THP-1 acute monocytic leukemia cells, human Raji Burkitt's lymphoma cells or mouse NIH-3T3 fibroblasts.

For immunoprecipitation, 2.5 µg to 5 µg of the antibody immunoprecipitates Bcl-10 from a RIPA lysate of 5×10^5 human Raji Burkitt's lymphoma cells or human THP-1 acute monocytic leukemia cells.

A minimum working dilution of 1:8,000 is determined by indirect immunoperoxidase staining of trypsin-digested, formalin-fixed, paraffin-embedded tissue sections of rat spleen.

Note: In order to obtain best results in different techniques and preparations we recommend determining optimal working dilution by titration test.

References

1. Willis, T.G., et al., *Cell*, **96**, 35-45 (1999).
2. Koseki, T., et al., *J. Biol. Chem.*, **274**, 9955-9961 (1999).
3. Yan, M., et al., *J. Biol. Chem.*, **274**, 10287-10292 (1999).
4. Thome, M., et al., *J. Biol. Chem.*, **274**, 9962-9968 (1999).
5. Zhang, Q., et al., *Nat. Genet.*, **22**, 63-68 (1999).
6. Srinivasula, S.M., et al., *J. Biol. Chem.*, **274**, 17946-17954 (1999).
7. Yoneda, T., et al., *J. Biol. Chem.*, **275**, 11114-11120 (2000).
8. Ruland, J., et al., *Cell*, **104**, 33-42 (2001).
9. Du, M.Q., et al., *Blood*, **95**, 3885-3890 (2000).
10. Ye, H., et al., *Am. J. Pathol.*, **157**, 1147-1154 (2000).
11. Bertin, J., et al., *J. Biol. Chem.*, **275**, 41082-41086 (2000).
12. Guiet, C., and Vito, P., *J. Cell Biol.*, **148**, 1131-1140 (2000).

lpg 6/01

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.