


**BIOSCOT®**

**Anti-K**

Monoclonal Human IgM  
Blood Grouping Reagent



For Slide, Tube, Microplate, Gel & Column Techniques.

Cell Line:MS-56

**INTENDED USE**

BIOSCOT Anti-K (cell line MS-56) monoclonal human IgM blood grouping reagent is used to ensure the immunological compatibility of blood and blood components intended for transfusion. This qualitative reagent will detect the presence or absence of the K (KEL1) antigen on the surface of human red blood cells when tested according to the slide, tube, microplate, gel and column techniques. The reagent is designed for in vitro diagnostic, professional use by operators trained in serological techniques.

**PRINCIPLE OF THE REAGENT**

When used by the recommended techniques this reagent will cause agglutination (clumping) of red cells carrying the specific antigen (positive test). Lack of agglutination of the red cells demonstrates the absence of the specific antigen (negative test).

The reagent is intended for use by non-automated, manual techniques.

The reagent has been characterised by the procedures recommended in these instructions for use, its suitability for use in other techniques must be determined by the user.

**PRECAUTIONS**

1. All blood products should be treated as potentially infectious. The human donors or the cell lines used to produce these reagents have been tested and found to be negative for Anti-HIV, Anti-HCV, HBsAg, EBV and Mouse Antibody Production (MAP) viruses. No known tests can guarantee that any product derived from human blood is free from infectious agents. Care must be taken in the use and disposal of each container and its contents.
2. The reagent contains 0.1% (w/v) sodium azide. Sodium azide may be toxic if ingested and may react with lead or copper plumbing to form highly explosive salts. On disposal, flush with large quantities of water.
3. The product should be clear. Turbidity may indicate bacterial contamination. This reagent should not be used if a precipitate, fibrin gel or particles are present.
4. The bovine material is obtained from USDA approved sources or from sources for which origin information is available. The donor animals have been inspected and certified disease free and are deemed to have low TSE (Transmissible Spongiform Encephalopathy) risk.
5. The product should be disposed of either by overnight immersion in disinfectants at appropriate concentrations or by autoclaving.

**CONTROLS**

It is recommended that a positive control and a negative control should be tested in parallel with each batch of tests. Tests must be considered invalid if controls do not show the expected reactions.

It is not required to use a reagent control in parallel with all tests using this reagent. Only in typing the red cells of patients known to have auto-antibodies or protein abnormalities is the use of a reagent control such as BIOSCOT Monoclonal Control (Product code: TT) recommended. This should be tested in parallel with the reagent.

**STORAGE**

Store the opened / unopened product at 2-8°C until the expiry date detailed on the product label.

Failure to store the product at the correct temperature, for example, storage at higher temperature or repeated freezing and thawing, may result in accelerated loss of reagent activity.

**SPECIMEN COLLECTION**

No special preparation of the patient/donor is required prior to specimen collection. Blood should be collected by an approved phlebotomy technique into tubes containing EDTA or CPD. The specimen should be tested as soon as possible following collection. Samples that cannot be tested within 24 hours of collection should be stored at 2-8°C. Testing should be carried out within 14 days of collection. Specimens displaying gross haemolysis or microbial contamination should not be tested with this reagent. Failure to store the specimens in the correct conditions may result in false positive or false negative results.

**MATERIALS PROVIDED**

Product code TR Anti-K blood grouping reagent contains antibody from cell line MS-56. The reagent is composed of monoclonal human IgM antibodies in a buffer solution containing macromolecular chemical potentiators. The reagents contain 0.1% (w/v) sodium azide and bovine material. The product is supplied filtered to 0.22 µm. The reagent has been optimised for use by the recommended techniques without further dilution or additions.

Contents:

1 x reagent vial **REF** TR-5ML-B

10 x reagent vials for **REF** TR-10X5ML-B

1 x information sheet

**MATERIALS REQUIRED BUT NOT PROVIDED**

**Slide Technique:**

- Microscope slide
- Isotonic saline or compatible plasma/serum
- Timer

**Tube Technique:**

- Test tube
- Isotonic saline
- 37°C Incubator
- Timer
- Centrifuge (1000 rcf)

**Microplate Technique:**

- U well microplate
- Isotonic saline
- Timer
- Centrifuge (100 rcf)
- Microplate shaker
- Microplate reader (optional)

**Bio-Rad ID Gel Technique:**

- Bio-Rad ID-Card "NaCl, enzyme test and cold agglutinins"
- Isotonic saline, phosphate buffered saline (PBS) or ID-diluent 2
- Micropipettes capable of delivering 10, 25 and 50 µL
- Timer
- Centrifuge suitable for Bio-Rad ID-Cards
- Reader (optional)

**Ortho BioVue Column Technique:**

- Ortho BioVue® Neutral Cassettes
- Isotonic saline, phosphate buffered saline (PBS) or Ortho® 0.8% Red Cell Diluent
- Micropipettes capable of delivering 10, 40 and 50 µL
- Timer
- Centrifuge suitable for Ortho BioVue Cassettes
- Reader (optional)

## **RECOMMENDED TECHNIQUE**

### **1. SLIDE TECHNIQUE**

- 1.1. Prepare a 35-50% suspension of test red cells in autologous (or compatible) plasma, serum or in isotonic saline.
- 1.2. Add 1 drop (40-50 µL) of Anti-K reagent to a clean, labelled microscope slide.
- 1.3. Add 1 drop (40-50 µL) of the suspension of test red cells.
- 1.4. Mix the antiserum and cells over an area about 2 cm in diameter by gently and continuously rocking the slide.
- 1.5. Read macroscopically after 2 minutes. Do not confuse any drying of the mixture with agglutination.

### **2. TUBE TECHNIQUE**

- 2.1. Prepare a 3-5% suspension of test red cells in isotonic saline.
- 2.2. Add 1 drop (40-50 µL) of Anti-K reagent to an appropriately labelled test tube
- 2.3. Add 1 drop (40-50 µL) of the suspension of test red cells.
- 2.4. Mix and centrifuge at 1000 rcf for 20 seconds.
- 2.5. Gently agitate the tube to dislodge the red cells and examine macroscopically for agglutination.
- 2.6. Incubate all negative or weakly positive tests at 37°C for 5 minutes and repeat steps 2.4 and 2.5. This may enhance the reaction strength in typing red cells of weak or variant phenotypes.

### **3. MICROPLATE TECHNIQUE**

- 3.1. Prepare a 3-5% suspension of test red cells in isotonic saline.
- 3.2. Add 1 drop (40-50 µL) of Anti-K reagent to the appropriate test wells of a microplate.
- 3.3. Add an equal volume (40-50µL) of the test cell suspension to the appropriate test wells .
- 3.4. Mix the contents of each well using manual means or a microplate shaker.
- 3.5. Incubate the microplate at ambient temperature for 15-20 minutes.
- 3.6. Centrifuge the microplate at 100 rcf for 40 seconds.
- 3.7. Resuspend the red cells by manual means or using the microplate shaker.
- 3.8. Read tests macroscopically or with a reader. The use of a plate reader must be validated by the user.

### **4. BIO-RAD ID GEL TECHNIQUE**

- 4.1. Prepare a 3-5% suspension of test red cells in isotonic saline/PBS or a 0.8% suspension in ID-Diluent 2.
- 4.2. Add 10 µL of 3-5% or 50 µL of 0.8% suspension of test red cells to the appropriate microtube of the ID-card
- 4.3. Add 25 µL of Anti-K reagent to the appropriate microtube.
- 4.4. Mix gently and centrifuge the ID-card at the speed and time recommended in the ID-card manufacturer's instructions.
- 4.5. Read macroscopically or with a reader. The use of a reader must be validated by the user.

### **5. ORTHO BIOVUE COLUMN TECHNIQUE**

- 5.1. Prepare a 3-5% suspension of test red cells in isotonic saline/PBS or a 0.8% suspension in Ortho 0.8% Red Cell Diluent.
- 5.2. Add 10 µL of 3-5% or 50 µL of 0.8% suspension of test red cells to the appropriate Ortho BioVue cassette reaction chamber.
- 5.3. Add 40 µL of Anti-K reagent to the appropriate reaction chamber.
- 5.4. Mix gently and centrifuge the cassette at the speed and time recommended in the cassette manufacturer's instructions.
- 5.5. Read macroscopically or with a reader. The use of a reader must be validated by the user.

## **LIMITATIONS**

Red cells that have a positive direct antiglobulin test (DAT) may produce false positive results. The use of BIOSCOT Monoclonal Control reagent (product code TT) is recommended for detection of such potentially false positive results.

Rigid polystyrene microplates are generally more suitable than those made from PVC. Each batch of microplates should be evaluated in the user's system prior to acceptance as suitable for routine usage.

Antigen variant cells may produce unexpected positive or negative reactions with samples previously typed with blood grouping reagents of polyclonal or other cell line-derived monoclonal sources.

Incorrect handling or storage of the ID-cards or Ortho BioVue cassettes could lead to incorrect results. Cards and cassettes should be stored and handled according to the manufacturer's instructions.

False positive and false negative results may occur through contamination of test materials or any deviation from the recommended techniques.

## **PERFORMANCE CHARACTERISTICS**

Anti-K (cell line MS-56) monoclonal blood grouping reagent TR has been tested by each of the recommended techniques with donor, clinical and neonatal specimens. The sample population represented all major phenotypes. The total number of tests (n), and the calculated sensitivity and specificity for each technique are displayed below:

| TECHNIQUE      | Anti-K Product Code TR |     |             |     |
|----------------|------------------------|-----|-------------|-----|
|                | SENSITIVITY            |     | SPECIFICITY |     |
|                | n                      | %   | n           | %   |
| Slide          | 10                     | 100 | 109         | 100 |
| Tube           | 24                     | 100 | 260         | 100 |
| Microplate     | 10                     | 100 | 109         | 100 |
| Bio-Rad ID Gel | 18                     | 100 | 82          | 100 |
| Ortho BioVue   | 16                     | 100 | 84          | 100 |

**Diagnostic Sensitivity:** The probability that the device gives a positive result in the presence of the target marker.

**Diagnostic Specificity:** The probability that the device gives a negative result in the absence of the target marker.

## **ANALYTICAL PERFORMANCE**

This blood grouping reagent(s) exhibited unequivocal positive or negative results by all recommended techniques. Performance was found to be acceptable in terms of repeatability, reproducibility, and robustness.

## **FURTHER INFORMATION**

For technical assistance contact: [SigmaAldrich.com/techservice](http://SigmaAldrich.com/techservice)  
Any serious incident that has occurred in relation to this reagent must be reported to Millipore (UK) Ltd and the competent authority of the Member State in which the user and/or the patient is established.

The summary of safety and performance (SSP) for this device is available in the European database on medical devices (Eudamed) at <https://ec.europa.eu/tools/eudamed>, where it is linked to the Basic UDI-4053252TRBTRVK.

## **BIBLIOGRAPHY**

1. Guidelines for the Blood Transfusion Services in the United Kingdom. 8th Edition. The Stationary Office, 2013
2. Issitt, P.D. and Anstee, D.J. Applied Blood Group Serology, 4th Edition, Montgomery Scientific Publications, 1998.
3. AABB Technical Manual 20<sup>th</sup> Edition, 2020.

**SUMMARY OF CHANGES FROM PI16/EN 2023-12 to PI16/EN 2025-10**

1. Update to Basic UDI number in further information section.
2. Changed heading for Summary of changes to Summary of changes from versions.
3. Date within cliché code updated to Oct 2025.



Millipore (UK) Ltd  
Fleming Road  
Kirkton Campus  
Livingston, EH54 7BN  
United Kingdom  
Tel: +44 (0)1506 404000

[www.sigmaaldrich.com](http://www.sigmaaldrich.com)

