

3-AMINO-1,2,4-TRIAZOLE
Sigma Prod. No. A8056**CAS NUMBER:** 61-82-5**SYNONYMS:** 3-Amino-1*H*-1,2,4-triazole;
1*H*-1,2,4-Triazol-3-amine; Aminotriazole; ATA; ENT
25445; Amizol; Cytrol; Weedazol; 3 A-T**PHYSICAL PROPERTIES:**

Appearance: White to light yellow powder

Melting Point: 150-153°C¹Molecular Formula: C₂H₄N₄

Molecular Weight: 84.08

STABILITY / STORAGE AS SUPPLIED:

Store desiccated at -20°C.

SOLUBILITY / SOLUTION STABILITY:

According to the Merck Index, 12th ed., p. 80, No. 506, the product is soluble in water, methanol, ethanol and chloroform, and sparingly soluble in ethyl acetate. It is insoluble in ether and acetone. Aqueous solutions are neutral. Sigma tests the solubility in hot water at 50 mg/mL, obtaining clear and faint yellow-tan solutions. A paper reports that an aqueous solution of radiolabeled product was stored refrigerated; the solution's shelf life is not given.²

USAGE / APPLICATIONS:


3-Amino-1,2,4-triazole is an inhibitor of mitochondrial and chloroplast function. Commercial grade 3-amino-1,2,4-triazole (which generally contains catalase anti-inhibitory impurities) is used as a herbicide and cotton defoliant.³

Enzyme or Reaction Affected:

Inhibits cyclizations leading to β -carotene α -Carotene accumulates.³
Causes chlorosis
Disorganizes chloroplast membranes
Blocks formation of 18S Fraction I protein
Blocks formation of chloroplast DNA
Blocks formation of 70S ribosomes
Inhibits riboflavin biogenesis
Catalase
Peroxidase
 δ -Aminolevulinatase synthetase
Methane monooxygenase
Imidazolglycerol-phosphate dehydratase

Comments:

Concentration = 1×10^{-4} M³
Concentration = 1×10^{-4} M³
³
³
³
³
Inhibits at 2×10^{-2} M in the presence of H₂O₂.³
No action.³
³
Inhibits 75% at 10 μ M.⁴
K_i = 30 μ M.⁴



3-AMINO-1,2,4-TRIAZOLE
Sigma Prod. No. A8056

REFERENCES:

1. *Aldrich Chemical Co. Catalog*, p. 97 (1996-1997).
2. M. C. DeTraglia, et al., *Anal. Biochem.*, 99, 464 (1979).
3. R. M. C. Dawson, et al., *Data for Biochemical Research*, 3rd ed., 302-303 (1993).
4. H. Zollner, *Handbook of Enzyme Inhibitors*, Part B, 563 (1993).

Sigma warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.