

3050 Spruce Street Saint Louis, Missouri 63103 USA Telephone 800-325-5832 • (314) 771-5765 Fax (314) 286-7828 email: techserv@sial.com sigma-aldrich.com

ProductInformation

Ethyl methanesulfonate

Product Number **M0880** Store at Room Temperature

Product Description

Molecular Formula: C₃H₈O₃S Molecular Weight: 124.2 CAS Number: 62-50-0 Boiling Point: 213-213.5 °C¹ Molarity: 9.7 (neat liquid) Density: 1.206 g/ml (20° C) Synonym: EMS

Ethyl methanesulfonate has been used as a mutagen for both mammalian cells² and plants.³ The mutagenicity¹ and carcinogenicity of this product has been reviewed.^{4,5,6,7} A comparative study of mutagenicity with ethyl and methyl methanesulfonate was reviewed.⁸

Precautions and Disclaimer

For Laboratory Use Only. Not for drug, household or other uses.

Preparation Instructions

EMS is miscible in chloroform (200 mg/ml).

EMS can be diluted to a concentration of 1% in water or serum-free medium and then added to a culture at a final concentration in the range of 200-300 μ g/ml.⁷

Storage/Stability

Reaction rate constants and half-lives for aqueous solutions of this product were calculated at increasing temperatures and at several pH's. The rate of hydrolysis is unaffected by varying (acidic) pH.

Temperature (°C)	Half-life in hours ³
0	1716.0
5	796.6
10	378.7
15	185.7
20	93.2
25	48.5
30	7.8
45	4.6
50	2.7.

References

- 1. IARC Monographs, 7, 245-252 (1974).
- 2. Thompson, L. H., Mutant Isolation. Meth. Enzymol., **5**8, 308-322 (1979).
- 3. Frose-Gertzen, E. E., et al., Nature, **198**, 447-8 (1963).
- 4. Merck Index, 11th ed., Entry# 3782.
- Alderson, T., Chemically Induced Delayed Germinal Mutation in Drosophila. Nature, 207(993), 164-167 (1965).
- Jenkins, J. B., The Induction of Mosaic and Complete Dumpy Mutants in Drosophila Melanogaster with Ethyl Methanesulfonate. Mutat. Res., 4(1), 90-92 (1967).
- Schalet, A. P., Interspecific Comparison of Ethyl Methanesulfonate-induced Mutation Rates in Relation to Genome Size. Mutat. Res., 49(3), 313-340 (1978).
- 8. Kondo, S., Environ. Sci. Res. 24, 743-785 (1981).

FEB/RXR 10/08

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.