

3050 Spruce Street Saint Louis, Missouri 63103 USA Telephone 800-325-5832 • (314) 771-5765 Fax (314) 286-7828 email: techserv@sial.com sigma-aldrich.com

# ProductInformation

### Betaine

Product Number **B 2629** Storage Temperature 2-8 °C

# **Product Description**

Molecular Formula:  $C_5H_{11}NO_2$ Molecular Weight: 117.1 CAS Number: 107-43-7 pK<sub>a</sub>: 1.83<sup>1</sup> Melting Point: 310 °C (with decomposition) Synonyms: oxyneurine, glycine betaine, (carboxymethyl)trimethylammonium inner salt, 1-carboxy-N,N,N-trimethylmethanammonium inner salt<sup>2</sup>

Betaine is a reagent that is used in soldering, resin curing fluxes, and organic synthesis.<sup>1</sup> Betaine compounds occur naturally in plants and animals, such as in sugar beets and marine algae. A review of the role of glycine betaine in protecting plants from environmental stress has been published.<sup>3</sup> The role of betaine and other organic osmolytes in maintaining osmolarity in mammalian cells has been reviewed.<sup>4</sup>

Betaine has been reported to diminish the base pair composition dependence of DNA thermal melting transitions, with a concentration of 5.2 M betaine being the value at which AT and GC base pairs are equally stable.<sup>5</sup> Betaine has subsequently found application in the polymerase chain reaction (PCR).<sup>6,7</sup> The usefulness of betaine in PCR has resulted from its ability to enhance DNA amplification by diminishing the formation of secondary structure in GC-rich DNA regions.<sup>8</sup> A protocol for the production of long-chain cDNAs that uses betaine and trehalose has been published.<sup>9</sup>

The permanent positive charge on betaine makes it a target of study for mass spectrometry (MS) analysis by such methods as electrospray ionization (ESI) MS.<sup>10</sup> A protocol for the analysis of betaine and other choline related compounds from tissues that combines HPLC with isotope dilution ESI-MS has been described.<sup>10</sup>

# **Precautions and Disclaimer**

For Laboratory Use Only. Not for drug, household or other uses.

# **Preparation Instructions**

This product is soluble in water (50 mg/ml), yielding a clear, colorless solution. 5 M solutions of betaine in water are available (Product No. B 0300). This product is soluble in also methanol (55 g/100 g, w/w) and ethanol (8.7 g/100 g).<sup>1</sup>

### References

- Data for Biochemical Research, 3rd ed., Dawson, R. M. C., et al., Oxford University Press (New York, NY: 1986), pp. 8-9.
- 2. The Merck Index, 12th ed., Entry# 1225.
- Sakamoto, A., and Murata, N., The role of glycine betaine in the protection of plants from stress: clues from transgenic plants. Plant Cell. Environ., 25(2), 163-171 (2002).
- Lang, F., et al., The diversity of volume regulatory mechanisms. Cell. Physiol. Biochem., 8(1-2), 1-45 (1998).
- Rees, W. A., et al., Betaine can eliminate the base pair composition dependence of DNA melting. Biochemistry, **32(1)**, 137-144 (1993).
- 6. Papp, A. C., et al., Strategies for amplification of trinucleotide repeats: Optimization of fragile X and androgen receptor PCR. Mol. Diagn., **1(1)**, 59-64 (1996).
- Weissensteiner, T., and Lanchbury, J. S., Strategy for controlling preferential amplification and avoiding false negatives in PCR typing. Biotechniques, **21(6)**, 1102-1108 (1996).
- Henke, W., et al., Betaine improves the PCR amplification of GC-rich DNA sequences. Nucleic Acids Res., 25(19), 3957-3958 (1997).
- Spiess, A. N., and Ivell, R., A highly efficient method for long-chain cDNA synthesis using trehalose and betaine. Anal. Biochem., **301(2)**, 168-174 (2002).

- Wood, K. V., et al., Characterization of betaines using electrospray MS/MS. Phytochemistry, 59(7), 759-765 (2002).
- Koc, H., et al., Quantitation of choline and its metabolites in tissues and foods by liquid chromatography/electrospray ionization-isotope dilution mass spectrometry. Anal. Chem., **74(18)**, 4734-4740 (2002).

The purchase of this product does not include a license to practice the claims of U.S. Patent No. 5,545,539, DE4411588, or DE4411594. The practice of the claims of these patents may require a license from the patent owners.

GCY/RXR 11/03

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.