www.sigmaaldrich.com



58030 Iodonitrotetrazolium chloride

(p-Iodonitrotetrazolium Violet, 2-(4-Iodophenyl)-3-(4-nitrophenyl)-5-phenyl-2*H*-tetrazolium chloride, INT)

CAS number: 146-68-9

Product Description:



Solubility:

INT dissolves in water to give a clear light yellow to light green solution at 4 mg/mL, although gentle heat and sonication are necessary for this concentration. INT dissolves up to 8 mg/mL in ethanol, and up to 30 mg/mL in ethylene glycol monomethyl ether. Aqueous solutions are ³ reported to be stable for several weeks when stored at 2-8°C in the dark.⁴ A 50 mg/mL solution in hot methanol/water (1:1) is clear and very deep yellow.²

Applications:

This tetrazolium dye precursor forms a purple formazan dye on reduction and has been used in a variety of applications. It is considered to have higher reactivity than some tetrazolium compounds, at least with respect to succinate dehydrogenase, with optimal results obtained using a concentration of 0.8 mM INT. The formazan (available as I7375) is ⁵ comparatively insoluble in water, so can be detected visually or extracted into ethanol or organic solvent and determined spectrophotometrically. Altman reviewed properties of a number of different tetrazolium salts.¹

For use in measurement of respiratory electron transport activity: *Marine Biology*, 30, 27-36 (1975). Used in a colorimetric assay to measure *serum* lactate dehydrogenase: *Clinical Chemistry*, 19, 766 (1973).

INT is used as an electron acceptor for the colorimetric assays of: lactate dehydrogenase^{1,11}, xanthine dehydrogenase¹⁰, lactyl-CoA dehydrogenase⁹, succinate dehydrogenase⁸, BBM II ketolisomerase⁷, histidinol dehydrogenase⁷ and diverse other hydrolases.

Active bacterial biomass was estimated using INT stain for counting¹² and cells have been stained with INT.



References:

- 1. Altman, F.P., Histochem. J., 8, 47-485 (1976). "Tetrazolium Salts: A Consumers Guide" (TLC chromatography systems for different salts, using UV detection, 300-400 nm).
- 2. Sigma-Aldrich quality control.
- 3. Green, F.J., The Sigma-Aldrich Handbook of Stains, Dyes and Indicators, p. 410 (1990).
- 4. Clin. Chem. Acta, 12, 210 (1965).
- 5. Green, J.D., and Narahara, H.T., J. Histochem. Cytochem., 28, 408 (1980).
- 6. A.L. Babson, S.R. Babson, Clin. Chem. 19, 766 (1973)
- 7. R.G. Martin, BBM II ketolisomerase and histidinol dehydrogenase Meth. Enzymol. 17B, 3 (1971)
- 8. D.J. Morré, succinate dehydrogenase Meth. Enzymol. 22, 130 (1971)
- 9. R.L. Baldwin, W.A. Wood, lactyl-CoA dehydrogenase Meth. Enzymol. 9, 683 (1966)
- 10. R. Fried, xanthine dehydrogenase Anal. Biochem. 16, 427 (1966)
- 11. M.N. Nachlas, Electron acceptor for the colorimetric assay of: lactate dehydrogenase Anal. Biochem. 1, 317 (1960)
- 12. Stamatiadis s., et al., Use of Staining inhibitors to separate fungal and bacterial activity in soil. Soil Bol. Biochem. 22, 81-88 (1990).
- 13. Suh Y.-A., ARNOLD R.S., Cell transformation by the superoxide-generating oxidase Mox1, Nature 401, 79 - 82 (1999)

Precautions and Disclaimer:

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

The vibrant M and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. Detailed information on trademarks is available via publicly accessible resources. © 2018 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the US and Canada.

