Instructions for Use

Automated Wright Stain Set

Procedure No. WSHT

Intended Use

Wright Stain Set is specifically designed for use with Hema-Tek® Automated Slide Stainers. Wright Stain reagents are for "In Vitro Diagnostic Use." For professional use only. The data obtained from this automated, qualitative procedure demonstrates cellular elements in blood and bone marrow films of human specimens. This data when reviewed in conjunction with other diagnostic tests and information may be used as an aid to diagnosis of viral and bacterial infections, allergenic reactions, and other pathological states.

When blood films are stained as described in the procedure, the white blood cell nucleus and cytoplasm take on characteristic blue or pink coloration. The purified dyes in the product eliminate inconsistent staining and yield lot-to-lot reproducible response.

Reagents

Thiazine Dye Solution (Cat. No. T3272-250ML)
Methanolic solution of a thiazine dye, 0.1% w/v, C.I. 52010.

Xanthene Dye Solution (Cat. No. X2251-475ML)
Buffered solution of a xanthene dye, 0.05% w/v, C.I. 45380, with surfactant, pH 6.5. Sodium azide, 0.02%, added as preservative

Rinse Solution No 2 (Cat. No. RS2-900ML)
Ethanol solution, 18% w/v, with surfactant. Sodium azide, 0.02%, added as preservative

Special Materials Required but Not Provided

- Hema-Tek® Stainer
- Microscope
- Microscope slides

Storage and Stability

Store reagents at room temperature (18-26°C). Reagents are stable until expiration date shown on the label.

Deterioration

Discard if turbidity develops.

Preparation

Wright Stain Set is supplied ready for use.

Precautions

These IVDs are intended for in vitro diagnostic use in a clinical laboratory environment. These IVDs are for professional use by qualified personnel only. Sigma-Aldrich IVDs may be operated by laboratory personnel who are trained to handle human specimens that can be infectious, use microscopes and other laboratory equipment and have color perception and visual acuity to distinguish colors and other objects under a microscope.

Normal precautions exercised in handling laboratory reagents should be followed. Dispose of waste observing all local, state, provincial or national regulations.

Procedure

Specimen Collection

No known test method can offer complete assurance that blood samples or tissue will not transmit infection. Therefore, all blood derivatives or tissue specimens should be considered potentially infectious.

Fresh whole blood films or fresh films from blood using EDTA as an anticoagulant must be used.

Notes

These settings have given satisfactory results in our laboratory. Individual preferences may dictate pump adjustments. If stain is too light, increase settings on Pump 1 until satisfactory color is obtained. For optimum results, it is recommended that the setting for Pump 2 does not exceed the setting for Pump 1.

Positive control slides should be included in each run.

Procedure

- General instrument set-up and use is as described in the Hema-Tek® procedure manual.
- Thiazine Dye Solution: Remove cap, leaving inner seal on bottle to minimize moisture pick-up. Insert cannula from Pump 1 directly through seal.
Pump Setting: -3 (maximum counterclockwise position of volume control knob).
- Xanthene Dye Solution: Remove cap, leaving inner seal on bottle. Insert cannula from Pump 2 directly through seal.
Pump Setting: -3 (maximum counterclockwise position of volume control knob).
- Rinse Solution 2: Remove cap, leaving inner seal on bottle. Insert cannula from Pump 3 directly through seal.
Pump Setting: -3 (maximum counterclockwise position of volume control knob).

Performance Characteristics

Nuclei will be varying shades of purple. Cytoplasmic staining will be varying shades of blue to light pink. Fine reddish to lilac granules may be present in cytoplasm of some cell types. Basophils will demonstrate dark blue to black granules in the cytoplasm. Eosinophils will demonstrate bright orange granules in the cytoplasm. Red blood cells should be pink to orange.¹

If observed results vary from expected results, please contact Sigma-Aldrich Technical Service for assistance.

Analytical Performance Characteristics

The analytical performance results for the given tests conducted on all target structures, confirm 100% sensitivity, specificity, and repeatability.

Cat. No	Product Description	Target	Intra-assay Specificity	Intra-assay Sensitivity	Inter-assay Specificity	Inter-assay Sensitivity
T3272	Thiazine Dye Solution	Nuclei	3 of 3	3 of 3	3 of 3	3 of 3
		Cytoplasm	3 of 3	3 of 3	3 of 3	3 of 3
		Granules	3 of 3	3 of 3	3 of 3	3 of 3
X2251	Xanthene Dye Solution	Nuclei	3 of 3	3 of 3	3 of 3	3 of 3
		Cytoplasm	3 of 3	3 of 3	3 of 3	3 of 3
		Granules	3 of 3	3 of 3	3 of 3	3 of 3
RS2	Rinse Solution No 2	Nuclei	3 of 3	3 of 3	3 of 3	3 of 3
		Cytoplasm	3 of 3	3 of 3	3 of 3	3 of 3
		Granules	3 of 3	3 of 3	3 of 3	3 of 3

Warnings and Hazards

Refer to Safety Data Sheet and product labeling for any updated risk, hazard or safety information.

WSHT:

H225 Highly flammable liquid and vapor.

H301 + H311 + H331 Toxic if swallowed, in contact with skin or if inhaled.

H319 Causes serious eye irritation.

H370 Causes damage to organs (Eyes, Central nervous system).

H412 Harmful to aquatic life with long lasting effects.

P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

P273 Avoid release to the environment.

P280 Wear protective gloves/ protective clothing/ eye protection/ face protection/ hearing protection.

P301 + P310 IF SWALLOWED: Immediately call a POISON CENTER/ doctor.

P303 + P361 + P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water.

P304 + P340 + P311 IF INHALED: Remove person to fresh air and keep comfortable for breathing. Call a POISON CENTER/ doctor.

If during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

Symbol Definitions

Symbols as defined in EN ISO 15223-1:2021

	Manufacturer		Catalogue Number
	Consult Instructions for Use		Batch Code
	Authorized Representative in the European Community/ European Union		European Union Declaration of Conformity (defined in IVDR 2017/746)
	Use-by Date		In vitro diagnostic medical device
	Temperature Limit		Caution
	Date of Manufacture		Importer
	Indicates the Authorised Representative in Switzerland		

References

- Hematology: Principles and Procedures, Sixth Edition, Brown AB, Lea & Febiger, Philadelphia 1993 p100-102

Contact Information

To place an order, please visit our web site at [SigmaAldrich.com](https://www.sigmaaldrich.com). For Technical Service, please visit the tech service page on our web site at [SigmaAldrich.com/techservice](https://www.sigmaaldrich.com/techservice).

Revision History

Rev. 3.0	2014
Rev. 4.0	2016
Rev. 5.0	2023
	Transferred to new template with current branding. Specified for professional use in intended use and precautions. Removed references to specific stainers. Moved aid to diagnosis statement to intended use. Revised intended use to align with IVDR guidelines. Updated Material Safety Data Sheet to Safety Data Sheet. Updated contact information. Removed instruction to follow CLSI for specimen collection. Removed EN 980 and changed to EN ISO 15223-1:2021 for symbols. Added adverse event contact information. Updated literature references. Added a revision history table. Added Warnings and Hazards. Added CH-REP information.
Rev. 6.0	2025
	Updated RS2-900ML name. Added CH-REP importer information.

Sigma-Aldrich, Inc.,
3050 Spruce Street,
St. Louis, MO 63103 USA
an affiliate of Merck KGaA,
Darmstadt, Germany
+1(314) 771-5765

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Merck KGaA,
64271
Darmstadt,
Germany

MDSS CH GmbH
Laurenzenvorstadt 61
5000 Aarau
Switzerland

The Initial M and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2023 Merck KGaA, Darmstadt, Germany and/or its affiliates. All rights reserved.

The life science business of Merck operates as MilliporeSigma in the U.S. and Canada.

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MDSS GmbH