



Product Information

SIGMA QUALITY CONTROL TEST PROCEDURE

Enzymatic Assay of PECTINASE (EC 3.2.1.15)

PRINCIPLE:

Polygalacturonic Acid + H₂O $\xrightarrow{\text{Pectinase}}$ Galacturonic Acid

I₂ + Galacturonic Acid \longrightarrow Oxidation products

I₂ + 2Na₂S₂O₃ \longrightarrow 2NaI + Na₂S₄O₆ (Titration of Excess Iodine)

CONDITIONS: T = 25°C, pH = 4.0

METHOD: Titrimetric

REAGENTS:

- A. 0.5% (w/v) Polygalacturonic Acid Solution (Polygalacturonic Acid)
(Prepare 100 ml in deionized water using Polygalacturonic Acid, Sigma Prod. No. P-3889.
Adjust to pH 4.0 at 25°C with 1 M NaOH.)
- B. 50 mM Iodine with 200 mM Potassium Iodide (I₂/KI)
(Prepare 500 ml in deionized water using Iodine,
Sigma Prod. No. I-3380 and Potassium Iodide, Sigma Prod. No. P-8256. Initially dissolve in
50 ml deionized water, allow to set at room temperature for 30 minutes and then QS to
500 ml with deionized water.)
- C. 1 M Sodium Carbonate (Na₂CO₃)
(Prepare 100 ml in deionized water using Sodium Carbonate, Anhydrous, Sigma Prod. No.
S-1641.)
- D. 2.0 N Sulfuric Acid (H₂SO₄)
(Prepare 200 ml in deionized water by diluting Sulfuric Acid, ACS Reagent, Sigma Prod. No.
S-1526 (1:17.8).)
- E. 100 mM Sodium Thiosulfate, Standardized (Na₂S₂O₃)
(Prepare 100 ml in deionized water using Sodium Thiosulfate, Pentahydrate, Sigma Prod.
No. S-8503. Standardize according to the ACS Reagent Procedure.)

**Enzymatic Assay of PECTINASE
(EC 3.2.1.15)**

REAGENTS: (continued)

- F. 1.0% (w/v) Starch Indicator
(Prepare 100 ml in deionized water using Starch Potato, Soluble, Sigma Prod. No. S-2630.
PREPARE FRESH.)
- G. Pectinase Enzyme Solution
(Immediately before use, prepare a solution containing 100 units/ml of Pectinase in cold deionized water.)

PROCEDURE:

Pipette (in milliliters) the following reagents into 50 ml Erlenmeyer flasks:

	<u>Test</u>	<u>Blank</u>
Reagent A (Polygalacturonic Acid)	4.90	5.00

Equilibrate to 25°C and then add the following reagent:

Reagent G (Pectinase)	0.10	-----
-----------------------	------	-------

Mix by swirling and incubate at 25°C for exactly 5.0 minutes. At the end of 5 minutes incubation, add:

Reagent B (I ₂ /KI)	5.0	5.0
Reagent C (Na ₂ CO ₃)	1.0	1.0

Mix by swirling and store in the dark for 20 minutes. Then add:

Reagent D (H ₂ SO ₄)	2.0	2.0
---	-----	-----

Mix by swirling and titrate the Test and Blank with Reagent E (Na₂S₂O₃) until it is light yellow. Then add 1 drop of Reagent F (Indicator) and continue to titrate with Reagent E (Na₂S₂O₃) until solutions are colorless. (Record the volume of Reagent E used.)

**Enzymatic Assay of PECTINASE
(EC 3.2.1.15)**

CALCULATIONS:

Units/ml enzyme =

$$\frac{(1)(100)(\text{ml of Reagent E for Blank} - \text{ml of Reagent E for Test})(\text{df})}{(5)(0.1)(2)}$$

1 = One μmole galacturonic acid is oxidized by 1 microequivalent of I_2

100 = microequivalents of S_2O_3 /ml of Reagent E

df = Dilution factor

5 = Time of reaction in minutes

0.1 = Volume (in milliliter) of enzyme used

2 = microequivalents of S_2O_3 oxidized per microequivalent of I_2 reduced

$$\text{Units/mg solid} = \frac{\text{units/ml enzyme}}{\text{mg solid/ml enzyme}}$$

$$\text{Units/mg protein} = \frac{\text{units/ml enzyme}}{\text{mg protein/ml enzyme}}$$

UNIT DEFINITION:

One unit will liberate 1.0 μmole of galacturonic acid from polygalacturonic acid per minute at pH 4.0 at 25°C.

FINAL ASSAY CONCENTRATIONS:

In a 5.00 ml reaction mix, the final concentrations are 0.49% (w/v) polygalacturonic acid and 10 units pectinase.

REFERENCE:

Kertesz, Z. I. (1955) *Methods in Enzymology*, Volume I, 162-164.

NOTES:

1. The standardization procedure is described in Reagent Chemicals, (1981) 6th ed., American Chemical Society Specifications, 551-552.
2. This assay is based on the cited reference.

**Enzymatic Assay of PECTINASE
(EC 3.2.1.15)**

NOTES: (continued)

3. Where Sigma Product or Stock numbers are specified, equivalent reagents may be substituted.

Sigma warrants that the above procedure information is currently utilized at Sigma and that all Sigma-Aldrich, Inc. products conform to the information in this and other Sigma-Aldrich, Inc. publications. Purchaser must determine the suitability of the information and product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.