

1.00482.1000

## Microscopy

### Periodic acid solution 0.5%

for the PAS reaction for the detection of aldehyde and mucosubstances in microscopy

#### For professional use only

**IVD** In Vitro Diagnostic Medical Device



#### Intended purpose

The PAS (periodic acid Schiff) reaction is one of the most frequently used chemical methods for histology.

This "Periodic acid solution 0.5% - for the PAS reaction for the detection of aldehyde and mucosubstances in microscopy" is used for human-medical cell diagnosis and serves the histological investigation of sample material of human origin. It is a ready-to-use solution that when used together with other in vitro diagnostic products from our portfolio makes target structures evaluable for diagnostic purposes (by sample preparation, fixing, embedding, staining, counterstaining, mounting) in human-histological specimen material, for example histological sections of e.g. the intestine or liver.

This ready-to-use Periodic acid solution 0.5% can be used to stain mucopolysaccharides in histological tissue specimens.

Unstained structures are relatively low in contrast and are extremely difficult to distinguish under the light microscope. The images created using the staining solutions help the authorized and qualified investigator to better define the form and structure in such cases. Further examinations may be necessary to reach a definitive diagnosis.

#### Principle

In the PAS reaction, the histological specimen material is first treated with periodic acid, resulting in the oxidation of the 1,2-glycols into aldehyde groups. The addition of Schiff's reagent (fuchsin-sulfuric acid) in the second step causes the aldehydes to react to form a brilliant red color. In the end result, the PAS reaction yields a specific color reaction with unsubstituted polysaccharides, neutral mucopolysaccharides, muco- and glycoproteins, and glyco- and phospholipids.

The PAS reaction can be further combined with the Alcian blue staining method to detect mucosubstances (glycosaminoglycans).

#### Sample material

Sections of formalin fixed, paraffin embedded tissue (3 - 4 µm thick paraffin sections) are used as starting material.

#### Reagents

Cat. No. 1.00482.1000

Periodic acid solution 0.5% 1 l  
for the PAS reaction for the detection of aldehyde and mucosubstances in microscopy

#### Also required (PAS staining):

Cat. No. 105174 Hematoxylin solution modified acc. to Gill III for microscopy 500 ml, 1 l, 2.5 l

Cat. No. 102572 Schiff's reagent Intense for the detection of aldehyde and mucosubstances in microscopy 1 l

or  
Cat. No. 109033 Schiff's reagent for microscopy 500 ml, 2.5 l

#### Optional (see "PAS staining - Procedure", footnote):

Cat. No. 105175 Hematoxylin solution modified acc. to Gill II for microscopy 500 ml, 2.5 l

#### Alternatively:

Cat. No. 101646 PAS staining kit for detection of aldehyde and mucosubstances 2x 500 ml

#### Sample preparation

The sampling must be performed by qualified personnel.

All samples must be treated using state-of-the-art technology.

All samples must be clearly labeled.

Suitable instruments must be used for taking samples and their preparation. Follow the manufacturer's instructions for application / use.

When using the corresponding auxiliary reagents, the corresponding instructions for use must be observed.

Deparaffinize and rehydrate sections in the conventional manner.

#### Reagent preparation

The Periodic acid solution 0.5% used for staining is ready-to-use, dilution of the solution is not necessary and merely produces a deterioration of the staining result and its stability.

#### Schiff's reagent

The stain can be carried out either with the ready-to-use solution provided in the PAS staining kit (Cat. No. 101646) or with a separately available Schiff's reagent Intense (Cat. No. 102572) or Schiff's reagent (Cat. No. 109033).

#### Procedure

When Schiff's reagent (Cat. No. 109033) is used, the PAS staining and the Alcian blue PAS staining procedure should be performed as per the instructions for use of the PAS staining kit (Cat. No. 101646).

#### PAS staining with Schiff's reagent Intense

##### Staining in the staining cell

Deparaffinize histological slides in the conventional manner and rehydrate in a descending alcohol series.

The slides should be allowed to drip off well after the individual staining steps, as a measure to avoid any unnecessary cross-contamination of solutions.

The stated times should be adhered to in order to guarantee an optimal staining result.

Slide with histological specimen	
Distilled water	rinse
Periodic acid solution 0.5%	5 min
Running tap water	3 min
Distilled water	rinse
Schiff's reagent Intense	10 min
Running tap water	3 min
Distilled water	rinse
Hematoxylin solution modified acc. to Gill III*	2 min
Running tap water	3 min
Ethanol 70 %	1 min
Ethanol 70 %	1 min
Ethanol 96 %	1 min
Ethanol 96 %	1 min
Ethanol 100 %	1 min
Ethanol 100 %	1 min
Xylene or Neo-Clear®	5 min
Xylene or Neo-Clear®	5 min
Mount the Neo-Clear®-wet slides with Neo-Mount® or the xylene-wet slides with e.g. Entellan® new and cover glass.	

\* To further enhance the brilliance and contrast of the PAS-positive structures, it is recommended to use hematoxylin solution modified according to Gill II (Cat. No. 105175).

After dehydration (ascending alcohol series) and clearing with xylene or Neo-Clear®, histological samples can be mounted with water-free mounting agents (e.g. Neo-Mount®, Entellan®, DPX new or Entellan® new) and a cover glass and can then be stored.

## Result

Nuclei	blue
Polysaccharides, glycogen, neutral mucopolysaccharides, muco- and glycoproteins, glyco- and phospholipids, basal membrane, collagen	purple

## Alcian blue PAS staining

When performing the Alcian blue PAS staining procedure, conventional Schiff's reagent (Cat. No. 109033) or PAS staining kit (Cat. No. 101646) should be used, since the intensive stain produced by Schiff's reagent Intense masks the typical blue stain of the Alcian blue dye and in certain tissues results in a mixed stain (violet).

## Technical notes

The microscope used should meet the requirements of a medical diagnostic laboratory.

When using histoprocessor systems or automatic staining systems, please follow the instructions for use supplied by the supplier of the system and software.

## Diagnostics

Diagnoses are to be made only by authorized and qualified personnel. Valid nomenclatures must be used.

This method can be supplementarily used in human diagnostics. Further tests must be selected and implemented according to recognized methods.

Suitable controls (e.g. ISOSLIDE® PAS, Cat. No. 1.00408.0001, ISOSLIDE® Alcian blue pH 2.5, Cat. No. 1.00425.0001) should be conducted with each application in order to avoid an incorrect result.

## Storage

Store the Periodic acid solution 0.5% - for the PAS reaction for the detection of aldehyde and mucosubstances in microscopy at +15 °C to +25 °C.

## Shelf-life

The Periodic acid solution 0.5% - for the PAS reaction for the detection of aldehyde and mucosubstances in microscopy can be used up to the stated expiry date.

After first opening of the bottle, the contents can be used up to the stated expiry date when stored at +15 °C to +25 °C.

The bottles must be kept tightly closed at all times.

## Capacity

2500 - 3000 applications / 500 ml

## Additional instructions

### For professional use only.

In order to avoid errors, the application must be carried out by qualified personnel only.

National guidelines for work safety and quality assurance must be followed. Microscopes equipped according to the standard must be used.

## Protection against infection

Effective measures must be taken to protect against infection in line with laboratory guidelines.

## Instructions for disposal

The package must be disposed of in accordance with the current disposal guidelines.

Used solutions and solutions that are past their shelf-life must be disposed of as special waste in accordance with local guidelines. Information on disposal can be obtained under the Quick Link "Hints for Disposal of Microscopy Products" at [www.microscopy-products.com](http://www.microscopy-products.com). Within the EU the currently applicable REGULATION (EC) No 1272/2008 on classification, labelling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No 1907/2006 applies.

## Auxiliary reagents

Cat. No. 100408	ISOSLIDE® PAS Control Slides with reference tissue for the detection of polysaccharides in histological tissue	25 tests
Cat. No. 100425	ISOSLIDE® Alcian blue pH 2.5 Control Slides with reference tissue for the detection of acid mucosubstances in histological tissue	25 tests
Cat. No. 100496	Formaldehyde solution 4%, buffered, pH 6.9 (approx. 10% Formalin solution) for histology	350 ml and 700 ml (in bottle with wide neck), 5 l, 10 l, 10 l Titripac®

Cat. No. 100579	DPX new non-aqueous mounting medium for microscopy	500 ml
Cat. No. 100974	Ethanol denatured with about 1 % methyl ethyl ketone for analysis EMSURE®	1 l, 2.5 l
Cat. No. 101646	PAS staining kit for detection of aldehyde and mucosubstances	2x 500 ml
Cat. No. 102572	Schiff's reagent Intense for the detection of aldehyde and mucosubstances in microscopy	1 l
Cat. No. 103699	Immersion oil Type N acc. to ISO 8036 for microscopy	100-ml dropping bottle
Cat. No. 104699	Immersion oil for microscopy	100-ml dropping bottle, 100 ml, 500 ml
Cat. No. 105174	Hematoxylin solution modified acc. to Gill III for microscopy	500 ml, 1 l, 2.5 l
Cat. No. 105175	Hematoxylin solution modified acc. to Gill II for microscopy	500 ml, 2.5 l
Cat. No. 107961	Entellan® new rapid mounting medium for microscopy	100 ml, 500 ml, 1 l
Cat. No. 108298	Xylene (isomeric mixture) for histology	4 l
Cat. No. 109016	Neo-Mount® anhydrous mounting medium for microscopy	100-ml dropping bottle, 500 ml
Cat. No. 109033	Schiff's reagent for microscopy	500 ml, 2.5 l
Cat. No. 109843	Neo-Clear® (xylene substitute) for microscopy	5 l
Cat. No. 111609	Histosec® pastilles solidification point 56-58°C embedding agent for histology	1 kg, 10 kg (4x 2.5 kg), 25 kg
Cat. No. 115161	Histosec® pastilles (without DMSO) solidification point 56-58°C embedding agent for histology	10 kg (4x 2.5 kg), 25 kg

## Hazard classification

Cat. No. 1.00482.1000

Please observe the hazard classification printed on the label and the information given in the safety data sheet.

The safety data sheet is available on the website and on request.

## Main components of the product

Cat. No. 1.00482.1000

H<sub>5</sub>IO<sub>6</sub> 5 g/l

## Other IVD products

Cat. No. 101647	Alcian blue solution pH 2.5 for microscopy	500 ml
Cat. No. 102439	Eosin Y-solution 0.5%, alcoholic for microscopy	500 ml, 2.5 l
Cat. No. 102472	ISOSLIDE® Warthin-Starry Control Slides with reference tissue for the detection of Helicobacter pylori and Spirochetes in histological tissue	25 tests
Cat. No. 102560	ISOSLIDE® AFB Control Slides with reference tissue for the detection of acid-fast bacteria in histological tissue	25 tests
Cat. No. 102561	ISOSLIDE® Congo Red Control Slides with reference tissue for the detection of amyloid structures in histological tissue	25 tests
Cat. No. 103999	Formaldehyde solution min. 37% free from acid stabilized with about 10% methanol and calcium carbonate for histology	1 l, 2.5 l, 25 l
Cat. No. 107164	Paraffin pastilles solidification point about 56-58°C for histology	10 kg (4x 2.5 kg)
Cat. No. 109204	Giemsa's azur eosin methylene blue solution for microscopy	100 ml, 500 ml, 1 l, 2.5 l
Cat. No. 117081	Eosin Y solution 1%, alcoholic for microscopy	1 l

## General remark

If during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorised representative and to your national authority.

## Literature

1. Romeis - Mikroskopische Technik, Editors: Maria Mulisch, Ulrich Welsch, 2015, Springer Spektrum, 19. Auflage
2. Welsch Sobotta - Lehrbuch Histologie, Editor: Ulrich Welsch, 2006, ELSEVIER Urban&Fischer, 2. Auflage
3. Theory and Practice of Histological Techniques, John D Bancroft, Marilyn Gamble, 2008, Churchill Livingstone ELSEVIER, 6th Edition
4. Laboratory Manual of Histochemistry, Linda L. Vacca, 1985, Raven Press
5. Basiswissen Histologie und Zytologie, Karl Heinz Stein, Hellmut Flenker, 1998, uZv, 2. Auflage
6. Histological & Histochemical Methods: Theory & Practice, J. A. Kiernan, 1990, Pergamon Press, 2nd Edition
7. Histological and Histochemical Methods, Theory and practice, J. A. Kiernan, 2015, Scion Publishing Ltd, 5th Edition
8. Conn's Biological Stains, R.W. Horobin, J.A. Kiernan, 2002, Biological Stain Commission Publication, 10th Edition



Consult instructions  
for use



Manufacturer



Catalog number



Batch code



Caution, consult  
accompanying documents



Use by  
YYYY-MM-DD



Temperature  
limitation

Status: 2021-Jan-25

Merck KGaA, 64271 Darmstadt, Germany,  
Tel. +49(0)6151 72-2440  
[www.microscopy-products.com](http://www.microscopy-products.com)

EMD Millipore Corporation, 400 Summit Drive  
Burlington MA 01803, USA, Tel. +1-978-715-4321  
Sigma-Aldrich Canada Co. or Millipore (Canada) Ltd.  
2149 Winston Park, Dr. Oakville, Ontario, L6H 6J8  
Phone: +1 800-565-1400

